TOXOREAGENT

TOXOREAGENT

REF RST7001

60 Tests

Art.-Nr. 487001

**Gebrauchsanweisung /
Instructions for Use /
Notice Technique**

**Nur zur Verwendung durch Fachpersonal /
For professional use only /
Exclusivement à un usage professionnel**

CE 2797

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TOXOREAGENT

Agglutination assay for the detection of *Toxoplasma gondii* antibodies in serum

Introduction

A latex reagent agglutination test for the semi-quantitative determination of anti-*Toxoplasma gondii* antibodies in human and animal serum samples.

FOR IN VITRO DIAGNOSTIC USE ONLY

Note that the CE-mark affixed to this product applies only to the use of the product with human samples, not with animal samples.

Important note for use of these kit instructions

Any assay relevant changes to the kit instructions for use (IFU) will lead to a change of the version number at the bottom on the last page of the IFU. All changes made will be identified on a separate sheet added to the IFU for a period of three months from the date of a change of the version. Please ensure that the latest version of the IFU is used for the assay procedure.

Contents

2. **AG/BEADS** Latex reagent or latex particles

1 x 12 mL, ready to use. Contains polystyrene latex particles coated with inactivated *Toxoplasma gondii* antigen and 0.1% proclin as preservative.

3. **CONTROL+** Positive Control

1 x 0.5 mL, ready to use; contains human serum with anti-*Toxoplasma gondii* antibodies and 0.1 % proclin as preservative. See titre concentration on QC document supplied in the kit.

4. **DIL** Diluent

1 x 50 mL. Contains 0.1% proclin as preservative.

5. Instructions for use

Further abbreviations

1. **RTU** ready to use

Stability and storage

- Store unopened at 2-8 °C in an upright position until the expiry date shown on the pack label. Once opened, the latex reagent and diluent should be stored at 2-8 °C with the lids tightly closed and may be used until the expiry date given on the label.
- The positive control should be stored at 2-8 °C with the lid tightly closed and may be used until the expiry date indicated on the label.
- **Do not freeze any of the reagents.**
- Serum samples may be stored at 2-8°C for up to 48 hours prior to testing. If longer storage is required,

store at -20°C for up to 1 year. Thawed samples must be mixed well prior to testing.

Warnings and precautions

- For *in vitro* diagnostic use only. Observe approved biohazard precautions and aseptic techniques when dealing with specimens of human origin.
- The positive control contains human serum. All sera, plasma and buffers containing human biological material were found to be negative for HIV, HCV and HBsAg. However, precautions such as the use of latex gloves must be taken.
- Sterilise all biohazard waste before disposal.
- Agglutination tests are sensitive to the effects of heat, direct sunlight and vibration. Keep away from such sources during test incubation periods.
- Some kit components contain bovine serum albumin. An infection risk is dominant after ingestion, inhalation or embrocation with reagents.
- The assay is based on latex particles. Latex may cause allergic reactions with skin. Wear gloves, coats, protection glasses etc. to reduce the risk.
- Latex particles must be resuspended thoroughly before use. In case of insufficiently resuspended Latex reagent, reactions may be hard to read due to latex button size.
- Use a separate disposable tip for each sample to prevent cross contamination.
- Replace caps on all reagents immediately after use.
- Do not use damaged or contaminated kit components.
- Kit components are matched and should not be interchanged between batches.
- For further safety information see Material Safety Data Sheet, available on request or from the Mast website.

Materials required but not provided

Standard microbiological supplies and equipment such as micropipettes, small dilution tubes and rigid U-shaped bottom microwell plates.

Sample Material

- a) Serum
- b) The serum or plasma samples can be stored at 2–8 °C for up to 3 days. They should be aliquoted and kept at -20°C for a longer storage.

The samples should not be frozen and thawed repeatedly.

After thawing samples should be briefly vortexed carefully before being used in the assay.

Lipemic, icteric, hemolytic or bacterially contaminated samples can cause false results.

Test Procedure

A. Sample preparation

1. Use fresh serum samples obtained by centrifugation of clotted blood. Serum samples may be stored at 2-8 °C for up to 48 hours or frozen for longer term storage.
2. Dilute all serum samples 1:8 by adding 50 µL of serum sample to a small tube and then add 350 µL of diluent. Mix well before use. Heat inactivation of specimens is not required prior to testing.

B. Positive Control

Note that the positive control is equivalent to a 1:8 dilution and requires no further dilution prior to use.

C. Procedure

1. Allow the TOXOREAGENT components to equilibrate to room temperature before use.
2. Gently shake the latex reagent to suspend the particles. Avoid frothing.
- Note: Make sure that all latex particles are thoroughly resuspended before use.**
3. Place 25 µL of diluent into wells 1-8 of row 1 of a rigid U-shaped bottom microwell plate using a micropipette. Repeat this operation for as many serum samples that are to be tested. All wells of column 9 are to be used as a control.
4. Dispense 25 µL of the reconstituted positive control to well 1 of the first row then add 25 µL of 1:8 diluted serum samples to well 1 of subsequent rows according to a pre-planned template.
5. Mix the contents of all wells in column 1 by pipetting up and down a few times and transfer 25 µL to the corresponding wells of column 2.
6. Continue mixing and transferring 25 µL amounts as in step 5 until you reach column 8. Finally mix and discard 25 µL from all wells of column 8. This will make a series of doubling dilutions from 1:16 in column 1 to 1:2048 in column 8. See diagram below.

a	1	2	3	4	5	6	7	8
b	1:16	1:32	1:64	1:128	1:256	1:512	1:1024	1:2048
c	25	25	25	25	25	25	25	25
d	25							
e	→	25 → dis- card						
f	Add 25µL latex reagent to all wells.							
	Mix contents of wells.							
	Leave tray at room temperature for at least 12 hours.							
	Read agglutination patterns.							

Key to Table

a – Well number.

b – Serum dilution.

c – µL diluent.

d – µL of positive control or 1:8 diluted serum samples.

e – Mix diluted serum + diluent and transfer.

f – Add latex reagent to all wells.

7. Add 25 µL of latex reagent to all wells. **Ensure that the latex reagent is fully resuspended before use.**
8. Mix the contents of the wells by gently tapping all four sides of the plate or by using a plate shaker set on high speed for at least 30 seconds.
9. Cover the plate to avoid evaporation and incubate at room temperature (15-25 °C) for at least 12 hours, **on a horizontal surface, in an undisturbed position that is free from vibration and away from heat and direct sunlight.**
10. Read the wells for agglutination patterns.

Reading and interpretation of results

Read the agglutination patterns according to the following criteria:

Result	Pattern	Interpretation
Intense agglutination with outer edge irregular and showing occasional tearing.		3+
Agglutination spread throughout the well.		2+
Agglutination is clearly present but does not cover the complete surface of the well.		1+
Weak agglutination pattern		0,5+
Negative reaction. A small and distinct circular sedimentation at the centre.		-

The titre is counted as the highest dilution of the test serum showing significant agglutination of 1+ or above.

The assay results can be considered as valid if the positive control matches the criteria given in the QC document.

For human serum

1. An antibody titre (serum dilution factor) of 1:32 or above indicates the presence of toxoplasma antibodies. It is suggested that for all positive samples the patient should be monitored every 2 weeks to assess for antibody levels. A rising titre level or maintenance of a high level indicates toxoplasmosis. For the diagnosis of ophthalmic toxoplasmosis even a low titre level is significant.

2. An antibody titre of 1:16 should be reported as positive but indicates the presence of a borderline level of toxoplasma antibodies. The probability of toxoplasmosis is low. Re-examination of the patient may be required by taking further blood samples to assess for a possible rise in titre value. If no rise in antibody titre is observed in subsequent samples the patient may be considered to be negative to toxoplasmosis.
3. Antibody titres of less than 1:16 may indicate the absence of clinically significant levels toxoplasma antibodies. In case of any doubt the test shall be repeated using another method.
4. If a negative result is obtained from a pregnant woman it is advisable to repeat the test during the period of pregnancy to detect a possible seroconversion.
In case of doubt it is recommended to request another serum sample and to consult the physician.

For pigs and cats

1. An antibody titre (serum dilution factor) of 1:64 or above indicates the presence of toxoplasma antibodies.
2. It is suggested that for all positive samples further samples should be taken and monitored every 2 weeks to assess for antibody levels.
3. An antibody titre of 1:32 indicates the presence of a borderline level of toxoplasma antibodies. The probability of toxoplasmosis is low. Re-examination of the animal by taking further blood samples may be required to assess for a possible rise in titre value. If no rise in antibody titre is observed in subsequent samples the animal may be considered to be negative to toxoplasmosis.
4. An antibody titre of 1:16 or less indicates the absence of clinically significant levels toxoplasma antibodies.

Limitations of use

- TOXOREAGENT detects IgG and IgM antibodies. The assay cannot differentiate between the two immuno globulin classes.
- TOXOREAGENT is useful for the serological diagnosis of toxoplasmosis. Human serum samples with titres of 1:16 are counted as weak positive or borderline and should be retested after approximately 2 weeks.
- When testing samples from immunocompromised, congenitally infected or renal dialysis patients serological test results should be considered along with clinical history and other aspects of patient management to be considered diagnostically significant. It is advisable to confirm the patient's anti-toxoplasma status with an alternative test.
- False positive reactions have been noted in patient samples positive for antinuclear antibodies, cytomegalovirus IgM antibodies, or both HBs and HBe antigens. Results from any serological test should always be considered along with the patient's clinical symptoms and history.

- Each diagnostic assay may give false-positive or false-negative results. In case of any doubt it is advisable to confirm the patient's anti-toxoplasma status with an alternative test.
- Serum samples only should be used in this test. Plasma samples should not be used as they may give rise to false positive results.
- Samples showing visibly lipaemic, haematic, icteric or other contaminations should be avoided.
- Samples up to the following concentrations were tested not showing any strange results, but nevertheless they should not be used:

Hemoglobin:	560 mg/dL
Bilirubin:	24 mg/dL
Triglycerides:	1000 mg/dL
- Occasionally a prozone effect may be observed with strongly positive sera. Please note that all samples should be fully titrated out. Screening is not recommended.

Quality control and reference material

It is recommended that quality control should be performed with the positive control provided to verify that the latex reagent is working correctly, and should be repeated at regular intervals. The antibody titre of the positive control is shown in the QC document supplied with the kit. If the positive control result falls outside of this range the test run should be considered invalid.

Also periodically check that latex reagents give correct results with known reference serum samples. Check for signs of deterioration. Do not use reagents if they are contaminated or cloudy.

Control materials used in the kit were tested against International Standards:

NIBSC 01/600:	1:2048
TOXM:	1:2048
NIBSC 09/B588:	1:1024

Performance characteristics

In a study of 146 human sera samples which were tested semi-quantitatively, the sensitivity of TOXOREAGENT compared to a Toxoplasma IgG EIA was found to be 93.8 % and the specificity was found to be 94.1 %. Positive results taken as $\geq 1:16$ for TOXOREAGENT.

All IgM positive samples in the tested panel gave positive reactions in TOXOREAGENT, too.

References

You can find the references at the end of the instruction for use.

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1. Rodrigues IMX, Castro AM, Gomes MBF, Amaral WN, Avelino MM; Congenital toxoplasmosis: evaluation of serological methods for the detection of anti-Toxoplasma gondii IgM and IgA antibodies, Mem. Inst. Oswaldo Cruz, Rio de Janeiro, 2009, Vol. 104 (3), 434-440
1. Kobayahshi A, Hirai N, Suzuki Y, Nishikawa H, Watanabe N. Evaluation of a commercial Toxoplasma latex agglutination test. Jap J Parasitol. 1977; 26: 175-180.
2. Balfour AH, Fleck DG, Hughes PA, Sharp D. Comparative study of three tests (dye test, indirect haemagglutination, latex agglutination test) for the detection of antibodies to *Toxoplasma gondii* in human sera. J Clin Pathol. 1982; 35: 228-232.
3. Hollimann, RE. Investigation of HIV-positive patients for Toxoplasmosis using the latex agglutination test. Serodiagn Immunother Inf Dis 1990; 4: 249-253.
4. <http://www.cdc.gov/parasites/toxoplasmosis/index.html>
5. <http://ecdc.europa.eu/en/healthtopics/toxoplasmosis/pages/index.aspx>
6. <http://www.rki.de/DE/Content/InfAZ/T/Toxoplasmose/Toxoplasmose.html>
7. <http://www.bfr.bund.de/en/toxoplasms-54399.html>
BgVV-Hefte Herausgegeben von M. Hartung, Bericht über die epidemiologische Situation der Zoonosen in Deutschland für 2001, 249 – 251

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Icons are used according to DIN EN ISO 15223-1

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Mast Diagnostica GmbH,
Mast Group Ltd.