

pecification shee

## Poly-HRP

## Streptavidin conjugated to Poly-Horseradish Peroxidase

Art.no. M2051

**Composition** Streptomyces avidinii.

**Reactivity** The HRP-conjugated streptavidin binds to biotin.

Conjugation Conjugate concentration: 1.0 mg/ml

Molecular E/P ratio : 20.0

Specificity The HRP-conjugated streptavidin was found specific for biotin in ELISA

and ELISPOT techniques.

Packing Each vial contains 0.1 ml Poly-HRP-conjugated streptavidin.

**Preservative** 5-Nitro-5-Bromo-1,3-Dioxan (600 ppm).

Storage and stability Poly-HRP-conjugated streptavidin should be stored below -18°C.

Repeated thawing and freezing may reduce the antibody activity and should be avoided. Do not use this product beyond the stated

expiration date.

Application Detection of biotinylated reagents in ELISA, ELISPOT and

immunostaining techniques.

**Directions for use** 1. Spin small volumes before use.

2. Store the stock solution below -18°C.

3. Before use: dilute the required quantity of the poly-HRP with dilution buffer and immediately restore the stock solution below -18°C.

4. A concentration of 1% casein in the dilution buffer is recommended to exclude an increased background, this casein solution is available as 5% stock solution (M2052).

5. For optimal results with ELISA techniques HPE buffer (M1940) is recommended, for ELISPOT the PeliSPOT (M2540) buffer is

recommended.

Reagents

Plesmanlaan 12

1066 CX Amsterdam

The Netherland