



Mast House, Derby Road, Bootle, Merseyside L20 1EA, United Kingdom.
Tel: +44 (0)151 933 7277 Fax: +44 (0)151 944 1332
www.mastgrp.com

MAST ID™ NITROCEFİN DISCS

**For the rapid detection of β -lactamase enzymes in
Neisseria gonorrhoeae, *Moraxella catarrhalis*, *Staphylococcus* spp.,
Haemophilus influenzae and anaerobic bacteria.**

The ability of certain bacteria to produce enzymes that inactivate β -lactam antibiotics i.e. penicillins and cephalosporins has long been recognised. Several clinical tests have been developed to detect β -lactamase enzymes. They include the iodometric method, the acidometric method and the use of a variety of chromogenic substrates. Nitrocefın is a chromogenic cephalosporin that has been found to be effective in detection of all known β -lactamase enzymes. β -lactamase enzymes hydrolyse the amide bond in the β -lactam ring of Nitrocefın resulting in a distinctive colour change from yellow to red.

Rapid β -lactamase tests can yield clinically relevant information earlier than an MIC test or disc diffusion test. MAST ID™ Nitrocefın discs are intended for use in the rapid testing of isolated colonies of *Neisseria gonorrhoeae*, *Moraxella catarrhalis*, *Staphylococcus* spp., *Haemophilus influenzae*¹ and anaerobic bacteria, notably *Bacteroides* spp.¹². This is because organisms within a taxonomic group or even a single strain may produce a diversity of β -lactamase enzymes with different substrate specificities. Nitrocefın is the only reliable test for detecting β -lactamase producing *Enterococcus* spp.

- | | |
|---------------------|-----------------------------------------------------------------------|
| Easy to use | • Smear colony onto surface of moistened disc |
| Convenient | • Tests can be performed using either clean Petri dish or glass slide |
| Cost effective | • No additional reagents required |
| Definable end point | • Colour change from yellow to red |