# MASTAZYME™ ANA PROFIL 8

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<b>Contents</b>	<b>Page</b>
1. Intended Use	13
2. Introduction	13
3. Principle of the Test	16
4. Kit Contents	17
5. Materials Required but not Provided	17
6. Warnings and Precautions	18
7. Storage and Stability	18
8. Specimen Collection and Handling	19
9. Assay Procedure	19
10. Results and Interpretation	21
11. Assay Performance	21
12. References	21

## 1. Intended Use

MASTAZYME™ ANA PROFIL 8 has been designed for the detection of specific antibodies against dsDNA, SS-A, SS-B, Sm, Sm/RNP, Scl-70, Jo-1 and Centromere in serum and plasma.

This assay is intended for *in vitro* diagnostic use only.

All laboratory test results should be interpreted in conjunction with other clinical data. The clinical judgement and further tests have to be taken into account additionally.

## 2. Introduction

In general ENA's ("extractable nuclear antigens") are considered a subgroup of ANA ("anti-nuclear antibodies"). Extractable nuclear antigens can be isolated by various biochemical methods from the nucleus of eukaryotic cells and most of them are well characterized. Their clinical relevance and their role in pathogenesis is widely accepted as well as the fact that for many diseases more than one autoantibody can be detected by diagnostic methods.

### Double stranded DNA (dsDNA) and single stranded DNA (ssDNA)

#### Antigen

Antibodies directed against double stranded DNA molecules can bind to single stranded DNA and RNA as well. The majority of ANA positive antibodies leading to homogeneous patterns in immunofluorescence assays are of anti-dsDNA specificity, some of them are also specific for DNA binding proteins only. The molecular structure of DNA allows binding of three groups of autoantibodies: antibodies directed to phosphoribose chains, antibodies to purine and/or pyrimidine bases and antibodies specific for respective conformational epitopes of the molecules. Although molecular components of DNA in prokaryotic and eukaryotic cells do not vary, human anti-dsDNA antibodies do not bind to DNA molecules of various species with comparable specificities and sensitivities.

The affinity and avidity of antibodies directed against sugar phosphates serving as backbone structures of DNA is highly dependent on negatively charged moieties of the phosphodiester chain. Therefore other autoantibodies capable of recognizing phospholipids (e.g. Cardiolipin) and chemical bondages of phosphodiester origin do bind to dsDNA and ssDNA antigens, too. Anti-EBV antibodies of the EBNA-1 type also show a strong cross-reactivity in all dsDNA and ssDNA ELISA. Purine and pyrimidine bases are forming the inner part of the native DNA molecule. Specific autoantibodies have access to bind to these purine and pyrimidine bases after melting native DNA into single-stranded molecules. Autoantibodies directed against adenosine, guanosine, cytosine and thymidine only bind to ssDNA antigens! GC rich epitopes in ssDNA are recognized by anti-ssDNA antibodies with sensitivity above average. To avoid any nonspecific cross-reactivities dsDNA molecules coated to the solid phase in the MASTAZYME™ Anti-dsDNA is free of ssDNA epitopes!

In contrary to anti-ssDNA binding activities the helical structure of dsDNA antigen as well as negatively charged moieties of the molecule are crucial for specific binding of anti-dsDNA antibodies.

#### Clinical relevance

Anti-dsDNA antibodies are highly predictive for systemic lupus erythematosus. Antibody titers correlate well with disease activity e.g. in cases of lupus nephritis. An increase in antibody titers indicates the development of a more severe state of the disease. The diagnosis of high avidity IgG antibodies is more predictive for the course of SLE than IgM type antibodies. Anti-dsDNA IgM antibodies are also found within healthy control groups. An intense monitoring of anti-dsDNA-positive patients is recommended, to detect an early rise of antibody titers and to start treatment at the earliest time possible. An increase of anti-dsDNA autoantibodies can be detected about 10 weeks before clinical symptoms occur.

Cases obviously not linked to an autoimmune pathogenesis express anti-ssDNA autoantibodies which correlate well in terms of coincidence with disease activity and clinical symptoms. For SLE diagnosis anti-ssDNA results are only of low specificity, in about 40 % of patients. Furthermore in cases of chronic liver diseases e.g. autoimmune hepatitis (AIH) and chronic biliary cirrhosis anti-ssDNA titers are detected above statistically normal prevalences.

In other liver diseases not directly associated with an autoimmune pathogenesis anti-ssDNA autoantibodies are present in serum. Anti-ssDNA is also present in metabolic diseases with a potential to develop vascular problems.

Testing for anti-dsDNA and anti-ssDNA is generally recommended

- when a specimen reacts homogenous in the nuclei of HEp-2 cells or in nuclei of liver tissue sections
- when clinical data are related to SLE symptoms
- when clinical data suggest a drug-induced form of SLE (the determination of anti-histone antibodies is recommended here as well).

## **SS-A**

### Antigen

SS-A ("soluble substance A") and the so-called Ro protein (Robert) are identical. SS-A consists of at least 4 short RNA molecules and two proteins of 60 kDa and 52 kDa. Although SS-A and calreticulin – a calcium binding protein – show some homology both proteins are clearly distinct from each other in their amino acid sequences. Regulatory proteins are usually conserved among animal species, so human autoantibodies bind easily to SS-A proteins isolated from non-human origin. Anti-SS-A antibodies bind to several domains of the native 60 kDa protein.

### Clinical relevance

A strong connection between anti-SS-A antibodies and lupus or lupus-like diseases is well accepted and confirmed in the literature (1). A high percentage of anti-SS-A is found within a group of SLE patients with skin involvement and among pregnant SLE patients with a newborn's risk of congenital heart block. The expression of both SS-A and SS-B indicates in most cases a mild progression of the disease. Sometimes SS-A can be the only autoantibodies found in so called „ANA negative“ and subacute cutaneous lupus patients. The autoantibody concentrations of anti-SS-A and anti-dsDNA correlate well with disease activity.

There is a high incidence of anti-SS-A and anti-SS-B for Sjögren's syndrome; antibodies reach high titers within this group of patients. Only a few patients suffering from scleroderma and mixed connective tissue disease express anti-SS-A and anti-SS-B.

Only IgG antibodies are of diagnostic and clinical interest to date, although specific antibodies of the IgM and IgA class can be detected in patient sera.

## **SS-B**

### Antigen

SS-B ("soluble substance B") and the so-called La protein (Lane) are identical. Despite almost identical names and quite similar clinical and diagnostic functions there is no homology between SS-A and SS-B regarding biochemical and molecular structures. SS-B is defined as a 47 kDa phosphoprotein associated with various RNAs like rRNA and tRNA. Even RNA molecules usually specific for SS-A proteins can bind to native SS-B. SS-B is involved in RNA transport from the nucleus to the cytoplasm. The protein is further part of the regulatory cascade of RNA polymerase III reactions. An ATPase activity and the ability to melt DNA/RNA hybrids support the view of the physiological functions of SS-B.

### Clinical relevance

In many patient sera both antibodies to SS-A and SS-B are often co-expressed. Both antibodies are generally found exclusively in Sjögren's syndrome with an underlying SLE. Up to 15 % of SLE patients express anti-SS-B antibodies. SLE patients positive for SS-A and SS-B usually develop a less severe form of systemic lupus erythematosus.

Isolated SS-B antibody responses are limited to an early form of the Sjögren's syndrome. SS-A and SS-B proteins are both expressed on heart muscle fibres increasing the risk of congenital heart block in the newborns after binding the respective autoantibodies. Both IgG antibodies pass through the placenta during pregnancy.

## **Sm and Sm/RNP**

### Antigen

The Sm antigen (Smith) consists of at least 9 different polypeptides which, together with nRNP proteins, are part of the nuclear splicing complex. Various snRNAs bind to a core complex of the Sm protein built of clearly defined proteins B/B', D and E. Antibodies directed against Sm in patient sera mainly bind to proteins B and D. The Sm antigen coated to the solid phase in the MASTAZYME™ Anti-Sm assay contains B and D proteins. Furthermore the Sm antigen does not contain any RNP moieties!

The Sm/RNP antigen, however, consists of various snRNPs (small nuclear ribonucleoprotein particles) in addition to the Sm part. These snRNPs and the proteins B/B', D, E, F or G plus a mRNA molecule are part of the spliceosome. Anti-RNP antibodies bind to a specific domain of the Sm/RNP protein complex on the solid phase which is different from the binding site of anti-Sm.

### Clinical relevance

Anti-Sm antibodies are considered as a marker for HLA DR4 associated SLE and in a few cases they are indicative of MCTD. Nevertheless anti-Sm is highly predictive for SLE. Although "only" 35 % of all SLE patients express anti-Sm positive results are highly specific for SLE. A correlation between antibody titers and disease activity has been observed. The clinical relevance of this data, however, is still under discussion.

Anti-snRNP antibodies are found in SLE but they are also found in a higher frequency in MCTD. Patients suffering from rheumatic arthritis, scleroderma, Sjögren's syndrome and progressive systemic sclerosis can be found positive for anti-snRNP autoantibodies.

## **Scl-70**

### Antigen

Anti-Scl-70 antibodies bind to eukaryotic topoisomerase I (topo I) an enzyme located in the nucleus where it is involved in relaxation of supercoiled genomic DNA during replication. In SDS polyacrylamide gels topo I migrates at a size of 100 kDa but smaller fragments even with proteolytic activity are also present, for example the Scl-70 antigen which is a 70 kDa fragment of the native 100 kDa topo I. Like other proteins involved in cellular regulation topo I is highly conserved among species leading to cross-reactions. An intact antigen conformation is most important for the detection of anti-Scl-70 autoantibodies.

### Clinical relevance

Anti-Scl-70 is detectable in 20–30 % of scleroderma patients where the antibody is considered as a disease defining marker. In Scl-70 positive scleroderma cases other diseases like SLE or Sjögren's syndrome may also be present and must be excluded. The prevalence of anti-Scl-70 is low among MCTD patients but even there the presence of the antibody is highly predictive for mixed connective tissue diseases. Antibody titers usually do not correlate with disease activity in scleroderma.

## **Jo-1**

### Antigen

Histidyl-tRNA synthetase is the corresponding antigen for anti-Jo-1 antibodies. The enzyme is found in prokaryotic and eukaryotic organisms but only antigens derived from higher eukaryotes or the human antigen itself are detected by anti-Jo-1 autoantibodies. Immunofluorescence assays clearly show the distribution of the antigen in the cell, concentrated in the perinuclear cytoplasm according to its physiological function where it is involved in protein biosynthesis.

### Clinical relevance

There is a strong correlation between anti-Jo-1 autoantibodies and inflammatory muscle diseases. Jo-1 is considered as a diagnostic marker for myositis. About 54 % of Jo-1 positives are found within patients suffering from primary myositis, 40 % were connected with dermatomyositis and 6% of Jo-1 positive patients showed a myositis together with other connective tissue diseases. Multi-systemic disorders defined by Jo-1 can be separated from other symptoms connected with an anti-synthetase syndrome.

### **Centromere**

#### Antigen

CENP-A (16 kDa), CENP-B (80 kDa) and CENP-C (140 kDa) are the main target antigens for anti-centromere antibodies with CENP-B. Further antigens like CENP-D and CENP-F are described but CENP-B is of most importance for diagnostic. Antibodies are directed against the inner and outer part of the kinetochor but do not bind to dsDNA molecules.

#### Clinical Relevance

Most frequently antibodies directed against kinetochor protein are associated with CREST syndrome and may be indicating a more benign disease progress. They are less common in within the other collagenosises but approximately 10–30 % of Raynaud's patients are found positive. Even healthy blood donors may show anti-centromere antibodies especially in the group of female donors (0.08 %).

Antibodies directed to CENP proteins are sometimes found together with AMA-M2 antibodies but never together with Scl-70 antibodies.

### **3. Principle of the Test**

The principle of the test reaction can be described in four stages.

#### **3.1 Serum incubation**

Specific antibodies bind to the antigens on the solid phase to form a stable immune complex. After a 30 minutes incubation at room temperature the wells are washed with prediluted washing buffer to remove all non-reactive serum components.

#### **3.2 Conjugate incubation**

The anti-human-IgG horseradish peroxidase conjugate is added to all wells. The conjugate binds to IgG antibodies on the solid phase antigen to form a stable sandwich. After a 30 minutes incubation at room temperature the excess conjugate is removed by washing all wells with washing buffer.

#### **3.3 Substrate reaction and stopping**

The TMB substrate is dispensed into each well and the peroxidase enzyme/substrate reaction forms a stable blue chromogen. The reaction and subsequently the colour development is stopped after 15 minutes incubation at room temperature by adding 1 N H<sub>2</sub>SO<sub>4</sub> to the wells. The change in pH also causes the chromogen to change colour from blue to yellow.

#### **3.4 Reading and interpretation**

The intensity of the colour is read in a microtiter plate reader at 450 nm (recommended reference wavelength for bichromatic measurement: 600–690 nm). The intensity of the colour (OD) is directly proportional to the concentration of the specific antibody in the patient sample.

#### 4. Kit Contents

The kit contains sufficient reagents for 12 x 8 = 96 determinations. The strips and solutions have to be stored at 2–8 °C. The expiry date is mentioned on the labels.

12 x	Microtiter strips	single strips each with 8 break-apart wells coated with a purified antigen per well of: dsDNA, SS-A, SS-B, Sm, Sm/RNP, Scl-70, Jo-1 and Centromere
1 x	Frame holder	
8 x 1.5 mL	Cut-off Calibrator	human pool-plasma containing antibodies against (antigen specific) the above listed antigens, diluted in buffer, ready to use
8 x 1.5 mL	Positive Control	human plasma containing antibodies against the (antigen specific) above listed antigens, diluted in buffer, ready to use
1 x 3,5 mL	Negative Control	human plasma, diluted in buffer, ready to use
2 x 60 mL	Sample diluent	PBS solution contains preservative, ready to use
1 x 12 mL	Enzyme conjugate solution	HRP-labelled goat anti-human-IgG, ready to use
1 x 12 mL	TMB substrate	3,3',5,5' Tetramethylbenzidine, ready to use
1 x 12 mL	Stopping solution	1 N sulfuric acid, ready to use
2 x 50 mL	Washing buffer 10 x concentrated	PBS/Tween buffer solution 10 x concentrated to be diluted 1:10 prior to use; the concentrate should be warmed up to 37 °C for 15 minutes to avoid any crystals

#### 5. Materials Required but not Provided

- 5 µL-, 100 µL- and 500 µL micro- and multichannel pipettes
- Microtiter plate reader with a 450 nm filter (reference filter 600–690 nm)
- Microtiter Plate Washer (in case of manual washing: wash bottle)
- Reagent tubes for the serum dilution
- Measuring cylinder
- Distilled water or water of higher quality

## 6. Warning and Precautions

- For *in vitro* diagnostic use only! Do not ingest or swallow! Laboratory safety precautions should be followed. Do not eat, drink or smoke in the laboratory.
- All sera and plasma or buffers based upon have been tested to HBsAg, HIV and HCV respectively with generally accepted methods and were found negative. Nevertheless, precautions like the use of latex gloves have to be taken.
- Serum and reagent spills have to be wiped off with a disinfecting solution (e.g. sodium hypochlorite 5 %) and have to be disposed of properly.
- All reagents have to be brought to room temperature (18–24 °C) before performing the test.
- Before pipetting all reagents should be mixed thoroughly by gentle tilting or swinging. Vigorous shaking with formation of foam should be avoided.
- It is important to pipet with constant intervals so that all the wells of the microtiter plate have the same conditions.
- When removing reagents out of the bottles care has to be taken that the stoppers are not contaminated. Further a possible mix-up has to be avoided. The content of the bottles is usually sensitive to oxidation so that they should be opened only for a short time.
- In order to avoid a carry-over or a cross-contamination separate disposable pipet tips have to be used.
- No reagents from different kit lots should be used and they should not be mixed with one another.
- All reagents have to be used within shelf life.
- In accordance with a Good Laboratory Practice (GLP) or following ISO 9001 all laboratory devices employed should be regularly checked regarding the accuracy and precision. This refers e.g. to microliter pipets and washing or reading (ELISA Reader) instrumentation.
- The contact of certain reagents especially the stopping solution and the substrate with skin, eye and mucosa has to be avoided because possible irritations and acid burns could arise and there exists a danger of intoxication.

## 7. Storage and Stability

Store all reagents at 2–8 °C.

The expiry date of each reagent is printed on the individual labels. Do not use any reagents after the expiry date has been exceeded.

The diluted washing buffer is stable for up to 4 weeks when stored at 2–8 °C.

The opened kit should be used within three months.

## 8. Specimen Collection and Handling

Both serum or plasma (EDTA, heparin) can be used for the determination. Serum is separated from the blood which is aseptically drawn by venipuncture after clotting and centrifugation. The serum or plasma samples can be stored at 2–8 °C for up to 3 days. They should be kept at -20 °C for a longer storage. The samples should not be frozen and thawed repeatedly. Lipemic, hemolytic or bacterially contaminated samples can cause false positive or false negative results.

Patient sera must be prediluted 1:101 in sample diluent (e.g. 10 µL serum + 1000 µL sample diluent) prior to testing.

## 9. Assay Procedure

### 9.1. Preparation of Reagents

**Allow all kit components and specimens to reach room temperature (RT, 18–24 °C) prior to use and mix well.**

**Washing buffer:** Dissolve any crystals which may be in the bottle by warming to 37 °C and then mix well.

Dilute the concentrated washing buffer 1:10 with distilled water (e.g. 50 mL buffer concentrate + 450 mL distilled water). Mix thoroughly.

- Strictly follow the instructions for reliable test performance. Any changes or modifications are within the responsibility of the user.
- All reagents and samples must be brought to room temperature before use but should not be left at this temperature for longer than necessary.
- A standard curve should be established with each assay.
- Put the unused microtiter strips back in the plastic bag and store them dry at 2–8 °C.

## 9.2. Assay Steps

Prepare a sufficient amount of microtiter wells for calibrators, controls and samples.

**Note: Other incubation conditions might be possible. In case of modifications of the recommended test procedure (e.g. incubation temperature 37 °C instead of RT) the user has to validate assay performance.**

		1	2	3	4	5	6	7	8	9	10	11	12
<b>A</b>	dsDNA	CO	Pos	S	S								
<b>B</b>	SS-A	CO	Pos	S	S								
<b>C</b>	SS-B	CO	Pos	S	S								
<b>D</b>	Sm	CO	Pos	S	S								
<b>E</b>	Sm/RNP	CO	Pos	S	S								
<b>F</b>	Scl-70	CO	Pos	S	S								
<b>G</b>	Jo-1	CO	Pos	S	S								
<b>H</b>	CENP	CO	Pos	S	S								

CO: Cut-off calibrator, Pos: positive control, S = Sample

1. Pipette 100 µL each of the diluted (1:101) samples and the ready to use cut-off calibrator and controls into the appropriate wells e.g.
2. Incubate the plate at room temperature for 30 minutes.
3. Discard the contents of the microwells and wash 3 times with 300 µL of diluted washing buffer. Afterwards remove residues of the washing solution by gentle tapping of the microtiter plate on a paper towel.
4. Pipette 100 µL of enzyme conjugate solution into each well.
5. Incubate the strips for 30 minutes at room temperature.
6. Discard the contents of the microwells and wash 3 times with 300 µL of diluted washing buffer. Afterwards remove residues of the washing solution by gentle tapping of the microtiter plate on a paper towel.
7. Dispense 100 µL of TMB substrate into each well.
8. Incubate for 15 minutes in the dark (e.g. drawer) at room temperature.
9. Add 100 µL of stopping solution to each well.
10. After thorough mixing and wiping the bottom of the plate, read the optical density at 450 nm and calculate the results. Blank against air. A bichromatic measurement using a reference wavelength of 600–690 nm is recommended.

**The developed colour is stable for at least 30 minutes. Read optical densities during this time.**

## 10. Results and Interpretation

### Qualitative Calculation

The calculated OD values for patient sera as mentioned above are compared with the value for the cut-off calibrator. If the value of the sample is higher, then it should be read as positive.

Normal range:

Negative: OD-sample < OD cut-off calibrator

Positive: OD-sample > OD cut-off calibrator

The positive calibrator must show an absorption value at least double the value received by the cut-off calibrator.

### Interpretation of results

No single laboratory result should be used for a diagnosis only. The results should be interpreted together with other laboratory parameters and other clinical findings.

As screening assay the ELISA detects antibodies directed against the antigens dsDNA, SS-A (Ro), SS-B (La), Sm, Sm/RNP, Scl-70, Jo-1 or Centromere. The assay is neither suitable for antibody typing nor for the identification of antibody specificity.

## 11. Assay Performance

The assay characteristics of the MASTAZYME ANA PROFIL 8 have been established and assessed according the European IVD directive. Detailed validation data can be provided on special request.

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Laboratoriumspräparate GmbH
**MAST GROUP Ltd.**