# MASTAZYME™ AMA-M2

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	<b>Contents</b>	<b>Page</b>
1.	Intended Use	11
2.	Introduction	11
3.	Principle of the Test	12
4.	Kit Contents	12
5.	Materials Required but not Provided	13
6.	Warning and Precautions	13
7.	Storage and Stability	14
8.	Specimen Collection and Handling	14
9.	Assay Procedure	14
10.	Results and Interpretation	15
11.	Assay Performance	16
12.	References	16

## 1. Intended Use

MASTAZYME™ AMA-M2 has been designed for the detection of AMA-M2.

This assay is intended for *in vitro* diagnostic use only.

All laboratory test results should be interpreted in conjunction with other clinical data. The clinical judgement and further tests have to be taken into account additionally.

## 2. Introduction

### Mitochondrial Proteins

#### Antigen

Anti-mitochondrial antibodies (AMA) are a heterogeneous group of autoantibodies directed against various proteins that are located in the outer and inner membrane of mitochondria. Expression of their distinct target antigens takes place in the tissues. Non-organ specific subtypes of antibodies have been described besides those antibodies which exhibit relatively high organ specificity.

Using fluorescence techniques with tissue sections from different organs, i.e. liver, stomach, kidney, heart and pancreas, the mitochondrial autoantibodies of different specificity can be discriminated. For some of the AMA subtypes their specific target protein remains unclear. HEp 2 cell monolayers for indirect immunofluorescence AMA-M2 are characterized as a fine-speckled cytoplasmic, perinuclear condensed fluorescence pattern. For differential diagnosis of the primary biliary cirrhosis (PBC) determination of AMA-M2 by ELISA is recommended because of its high sensitivity and specificity.

Subtypes	Autoantigen	Mitochondrial localisation	Clinical Relevance
M1	Cardiolipin	inner membrane	Lues II
M2	Proteins of the $\alpha$ -keto acid dehydrogenase complex	inner membrane	Primary biliary cirrhosis (PBC)
M3		outer membrane	Drug-induced LE
M4	Associated with sulfite oxidase	outer membrane	PBC
M5		outer membrane	SLE and undifferentiated collagenoses, autoimmune haemolytic anaemia
M6		outer membrane	Drug induced hepatitis
M7	Sarcosin dehydrogenase	inner membrane	Cardiomyopathy, myocarditis
M8		outer membrane	PBC
M9	Associated with glycogen phosphorylase	outer membrane (probably also cytoplasmic)	PBC

#### Clinical relevance

Specific anti-mitochondrial antibodies have been described for the primary biliary cirrhosis (PBC) as subtypes M2, M4, M8 and M9. Other AMA subtypes are related to other diseases like collagenoses (AMA-M5) and drug induced LE and hepatitis (AMA-M3 and AMA-M6). The heterogeneously reacting specific anti-mitochondrial antibodies of the M2 subtype are directed against three related proteins of the  $\alpha$ -keto acid dehydrogenase complex which is located at the inside of the mitochondrial membrane. The major epitope recognized is located on the E2 subunit and the protein X of the pyruvate dehydrogenase complex (PDC). Additionally AMA-M2 recognize the (E1 und E3) subunits of the same complex

and the E2 subunit of several other multi enzyme complexes such as the 2-oxoglutarate dehydrogenase complex (OGDC) and the branched chain 2-oxo acid dehydrogenase complex (BCOADC).

In patients with other autoimmune diseases determination of AMA antibodies allows an early screening for the occurrence of subtype M2 and M9 antibodies which may be related with the development and / or association of PBC. Profiling the AMA subtypes allows an immunological and prognostic classification of the primary biliary cirrhosis. Early cases of symptomatic PBC often exhibit only AMA-M2 subtype (sometimes in combination with AMA-M9), whereas progressive cases and mixed syndromes with chronic acute hepatitis (CAH) are related with the occurrence of AMA-M2, -M4 and -M8 subtypes.

### 3. Principle of the Test

The principle of the test reaction can be described in four stages.

#### 3.1 Serum incubation

Specific antibodies bind to the antigens on the solid phase to form a stable immune complex. After a 30 minutes incubation at room temperature the wells are washed with pre-diluted wash buffer to remove all non-reactive serum components.

#### 3.2 Conjugate incubation

The anti-human-IgG horseradish peroxidase conjugate is added to all wells. The conjugate binds to IgG antibodies on the solid phase antigen to form a stable sandwich. After a 30 minutes incubation at room temperature the excess conjugate is removed by washing all wells with washing buffer.

#### 3.3 Substrate reaction and stopping

The TMB substrate is dispensed into each well and the peroxidase enzyme/substrate reaction forms a stable blue chromogen. The reaction and subsequently the colour development is stopped after 15 minutes incubation at room temperature by adding 1 N H<sub>2</sub>SO<sub>4</sub> to the wells. The change in pH also causes the chromogen to change colour from blue to yellow.

#### 3.4 Reading and interpretation

The intensity of the colour is read in a microtiter plate reader at 450 nm (recommended reference wavelength for bichromatic measurement: 600–690 nm). The intensity of the colour (OD) is directly proportional to the concentration of the specific antibody in the patient sample.

### 4. Kit Contents

The kits contain sufficient reagents for 12 x 8 = 96 determinations. The strips and solutions have to be stored at 2–8 °C. The expiry date is mentioned on the labels.

12	Microtiter strips	single strips each with 8 break-apart wells coated with pyruvate dehydrogenase
1 x	Frame holder	
6 x 1,5 mL	Calibrators 1–6	human plasma containing antibodies to mitochondrial antigen (AMA-M2), diluted in buffer, ready to use

	AMA-M2 concentration [U/mL]
Calibrator 1	1
Calibrator 2 (Cut-off)	15
Calibrator 3	25
Calibrator 4	50
Calibrator 5	100
Calibrator 6	200

1 x 150 µL	Positive Control	human plasma containing antibodies to mitochondrial antigen (AMA-M2), neat, to be diluted 1:101
1 x 150 µL	Negative Control	human plasma free of AMA-M2, neat, to be diluted 1:101
2 x 60 mL	Sample diluent	PBS solution contains preservative, ready to use
1 x 12 mL	Enzyme conjugate solution	HRP-labelled goat anti-human-IgG, ready to use
1 x 12 mL	TMB substrate	3,3',5,5' Tetramethylbenzidine, ready to use
1 x 12 mL	Stopping solution	1 N sulfuric acid, ready to use
2 x 50 mL	Washing buffer 10 x concentrated	PBS/Tween buffer solution 10x concentrated to be diluted 1:10 prior to use; the concentrate should be warmed up to 37 °C for a short time to avoid any crystals

## 5. Materials required but not provided

- 5 µL-, 100 µL- and 500 µL micro- and multichannel pipettes
- Microtiter plate reader with a 450 nm filter (reference filter 600–690 nm)
- Microtiter Plate Washer (in case of manual washing: wash bottle)
- Reagent tubes for the serum dilution
- Measuring cylinder
- Distilled water or water of higher quality

## 6. Warning and Precautions

- For *in vitro* diagnostic use only! Do not ingest or swallow! Laboratory safety precautions should be followed. Do not eat, drink or smoke in the laboratory.
- All sera and plasma or buffers based upon have been tested to HBsAg, HIV and HCV respectively with generally accepted methods and were found negative. Nevertheless, precautions like the use of latex gloves have to be taken.
- Serum and reagent spills have to be wiped off with a disinfecting solution (e.g. sodium hypochlorite 5 %) and have to be disposed of properly.
- All reagents have to be brought to room temperature (18–24 °C) before performing the test.
- Before pipetting all reagents should be mixed thoroughly by gentle tilting or swinging. Vigorous shaking with formation of foam should be avoided.
- It is important to pipet with constant intervals so that all the wells of the microtiter plate have the same conditions.
- When removing reagents out of the bottles care has to be taken that the stoppers are not contaminated. Further a possible mix-up has to be avoided. The content of the bottles is usually sensitive to oxidation so that they should be opened only for a short time.
- In order to avoid a carry-over or a cross-contamination separate disposable pipet tips have to be used.
- No reagents from different kit lots should be used and they should not be mixed with one another.
- All reagents have to be used within shelf life.

- In accordance with a Good Laboratory Practice (GLP) or following ISO 9001 all laboratory devices employed should be regularly checked regarding the accuracy and precision. This refers e.g. to microliter pipets and washing or reading (ELISA Reader) instrumentation.
- The contact of certain reagents especially the stopping solution and the substrate with skin, eye and mucosa has to be avoided because possible irritations and acid burns could arise and there exists a danger of intoxication.

## 7. Storage and Stability

Store all reagents at 2–8 °C.

The expiry date of each reagent is printed on the individual labels. Do not use any reagents after the expiry date has been exceeded.

The diluted washing buffer is stable for up to 4 weeks when stored at 2–8 °C.

The opened kit should be used within three months.

## 8. Specimen Collection and Handling

Both serum and plasma (EDTA, heparin) can be used for the determination. Serum is separated from the blood which is aseptically drawn by venipuncture after clotting and centrifugation. The serum or plasma samples can be stored at 2–8 °C for up to 3 days. They should be kept at -20 °C for a longer storage. The samples should not be frozen and thawed repeatedly. Lipemic, haemolytic or bacterially contaminated samples can cause false positive or false negative results.

Patient sera must be pre-diluted 1:101 in sample diluent (e.g. 5 µL serum + 500 µL sample diluent) prior to testing.

## 9. Assay Procedure

Dilute controls and patient samples 1:101 with sample diluent before testing. For example pipette 10 µL of sample to 1,000 µL of sample diluent in a suitable tube.

**Controls are not pre-diluted and must therefore be diluted 1:101 before use!**

### 9.1. Preparation of Reagents

**Allow all kit components and specimens to reach room temperature (RT, 18–24 °C) prior to use and mix well.**

**Washing buffer:** Dissolve any crystals which may be in the bottle by warming to 37 °C and then mix well.

Dilute the concentrated washing buffer 1:10 with distilled water (e.g. 50 mL buffer concentrate + 450 mL distilled water). Mix thoroughly.

**Controls:** Dilute controls (1:101) as for samples with sample diluent

- Strictly follow the instructions for reliable test performance. Any changes or modifications are within the responsibility of the user.
- All reagents and samples must be brought to room temperature before use, but should not be left at this temperature for longer than necessary.
- A standard curve should be established with each assay.
- Put the unused microtiter strips back in the plastic bag and store them dry at 2–8 °C.

## 9.2. Assay Steps

Prepare a sufficient amount of microtiter wells for calibrators, controls and samples.

**Note: Other incubation conditions might be possible. In case of modifications of the recommended test procedure (e.g. incubation temperature 37 °C instead of RT) the user has to validate assay performance.**

1. Pipette 100 µL each of the diluted (1:101) samples, ready to use calibrators and diluted (1:101) controls into the appropriate wells.
2. Incubate the plate at room temperature for 30 minutes.
3. Discard the contents of the microwells and wash 3 times with 300 µL of diluted washing buffer. Afterwards remove residues of the washing solution by gentle tapping of the microtiter plate on a paper towel.
4. Pipette 100 µL of enzyme conjugate solution into each well.
5. Incubate the strips for 30 minutes at room temperature.
6. Discard the contents of the microwells and wash 3 times with 300 µL of diluted washing buffer. Afterwards remove residues of the washing solution by gentle tapping of the microtiter plate on a paper towel.
7. Dispense 100 µL of TMB substrate into each well.
8. Incubate for 15 minutes in the dark (e.g. drawer) at room temperature.
9. Add 100 µL of stopping solution to each well.
10. After thorough mixing and wiping the bottom of the plate, read the optical density at 450 nm and calculate the results. Blank against air. A bichromatic measurement using a reference wavelength of 600–690 nm is recommended.

**The developed colour is stable for at least 30 minutes. Read optical densities during this time.**

## 10. Results and Interpretation

### Qualitative

Calibrator 2 or calibrator 3 can be used for qualitative interpretation of the assays, depending on specificity requirements to include or eliminate borderline results. A quantitative interpretation is not possible by using calibrator 2 or 3 only.

### Quantitative

A standard curve is plotted by entering the mean absorbance of each calibrator on the Y-axis and the corresponding concentration on the X-axis using either log/log or linear/log graph paper. The concentration of the controls and specimen can then be read directly from the graph.

The calculation of the results can also be performed using a computer and a suitable software program. In this case the results are calculated using a 4 parameter fit algorithm.

### Suggested normal values

In a normal range study using serum samples from healthy blood donors and disease-state sera the following normal and elevated ranges have been established.

Negative:	< 15	U/mL
Borderline	15–25	U/mL
Positive:	> 25	U/mL

## **Recommendation**

Each laboratory should establish its individual normal ranges based on results obtained from the local population.

## **Interpretation of results**

MASTAZYME™ AMA-M2 detects antibodies directed against the E2 subunit of the mitochondrial dehydrogenase complex. Autoantibodies reacting with 70 kd, 50 kd and 42 kd antigens of pyruvate-dehydrogenase located in the inner mitochondrial membrane give positive results. There is a strong correlation between the prevalence of AMA-M2 and primary biliary cirrhosis; within 5 % of PBC patients no AMA can be detected with any of the widely used diagnostic tools.

No cross-reactivity was found between MASTAZYME™ AMA-M2 and sera derived from second stage syphilis patients but should be considered for differential diagnostic reasons.

Anti mitochondrial antibody concentrations do not correlate with stage of disease or prognosis. Environmental factors, genetical predispositions or chronic graft versus host disease may lead to elevated antibody levels.

Results are affected by several patient factors. Clinical diagnosis and disease prognosis should not be based on a single result. Other lab data and clinical observations must be taken into consideration before a conclusive diagnosis can be made.

## **11. Assay Performance**

The assay characteristics of the MASTAZYME™ AMA-M2 have been established and assessed according to the European IVD directive. Detailed validation data can be provided on special request.

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Laboratoriumspräparate GmbH
**MAST GROUP Ltd.**