*ENZYWELL*

**TOXOPLASMA IgG AVIDITY**

**REF 91098 (24 tests)**

Manufactured by: DIESSE Diagnostica Senese  
Via delle Rose 10  
53035 Monteriggioni (Siena) - Italy

**INDICE / INDEX**

1. UTILIZZAZIONE / INTENDED USE / UTILISATION
2. INTRODUZIONE / SUMMARY AND EXPLANATION OF TEST /INTERET CLINIQUE
3. PRINCIPIO DEL METODO / PRINCIPLE OF THE TEST / PRINCIPE
4. COMPOSIZIONE DEL KIT E PREPARAZIONE DEI REAGENTI / KIT CONTENTS AND REAGENT PREPARATION / CONTENU DU COFFRET ET PREPARATOIN DES REACTIFS
5. MODALITA' DI CONSERVAZIONE E STABILITA' DEI REAGENTI / STORAGE AND STABILITY OF REAGENTS / MODALITE DE CONSERVATION ET STABILITE DES REACTIFS
6. PRECAUZIONI / PRECAUTIONS / PRECAUTIONS D'UTILISATION
7. TIPO DI CAMPIONE E CONSERVAZIONE / TYPE AND STORAGE OF SAMPLE / TYPE D'ECHANTLLON
8. PROCEDIMENTO / TEST PROCEDURE / MODE OPERATOIRE
9. SCHEMA DEL SAGGIO / SCHEME OF TEST PROCEDURE / SCHEMA D'ANALYSE
10. VALIDAZIONE DEL TEST / VALIDATION OF THE TEST / VALIDATION DU TEST
11. INTERPRETAZIONE DEL TEST / INTERPRETATION OF RESULTS / INTERPRETATION DES RESULTATS
12. LIMITAZIONI DELLA PROCEDURA /LIMITATIONS OF THE PROCEDURE / LIMITES DE LA METHODE
13. SENSIBILITA' E SPECIFICITA' DIAGNOSTICA / DIAGNOSTIC SENSITIVITY AND SPECIFICITY / SENSIBLITE ET SPECIFICITE DIAGNOSTIQUE
14. PRECISIONE / PRECISION
15. GUIDA DEI PROBLEMI DI UTILIZZO / "TROUBLE SHOOTING"
16. BIBLIOGRAFIA / REFERENCES / BIBLIOGRAPHIE

## ENZYWELL TOXOPLASMA IgG AVIDITY KIT

**REF** 91098

(Italiano)

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(English)

**1. INTENDED USE:**

**IMMUNOENZYMATIC METHOD FOR THE DETERMINATION OF THE AVIDITY OF ANTI TOXOPLASMA GONDII IgG ANTIBODIES IN HUMAN SERUM, AS AN AID IN THE DIAGNOSIS OF TOXOPLASMA INFECTION.**

**2. INTRODUCTION:**

*Toxoplasma gondii* is a ubiquitous protozoa which causes infection in all species of mammals. In man, serological studies have shown that a high percentage of the adult population becomes infected by this parasite. The disease is particularly dangerous during pregnancy, being transmitted from mother to fetus and causing serious malformations.

For this reason it is important to determine the antibody titer before the beginning of pregnancy if possible; those subjects who do not present immunity should be followed up with frequent analyses to check for serum conversion.

Sometimes it is not sufficient to test for IgM antibodies, as the presence of this class may be due to the persistence of IgM or asymptomatic reinfection by *Toxoplasma*, without risk for the fetus.

For this reason it is useful to assay the avidity of the IgG antibodies, in as far as it has been demonstrated that avidity is low in the acute phase of infection and increases over a period of time. The presence of low avidity antibodies is therefore an indication of a recent or current infection. The avidity of IgG antibodies in the serum can be assayed with this ELISA kit.

**3. PRINCIPLE OF THE TEST:**

The test is based on the ELISA technique (Enzyme linked Immunosorbent Assay).

The purified and inactivated *Toxoplasma gondii* is bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with dilute human serum. The sample is placed in two wells (or in four wells if the test is performed in duplicate) and left to incubate.

After washing of the solid phase, the Avidity Buffer is placed in one well and the Wash Buffer in the other. The avidity buffer will cause dissociation of the antigen/antibody bond, proportional to the avidity of the antibodies. After washings to eliminate the proteins that have not reacted, incubation is performed with the conjugate composed of anti-human IgG monoclonal antibodies labelled with peroxidase HRP.

The unbound conjugate is eliminated and the peroxidase substrate specific for the enzyme is added.

The colour which develops is proportional to the concentration of specific antibodies present in the serum sample still bound to the antigen. The ratio between the OD of the well containing the Avidity buffer and that containing the Wash Buffer is expressed as a percentage and gives the degree of avidity of the serum sample for the *Toxoplasma* antigen.

**4. KIT CONTENTS AND REAGENT PREPARATION:**

- Reagents are sufficient for 48 determinations.

- **Bring reagents to room temperature before use.**

**MT PLATE** MICROPLATE COATED WITH *TOXOPLASMA GONDII* (COLOURED IN YELLOW). 3x16 wells.

Use: open the package at the opposite end from the code (T, followed by the lot number) which is useful for identification purposes, remove the support and wells to be used from the foil package, and place the unused wells in the polythene bag with the silica gel, expel the air and seal by pressing the closure.

**CONJ** CONJUGATE. 1x12 mL.

Contents: monoclonal antibodies labelled with Peroxidase, in phosphate buffer solution containing phenol 0.05% and Bronidox 0.02%. Ready for use without further dilution.

**CONTROL HIGH** HIGH AVIDITY CONTROL (2.5 mL)

Contents: Human serum containing IgG anti-*Toxoplasma* with high avidity, in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

**CONTROL LOW** LOW AVIDITY CONTROL (2.5 mL)

**Contents:** Human serum containing IgG anti-Toxoplasma with low avidity in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

**AVIDITY BUF** AVIDITY BUFFER (1 x 10 mL)

**Contents:** Urea in saline buffered solution containing phenol 0.05% and Bronidox 0.02%, ready for use. If crystals are present, they should be dissolved at 37°C before use.

**WASH BUF 10x** WASH BUFFER 10X. 1 x 100 mL. **INTERCHANGEABLE BETWEEN LOTS**

**Contents:** Phosphate buffered saline, concentrated 10 times; contains Brij 0.5% w/v.

**Preparation:** dilute 1:10 with distilled water in order to obtain the washing buffer ready for use. If crystals are present, they should be dissolved at 37°C before dilution.

**SAMP DIL** DILUENT 2 (PF93611). 1 x 100 mL. **INTERCHANGEABLE BETWEEN LOTS**

To be used to dilute samples.

**Contents:** Proteic solution in phosphate buffer with sodium azide 0,09% containing methyl orange as dye.

**SUBS TMB** SUBSTRATE. 12 mL. Ready for use. **INTERCHANGEABLE BETWEEN LOTS**

**Contents:** Tetramethylbenzidine and hydrogen peroxide stabilised in citrate buffer, pH 3.8.

**H<sub>2</sub>SO<sub>4</sub> 0.3 M** STOP SOLUTION. 1x16 mL. **INTERCHANGEABLE BETWEEN LOTS**

H<sub>2</sub>SO<sub>4</sub> 0.3 mol/L, in solution ready for use.

ADHESIVE FILMS (2).

POLYTHENE BAG (1).

#### **MATERIALS REQUIRED BUT NOT PROVIDED.**

- Incubator at 37°C
- Microplate reader, wavelength 450 nm or 450/620 nm and 405 nm, with OD linearity up to 2,000 (at least).
- Microplate washer (preferable) able to dispense volumes in the range 225-375 µL
- Distilled or deionized water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 10, 100, 1000 µl solution
- Disposable gloves
- Timer
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials
- Absorbent tissue.

#### **5. STORAGE AND STABILITY OF REAGENTS**

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

#### **Reagents have a limited stability after opening and/or preparation**

REAGENT	CONDITIONS
Microplate	5 weeks at 2/8°C, polythene bag
Control sera	5 weeks at 2/8°C
Conjugate	5 weeks at 2/8°C
Substrate	up to the expiry date at 2/8°C, 1 week at 15-30°C; store in the dark
Sample Diluent	up to the expiry date at 2/8°C
Wash Buffer	2 weeks at 2/8°C, 5 days at 15/30°C.
Avidity Buffer	up to the expiry date at 2/8°C
Stop Solution	up to the expiry date at 2/8°C

#### **6. PRECAUTIONS**

*FOR IN VITRO DIAGNOSTIC USE ONLY.*

*This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.*

### Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
2. The following reagents contain low concentrations of harmful or irritant substances:
  - a) The Wash Buffer contains detergents
  - b) The conjugate contains phenol
  - c) The substrate is acid
  - d) Calibrators contain Sodium Azide 0.09% as preservative which may react with copper and lead in plumbing to form potentially explosive metal azides. On disposal, flush with large volumes of water to prevent azide buildup. If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
3. Non-disposable apparatus should be sterilized after use. The preferred method is to autoclave for 1 h at 121°C; disposables should be autoclaved or incinerated.
4. Sulphuric acid required for the Stop Solution and hydrochloric acid used for washing glassware are corrosive and should be handled with appropriate care. If they come into contact with the skin or eyes, wash thoroughly with water.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

### Analytical precautions

1. Allow all reagents and samples to come to room temperature (18-30°C) before use. Immediately after use return reagents to the recommended storage temperature. **It is important to work at the correct temperature. Check that the thermostat does not go below 35°C or over 39°C.** Open the envelope containing the strips after at least ½ hr at room temperature.
2. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
3. Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.
4. Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or high quality deionized water.
5. Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps.
6. Do not allow wells to become dry during the assay procedure.
7. Care must be taken not to cross-contaminate reagents. It is important that pipettes are dedicated for exclusive use with the various reagents.
8. Care should be taken to avoid touching or splashing the rim of the well with conjugate. Do not "blow-out" from microplates.
9. Enzyme immunoassays can occasionally exhibit an "edge effect" which must be minimized by increasing the humidity during incubation steps. Plates must be covered with their covers and incubated at 37°C either in a water bath with a rack or float to support the plates if necessary, or in an incubator. Alternatively, plates can be incubated in an approved analyzer. See the appropriate operating manual for further details. CO<sub>2</sub> incubators must not be used.
10. Ensure that the bottom of the plate is clean and dry, and that no bubbles are present on the surface of the liquid before reading the plate.
11. Use of highly hemolyzed samples, incompletely clotted sera, or samples with microbial contamination may give rise to erroneous results.
12. The use of the kit with automate4d equipment has to be validated by the users.
13. For each instrument used, read the manufacturer's instructions manual carefully to obtain additional information on the following points:
  - installation and particular requisites
  - operating principles, instructions, precautions and risks
  - manufacturer's specifications and instrument performance
- servicing and maintenance.

### 7. TYPE AND STORAGE OF SAMPLE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times. Defrosted samples must be carefully shaken before performing the test. Heat inactivation can lead to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric or contaminated samples should be avoided. If a new sample cannot be obtained, such samples should be clarified by filtration (0.45 µm) or centrifugation (3000 rpm x 10').

**The test cannot be performed on human plasma.**

### **8. TEST PROCEDURE:**

- Prepare the required number of strips.
- Prepare the washing buffer by diluting the Wash Buffer 10x (100 mL + 900 mL H<sub>2</sub>O).

Leave one well empty for the blank. Dilute samples 1:101 distributing 10 µL of serum into 1 mL of diluent. Dispense 100 µL of each diluted sample in at least two wells, and 100 µL of High Avidity Control and 100 µL of Low Avidity Control in duplicate. Wells are covered with protective film and incubated for 60 minutes at 37°C. After washing four times for 30 seconds, dispense Wash Buffer in one well of the serum samples and of the controls and Avidity Buffer in the other, and incubate for 10' at 37°C. Aspirate the buffer and perform 4 washings for 30 seconds each (300 µL). Add 100 µL of conjugate to each well and incubate again for 45 minutes at 37°C, covering the wells with the protective film. The plate is washed again 4 times, as described above. Finally, the substrate is distributed, 100 µL/well, including the blank.

After 15 minutes at room temperature the enzymatic reaction is stopped with 100 µL of Stop Solution.

The adsorbance (O.D.) is read at 450 nm or at 450/620 nm within 30 min.

<b>9. TEST PROCEDURE FOR TOXOPLASMA IgG AVIDITY KIT</b>
---

STEP 1                    Place 100 µL of diluted sample/controls in duplicate in the wells of the strips.

-

Incubate for 60 min. at 37°C

-

Wash 4 times (30 seconds each; 300 µL)

STEP 2                    For each sample and control, place 100 µL of Wash Buffer in one well and 100 µL of Avidity buffer in the other; incubate for 10' at 37°C

-

Wash four times (30 seconds each; 300 µL)

STEP 3                    Add 100 µL of conjugate to each well

-

Incubate for 45 min. at 37°C

-

Wash 4 times

STEP 4                    Add 100 µL of Substrate to each well including the blank

-

Incubate for 15 min. at R.T.

STEP 5                    Add 100 µL of Stop Solution

-

Read absorbance at 450 nm within 30 min

### **10. TEST VALIDATION**

The OD Blank must be ≤ 0.15.

The high and low avidity controls in the wells incubated with the Wash Buffer must give an OD value > 0.500. The percentage of avidity of the low control must be less than 35%, and the high avidity control must be higher than 50% compared to that found with the Wash Buffer (calculation performed as reported in the following paragraph).

### **11. INTERPRETATION OF THE RESULTS**

The percentage of avidity of the samples is expressed and calculated using the ratio between the OD found in the wells containing Avidity buffer and those with Wash Buffer, subtracting the value of the test blank.

$$\frac{\text{O.D. with Avidity Buffer} - \text{O.D. blank}}{\text{O.D. with Wash Buffer} - \text{O.D. blank}} \times 100 = \% \text{ of avidity}$$

A ratio over 35% indicates the presence of high avidity IgG antibodies; when it is lower than 30%, this indicates the presence of low avidity IgG. If the percentage is between 30 and 35%, there is a medium degree of avidity (borderline).

In this case the interpretation of the test is doubtful.

The calculation can only be performed on samples with give an OD  $\geq 0.40$  in the well washed with Wash Buffer.

#### Warning

If the OD of the reference well is higher than 2000, read the test at 405 nm. If the OD remains  $> 2000$  repeat the test diluting 1/400.

If the OD of the reference well is  $\leq 0.4$  repeat the test diluting the sample 1/ 25.

### **12. LIMITATIONS**

A result of high avidity does not exclude the possibility of a recent infection. On the other hand, given the high specificity of the test, a positive result strongly indicates an infection within the previous 3 months.

### **13. DIAGNOSTIC SENSITIVITY AND SPECIFICITY**

In a clinical trial, 81 samples from patients who had undergone an infection during the previous three months were tested. The method being tested was compared with 4 other commercial kits. The Diesse kit showed the highest specificity (84%). The other 4 methods gave values of 50%, 33%, 58%, 69% respectively. The sensitivity of the Diesse kit was 64%.

Diesse showed 100% specificity in the control population (49 samples) while the other methods gave many false positive results.

### **14. PRECISION**

Table 1 Precision "in run" (n=6)

Sample	% average avidity	S.D..	CV%
Low avidity	19	0.7	4
High avidity	78	2	2

Table 2 Precision between runs

Sample	% average avidity	S.D.	CV%
S1	79	4	5
S2	52	1	2
S3	14	1	4

### **15. TROUBLE SHOOTING GUIDE**

PROBLEM	POSSIBLE SOURCE	TEST OR ACTION
Invalid run (all negative)	One or more reagents not added or added in wrong sequence	Recheck procedure Check for unused solutions. Repeat test.
	Unreactive plate	Check the code on the package containing the plate (see package insert point 4 for correct code).
		Check for moisture in unused plate. (Silica gel dessiccant must be pale yellow ).Repeat test
Invalid run (all positive)	Contamination of substrate	Take new aliquot of substrate.
	Inadequate washing	Ensure that wash apparatus works well
Poor precision	Incomplete washing of wells	Ensure that wash apparatus works well
	Inadequate aspiration of wells	Ensure that wash apparatus works well
	Pipetting error	Check pipette function
	Reagent addition too slow	Avoid drying of the plate after washing step. Add reagents immediately
	Presence of bubbles	Avoid air bubbles during pipetting.
	Optical pathway not clean	Check instrument light source and detector for dirt. Wipe bottom of plate with soft tissue.
Inadequate Color development	Incorrect incubation times or temperature	Check for temperature control and time monitoring
		Adhere to recommended instruction for use.
	Inadequate volume of substrate added to the plate	Check pipette function.

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