ENZYWELL

RUBELLA IgM

REF 91031 (96 tests)

Manufactured by: DIESSE Diagnostica Senese
Via delle Rose 10
53035 Monteriggioni (Siena) - Italy

INDICE / INDEX

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(English)

1. INTENDED USE

IMMUNOENZYMATIC CAPTURE METHOD FOR THE QUALITATIVE DETERMINATION OF IgM-CLASS ANTIBODIES TO RUBELLA VIRUS IN HUMAN SERUM, AS AN AID IN THE DIAGNOSIS OF RUBELLA VIRUS INFECTION.

2. SUMMARY AND EXPLANATION OF THE TEST

Rubella virus causes a very common infectious disease in man, accompanied by a typical diffused rash. The course of the illness is generally short and without complications. However, it is harmful for the fetus if contracted during pregnancy, because the virus is teratogenic and may cause very serious malformations especially during the first month of pregnancy. It is therefore essential to determine the patient's immune status before the start of pregnancy, if possible, in order to perform vaccination or, if this is impossible, to monitor the course of the pregnancy to observe whether seroconversion occurs. Several methods are employed for the assay of specific immunoglobulins (1). ELISA is one of these, and has the advantage of being completely automated. The serological diagnosis of Rubella may be performed with the assay of specific IgM because these antibodies, together with the clinical symptoms, are of considerable diagnostic relevance.

3. PRINCIPLE OF THE TEST

The test for the assay of anti Rubella Virus IgM is based on the principle of the capture of these immunoglobulins and the subsequent identification of those which are specific, due to their ability to bind an antigen conjugated to peroxidase. The capture is performed using monoclonal antibodies bound to the solid phase (microtiter wells) (2,3,4). The antigen is made of purified and inactivated Rubella Virus labelled with peroxidase bound to specific anti Rubella Virus monoclonal antibodies.

4. KIT CONTENTS AND REAGENT PREPARATION

- Reagents are sufficient for 96 determinations.
- **Bring to room temperature before use.**

MT PLATE *MICROPLATE. 12x8 wells coated with anti-human IgM monoclonal antibodies.*

Use: open the package at the opposite end from the code (MM followed by the lot number) which is useful for identification purposes, remove the support and strips to be used from the foil package, and place the unused strips in the polythene bag with the silica gel, expel the air and seal by pressing the closure.

CONTROL + *POSITIVE CONTROL (1 x 1.6 mL)*

Contents: Diluted human serum in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

Colour: the colour is proportional to the relative antibody titer.

CONTROL CUT OFF *CUT OFF CONTROL (1 x 2.5 mL)*

Contents: Diluted human serum in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

Colour: the colour is proportional to the relative antibody titer.

Ag *ANTIGEN. Freeze-dried powder x 6 vials.*

Contents: Purified Rubella virus in Phosphate buffer with mouse ascitic fluid and lactose.

Preparation: reconstitute with the conjugate volume shown on the label, mixing by inversion.

CONJ CONJUGATE.

Contents: monoclonal antibodies labelled with Peroxidase, in phosphate buffer solution containing phenol 0.05% and Bronidox 0,02%.

Preparation: ready for use.

The immunocomplex should be prepared about 45 min. before use.

CONTROL - IgM NEGATIVE CONTROL (PF93900) (1 x 1.6 mL) INTERCHANGEABLE BETWEEN LOTS

Contents: Diluted human serum in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

WASH BUF 10x WASH BUFFER 10X (PF93603). 1 x 100 mL. INTERCHANGEABLE BETWEEN LOTS

Contents: Phosphate buffered saline, concentrated 10 times; contains Brij 0.5% .

Preparation: dilute the required volume 1:10 with distilled water in order to obtain the washing buffer ready for use. If crystals are present, they should be dissolved at 37°C before dilution.

SAMP DIL DILUENT 2 (PF93611). 1 x 100 mL. INTERCHANGEABLE BETWEEN LOTS

To be used to dilute samples.

Contents: Proteic solution in phosphate buffer with sodium azide 0,09% containing methyl orange as dye.

SUBS TMB SUBSTRATE (PF93619). 12 mL. Ready for use. INTERCHANGEABLE BETWEEN LOTS

Contents: Tetramethylbenzidine 0.26 mg/mL and hydrogen peroxide 0.01% stabilised in citrate buffer 0.05 mol/L (pH 3.8).

H₂SO₄ 0.3 M STOP SOLUTION (PF93602). 1x16 mL. INTERCHANGEABLE BETWEEN LOTS

H₂SO₄ 0.3 mol/L, in solution ready for use.

ADHESIVE FILMS (2).

POLYTHENE BAG (1).

MATERIALS REQUIRED BUT NOT PROVIDED.

- Incubator at 37°C
- Microplate reader (wave length 450 or at 450/620 nm, with linearity up to OD >= 2000)
- Microplate washer (preferable) able to dispense volumes in the range 225-375 µl
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 10, 100, 1000 µl solution
- Disposable gloves
- Timer
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials
- Absorbent tissue.

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

Reagents have a limited stability after opening and/or preparation

REAGENT	CONDITIONS
Microplate	5 weeks at 2/8°C, polythene bag
Control sera	5 weeks at 2/8°C
Conjugate	5 weeks at 2/8°C
Reconstituted antigen	must be used the same day, it cannot be frozen
Substrate	up the expiry date at 2/8°C, 1 week at 15-30°C in the dark
Sample Diluent	up the expiry date at 2/8°C
Wash Buffer	2 weeks at 2/8°C, 5 days at 15/30°C.
Stop Solution	up the expiry date at 2/8°C

6. PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY.

Caution:

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
2. The following reagents contain low concentrations of harmful or irritant substances:
 - a) The Wash Buffer contains detergents
 - b) The conjugate contains phenol
 - c) The substrate is acid
 - d) The controls contain 0.09% Sodium Azide which can react with lead and copper in plumbing forming highly explosive deposits of metal azides; dilute with large amounts of water to eliminate.

If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
3. Non-disposable apparatus should be sterilized after use. The preferred method is to autoclave for 1 h at 121°C; disposables should be autoclaved or incinerated.
4. Sulphuric acid required for the Stop Solution and hydrochloric acid used for washing glassware are corrosive and should be handled with appropriate care. If they come into contact with the skin or eyes, wash thoroughly with water.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

1. Allow all reagents and samples to come to room temperature (18-30°C) before use. Immediately after use return reagents to the recommended storage temperature. **It is important to work at the correct temperature. Check that the thermostat does not go below 35°C or over 39°C.** Open the envelope containing the strips after at least ½ hr at room temperature.
2. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
3. Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.
4. Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or high quality deionized water.
5. Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps.
6. Do not allow wells to become dry during the assay procedure.
7. Care must be taken not to cross-contaminate reagents. It is important that pipettes are dedicated for exclusive use with the various reagents.
8. Care should be taken to avoid touching or splashing the rim of the well with conjugate. Do not "blow-out" from microplates.
9. Enzyme immunoassays can occasionally exhibit an "edge effect" which must be minimized by increasing the humidity during incubation steps. Plates must be covered with their covers and incubated at 37°C either in a water bath with a rack or float to support the plates if necessary, or in an incubator. Alternatively, plates can be incubated in an approved analyzer. See the appropriate operating manual for further details. CO₂ incubators must not be used.
10. Ensure that the bottom of the plate is clean and dry, and that no bubbles are present on the surface of the liquid before reading the plate.
11. Use of highly hemolyzed samples, incompletely clotted sera, or samples with microbial contamination may give rise to erroneous results.
12. The use of the kit with automated equipment has to be validated by the user.
13. For each instrument used, read the manufacturer's instructions manual carefully to obtain additional information on the following points:
 - installation and particular requisites
 - operating principles, instructions, precautions and risks

- manufacturer's specifications and instrument performance
- servicing and maintenance.

7. TYPE AND STORAGE OF SAMPLE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times. Defrosted samples must be carefully mixed before the test. Heat inactivation can lead to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric or contaminated samples should be avoided. If a new sample cannot be obtained, such samples should be clarified by filtration (0.45 µm) or centrifugation (3000 rpm x 10').

The test cannot be performed on human plasma.

8. TEST PROCEDURE

- Prepare the required number of strips.
- Prepare the washing buffer by diluting the Wash Buffer 10x (100 mL + 900 mL H₂O).
- Prepare the immunocomplex by adding the Conjugate to the antigen (volume shown on the label).

Dilute samples 1:101 distributing 10 µL of serum into 1 mL of diluent. Dispense 100 µL of each diluted sample per well (duplicate testing is recommended). Place UNDILUTED controls in a strip (100 µL in each well). The minimum requisite is 1 negative control, 2 cut-off and 1 positive control. Leave one well for the blank, performed using 100 µL of the substrate mixture.

Wells are covered with protective film and incubated for 45 minutes at 37°C. After washing four times for 30 seconds (300 µl) , add 100 µL of the immunocomplex (Antigen-anti Rubella Virus monoclonal antibodies labelled with peroxidase) to each well and incubate again for 45 minutes at 37°C, covering the wells with the protective film. The plate is washed again 4 times, as described above. Finally, the substrate is distributed, 100 µL/well.

After 15 minutes at room temperature the enzymatic reaction is stopped with 100 µL of Stop Solution.

The adsorbance (O.D.) is read at 450 nm or 450/620 nm within 30 min.

9. TEST PROCEDURE FOR RUBELLA IgM
--

STEP 1 Place 100 µL of diluted sample//controls in the wells of the strips.

-

Incubate for 45 min. at 37°C

-

Wash 4 times (300 µl)

STEP 2 Add 100 µL of immunocomplex to each well

-

Incubate for 45 min. at 37°C

-

Wash 4 times (300 µl)

STEP 3 Add 100 µL of Substrate to each well

-

Incubate for 15 min. at R.T.

STEP 4 Add 100 µL of Stop Solution

-

Read absorbance at 450 nm within 30 min

10. TEST VALIDATION

Subtract the value of the blank (<= 0.150) from all the other readings. The O.D. values of the control Cut-off serum tested in triplicate must be within 25% of the mean value. Disregard any abnormal value and recalculate the mean. The Positive control must have an O.D. at least 1.5 times that of the Cut-off serum. The ratio between Negative Control and Cut-off must be <= 0.6. The OD of the cut-off must be >= 0.2 at 450 nm and >= 0.16 at 450/620 nm.

11. INTERPRETATION OF THE RESULTS

Qualitative results

If the adsorbance of the sample is higher than that of the Cut-off, the sample is positive for the presence of specific IgM.

Calculate the ratio between the O.D. value of the sample and that of the Cut-off (INDEX).

The sample is considered:

Positive: if the ratio is > 1.2 .

Doubtful: $\pm 20\%$ of the Cut-off.

Negative: if the ratio is < 0.8 .

If the result is doubtful, repeat the test. If it remains doubtful, collect a new serum sample.

12. LIMITATIONS OF THE PROCEDURE

Specimens collected earlier than 5 days after onset of symptoms may not contain sufficient IgM antibody to be positive in the test.

Since maternal IgM antibody does not cross the placenta, virus-specific IgM antibody in the new-born can be assumed to be a result of congenital infection. Because of the possibility of contamination of cord blood with maternal IgM, it is prudent to confirm positive viral IgM antibody results on cord blood samples by testing a follow-up specimen from the infant.

Rubella IgM antibody responses recently have been reported in a small proportion of patients with infectious mononucleosis; some of these showed narrow specificity for certain strains of Rubella virus. Polyclonal stimulation of B memory cells by EBV was considered in some cases to be a possible mechanism for heterotypic IgM antibody responses.

The results must always be interpreted together with the clinical picture and data from other diagnostic investigations.

Once reconstituted, the antigen must be used the same day.

13. ANALYTICAL SPECIFICITY

65 serum samples containing potential interfering substances were test:

- Antinuclear antibodies (ANA) (n=11)
- Rheumatoid factor (n=8)
- Heterophile antibodies (n=5)
- Cytomegalovirus igM (23)
- Bilirubin (n=8)
- Triglycerides (n=10).

No interference was noted in any case.

14. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

In a clinical trial performed in a hospital laboratory, 294 samples were analysed, 27 of which came from the serum collection. The samples were analysed with the Diesse kit and with another commercial immunoenzymatic method: the comparison showed that Diesse found 1 false negative and 1 false positive result.

Results are summarized in the following table:

	REFERENCE METHOD	
	+	-
DIESSE +	36	1
DIESSE -	1	256

The Enzywell kit has a diagnostic sensitivity of 97.3% and a diagnostic specificity of 99.6%.

15. PRECISION

“In run” Precision:

Sample	RVM 1 (Negative < Cut Off)	RVM 2 (Positive > Cut Off)	RVM 3 (Positive)	Cut Off	Positive Control
<i>n</i> (replicates)	24	24	24	12	12
<i>O.D.</i>	0.281	0.941	1,596	0.501	1,799
<i>CV%</i>	3	4	2	8	1

“Between run” Precision:

<i>Sample</i>	<i>INDEX</i>	
	<i>Average</i>	<i>CV%</i>
<i>Pos.Con.</i>	3,8	6
<i>RVM1</i>	0,4	8
<i>RVM2</i>	1,5	10
<i>RVM3</i>	2,6	8

Precision between lots:

<i>Sample</i>	<i>INDEX</i>			<i>Media</i>	<i>CV%</i>
	<i>Lot n. 136</i>	<i>Lot n. 137</i>	<i>Lot n. 138</i>		
<i>Positive Control</i>	5.25	4.46	5.36	5.25	9.8
<i>RVM1</i>	0.65	0.54	0.61	0.6	9.3
<i>RVM2</i>	2.57	2.17	2.50	2.41	8.9
<i>RVM3</i>	3.95	4.46	3.95	4.12	7.14

16. TROUBLE SHOOTING GUIDE

PROBLEM	POSSIBLE SOURCE	TEST OR ACTION
Invalid run (all negative)	One or more reagents not added or added in wrong sequence	Recheck procedure Check for unused solutions. Repeat test.
	Unreactive plate	Check the code on the package containing the plate (see package insert paragraph 4 for correct code).
Invalid run (all positive)		Check for moisture in unused plate. (Silica gel dessiccant must be pale yellow).Repeat test
	Contamination of substrate	Take new aliquot of substrate.
Poor precision	Inadequate washing	Ensure that wash apparatus works well
	Incomplete washing of wells	Ensure that wash apparatus works well
	Inadequate aspiration of wells	Ensure that wash apparatus works well
	Pipetting error	Check pipette function
	Reagent addition too slow	Avoid drying of the plate after washing step. Add reagents immediately
	Presence of bubbles	Avoid air bubbles during pipetting.
Inadequate Color development	Optical pathway not clean	Check instrument light source and detector for dirt. Wipe bottom of plate with soft tissue.
	Incorrect incubation times or temperature	Check for temperature control and time monitoring
		Adhere to recommended instruction for use.
	Inadequate volume of substrate added to the plate	Check pipette function.

17. REFERENCES

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. S. Chantler and C.J. Evans: Selection and performance of monoclonal and polyclonal antibodies in an IgM antibody capture enzyme immunoassay for rubella. J. Immunol. Meth. 87: 109 (1986).
3. R.S. Tedder and J.L. Yao: The production of monoclonal antibodies to rubella haemagglutinin and their use in antibody-capture assays for rubella-specific IgM. J. Hyg. Camb. 88: 335 (1982).
4. B. Forghani, C. Myoraku, N. Schmidt: Use of monoclonal antibodies to human immunoglobulin M in "capture" assays for measles and Rubella immunoglobulin M. J. Clin. Microbiol. 18, 652 (1983).

DIESSE Diagnostica Senese
Via delle Rose 10
53035 Monteriggioni (Siena) Italy
Tel. 0577-587111