

Mast House, Derby Road, Bootle, Merseyside L20 1EA, United Kingdom. Tel: +44 (0)151 933 7277 Fax: +44 (0)151 944 1332 www.mastgrp.com

# ENZYWELL HELICOBACTER PYLORI IgG EIA KIT

### Principle of the test

The DIESSE ENZYWELL HELICOBACTER PYLORI IgG kit utilises a stepwise EIA method for the determination of IgG class antibodies to HELICOBACTER PYLORI in human serum.

The antigen, HELICOBACTER PYLORI, is bound to the microtitre wells.

#### Features and benefits of the kit

The DIESSE ENZYWELL HELICOBACTER PYLORI IgG kit has been designed to be "user-friendly" and to minimise waste.

- 1. Short incubation time Results available within two hours.
- 2. 12 x 8 breakwell presentation This gives good flexibility to meet varying specimen throughput.
- **3. Sensitivity 96%, Specificity 99%** Meets MDA requirements, helps to ensure reliable results.
- 4. Ready for use reagents Calibrator, conjugate, substrate and stop solution, are all liquid and ready for use. This helps to reduce preparation time and minimise wastage of reagents.

#### **ORDERING INFORMATION**

PRODUCT CODE	91060
PACK SIZE	12 x 8 tests



## **TEST PROCEDURE FOR HELICOBACTER PYLORI IgG**

<u>STEP1</u>	Place 100μL of diluted sample/blank/calibrators in the wells of the strips ↓
	(Incubate for 45 mins at 37°C)
	Wash 4 times ↓
STEP2	Add 100µL of Conjugate to each well
	↓ (Incubate for 45 mins at 37° C)
	Wash 4 times
	$\downarrow$
STEP3	Add 100µL of Substrate to each well $\downarrow$
	(Incubate for 15 minutes at room temperature) $\downarrow$
<u>STEP4</u>	Add 100µL of Stop Solution $\downarrow$
	Read absorbance at 450nm within 15 minutes $\downarrow$
STEP5	Interpretation of results
	(O.D. sample > O.D. cutoff = positive result)
	For more details see pack insert

REFERENCES

Marshall B.J. and Warren J.R.: Unidentified curved bacilli in the stomach of patients with gastritis and peptic ulceration. Lancet i, 1311 (1984).
Jones D.M., Lessels A.M., Eldridge J.: Campylobacter-like organisms on the gastric mucosa: culture, histological and serological studies. J. Clin. Pathol. 37: 1002 (1984).
Blaser M.J.: H. pylori and the pathogenesis of gastroduodenal inflammation. J. Inf. Dis. 161: 626 (1990).
Valle J., SeppJ., K., Sipponen P., Kasunen T.U.: Disappearance of gastritis after eradication of H. pylori: a morphometric study. Scand. J. Gastroenterology 26: 1057 (1991).

G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
Laheij R. J. F. et al Evaluation of Commercially Available Helicobacter pylori Serology Kits: a Review. J Clin Micro 36:10: (1998)

FTL 07/02 V1.0