

Mast House, Derby Road, Bootle, Merseyside L20 1EA, United Kingdom. Tel: +44 (0)151 933 7277 Fax: +44 (0)151 944 1332 www.mastgrp.com



# ENZYWELL CYTOMEGALOVIRUS IgM EIA KIT

## Principal of the test

The DIESSE ENZYWELL CYTOMEGALOVIRUS IgM KIT is a stepwise EIA capture method for the determination of IgM class antibodies to Cytomegalovirus in human serum.

The capture is performed using monoclonal antibodies bound to the microtitre wells

The antigen is composed of purified Cytomegalovirus virus.

### Features and benefits of the kit

The DIESSE ENZYWELL CYTOMEGALOVIRUS IgM kit has been designed to be "user-friendly" and to minimise waste

- 1. IgM capture system is integral within test method
  Saves time
- 2. Short incubation time

Results available within two hours

- 2. 12 x 8 breakwell test presentation
  Good flexibility to meet varying specimen throughput
- 4. 96.4% Sensitivity, 100% Specificity

Helps to ensure valid results

5. Ready for use reagents

Positive, negative and cut-off controls, conjugate and stop solution are all liquid and ready for use.

Reduces preparation time and minimises wastage

### ORDERING INFORMATION

PRODUCT CODE 91011

PACK SIZE 12 x 8 tests



Mast House, Derby Road, Bootle, Merseyside L20 1EA, United Kingdom. Tel: +44 (0)151 933 7277 Fax: +44 (0)151 944 1332 www.mastgrp.com

# TEST PROCEDURE FOR CYTOMEGALOVIRUS IGM EIA **KIT**

STEP 1 Place 100µL of diluted sample/controls in the wells of the strips

Incubate for 45 mins at 37 C

Wash 4 times

STEP2 Add 100µL of Antigen/conjugate to each well

Incubate for 45 mins at 37 C

Wash 4 times

STEP3 Add 100µL of TMB Substrate to each well

Incubate for 15 minutes at room temperature

STEP4 Add 100µL of Stop Solution

Read absorbance at 450nm within 15 minutes

STEP5 Interpretation of results

O.D. sample > O.D. cutoff = positive result

For more details see pack insert

### REFERENCES

G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
H.O. Kangro: Cytomegalovirus serology: does it give the answer? Serodiagnosis and Immunotherapy 1: 91 (1987).
Grint P.C.A. et al.: Screening tests for antibodies to cytomegalovirus: an evaluation of five commercial products. J. Clin. Pathol. 38: 1059 (1985). 3.

Van Loon A.M. et al.: Direct enzyme-linked immunosorbent assay that uses peroxidase-labelled antigen for dtermination of immunoglobulin M antibody to cytomegalovirus. J. Clin. Microbiol. 13: 416

M. Musiani et al.: Rapid detection of antibodies against cytomegalovirus induced immediate early and early antigens by an enzyme linmed immunosorbent

<sup>5.</sup> assay. J. Clin. Pathol. 37: 122 (1984).

F. de Ory et al.: Serological diagnosis of cytomegalovirus infections: comparison of six commercial methods of ELISA. Serological methods of ELISA.

<sup>6.</sup> 423 (1988)

P. Dal Monte et al.: Citomegalovirus umano. Diagnosis 2: 67 (1990).