

Mast House, Derby Road, Bootle, Merseyside L20 1EA, United Kingdom. Tel: +44 (0)151 933 7277 Fax: +44 (0)151 944 1332 www.mastgrp.com

# ENZYWELL CYTOMEGALOVIRUS IgG EIA KIT

# Principle of the test

The DIESSE ENZYWELL CYTOMEGALOVIRUS IgG kit utilises a stepwise EIA method for the determination of IgG class antibodies to Cytomegalovirus IgG antibodies in human serum.

The antigen, Cytomegalovirus, is bound to the microtitre wells.

## Features and benefits of the kit

The DIESSE ENZYWELL CYTOMEGALOVIRUS IgG kit has been designed to be "user-friendly" and to minimise waste.

- 1. Short incubation time Results available within two hours.
- 2. 12 x 8 break-well presentation This gives good flexibility to meet varying specimen throughput.
- **3. Sensitivity 100%, Specificity 100%** Helps to ensure reliable results

#### 4. Ready for use reagents

Calibrators, control,conjugate, substrate and stop solution are all liquid and ready for use. This helps to reduce preparation time and minimise wastage of reagents.

5. Semi-quantitative May be used as an indication of

Immunity

### **ORDERING INFORMATION**

PRODUCT CODE	91010
PACK SIZE	12 x 8 tests



# **TEST PROCEDURE FOR CYTOMEGALOVIRUS IgG**

<u>STEP1</u>	Place 100μL of diluted sample/blank/calibrators in the wells of the strips ↓	
	(Incubate for 45 mins at 37°C)	
	Wash 4 times ↓	
STEP2	Add 100µL of Conjugate to each well ↓	
	(Incubate for 45 mins at 37° C)	
	Wash 4 times	
	$\downarrow$	
<u>STEP3</u>	Add 100μL of Substrate to each well ↓	
(Incubate for 15 minutes at room temperature) ↓		
<u>STEP4</u>	Add 100μL of Stop Solution ↓	
	Read absorbance at 450nm within 15 minutes $\downarrow$	
STEP5	Interpretation of results	
	(Ratio O.D. sample and O.D. cutoff $> 1.1 = positive result)$	
	For more details see pack insert	

REFERENCES

 G. B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
 H.O. Kangro: Cytomegalovirus serology: does it give the answer? Serodiagnosis and Immunotherapy 1: 91 (1987).
 Grint P.C.A. et al.: Screening tests for antibodies to cytomegalovirus: an evaluation of five commercial products. J. Clin. Pathol. 38: 1059 (1985).
 Van Loon A.M. et al.: Direct enzyme-linked immunosorbent assay that uses peroxidase-labelled antigen for dtermination of immunoglobulin M antibody to cytomegalovirus. J. Clin. Microbiol. 13: 416
 M. Musiani et al.: Rapid detection of antibodies against cytomegalovirus induced immediate early and early antigens by an enzyme linmed immunosorbent assay. J. Clin. Pathol. 37: 122 (1984).
 F. de Ory et al.: Serological diagnosis of cytomegalovirus infections: comparison of six commercial methods of ELISA. Serodiagnosis and Immunotherapy 2: 423 (1988).
 P. Dal Monte et al.: Citomegalovirus umano. Diagnosis 2: 67 (1990).