ENZYWELL

CYTOMEGALOVIRUS IgG AVIDITY

REF 91092 (24 tests)

Manufactured by: DIESSE Diagnostica Senese
Via delle Rose 10
53035 Monteriggioni (Siena) - Italy

INDICE / INDEX

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17. BIBLIOGRAFIA / REFERENCES

1. L. Grangeot-Keros, M.J. Mayaux, P. Lebon et al. Value of CMV IgG avidity index for the diagnosis of primary CMV infection in pregnant women. *J. Inf. Dis.* 175: 944 (1997).
2. N.K. Blackburn, T.G. Besselaar et al. Differentiation of primary cytomegalovirus infection from reactivation using the urea denaturation test for measuring antibody avidity. *J. Med. Virol.* 33: 6 (1991).
3. G.B. Wisdom. Enzyme-immunoassay. *Clin. Chem.* 22: 1243 (1976).

DIESSE Diagnostica Senese
Via delle Rose 10
53035 Monteriggioni (Siena) Italy.
Tel. 0577-587111

INSTRUCTIONS FOR USE

ENZYWELL CYTOMEGALOVIRUS IgG AVIDITY KIT

REF 91092

1. INTENDED USE:

IMMUNOENZYMATIC METHOD FOR THE DETERMINATION OF THE AVIDITY OF IgG-CLASS ANTIBODIES TO CYTOMEGALOVIRUS IN HUMAN SERUM

2. SUMMARY AND EXPLANATION OF THE TEST:

Cytomegalovirus is a herpes virus transmitted by close human contact. No symptoms of infection are apparent in the majority of cases. However, the virus is very dangerous and may be fatal in immunodepressed patients. Serum-negative female patients who become infected during pregnancy may transmit the disease to the fetus. In 95% of cases this occurs without symptoms, but some neonates may present jaundice, hepato-splenomegaly and retarded psycho-motorial development. For this reason it is of great importance to determine the immunitary state of the patient and check for serum conversion.

Sometimes it is not sufficient to test for IgM antibodies, as the presence of this class may be due to the persistence of IgM or asymptomatic reinfection by Cytomegalovirus, without risk for the fetus.

For this reason it is useful to assay the avidity of the IgG antibodies, in as far as it has been demonstrated that avidity is low in the acute phase of infection and increases over a period of time. The presence of low avidity antibodies is therefore an indication of a recent or current infection. The avidity of IgG antibodies in the serum can be assayed with this ELISA kit.

3. PRINCIPLE OF THE TEST:

The test is based on the ELISA technique (Enzyme linked Immunosorbent Assay) (1-3).

The partially purified and inactivated Cytomegalovirus antigen is bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with dilute human serum. The sample is placed in two wells (or in four wells if the test is performed in duplicate) and left to incubate.

After washing of the solid phase, the Avidity Buffer is placed in one well and the rinsing buffer in the other. The avidity buffer will cause dissociation of the antigen/antibody bond, proportional to the avidity of the antibodies. After washings to eliminate the proteins that have not reacted, incubation is performed with the conjugate composed of anti-human IgG monoclonal antibodies labelled with horse-radish peroxidase.

The unbound conjugate is eliminated and the peroxidase substrate specific for the enzyme is added.

The colour which develops is proportional to the concentration of specific antibodies present in the serum sample still bound to the antigen on the solid phase. The ratio between the OD of the well containing the Avidity buffer and that containing the Wash Buffer is expressed as a percentage and gives the degree of avidity of the serum sample for the Cytomegalovirus antigen.

4. KIT CONTENTS AND REAGENT PREPARATION:

Reagents are sufficient for 48 determinations.

⚠ Bring reagents to room temperature before use.

MT PLATE MICROPLATE COATED WITH CYTOMEGALOVIRUS. 3x16 wells.

Use: open the package at the opposite end from the code (C, followed by the lot number) which is useful for identification purposes, remove the support and wells to be used from the foil package, and place the unused wells in the polythene bag with the silica gel, expel the air and seal by pressing the closure.

CONJ CONJUGATE. 1x12 mL.

Contents: anti-human IgG monoclonal antibodies labelled with Peroxidase, in phosphate buffer solution. Ready for use without further dilution.

CONTROL HIGH HIGH AVIDITY CONTROL (2.5 mL)

Contents: Human serum containing low avidity IgG anti-CMV in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

CONTROL LOW *LOW AVIDITY (2.5 mL) CONTROL (2.5 mL)*

Contents: Human serum containing high avidity IgG anti-CMV in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

AVIDITY BUF *AVIDITY BUFFER (1 x 10 mL)*

Contents: Urea in buffered saline solution containing phenol (0.05%) and Bronidox (0.02%), ready for use. If crystals are present, they should be dissolved at 37°C before use.

WASH BUF 10x *WASH BUFFER 10X. 1 x 100 mL.* **INTERCHANGEABLE BETWEEN LOTS**

Contents: Phosphate buffered saline, concentrated 10 times; contains Brij 0.5% w/v.

Preparation: dilute 1:10 with distilled water in order to obtain the washing buffer ready for use. If crystals are present, they should be dissolved at 37°C before dilution.

SAMP DIL 50x *DILUENT 50X. 1 x 4.5 mL. For dilution of serum samples.* **INTERCHANGEABLE BETWEEN LOTS**

Contents: Proteic solution conc. 50 times, containing BSA 5% w/v and bronidox 0.02%.

Preparation: Dilute 1:50 in the washing buffer to obtain the diluent ready for use.

SUBS TMB *SUBSTRATE. 12 mL. Ready for use.* **INTERCHANGEABLE BETWEEN LOTS**

Contents: Tetramethylbenzidine and hydrogen peroxide stabilised in citrate buffer, pH 3.8.

H₂SO₄ 0.3 M *STOP SOLUTION. 1x16 mL.* **INTERCHANGEABLE BETWEEN LOTS**

H₂SO₄ 0.3 mol/L, in solution ready for use.

ADHESIVE FILMS (2).

POLYTHENE BAG (1).

MATERIALS REQUIRED BUT NOT PROVIDED.

- Microplate reader, wavelength 450 nm or 450/620 nm, with OD linearity up to 2,000 (at least).
- Microplate washer (preferable) able to dispense volumes in the range 225-375 µL
- Distilled or deionized water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 10, 100, 1000 µl solution
- Disposable gloves
- Timer
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials
- Absorbent tissue.

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

Reagents have a limited stability after opening and/or preparation

REAGENT	CONDITIONS
Microplate	5 weeks at 2/8°C, polythene bag
Control sera	5 weeks at 2/8°C
Conjugate	5 weeks at 2/8°C
Substrate	up to the expiry date at 2/8°C, 1 week at 15-30°C; store in the dark
Sample Diluent	ready for use, 2 weeks at 2/8°C
Wash Buffer	2 weeks at 2/8°C, 5 days at 15/30°C.
Avidity Buffer	up to the expiry date at 2/8°C
Stop Solution	up to the expiry date at 2/8°C

6. PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY.

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
2. The following reagents contain low concentrations of harmful or irritant substances:
 - a) The Wash Buffer contains detergents
 - b) The conjugate contains phenol
 - c) The substrate is acid
 - d) Calibrators contain Sodium Azide 0.09% as preservative which may react with copper and lead in plumbing to form potentially explosive metal azides. On disposal, flush with large volumes of water to prevent azide buildup.

If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
3. Non-disposable apparatus should be sterilized after use. The preferred method is to autoclave for 1 h at 121°C; disposables should be autoclaved or incinerated.
4. Sulphuric acid required for the Stop Solution and hydrochloric acid used for washing glassware are corrosive and should be handled with appropriate care. If they come into contact with the skin or eyes, wash thoroughly with water.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

1. Allow all reagents and samples to come to room temperature (18-30°C) before use. Immediately after use return reagents to the recommended storage temperature. **It is important to work at the correct temperature. Check that the thermostat does not go below 35°C or over 39°C.** Open the envelope containing the strips after at least ½ hr at room temperature.
2. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
3. Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.
4. Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or high quality deionized water.
5. Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps.
6. Do not allow wells to become dry during the assay procedure.
7. Care must be taken not to cross-contaminate reagents. It is important that pipettes are dedicated for exclusive use with the various reagents.
8. Care should be taken to avoid touching or splashing the rim of the well with conjugate. Do not "blow-out" from microplates.
9. Enzyme immunoassays can occasionally exhibit an "edge effect" which must be minimized by increasing the humidity during incubation steps. Plates must be covered with their covers and incubated at 37°C either in a water bath with a rack or float to support the plates if necessary, or in an incubator. Alternatively, plates can be incubated in an approved analyzer. See the appropriate operating manual for further details. CO₂ incubators must not be used.
10. Ensure that the bottom of the plate is clean and dry, and that no bubbles are present on the surface of the liquid before reading the plate.
11. Use of highly hemolyzed samples, incompletely clotted sera, or samples with microbial contamination may give rise to erroneous results.

~~16.12.~~ The use of the kit with automated instruments must be validated by the user.

~~17.13.~~ For each instrument used, read the manufacturer's instructions manual carefully to obtain additional information on the following points:

- installation and particular requisites
- operating principles, instructions, precautions and risks

- manufacturer's specifications and instrument performance
- servicing and maintenance.

7. TYPE AND STORAGE OF SAMPLE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C; avoid repeated freezing/thawing cycles.. Defrosted samples must be carefully shaken before performing the test. Heat inactivation can lead to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric or contaminated samples should be avoided. If a new sample cannot be obtained, such samples should be clarified by filtration (0.45 µm) or centrifugation (3000 rpm x 10').

The test cannot be performed on human plasma.

8. TEST PROCEDURE:

~~2.-~~ Prepare the required number of strips.

~~3.-~~ Prepare the washing buffer by diluting the Wash Buffer 10x (100 mL + 900 mL H₂O).

~~4.-~~ Prepare the diluent for samples by adding Diluent 50x to the washing buffer (add 4 mL to 200 mL).

Leave a well empty for the test blank. Dilute samples 1:301 distributing 5 µl of serum into 1.5 mL of diluent. Dispense 100 µL of each diluted sample in at least two wells, and 100 µl of High Avidity Control and 100 µl of Low Avidity Control in duplicate. Wells are covered with protective film and incubated for 60 minutes at 37°C. After washing four times for 30 seconds (300 µL), dispense Wash Buffer in one well of the serum samples and of the controls and Avidity Buffer in the other, and incubate for 10' at 37°C. Aspirate the buffer and perform 4 washings for 30 seconds each. Add 100 µl of conjugate to each well and incubate again for 45 minutes at 37°C, covering the wells with the protective film. The plate is washed again 4 times, as described above. Finally, the substrate is distributed, 100 µl/well., including the blank. After 15 minutes at room temperature the enzymatic reaction is stopped with 100 µl of Stop Solution. The adsorbance (O.D.) is read at 450 nm or at 450/620 nm within 30 min.

9. TEST PROCEDURE FOR CYTOMEGALOVIRUS IgG AVIDITY KIT	
--	--

STEP 1	Place 100 µl of diluted sample /controls in duplicate in the wells of the strips.
--------	---

-
Incubate for 60 min. at 37°C

-
Wash 4 times (30 seconds each, 300 µL)

STEP 2	For each sample/control , place 100 µl of Wash Buffer in one well and 100 µl of Avidity buffer in the other; incubate for 10' at 37°C
--------	---

-
Wash four times (30 seconds each, 300 µL)

STEP 3	Add 100 µl of conjugate to each well
--------	--------------------------------------

-
Incubate for 45 min. at 37°C

-
Wash 4 times

STEP 4	Add 100 µl of Substrate to each well including the blank
--------	--

-
Incubate for 15 min. at R.T.

STEP 5	Add 100 µl of Stop Solution
--------	-----------------------------

-
Read absorbance at 450 nm within 30 min

10. QUALITY CONTROL VALUES

The OD of the Blank must be ≤ 0.15.

The high and low avidity controls in the wells incubated with the Wash Buffer must give an OD value > 0.500. The percentage of avidity of the low control must be less than 35%, and the high avidity control must be higher than 50% compared to that found with the Wash Buffer. (Calculation performed as reported in the following paragraph).

11. INTERPRETATION OF THE RESULTS

The percentage of avidity of the samples is expressed and calculated using the ratio between the OD found in the wells containing Avidity buffer and those with Wash Buffer., subtracting the value of the blank.

$$\frac{\text{O.D. with Avidity Buffer} - \text{O.D. blank}}{\text{O.D. with Wash Buffer} - \text{O.D. blank}} \times 100 = \% \text{ of avidity}$$

O.D. with Wash Buffer – O.D. blank

A ratio over 40% indicates the presence of high avidity IgG antibodies; when it is lower than 30%, this indicates the presence of low avidity IgG. If the percentage is between 30 and 40%, there is a medium degree of avidity (borderline). The calculation can only be performed on samples with give an OD \geq 0.40 in the well washed with Wash Buffer.

Warning: In the case of O.D. $<$ 0.40, repeat the test diluting the sample 1/101.

12. LIMITATIONS OF THE TEST

A finding of high avidity (especially in transplanted subjects) does not exclude the possibility of a recent infection. In the case of newborn babies, low avidity is not a sign of an acute infection in progress.

13. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

In an external evaluation, 100 serum samples were tested with two commercial kits for the determination of anti-CMV IgG and IgM. The clinical picture of all patients was known.

The % avidity calculated was correlated with the stage of infection. Of the 100 samples tested, 14 were from patients with an acute infection, 15 from patients who had undergone transplant (7 of which with an acute infection), the rest came from patients with a prior infection.

The percentage of avidity in the 14 patients with acute infections varied from a minimum of 6 to a maximum, of 35-40 (advanced primary infections); in the cases of previous infections the avidity varied from a minimum of 50 to a maximum of 100 (\geq 70 in 61% of cases). The transplanted patients showed high avidity levels; the result was falsified by the administration of gammaglobulins, which is often given in these patients).

14. PRECISION**Table 1 Precision “in run” (n=6)**

Sample	average % avidity	S.D.	CV%
Low avidity	28	2	6
High avidity	92	3	3

Table 2 Precision “between runs”

Sample	average % avidity	S.D.	CV%
S1	72	1	15
S2	31	1	2
S3	8	2	28

15. TROUBLE SHOOTING GUIDE

PROBLEM	POSSIBLE SOURCE	TEST OR ACTION
Invalid run (all negative)	One or more reagents not added or added in wrong sequence	Recheck procedure Check for unused solutions. Repeat test.
	Unreactive plate	Check the code on the package containing the plate (see package insert point 4 for correct code).
		Check for moisture in unused plate. (Silica gel dessiccant must be pale yellow).Repeat test
Invalid run (all positive)	Contamination of substrate	Take new aliquot of substrate.
	Inadequate washing	Ensure that wash apparatus works well

Poor precision	Incomplete washing of wells	Ensure that wash apparatus works well
	Inadequate aspiration of wells	Ensure that wash apparatus works well
	Pipetting error	Check pipette function
	Reagent addition too slow	Avoid drying of the plate after washing step. Add reagents immediately
	Presence of bubbles	Avoid air bubbles during pipetting.
	Optical pathway not clean	Check instrument light source and detector for dirt. Wipe bottom of plate with soft tissue.
Inadequate Color development	Incorrect incubation times or temperature	Check for temperature control and time monitoring
		Adhere to recommended instruction for use.
	Inadequate volume of substrate added to the plate	Check pipette function.

16. REFERENCES

1. L. Grangeot-Keros, M.J. Mayaux, P. Lebon et al. Value of CMV IgG avidity index for the diagnosis of primary CMV infection in pregnant women. *J. Inf. Dis.* 175: 944 (1997).
2. N.K. Blackburn, T.G. Besselaar et al. Differentiation of primary cytomegalovirus infection from reactivation using the urea denaturation test for measuring antibody avidity. *J. Med. Virol.* 33: 6 (1991).
3. G.B. Wisdom. Enzyme-immunoassay. *Clin. Chem.* 22: 1243 (1976).

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