

IVD solutions through partnership

**Mast  
Group**



# **MASTDISCS<sup>®</sup>** *Combi*

**ESBL & AmpC Detection Disc Sets - D68C**

- Differentiate resistance enzyme types
- Simple comparative interpretation
- Low cost implementation
- Compliance with international standard methodologies (EUCAST & CLSI)

# EXTENDED SPECTRUM $\beta$ LACTAMASE DETECTION

Extended spectrum beta-lactamases (ESBL) are bacterial enzymes which confer resistance to penicillin and cephalosporin antibiotics. The emergence of ESBL producing pathogens has become increasingly significant in limiting the antibiotic treatment options, representing a serious complication for antibiotic management.

Mast Group Ltd. range of ESBL detection discs offers laboratories a simple, reliable and low cost means of identification and detection by double disc diffusion, using paired and combination disc sets.

The presence of an ESBL and/or AmpC is easily determined by zone size comparison when simultaneously tested with antibiotic and antibiotic plus inhibitor combinations.

## Visual of results

**D68C<sup>1</sup>**

**AmpC & ESBL Detection Set**

Confirmation of AmpC and/or ESBL production in Enterobacteriaceae

A CPD10

B CPD10 + ES $\beta$ L inhibitor

C CPD10 + AmpC inhibitor

D CPD10 + ES $\beta$ L inhibitor  
+ AmpC Inhibitor


Calculator programme is  
available to download from  
[www.mastgrp.com](http://www.mastgrp.com)

**ESBL Positive**




$Z_B - Z_A$  and  $Z_D - Z_C \geq 5\text{mm}$   
and the differences of each of  
 $Z_B$  and  $Z_D$  and  $Z_A$  and  $Z_C$  are  $< 4\text{mm}$

**AmpC Positive**




$Z_D - Z_B$  and  $Z_C - Z_A \geq 5\text{mm}$   
and the differences of each of  
 $Z_A$  and  $Z_B$  and  $Z_C$  -  $Z_D$  are  $< 4\text{mm}$

**AmpC and ESBL Positive**



$Z_D - Z_C \geq 5\text{mm}$  and the difference of  
 $Z_A$  and  $Z_B$  is  $< 4\text{mm}$

**AmpC and ESBL Negative**



All zones differ by  $\leq 2\text{mm}$

CPD – Cefpodoxime



# D69C<sup>2</sup>

## AmpC Detection Set

Confirmation of either chromosomal or plasmid acquired AmpC

A CPD10 + AmpC inducer

### AmpC Positive

B CPD10 + AmpC inducer  
+ ESβL inhibitor


C CPD10 + AmpC inducer  
+ ESβL inhibitor  
+ AmpC inhibitors



$Z_C - Z_A$  and  $Z_C - Z_B \geq 5\text{mm}$

CPD – Cefpodoxime

### AmpC Negative



All zones differ by  $\leq 3\text{mm}$

# D63C<sup>3</sup>

## Cefepime 30 & Cefepime 30/Clavulanic Acid 10

Confirmation of ESβL production in Enterobacteriaceae with chromosomal AmpC

CPM30

CPM30/CLAV10

### ESBL Positive




$Z_2 - Z_1 \geq 5\text{mm}$

CPM – Cefepime

CLAV – Clavulanic Acid

### ESBL Negative



All zones differ by  $\leq 2\text{mm}$

# D52C<sup>4</sup>

## Extended Spectrum βLactamase Set

Confirmation of ESβL production in Enterobacteriaceae with no chromosomal de-repressed or inducible AmpC

CAZ30

CAZ30/CLAV10


CTX30

CTX30/CLAV10

CPD30

CPD30/CLAV10

### ESBL Positive



$Z_4 - Z_1 \geq 5\text{mm}$  and/or  
 $Z_4 - Z_3 \geq 5\text{mm}$  and/or  $Z_6 - Z_5 \geq 5\text{mm}$


CAZ – Ceftazidime

CLAV – Clavulanic Acid

CTX – Cefotaxime

CPD – Cefpodoxime

### ESBL Negative



All zones differ by  $\leq 2\text{mm}$

# D67C<sup>5</sup>

## Extended Spectrum $\beta$ Lactamase Set (CPD10)

Confirmation of ES $\beta$ L production in Enterobacteriaceae with no chromosomal de-repressed or inducible AmpC


CAZ30  
CAZ30/CLAV10

CTX30  
CTX30/CLAV10

CPD10  
CPD10/CLAV1


CAZ - Ceftazidime  
CLAV - Clavulanic Acid  
CTX - Cefotaxime  
CPD - Cefpodoxime

### ESBL Positive



$Z_2 - Z_1 \geq 5\text{mm}$  and/or  $Z_4 - Z_3 \geq 5\text{mm}$  and/or  $Z_6 - Z_5 \geq 5\text{mm}$

### ESBL Negative



All zones differ by  $\leq 2\text{mm}$

# D62C<sup>6</sup>

## Cefotaxime 30 & Cefotaxime 30 /Clavulanic Acid 10


# D64C<sup>6</sup>

## Ceftazidime 30 & Ceftazidime 30/Clavulanic Acid 10

Confirmation of ES $\beta$ L production in Enterobacteriaceae with no chromosomal de-repressed or inducible AmpC when both disc sets used concurrently

CTX30  
CTX30/CLAV10  
CAZ30  
CAZ30/CLAV10


### ESBL Positive



$Z_2 - Z_1 \geq 5\text{mm}$  and/or  $Z_4 - Z_3 \geq 5\text{mm}$

CTX - Cefotaxime  
CAZ - Ceftazidime  
CLAV - Clavulanic Acid

### ESBL Negative



All zones differ by  $\leq 2\text{mm}$

# D66C<sup>7</sup>

## Cefpodoxime 10 & Cefpodoxime 10/Clavulanic Acid 1

Confirmation of ES $\beta$ L production in Enterobacteriaceae with no chromosomal de-repressed or inducible AmpC

CPD10  
CPD10/CLAV1


### ESBL Positive



$Z_2 - Z_1 \geq 5\text{mm}$

CPD - Cefpodoxime  
CLAV - Clavulanic Acid


### ESBL Negative



All zones differ by  $\leq 2\text{mm}$

# Flowchart for laboratory use of **MASTDISCS® Combi**

- D52C - Extended Spectrum  $\beta$  Lactamase Set (50 tests)
- D62C - Cefotaxime 30ug & Cefotaxime 30ug/Clavulanic Acid 10ug (150 tests)
- D63C - Cefepime 30ug & Cefepime 30ug/Clavulanic Acid 10ug (150 tests)
- D64C - Ceftazidime 30ug & Ceftazidime 30ug/Clavulanic Acid 10ug (150 tests)
- D66C - Cefpodoxime 10ug & Cefpodoxime 10ug/Clavulanic Acid 1ug (150 tests)
- D67C - Extended Spectrum  $\beta$  Lactamase Set (CPD10) (50 tests)
- D68C - AmpC & ESBL Detection Set (50 tests)
- D69C - AmpC Detection Set (50 tests)
- D71C - MAST® CAT-ID - Faropenem discs (250 tests)
- D73C - MAST® Carba plus Detection Set (50 tests)
- D74 - MAST® ICT Indirect Carbapenemase Test (25 tests)



\* Concentrations and combinations are dependent on International Standard used.  
 NB. All carbapenemase producing organisms should be sent to a reference laboratory for confirmation

## Ordering Information

Order Code	Product/Usage	Pack size	No. Tests
171682	<b>MASTDISCS® Combi AmpC &amp; ESBL Detection Set - D68C<sup>1</sup></b> Confirmation of AmpC and/or ESBL production in isolates of Enterobacteriaceae. When interpreted as 'further work required' use D69C to confirm AmpC production and D63C for confirmation of ESBL production when AmpC is also present.	4 × 50 discs	50 tests
171692	<b>MASTDISCS® Combi AmpC Detection Set -D69C<sup>2</sup></b> Confirmation of AmpC production in isolates of Enterobacteriaceae with either plasmid acquired or chromosomal AmpC. Can be used when interpreted as 'further work required' on D68C for confirmation of AmpC production.	3 × 50 discs	50 tests
171632	<b>MASTDISCS® Combi Cefepime ESBL ID Disc Set - D63C<sup>3</sup></b> Confirmation of ESBL production in isolates of Enterobacteriaceae with chromosomal AmpC e.g. <i>Enterobacter</i> spp., <i>Citrobacter freundii</i> , <i>Morganella morganii</i> , <i>Providencia</i> spp., <i>Hafnia alvei</i> , <i>Serratia</i> spp. Can be used when interpreted as 'further work required' on D68C for confirmation of ESBL production when AmpC is also present.	6 × 50 discs	150 tests
701055	<b>MASTDISCS® Combi ESBL ID Disc Set - D52C<sup>4</sup> (EUCAST &amp; CLSI)</b> Confirmation of ESBL production in isolates of Enterobacteriaceae with no chromosomal de-repressed or inducible AmpC e.g. <i>Escherichia coli</i> , <i>Klebsiella</i> spp., <i>Proteus mirabilis</i> , <i>Salmonella</i> spp., <i>Shigella</i> spp.	6 × 50 discs	50 tests
171672	<b>MASTDISCS® Combi ESBL ID Disc Set (CPD10) - D67C<sup>5</sup> (EUCAST &amp; CLSI)</b> Confirmation of ESBL production in isolates of Enterobacteriaceae with no chromosomal de-repressed or inducible AmpC e.g. <i>Escherichia coli</i> , <i>Klebsiella</i> spp., <i>Proteus mirabilis</i> , <i>Salmonella</i> spp., <i>Shigella</i> spp.	6 × 50 discs	50 tests
171622	<b>MASTDISCS® Combi Cefotaxime ESBL ID Disc Set D62C<sup>6</sup> (CLSI)</b> Confirmation of ESBL production in isolates of Enterobacteriaceae with no chromosomal de-repressed or inducible AmpC e.g. <i>Escherichia coli</i> , <i>Klebsiella</i> spp., <i>Proteus mirabilis</i> , <i>Salmonella</i> spp., <i>Shigella</i> spp.	6 × 50 discs	150 tests
171642	<b>MASTDISCS® Combi Cefepime ESBL ID Disc Set - D64C<sup>6</sup> (CLSI)</b> Confirmation of ESBL production in isolates of Enterobacteriaceae with no chromosomal de-repressed or inducible AmpC e.g. <i>Escherichia coli</i> , <i>Klebsiella</i> spp., <i>Proteus mirabilis</i> , <i>Salmonella</i> spp., <i>Shigella</i> spp. <b>D62C &amp; D64C must be used concurrently</b>	6 × 50 discs	150 tests
171662	<b>MASTDISCS® Combi Cefpodoxime ESBL ID Disc Set - D66C<sup>7</sup> (EUCAST)</b> Confirmation of ESBL production in isolates of Enterobacteriaceae with no chromosomal de-repressed or inducible AmpC e.g. <i>Escherichia coli</i> , <i>Klebsiella</i> spp., <i>Proteus mirabilis</i> , <i>Salmonella</i> spp., <i>Shigella</i> spp.	6 × 50 discs	150 tests

## Additional products from Mast Group Ltd.'s Range

Order Code	Product	Pack size	No. Tests
171712	<b>MASTDISCS® ID - Carbapenemase Activity Test (CAT) - D71C</b>	6 × 50 discs	250 tests
171722	<b>MASTDISCS® Combi AmpC, ESBL &amp; Carbapenemase Detection Set - D72C</b>	6 × 50 discs	50 tests
171732	<b>MASTDISCS® Combi Carba plus - D73C</b>	5 × 50 discs	50 tests
171742	<b>MAST®/ID - Indirect Carbapenemase Test (ICT) - D74</b>	25 devices	25 tests

v 5.0 CA 09-18 (MD2019-01-30)

**United Kingdom**  
Mast Group Ltd.  
Mast House, Derby Road, Bootle  
Merseyside L20 1EA

Tel: +44 (0)151 933 7277  
Fax: +44 (0)151 944 1332  
e-mail: sales@mastgrp.com

**Germany**  
Mast Diagnostica GmbH  
Feldstraße 20  
DE-23858 Reinfeld

Tel: +49 (0)4533 2007 0  
Fax: +49 (0)4533 2007 68  
e-mail: mast@mast-diagnostica.de

**France**  
Mast Diagnostic  
12 rue Jean-Jacques Mention  
CS91106, 80011 Amiens, CEDEX 1

Tél: +33 (0)322 80 80 67  
Fax: +33 (0)322 80 99 22  
e-mail: info@mast-diagnostic.fr

[www.mast-group.com](http://www.mast-group.com)