MASTABLOT™ TP

| | | |
|-------------------|---|---------|
| MASTABLOT™ TP IgG | REF 6653G08 UDI-DI 4250729700149 | 8 Tests |
| MASTABLOT™ TP IgM | REF 6653M08 UDI-DI 4250729700156 | 8 Tests |

Gebrauchsanweisung / Instructions for Use / Notice d'utilisation

**Nur zur Verwendung durch Fachpersonal /
For professional use only
Exclusivement pour un usage professionnel**

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MASTABLOT™ TP IgG / IgM

Intended Purpose

Qualitative indirect immunoblot for the detection of IgM / IgG antibodies against *T. pallidum* antigens in human serum and plasma as an aid in the clarification of an acute or past syphilis infection.

The assay is suitable for manual use and is intended for professional in-vitro diagnostic use only. All laboratory test results should be interpreted in conjunction with other clinical data. The clinical judgement and further tests have to be taken into account additionally.

Important Note for Use of these Kit Instructions

Any assay relevant changes to the kit instructions for use (IFU) will lead to a change of the version number at the bottom on the last page. All assay relevant changes made will be identified in the change history at the end of this document. Please ensure that the latest version of the IFU is used for the assay procedure.

Principle of the Test

Antigens specific for *T. pallidum* are coated onto nitrocellulose strips at defined positions. Each strip contains a cut-off-control and one reaction control. Patient serum or plasma specimens are applied together with the strips in the incubation chamber and incubated. If specific antibodies are present in the serum / plasma they will bind to the fixed antigens on the strips forming a stable antigen-antibody complex. Strips are then washed to remove any unbound materials and complexed antibodies are detected by the addition of an alkaline phosphatase labelled anti-human immunoglobulin conjugate.

After a further washing step to remove any unbound conjugate, strips are incubated with BCIP/NBT substrate to form a visible dark bluish precipitate. Adding water to this substrate reaction will stop the formation of the precipitating colour.

Kit Contents

1. **STRIPS** Nitrocellulose strips
8 Nitrocellulose strips coated with recombinant *T. pallidum* antigens (p15, p17, TmpA, p47); depending on the intended use of the kit the strips are specific for IgG or IgM detection.
2. **BUFFER** PBS
15 mL sample diluent. Green coloured solution, ready to use, contains proclin (< 0.1 %)
3. **CONJ G** **CONJ M** Conjugate
13 mL conjugate; according to kit specification:

Anti-human-IgG H+L chain (goat), labelled with alkaline Phosphatase.
/
Anti-human-IgM μ -chain (goat) labelled with alkaline Phosphatase.

Red coloured solution, ready for use, the buffer contains 0.02 % methylisothiazolone and 0.02 % bromnitrodioxane.
4. **SUBS** Substrat
13 mL BCIP/NBT solution, ready to use
5. **WASH** **CONC** Wash buffer,
15 mL Waschpufferkonzentrat 10 x concentrated, contains Proclin (< 0.1%)
6. Incubation Chamber
7. Interpretation / documentation sheet

Further Abbreviations

RTU Ready to use

Materials required but not provided

1. Sterile test tubes.
2. 5 μ L-, 100 μ L- and 500 μ L micropipettes or a multichannel pipette (optional).
3. Disposable pipette tips
4. Aqua dest. or HPLC grade water.
5. Measuring cylinder.
6. Forceps (optional).
7. Shaker.
8. Wash bottle.

Warnings

1. Comply with the general health and safety guidelines for working with potentially infectious materials. Wear appropriate protective clothing and use appropriate lab facilities.
2. Do not cross-contaminate reagents or interchange caps on bottles. Do not use contaminated reagents. Do not use reagents with damaged bottles or packaging. Use a separate pipette or pipette tips for each sample and reagent.
3. In the event that haemolysed or lipaemic serum must be used, heat inactive at 56 °C before use.
4. The kit reagents may contain bovine components. To protect against infections with prions or zoonotic pathogens, standard occupational safety regulations (GLP) must be observed.
5. Proclin is used as a preservative as marked. It may be toxic if ingested. Always dispose of by flushing to drain with plenty of water.

Precautions / Notes

1. Read instructions carefully before conducting the assay. Do not modify the procedure without prior validation.
2. Do not use reagents beyond the expiry date.
3. Do not mix reagents between different lots, as reagents have been calibrated for each batch.
4. Use disposable plasticware wherever possible. Reusable glassware should be washed thoroughly and rinsed free of detergents before use.
5. Always use a new pipette tip for each pipetting step.
6. Only use laboratory water of the highest quality (e.g. Aqua dest., Aqua bidest., Aqua iniectionabilia, HPLC grade).
7. Do not allow nitrocellulose strips to dry out during the assay procedures.
8. Do not expose slides to intense sunlight or similar adverse conditions while incubating.

Disposal

Kit reagents, samples and contaminated disposables should be disposed of in accordance with the relevant disposal guidelines and regulations for infectious materials.

Stability and Storage

MASTABLOT™ TP IgG / IgM can be used until the end of expiry date as indicated on kit label. All kits components and reagents should be stored at 2–8 °C and brought to room temperature before use.

Diluted wash buffer can be stored at 2–8 °C for up to 30 days.

Sample Material

- a) Serum, Plasma
- b) Samples can be stored at 2–8 °C for up to 3 days. Serum should be aliquoted immediately after sampling and kept at -20°C for longer storage. The samples should not be frozen and thawed repeatedly. After thawing samples should be briefly vortexed carefully before being used in the assay. Lipaemic, haemolytic, icteric or bacterially contaminated samples can cause erroneous results (Ref. 6, 7, 8).

Test Procedure

1. Optional: Inactivate all sera by heating for 30 min at 56 °C. Do not inactivate plasma.
2. Allow all materials to reach room temperature prior to use.
3. Dilute the wash buffer concentrate 1:10 with distilled water:
Example:
25 mL wash buffer + 225 mL of distilled water
4. **Important note:** Always use the IgG labelled strips for IgG determinations, and IgM-labelled strips for the IgM determinations. Only touch the strips at the label end. Always use clean disposable gloves when touching the nitrocellulose strips! Alternatively use forceps for touching strips. Strips must be soaked in the incubation chamber with serum diluent before use. The nitrocellulose side of the strips with the label on top should be upward. For soaking add 1.5 mL of serum diluent to the strip in the incubation chamber. Incubate the strips with the diluent for 10 min ± 2 min.
5. For serum incubation directly add the specimen to the diluent into the incubation chamber.

Identification of IgG

Sample dilution for IgG 1:100:

Add 15 µL of serum / plasma to the diluent.

1. Incubate IgG assays for 60 min at room temperature on a shaker.
2. Carefully discard the serum diluent solutions of each incubation chamber. The use of appropriate pipettes or a vacuum pump is helpful. Immediately add 1.5 mL of 1x wash buffer to each strip. Gently shake for 5 min, discard the wash buffer then repeat the washing step 2–3 times.
3. After the washing step add 1.5 mL of IgG conjugate to each strip. The conjugate should cover the strips completely.
4. Incubate the **assay 60 min** while gently shaking.
5. After the conjugate reaction wash the strips as described in step 2 (Identification of IgG).
6. Add 1.5 mL of BCIP/NBT substrate solution to each strip.

Gently shake the IgG reactions for 5–8 min at room temperature.

The substrate reaction can be stopped at any time by repeatedly adding tap water or distilled water. It is recommended to stop the reaction immediately after a bluish background colour starts to develop heavily.

Identification of IgM

Sample dilution for IgM 1:50:

Add 30 µL of serum / plasma to the diluent.

1. Incubate IgM assays for 120 min at room temperature on a shaker.
2. Carefully discard the serum diluent solutions of each incubation chamber. The use of appropriate pipettes or a vacuum pump is helpful.

Immediately add 1.5 mL of 1x wash buffer to each strip. Gently shake for 5 min, discard the wash buffer then repeat the washing step 2–3 times.

3. After the washing step add 1.5 mL of IgM conjugate to each strip.

The conjugate should cover the strips completely.

4. Incubate the assay 60 min while gently shaking.
5. After the conjugate reaction wash the strips as described in step 2 (Identification of IgM).
6. Add 1.5 mL of BCIP/NBT substrate solution to each strip.

Gently shake the IgM reactions for 8-10 min at room temperature.

The substrate reaction can be stopped at any time by repeatedly adding tap water or distilled water. It is recommended to stop the reaction immediately after a bluish background colour starts to develop heavily.

Interpretation of Results

Controls have to fit the lot specific reactions stated in the QC certificate. The figures in the table below are for orientation only and do not reflect the data of the current lot.

| Sample | Expected result |
|------------------|-----------------|
| Positive control | 3 + bis 4 + |
| Negative control | negative |

For a correct interpretation the results should be compared with positive and negative controls.

If further dilution on the samples or controls are requested on the basis of the screening titer, this can be achieved by serial dilution with PBS on screening titer.

Interpretation of Results

Test Validation

A reaction control band as well as a serum control band are coated onto each strip. The reaction control band must show a clear intense signal.

Each strip is coated with a specific cut-off defining control band, too, which is at the lower end of the strip.

The cut-off band is of low reactivity in both IgG and IgM assay versions.

Control band reactivities shall be as follows:

| | Expected result |
|---------------------------|--------------------|
| Serum control band IgG | 1+ to 2+ intensity |
| Serum control band IgM | 1+ to 2+ intensity |
| Reaction control band IgG | 2+ to 3+ intensity |
| Cut-off band IgG | Weak intensity (+) |
| Reaction control band IgM | 2+ to 3+ intensity |
| Cut-off band IgM | Weak intensity (+) |

As long as strips have not dried completely cut-off bands appear to be very weak. With drying the intensity of the cut-off band increases. It reaches its optimal visibility after the strips are completely dried.

A result is only valid if the serum, reaction and cut-off bands are visible. The intensity of each band should be within the specifications listed in the table above or the QC document.

Interpretation of Specimen Results

It is advisable to remove the strips while still wet to fix them onto the lamination. If necessary, apply clear tape over the labelled end to ensure adhesion to the page.

After completely drying the band on the strips can be easily read and interpreted. Drying can be fastened by using a hairdryer or a placing the strips into a drying incubator.

IgG positive reactivity:

IgG positive specimens must show at least **1 strong reactive antigen band**, which should belong to one of the antigen specificities: p15, p17, or p47, respectively. The intensity of that strongly reacting band must be far beyond the strength of the cut-off control band.

In most cases further bands are visible, too. These bands must have an intensity equal to or greater than the cut-off band.

The TmpA band may confirm a positive result. However, a TmpA band on its own is not sufficient to be considered a positive result.

Blots with **negative specimens** do not show any reactive antigen bands of the p15, p17, TmpA and p47 specificity. Strips with just 1 single band in the intensity of the cut-off band are considered as positive.

IgM positive reactivities:


IgM positive specimens must show at least **1 band** of either the p15, p17 or p47 antigen specificity. The intensity of the band must be equal to or greater than the cut-off control.

The TmpA band may confirm a positive result. However, a TmpA band on its own is not sufficient to be considered a positive result.

Blots with **negative specimens** do not show any reactive antigen bands of the p15, p17, TmpA and p47 specificity. Any antigen band reactivity less than the cut-off band must be read as negative.

IgG-Blot example

Abbreviation: CO = cut-off, SC = Serum Control, RC = reactive control

Correlation of the results to the MASTAFLUOR™ FTA-ABS

Low reactivity in MASTABLOT™ TP IgM does correlate with low reactivity in 19S IgM FTA-ABS (titre ≤ 1:20).

Limitations / Interferences

1. The MASTABLOT™ TP IgG / IgM test may occasionally give a false positive result under certain patient conditions or diseases e.g. pregnancy, leprosy and systemic lupus erythematosus.
2. The MASTABLOT™ TP IgG / IgM test has been demonstrated to be highly sensitive and specific, however results should be considered along with all serological tests, clinical history and other aspects of patient management to be considered diagnostically significant.

Performance Characteristics

Sensitivity and Specificity

MASTABLOT™ TP IgG and IgM were tested using various plasma and serum samples. The results were validated through comparisons with two comparison devices (MASTAFLUOR™ FTA-ABS IgG / IgM and SERODIA® TPPA). MASTABLOT™ TP IgG was tested with 364 samples (209 positives and 155 negatives). MASTABLOT™ TP IgM was tested with 365 samples (119 positives and 246 negatives).

The analytical performance was calculated as follows:

| | Formula | Value | |
|---------------------------------|---------------------------|--------|--------|
| | | IgG | IgM |
| Sensitivity (Sens) | $\frac{TP}{(TP + FN)}$ | > 99% | > 99% |
| Specificity (Spec) | $\frac{TN}{(TN + FP)}$ | 99% | > 99% |
| Positive predictive value (PPV) | $\frac{TP}{(TP + FP)}$ | 99% | > 99% |
| Negative predictive value (NPV) | $\frac{TN}{(TN + FN)}$ | > 99% | > 99% |
| Efficiency (Eff) | $\frac{(TP + TN)}{Total}$ | > 95% | > 95% |
| Positive likelihood ratio (LR+) | $\frac{Sens}{(1 - Spec)}$ | > 10 | > 10 |
| Negative likelihood ratio (NR+) | $\frac{(1 - Sens)}{Spec}$ | < 0.01 | < 0.01 |

Traceability

NIBSC QCRSYPHQC2 can be used to monitor the activity of MASTABLOT™ TP IgG.

The INSTAND e.V. RV 311, Pr. 62, 2019-05, is an IgM positive sample and can be used to monitor activity of MASTABLOT™ TP IgM.

Precision

No deviation was found when testing specific samples using multiple lots MASTABLOT™ TP IgG / IgM. Retesting within one lot showed identical results. Thus, precision of the assay was demonstrated.

Reporting serious incidents

All serious incidents that have occurred in connection with the device must be reported to the manufacturer and the competent authority of the Member State in which the user and / or patient is established.

References

You can find the references at the end of the instruction for use.

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Change history includes changes to significant aspects of the assay / IFU.

| Chapter | Description of change |
|-----------------------------|---|
| Warnings and Precautions | Separation into two different sections Inclusion of additional points for assay handling |
| Disposal | Separation from warnings and precautions |
| Performance Characteristics | Expansion of existing data |

**Mast Diagnostica GmbH,**

Feldstraße 20,
23858 Reinfeld,
Deutschland

Tel: +49 (0)4533 2007 0
Fax: +49 (0)4533 2007 68
email: mast@mast-diagnostica.de
Web: www.mast-group.com

Mast Group Ltd.

Mast House, Derby Road
Bootle, Merseyside, L20 1EA
United Kingdom

Tel: +44 (0)151 472 1444
Fax: +44 (0)151 944 1332
email: sales@mastgrp.com
Web: www.mast-group.com

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Icons are used according to DIN EN ISO 15223-1

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