MASTAZYME™ ANCA
PR3 (c-ANCA)
MPO (p-ANCA)

MASTAZYME™ Proteinase 3 (PR3)	REF 735013 UDI-DI 4250729700392	12 x 8 Tests
MASTAZYME™ Myeloperoxidase (MPO)	REF 735012 UDI-DI 4250729700385	12 x 8 Tests

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MPO (p-ANCA)

Intended purpose

The MASTAZYME™ PR3 / MPO is an enzyme immunoassay for the detection and semi-quantification of specific autoantibodies against proteinase 3 and myeloperoxidase respectively in human serum and plasma, used in the analysis of autoimmune diseases.

The assay is suitable for manual or automated use on open EIA processor systems and is intended for professional in vitro diagnostic use only. All laboratory test results should be interpreted in conjunction with other clinical data. The clinical judgement and further tests have to be taken into account.

Important Note for Use of these Kit Instructions

Any assay relevant changes to the kit instructions for use (IFU) will lead to a change of the version number at the bottom on the last page. All changes made will be identified on a separate sheet added to the IFU for a period of three months from the date of a change of the version. Please ensure that the latest version of the IFU is used for the assay procedure.

Principle of the Test

The ELISA test can be described in four stages:

Serum incubation

Specific antibodies bind to antigens on the solid phase to form a stable immune complex. After incubation for 30 minutes at room temperature the wells are washed to remove all non-reactive serum components.

Conjugate incubation

The anti-human-IgG horseradish peroxidase conjugate is added to all wells. The conjugate binds to IgG antibodies on the solid phase antigen to form a stable sandwich immune complex. After incubating for 30 minutes at room temperature the excess conjugate is removed by washing all wells with wash buffer.

Substrate reaction and stopping

The TMB substrate is dispensed into each well and the peroxidase enzyme/substrate reaction forms a stable blue chromogen. The reaction and consequent colour development is stopped after a 15 minute incubation at room temperature by adding 0.25 M H₂SO₄ to the wells. The shift in pH also causes the chromogen to change colour from blue to yellow.

Interpretation

The colour intensity is read in a microtiter plate reader at 450 nm (recommended reference wavelength for bichromatic measurement: 600–690 nm). The intensity of the colour (OD) is directly proportional to the

concentration of the specific antibody in the patient sample.

The results can be read from a calibration curve or with an electronic graphing package using a 4 parameter sigmoidal curve.

Kit Contents

- STRIPS** Microtiter strips
12 single strips in a holder, each with 8 break-apart wells coated with purified PR3 or MPO respectively.
- CONTROL+** Positive Control
1 x 1.5 mL stabilized human plasma, ready-to-use, can be integrated into the calibration curve, contains Proclin (< 0.1%), potentially infectious (see warnings)
- CONTROL-** Negative Control
1 x 1.5 mL stabilized human plasma, ready-to-use, can be integrated into the calibration curve, contains Proclin (< 0.1%), potentially infectious (see warnings)
- CAL 1-4** Calibrators 1-4
4 x 1.5 mL stabilized human plasma, ready-to-use, containing antibodies against the above mentioned antigens, contains Proclin (< 0.1%), potentially infectious (see warnings)

	Anti-PR3 / MPO concentration [U/mL]
Calibrator 1	1
Calibrator 2 (Cut-off)	5
Calibrator 3	25
Calibrator 4	100
Positive Control *	12.5
Negative Control *	2.5

* For batch specific concentration see QC document.

- DIL** Diluent
2 x 60 mL sample diluent, ready-to-use
- CONJ G** Enzyme conjugate
1 x 12 mL HRP-labelled goat anti-human-IgG, ready-to-use
- SUBS** TMB substrate
1 x 12 mL 3,3', 5,5'-Tetramethylbenzidine, ready-to-use
- STOP** Stopping solution
1 x 12 mL 0.25 M H₂SO₄ (sulfuric acid), ready-to-use
- WASH CONC** Wash buffer 10 x concentrated
2 x 50 mL Wash buffer solution, to be diluted 1:10 prior to use; the concentrate should be warmed up to 37 °C for 15 min to dissolve any crystals

Further Abbreviations

- RTU** ready to use

Materials Required but not Provided

1. 5 µL-, 100 µL- and 500 µL micropipettes and a multichannel pipette (optional).
2. Microtiter plate reader with a 450 nm filter (reference filter 600–690 nm).
3. Microtiter plate washer (if washing manually: wash bottle).
4. Reagent tubes for serum dilution.
5. Measuring cylinder.
6. Distilled or deionized water.

Warnings

1. Comply with the general health and safety guidelines for working with potentially infectious materials. Wear appropriate protective clothing and use appropriate lab facilities.
2. Do not cross-contaminate reagents or interchange caps on bottles. Do not use contaminated reagents. Do not use reagents with damaged bottles or packaging. Use a separate pipette or pipette tips for each sample and reagent.
3. The control/calibrator materials of human origin provided have been tested for antibodies to HIV, HCV and for HBsAg and found to be negative. However, they should be treated as potentially infectious materials and capable of transmitting diseases. No guarantee is given that the samples are free of infections or microbial contamination.
4. The kit reagents may contain bovine components. To protect against infections with prions or zoonotic pathogens, standard occupational safety regulations (GLP) must be observed.
5. Do not use haemolysed, lipemic or icteric samples (Lit. 5, 6, 7).
6. Certain reagents are irritants, especially the Stopping solution and TMB substrate. Avoid contact with skin, eyes and mucosal membranes. In case of accident, rinse with water and seek medical attention. Clean all equipment after use to avoid secondary contact incidents.
7. Proclin is used as a preservative as marked. It may be toxic if ingested. Always dispose of by flushing to drain with plenty of water.

Precautions

1. Read instructions carefully before conducting the assay. Do not modify the procedure without prior validation.
2. Do not use reagents beyond the expiry date.
3. Do not mix reagents between different lots, as reagents have been calibrated for each batch.
4. Use disposable plasticware wherever possible. Re-usable glassware should be washed thoroughly and rinsed free of detergents before use.
5. Always use a new pipette tip for each pipetting step.
6. All reagents should be brought to room temperature (18 to 24 °C) before starting the test procedure.

7. Before pipetting, all reagents should be mixed thoroughly. Vigorous shaking leading to formation of foam should be avoided.
8. It is important to pipette with constant time intervals, so that all the wells of the microtiter plate have the same reaction conditions.
9. When pipetting reagents out of the bottles care has to be taken that the stoppers are not contaminated. Furthermore, disposable pipette tips are strongly recommended to avoid cross contamination. The contents of the bottles are usually sensitive to oxidation so should be opened only for a short time.
10. Only use laboratory water of the highest quality (e.g. Aqua dest., Aqua bidest., Aqua iniectionabilia, HPLC grade).
11. In accordance with Good Laboratory Practice (GLP) all laboratory equipment should be checked regularly regarding accuracy and precision. This refers e.g. to microliter pipets and washing or ELISA reader instrumentation.

Disposal

Kit reagents, samples and contaminated disposables should be disposed of in accordance with the relevant disposal guidelines and regulations for infectious materials.

Storage and Stability

Store all reagents at 2–8 °C.

The expiry date of each reagent is printed on the individual labels. Do not use any reagents after the expiry date has been exceeded.

The diluted wash buffer is stable for up to 4 weeks when stored at 2–8 °C. However ensure that this is at room temperature before testing.

Kits should be used within three months of opening.

Sample Material

Both serum and plasma can be used for testing. Samples can be stored at 2–8 °C for up to 3 days, although for longer storage periods they should be aliquoted and kept at -20 °C.

Repeated freezing and thawing is contra-indicated.

Thawed samples should be mixed (vortex) before used in the test!

Lipemic, hemolytic, icteric or bacterially contaminated samples can cause false positive or false negative results (Lit. 5, 6, 7).

Test Procedure

1. Preparation of Reagents Samples

Allow all kit components and specimens to reach room temperature (RT, 18–24 °C) prior to use and mix well.

Sample preparation: Patient sera should be prediluted 1:101 in sample diluent (e.g. 5 µL serum + 500 µL sample diluent) prior to testing.

Washing buffer: Dissolve any crystals which may be in the bottle by warming to 37 °C and then mix well.

Dilute the concentrated washing buffer 1:10 with distilled or deionized water (e.g. 50 mL buffer concentrate + 450 mL distilled water). Mix thoroughly.

- Strictly follow the instructions. The user takes sole responsibility for any modifications to the test procedure.
- Reagents should not be left at room temperature for longer than necessary.
- A standard calibration curve should be established with each assay for quantitative results.
- Store any unused microtiter strips in the resealable aluminium bag provided and store under dry conditions at 2–8 °C.

2. Assay steps

Prepare a sufficient amount of microtiter wells for the calibrators, controls and samples.

Note: Other incubation conditions are possible. However, in case of modifications to the recommended test procedure (e.g. an incubation temperature of 37 °C instead of RT) assay performance should be verified by the user.

1. Pipette 100 µL each of the diluted (1:101) samples and the ready to use calibrators and controls into the appropriate wells.
2. Incubate the test strips at room temperature for 30 minutes.
3. Discard the contents of the microwells and wash 3 times with 300 µL of diluted wash buffer. Afterwards remove residual wash buffer by gentle tapping of the microtiter plate onto a paper towel.
4. Pipette 100 µL of enzyme conjugate solution into each well.
5. Incubate the test strips for 30 minutes at room temperature.
6. Perform the washing step according to point 3.
7. Pipette 100 µL of TMB substrate into each well.
8. Incubate for 15 minutes in the dark (e.g. drawer) at room temperature.
9. Add 100 µL of stopping solution to each well.
10. After carefully mixing read the optical density at 450 nm. Blank against air. A bichromatic measurement using a reference wavelength of 600–690 nm is recommended. Concentrations can be plotted with an electronic graphing package or by hand against the calibration curve. After adding the stopping solution the developed colour is stable for at least 30 minutes. Read optical densities during this time.

Interpretation of Results

Qualitative

Calibrator 2 is used for qualitative interpretation of the assays. A quantitative interpretation is not possible when only using the cut-off defining calibrator 2 alone.

Quantitative

The concentration in patient sample can be calculated using a computer and suitable software program. In this case the results are calculated using e.g. a 4 parameter fit algorithm to produce a sigmoidal curve.

Optional: A standard curve is plotted by entering the mean absorbance of each calibrator on the Y-axis and the corresponding concentration on the X-axis using either log/log or linear/log graph paper. The concentration of the controls and specimen can then be read directly from the graph.

Suggested normal values

In a normal range study using samples from healthy blood donors and disease-state samples the following normal and elevated ranges have been established:

Negative: < 5 U/mL

Positive: > 5 U/mL

Interpretation

Results can be affected by several patient factors. Clinical diagnosis and disease prognosis should not be based on a single test result. Other lab data and clinical observations must be taken into consideration before a conclusive diagnosis can be made.

c-ANCA: Proteinase 3

The solid phase of the assay is coated with highly purified human proteinase 3. The antigen is separated from other ANCA related serine proteases in the α-granules of neutrophil granulocytes (i. e. cathepsin G or elastase) which therefore eliminates cross-reactions due to antigen contamination. The MASTAZYME™ PR3 is ideal for confirmation or subspecification of sera which have reacted positive in a c-ANCA specific immunofluorescence assay.

Anti-PR3 titers correlate directly with the activity of Wegener's granulomatosis.

Some HIV positive sera may show unspecific cross-reaction patterns. Malignancies may also give unspecific positive results without any clinical data for vasculitis. Most of these unspecific results are limited to IFA and can be excluded by ELISA.

p-ANCA: Myeloperoxidase

The solid phase of the assay is coated with human myeloperoxidase derived from granulocytes. The use of native antigen and a special coating technique allows the antigen to maintain its stable protein conformation. For this reason no cross-reactivity with GS-ANA is observed.

Positive results are highly indicative for necrotizing forms of vasculitis, Churg-Strauss syndrome and Wegener's granulomatosis.

Thyroid peroxidase and myeloperoxidase enzymes physiologically show homologies in their amino acid sequence. No cross-reactions have been observed but

should be taken into consideration when unusual results are obtained.

Performance data

Clinical performance

To demonstrate clinical performance 44 (PR3), 27 (MPO) positive and 42 (PR3), 69 (MPO) negative clinical samples (serum and plasma) were tested. Based on these the following clinical performance parameters were calculated

	Formula	Value	
		PR3	MPO
Diagnostic Sensitivity	$\frac{TP}{(TP + FN)}$	97.7%	96.3%
Diagnostic Specificity	$\frac{TN}{(TN + FP)}$	94.2%	90.0%
Positive predictive value (PPV)	$\frac{TP}{(TP + FP)}$	93.5%	78.8%
Negative predictive value (NPV)	$\frac{TN}{(TN + FN)}$	98.0%	98.4%
Efficiency	$\frac{(TP + TN)}{total}$	95.8%	91.7%
Positive likelihood ratio (PLR)	$\frac{Sens}{(1 - Spec)}$	16.94	9.49
Negative likelihood ratio (NLR)	$\frac{(1 - Sens)}{Spec}$	0.02	0.04

(Abbreviations: TP: True positive, TN: True negative, FN: False negative, FP: False positive)

Precision

Repeatability was demonstrated on one batch with two different parameters over multiple measurement dates.

Reproduceability was demonstrated on multiple batches over multiple measurement dates.

Reporting serious incidents

All serious incidents that have occurred in connection with the device must be reported to the manufacturer and competent authority of the member state in which the user and / or patient is established.

References

You can find the references at the end of the instruction for use.

Referenzen / References:

1. Gross WL. et al., Antineutrophil cytoplasmic autoantibodies: immunobiological aspects, Klin.Wochenschr, 69, 558-566 (1991)
2. Gross WL et al., ANCA-assoziierte Vaskulitiden (Wegener-Granulomatose, Churg-Strauss-Syndrom, Mikroskopische Polyangiitis) 1.Systematik, Pathogenese und Klinik; Z.Rheumatol, 54, 279-290 (1995)
3. Seibold F et al., Clinical significance of antibodies against neutrophils in patients with inflammatory bowel disease and primary sclerosing cholangitis; Gut, 33, 657-662 (1992)
4. Seibold F et al., Neutrophil Autoantibodies: A genetic marker in primary sclerosing cholangitis and ulcerative colitis, Gastroenterol, 107, 532-536 (1994)
5. Thomas, L.: Haemolysis as influence & interference factor (2002) eJIFCC 13(4) 95-8.
6. Nikolac, N.: Lipemia: causes, interference mechanisms, detection and management (2014) Biochem 24(1): 57-67
7. Mainali, S. et al: Frequency of icteric interference in clinical chemistry laboratory tests and causes of severe icterus (2021) Pract. Lab. Med. e00259

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Change history includes changes to significant aspects of the assay / IFU.

Chapter	Description of change
Warnings and Precautions	Separation into two different sections Inclusion of additional points for assay handling
Disposal	Separation from warnings and precautions
References	Addition of References for lipemic, haemolytic and icteric samples

**Mast Diagnostica GmbH,**

Feldstraße 20,
23858 Reinfeld,
Deutschland

Tel: +49 (0)4533 2007 0

Fax: +49 (0)4533 2007 68

email: mast@mast-diagnostica.de

Web: www.mast-group.com

Mast Group Ltd.

Mast House, Derby Road
Bootle, Merseyside, L20 1EA
United Kingdom

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Icons are used according to DIN EN ISO 15223-1