MASTAZYME™ ANA Profil HJS

MASTAZYME™ ANA Profil HJS

REF 733018
UDI-DI 4250729700323

12 x 8 Tests

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	Sens	Spec	PPV	NPV	Eff.	PLR	NLR
His	96.7%	99.2%	96.7%	99.2%	98.7%	>10	<0.1
SS-A	97.6%	99.0%	98.8%	98.1%	98.4%	>10	<0.1
SS-B	91.0%	99.3%	96.8%	97.9%	97.7%	>10	<0.1
Sm	92.9%	>99%	96.3%	99.1%	98.8%	>10	0.1
Sm/RNP	96.3%	98.5%	98.1%	97.0%	97.5%	>10	<0.1
Sci-70	>99%	>99%	>99%	>99%	>99%	>10	<0.1
Jo-1	>99%	98.1%	85.7%	>99%	98.3%	>10	<0.1
CENP	90.0%	93.3%	94.6%	87.5%	91.3%	>10	0.1

(Abkürzungen: PPV: Positiver Prädiktiver Wert, NPV: Negativer Prädiktiver Wert, PLR: Positive Likelihood Ratio, NLR: Negative Likelihood ratio, His: Histon, CENP: Centromer)

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MASTAZYME™ ANA Profil HJS

Intended Purpose

The MASTAZYME™ ANA Profil HJS is an enzyme immunoassay for the qualitative detection of specific autoantibodies to Histones, SS-A, SS-B, Sm, Sm/RNP, Scl-70, Jo-1 and Centromer in human serum and plasma used in the analysis of autoimmune diseases.

The assay is suitable for manual or automated use on open EIA processor systems and is intended for professional in vitro diagnostic use only. All laboratory test results should be interpreted in conjunction with other clinical data. The clinical judgement and further tests have to be taken into account.

Important Note for Use of these Kit Instructions

Any assay relevant changes to the kit instructions for use (IFU) will lead to a change of the version number at the bottom on the last page. All changes made will be identified on a separate sheet added to the IFU for a period of three months from the date of a change of the version. Please ensure that the latest version of the IFU is used for the assay procedure.

Principle of the Test

The ELISA test can be described in four stages:

Serum incubation

Specific antibodies bind to antigens on the solid phase to form a stable immune complex. After incubation for 30 minutes at room temperature the wells are washed to remove all non-reactive serum components.

Conjugate incubation

The anti-human-IgG horseradish peroxidase conjugate is added to all wells. The conjugate binds to IgG antibodies on the solid phase antigen to form a stable sandwich immune complex. After incubating for 30 minutes at room temperature the excess conjugate is removed by washing all wells with wash buffer.

Substrate reaction and stopping

The TMB substrate is dispensed into each well and the peroxidase enzyme/substrate reaction forms a stable blue chromogen. The reaction and consequent colour development is stopped after a 15 minute incubation at room temperature by adding 0.25 M H₂SO₄ to the wells. The shift in pH also causes the chromogen to change colour from blue to yellow.

Interpretation

The colour intensity is read in a microtiter plate reader at 450 nm (recommended reference wavelength for bichromatic measurement: 600–690 nm). The intensity of the colour (OD) is directly proportional to the concentration of the specific antibody in the patient sample.

Kit Contents

- STRIPS** Microtiter strips
12 single strips in a holder, each with 8 break-apart wells coated with a purified antigen per well of: Histone, SS-A, SS-B, Sm, Sm/RNP, Scl-70, Jo-1 and centromere.
- CONTROL+** Positive Control
8 x 1.5 mL human plasma containing antibodies against the (antigen specific) above listed antigens, diluted in buffer, ready to use, contains Proclin (< 0.1%), potentially infectious (see warnings)
- CONTROL-** Negative Control
1 x 3.5 mL human plasma, diluted in buffer, ready to use, contains Proclin (< 0.1%), potentially infectious (see warnings)
- CALCO** Cut-off calibrator
8 x 1.5 mL human plasma containing antibodies against the above listed antigens, diluted in buffer, ready to use, contains Proclin (< 0.1%), potentially infectious (see warnings)
- DIL** Diluent
2 x 60 mL sample diluent, ready-to-use
- CONJ G** Enzyme conjugate
1 x 12 mL HRP-labelled goat anti-human-IgG, ready-to-use
- SUBS** TMB substrate
1 x 12 mL 3,3', 5,5'-Tetramethylbenzidine, ready-to-use
- STOP** Stopping solution
1 x 12 mL 0.25 M H₂SO₄ (sulfuric acid), ready-to-use
- WASH** **CONC** Wash buffer 10 x concentrated
2 x 50 mL Wash buffer solution, to be diluted 1:10 with distilled or deionized water prior to use; the concentrate should be warmed up to 37 °C for 15 min to dissolve any crystals.

Further Abbreviations

- RTU** ready to use

Materials Required but not Provided

- 5 µL-, 100 µL- and 500 µL micropipettes or a multichannel pipette (optional).
- Microtiter plate reader with a 450 nm filter (reference filter 600–690 nm).
- Microtiter plate washer (if washing manually: wash bottle).
- Reagent tubes for serum dilution.
- Measuring cylinder.
- Distilled or deionized water.

Warnings

1. Comply with the general health and safety guidelines for working with potentially infectious materials. Wear appropriate protective clothing and use appropriate lab facilities.
2. Do not cross-contaminate reagents or interchange caps on bottles. Do not use contaminated reagents. Do not use reagents with damaged bottles or packaging. Use a separate pipette or pipette tips for each sample and reagent.
3. The control/calibrator materials of human origin provided have been tested for antibodies to HIV, HCV and for HBsAg and found to be negative. However, they should be treated as potentially infectious materials and capable of transmitting diseases. No guarantee is given that the samples are free of infections or microbial contamination.
4. The kit reagents may contain bovine components. To protect against infections with prions or zoonotic pathogens, standard occupational safety regulations (GLP) must be observed.
5. Do not use haemolysed, lipaemic or icteric samples (Lit. 5, 6, 7).
6. Certain reagents are irritants, especially the Stopping solution and TMB substrate. Avoid contact with skin, eyes and mucosal membranes. In case of accident, rinse with water and seek medical attention. Clean all equipment after use to avoid secondary contact incidents.
7. Proclin is used as a preservative as marked. It may be toxic if ingested. Always dispose of by flushing to drain with plenty of water.

Precautions

1. Read instructions carefully before conducting the assay. Do not modify the procedure without prior validation.
2. Do not use reagents beyond the expiry date.
3. Do not mix reagents between different lots, as reagents have been calibrated for each batch.
4. Use disposable plasticware wherever possible. Re-usable glassware should be washed thoroughly and rinsed free of detergents before use.
5. Always use a new pipette tip for each pipetting step.
6. All reagents should be brought to room temperature (18 to 24 °C) before starting the test procedure.
7. Before pipetting, all reagents should be mixed thoroughly. Vigorous shaking leading to formation of foam should be avoided.
8. It is important to pipette with constant time intervals, so that all the wells of the microtiter plate have the same reaction conditions.
9. The contents of the bottles are usually sensitive to oxidation so should be opened only for a short time.
10. Only use laboratory water of the highest quality (e.g. Aqua dest., Aqua bidest., Aqua iniectionis, HPLC grade).

11. In accordance with Good Laboratory Practice (GLP) all laboratory equipment should be checked regularly regarding accuracy and precision. This refers e.g. to microliter pipets and washing or ELISA reader instrumentation.

Disposal

Kit reagents, samples and contaminated disposables should be disposed of in accordance with the relevant disposal guidelines and regulations for infectious materials.

Storage and Stability

Store all reagents at 2–8 °C.

The expiry date of each reagent is printed on the individual labels. Do not use any reagents after the expiry date has been exceeded.

The diluted wash buffer is stable for up to 4 weeks when stored at 2–8 °C. However, ensure that this is at room temperature before testing.

Kits should be used within three months of opening.

Sample Material

Both serum and plasma can be used for testing. Samples can be stored at 2–8 °C for up to 3 days, although for longer storage periods they should be aliquoted and kept at -20 °C.

Repeated freezing and thawing is contra-indicated.

Thawed samples should be mixed (vortex) before used in the test!

Lipaemic, haemolytic, icteric or bacterially contaminated samples can cause false positive or false negative results (Lit. 5, 6, 7).

Test Procedure

1. Preparation of Reagents Samples

Allow all kit components and specimens to reach room temperature (RT, 18–24 °C) prior to use and mix well.

Sample preparation: Patient sera should be prediluted 1:101 in sample diluent (e.g. 5 µL serum + 500 µL sample diluent) prior to testing.

Washing buffer: Dissolve any crystals which may be in the bottle by warming to 37 °C and then mix well.

Dilute the concentrated washing buffer 1:10 with distilled water (e.g. 50 mL buffer concentrate + 450 mL distilled water). Mix thoroughly.

- Strictly follow the instructions. The user takes sole responsibility for any modifications to the test procedure.
- Reagents should not be left at room temperature for longer than necessary.
- Store any unused microtiter strips in the resealable aluminium bag provided and store under dry conditions at 2–8 °C.

2. Assay steps

Prepare a sufficient amount of microtiter wells for the calibrators, controls and samples.

Note: Other incubation conditions are possible. However, in case of modifications to the recommended test procedure (e.g. an incubation temperature of 37 °C instead of RT) assay performance should be verified by the user.

		1	2	3	4	5	6	7	8	9	10	11	12
A	His-tone	CO	Pos	S	S								
B	SS-A	CO	Pos	S	S								
C	SS-B	CO	Pos	S	S								
D	Sm	CO	Pos	S	S								
E	Sm/RNP	CO	Pos	S	S								
F	Scl-70	CO	Pos	S	S								
G	Jo-1	CO	Pos	S	S								
H	CENP	CO	Pos	S	S								

CO: Antigen specific cut-off calibrator, Pos: positive control,
S = Sample

1. Pipette 100 µL each of the diluted (1:101) samples and the ready to use calibrators and controls into the appropriate wells.
2. Incubate the test strips at room temperature for 30 minutes.
3. Discard the contents of the microwells and wash 3 times with 300 µL of diluted wash buffer. Afterwards remove residual wash buffer by gentle tapping of the microtiter plate onto a paper towel.
4. Pipette 100 µL of enzyme conjugate solution into each well.
5. Incubate the test strips for 30 minutes at room temperature.
6. Perform the washing step according to point 3.
7. Pipette 100 µL of TMB substrate into each well.
8. Incubate for 15 minutes in the dark (e.g. drawer) at room temperature.
9. Add 100 µL of stopping solution to each well.
10. After carefully mixing read the optical density at 450 nm. Blank against air. A bichromatic measurement using a reference wavelength of 600–690 nm is recommended. Concentrations can be plotted with an electronic graphing package or by hand against the calibration curve. After adding the stopping solution the developed colour is stable for at least 30 minutes. Read optical densities during this time.

Interpretation of Results

Qualitative

The Cut-off calibrator is used for qualitative interpretation of the assays, depending on specificity requirements to include or eliminate borderline results. A quantitative interpretation is not possible when only using the Cut-off calibrator alone.

The calculated OD values for patient sera as mentioned above are compared with the value for the Cut-off calibrator. If the value of the sample is higher, then it should be read as positive, if lower, the result should be read as negative.

Suggested normal values

In a normal range study using serum samples from healthy blood donors and disease-state sera the following normal and elevated ranges have been established. It is recommended that each laboratory should establish its individual normal ranges based on results obtained from the local population:

Negative: OD-sample < OD cut-off calibrator

Positive: OD-sample > OD cut-off calibrator

Interpretation

Results can be affected by several patient factors. Clinical diagnosis and disease prognosis should not be based on a single test result. Other lab data and clinical observations must be taken into consideration before a conclusive diagnosis can be made.

The MASTAZYME™ ANA Profile HJS detects antibodies directed against the antigens histone, SS-A (Ro), SS-B (La), Sm, Sm/RNP, Scl-70, Jo-1 and/or centromere. A typing of the specific antibody is possible at the very antigen position.

References

You can find the references at the end of the instruction for use.

Performance Characteristics

Clinical performance

To demonstrate clinical performance clinical samples (serum and plasma) were tested. Based on these the following clinical performance parameters were calculated

	Formula
Diagnostic Sensitivity	$\frac{TP}{(TP + FN)}$
Diagnostic Specificity	$\frac{TN}{(TN + FP)}$
Positive predictive value (PPV)	$\frac{TP}{(TP + FP)}$
Negative predictive value (NPV)	$\frac{TN}{(TN + FN)}$
Efficiency	$\frac{(TP + TN)}{total}$
Positive likelihood ratio (PLR)	$\frac{Sens}{(1 - Spec)}$
Negative likelihood ratio (NLR)	$\frac{(1 - Sens)}{Spec}$

(Abbreviations: TP: True positive, TN: True negative, FN: False negative, FP: False positive)

	Sens	Spec	PPV	NPV	Eff.	PLR	NLR
His	96.7%	99.2%	96.7%	99.2%	98.7%	>10	<0.1
SS-A	97.6%	99.0%	98.8%	98.1%	98.4%	>10	<0.1
SS-B	91.0%	99.3%	96.8%	97.9%	97.7%	>10	<0.1
Sm	92.9%	>99%	96.3%	99.1%	98.8%	>10	0.1
Sm/RNP	96.3%	98.5%	98.1%	97.0%	97.5%	>10	<0.1
Scl-70	>99%	>99%	>99%	>99%	>99%	>10	<0.1
Jo-1	>99%	98.1%	85.7%	>99%	98.3%	>10	<0.1
CENP	90.0%	93.3%	94.6%	87.5%	91.3%	>10	0.1

(Abbreviations: PPV: positive predictive value, NPV: negative predictive value, PLR: positive likelihood ratio, NLR: negative likelihood ratio, His: Histone, CENP: Centromer)

Precision

Repeatability was demonstrated on one batch with at least two different parameters over multiple measurement dates.

Reproduceability was demonstrated on multiple batches over multiple measurement dates.

Reporting serious incidents

All serious incidents that have occurred in connection with the device must be reported to the manufacturer and competent authority of the member state in which the user and / or patient is established.

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Change history includes changes to significant aspects of the assay / IFU.

Chapter	Description of change
Warnings and Precautions	Separation into two different sections Inclusion of additional points for assay handling
Disposal	Separation from warnings and precautions
References	Addition of References for lipemic, haemolytic and icteric samples

**Mast Diagnostica GmbH,**

Feldstraße 20,
23858 Reinfeld,
Deutschland

Tel: +49 (0)4533 2007 0

Fax: +49 (0)4533 2007 68

email: mast@mast-diagnostics.de

Web: www.mast-group.com

Mast Group Ltd.

Mast House, Derby Road
Bootle, Merseyside, L20 1EA
United Kingdom

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Icons are used according to DIN EN ISO 15223-1