MASTAFLUOR™ SKMA

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1. Autoantibodies, Eds. JB Peter, Y Shoenfeld, 1996, Elsevier

Intended Use

For the detection of skeleton muscles autoantibodies on rat thigh muscle tissue sections. A sample is considered positive if a green fluorescence is visible on the respective cellular structures.

Contents

Content	50 Tests	Specifications
Slides	10 x 5 Wells	Coated with tissue section of rat thigh muscle

Stability and Storage

MASTAFUOR™ SKMA slides are stable until the end of shelf life as indicated on the label. The slides should be stored at 2–8 °C.

Serum and plasma specimens may be stored at 2–8 °C for up to 3 days prior to use or for longer storage at -20 °C. Repeated freezing and thawing of specimens should be avoided.

Warnings and Precautions

1. The reagents supplied in this kit are for *in vitro* diagnostic use only.
2. Read instructions carefully before conducting the assay. Do not modify the procedure.
3. Do not use kit beyond the expiry date.
4. Wear appropriate protective clothing and use appropriate facilities.
5. Do not mouth pipette.
6. Use disposable plastic ware where ever possible. Reusable glassware should be washed thoroughly and rinsed free of detergents before use.
7. Use distilled water or water of higher quality.
8. Do not allow wells to dry out during the assay procedures.
9. Do not expose slides to intense sunlight or similar adverse conditions while incubating.
10. Contaminated plastic ware should be disposed and incinerated. Contaminated glassware should either be autoclaved at 121 °C or decontaminated with a solution of 2.5 % sodium hypochlorite. All liquid waste should be decontaminated appropriately e.g. using sodium hypochlorite.
11. Microbial contaminated serum samples should not be used.
12. In the event that haemolysed or lipemic serum must be used, heat inactive at 56 °C before use

Test Procedure (valid only for assay specific MASTAFUOR reagents)

1. Allow all materials to reach room temperature (at least 20 °C) prior to use.
2. Reconstitute one sachet of PBS powder with 1 litre of distilled water or water of higher quality.
3. Dilute patient samples 1:20 with PBS.
Example: 10 µL serum + 190 µL of PBS.

4. Controls are supplied ready to use and must not be diluted.
5. Carefully remove the required number of slides from their sachets and mark accordingly. Hold the slide at the edge and do not touch the wells. Place the slides in a moist chamber.
6. Apply 20–25 µL of diluted specimens and ready to use controls to respective wells on the slides according to a prepared work format. Ensure that all the wells are covered and that serum does not escape from the wells.

Direct contact of a pipette tip with the slide surface may result in damage to the antigen substrate and should be avoided.

7. Cover the moist chamber and incubate at room temperature for 30 min.
8. After incubation wash slides carefully with PBS from a wash bottle, taking care not to direct the jet onto the test wells. This may be done by directing the jet of PBS along the centre of the slide, tilting the slide first towards wells 1–5 and then towards 6–10.
9. Immerse the slides in a staining dish or Coplin jar containing PBS and wash for 15 min. A change of the PBS buffer after 5 min increases washing strength.
10. Remove slides one at a time from the staining jar and drain off any excess PBS. Using a blotter, dry the area outside wells. Wipe off any excess PBS from the back of the side. Do not touch the well surface.
11. Immediately transfer the slide to a moist chamber and add 20–25 µL of FITC-conjugate to each well. Ensure that all the wells are covered by the conjugate.
12. Incubate slides at room temperature in a covered moist chamber in the dark for 30 min.
13. Switch on microscope 15 min prior to use.
14. After incubation wash the slides with PBS as in steps 9 and 10.
15. Optional counter stain: add 3–5 drops of Evans Blue to the PBS wash buffer in the last 5 min washing step. The counter stain increases the contrast of the immunofluorescent pattern.
Briefly wash the slides in PBS to remove all Evans Blue.
16. Apply one small drop of Mounting Medium to each well on a slide. Place a cover slip on the slide and read the results immediately. Take care not to trap air pockets under the cover slip.
Excess Mounting Medium on a slide may result in high background fluorescence due to light scattering. Remove excess Mounting Medium with a paper towel avoiding any direct movement of the cover slip.
17. Examine under a fluorescence microscope at a total magnification of 200x to 400x (or any other FITC suitable light source).

Interpretation of Specimen Results (valid only for assay specific MASTAFLUOR reagents)

Test Validation

The evaluation should always be performed with the positive and negative controls. Specimens are assessed by way of comparison with the controls.

Interpretation

The presence of anti-skeletal muscle antibodies is indicated by fine and continuous fluorescent cross lines vertical to the myofibrils.

This is caused by the binding of specific antibodies to the various regions of the sarcomere mainly in the regions of the I- but also of the A-band and Z-stripes.

To distinguish between the fluorescence of the I- and the A-band, a phase-contrast picture is advised.

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