MASTAFLUOR™ Combi-3

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Intended Use

MASTAFLUOR™ Combi-3 slides are coated with rat liver / kidney / stomach tissue sections for the immunofluorescence detection of specific autoantibodies in human serum or plasma.

Kit Contents

Content	100 Tests	Specifications
Slides	10 x 10 Wells	Coated with combined tissue sections of rat stomach, liver and kidney

Stability and Storage

MASTAFLUOR™ Combi-3 slides are stable until the end of shelf life as indicated on the labels. The slides should be stored at 2–8 °C.

Serum and plasma specimens may be stored at 2–8 °C for up to 3 days prior to use or for longer storage at -20 °C. Repeated freezing and thawing of specimens should be avoided.

Warnings and Precautions

1. The reagents supplied in this kit are for *in vitro* diagnostic use only.
2. Read instructions carefully before conducting the assay. Do not modify the procedure.
3. Do not use kit beyond the expiry date.
4. Wear appropriate protective clothing and use appropriate facilities.
5. Do not mouth pipette.
6. Use disposable plastic ware where ever possible. Re-usable glassware should be washed thoroughly and rinsed free of detergents before use.
7. Use distilled water or water of higher quality.
8. Do not allow wells to dry out during the assay procedures.
9. Do not expose slides to intense sunlight or similar adverse conditions while incubating.
10. Contaminated plastic ware should be disposed and incinerated. Contaminated glassware should either be autoclaved at 121 °C or decontaminated with a solution of 2.5 % sodium hypochlorite. All liquid waste should be decontaminated appropriately e.g. using sodium hypochlorite.
11. Microbial contaminated serum samples should not be used.
12. In the event that haemolysed or lipemic serum must be used, heat inactive at 56 °C before use.

Test Procedure (valid only for assay specific MASTAFLUOR reagents)

1. Allow all materials to reach room temperature (at least 20 °C) prior to use.
2. Reconstitute one sachet of PBS powder with 1 litre of distilled water or water of higher quality.
3. Dilute patient samples 1 : 20 with PBS.
Example: 10 µL serum + 190 µL of PBS.
4. Controls are supplied ready to use and must not be diluted.
5. Carefully remove the required number of slides from their sachets and mark accordingly. Hold the slide at the edge and do not touch the wells. Place the slides in a moist chamber.
6. Apply 20–25 µL of diluted specimens and ready to use controls to respective wells on the slides according to a prepared work format. Ensure that all the wells are covered and that serum does not escape from the wells.

Direct contact of a pipette tip with the slide surface may result in damage to the antigen substrate and should be avoided.
7. Cover the moist chamber and incubate at room temperature for 30 min.
8. After incubation wash slides carefully with PBS from a wash bottle, taking care not to direct the jet onto the test wells. This may be done by directing the jet of PBS along the centre of the slide, tilting the slide first towards wells 1–5 and then towards 6–10.
9. Immerse the slides in a staining dish or Coplin jar containing PBS and wash for 15 min. A change of the PBS buffer after 5 min increases washing strength.
10. Remove slides one at a time from the staining jar and drain off any excess PBS. Using a blotter, dry the area outside wells. Wipe off any excess PBS from the back of the side. Do not touch the well surface.
11. Immediately transfer the slide to a moist chamber and add 20–25 µL of FITC-conjugate to each well. Ensure that all the wells are covered by the conjugate.
12. Incubate slides at room temperature in a covered moist chamber in the dark for 30 minutes.
13. Switch on microscope 15 minutes prior to use.
14. After incubation wash the slides with PBS as in steps 9 and 10.
15. Optional counter stain: add 3–5 drops of Evans Blue to the PBS wash buffer in the last 5 min washing step. The counter stain increases the contrast of the immunofluorescent pattern.

Briefly wash the slides in PBS to remove all Evans Blue.

16. Apply one small drop of Mounting Medium to each well on a slide. Place a cover slip on the slide and read the results immediately. Take care not to trap air pockets under the cover slip.

Excess Mounting Medium on a slide may result in high background fluorescence due to light scattering. Remove excess Mounting Medium with a paper towel avoiding any direct movement of the cover slip.

17. Examine under a fluorescence microscope at a total magnification of 200x to 400x.

Interpretation of Results (valid for assay specific MASTAFLUOR reagents)

Test Validation

Specimen are assessed by way of comparison with the controls included in the test.

Additional specific controls can be ordered separately.

Interpretation of Specimen Results

Screening titre: 1:20

The evaluation should always be performed with positive and negative controls

The combined tissue section allows the differentiation of various antibodies within one test area and may thus be applied as a diagnostic test for the following autoimmune antibodies (in case of diverse antibodies it is advisable to look for further diagnostic identification).

ABBA The presence of ABBA is indicated by fluorescence of the brush border of the proximal tubules. These are heterophilic antibodies.

AMA The presence of anti-mitochondrial antibodies displays a fine granular cytoplasmatic fluorescence of the renal tubules. The distal tubules are richer in mitochondria and therefore display a more intense fluorescence in contrast to the proximal tubules.

ANA ANA-positive samples cause a fluorescence of the cell nuclei of the liver tissue. The sample should be assessed as negative if no or a faintly green staining of the nuclei is visible. Highly positive samples should be titrated to determine hidden fluorescence patterns or mixed reactions. For the differentiation of various patterns (e.g. homogenous, nucleolar, speckled, centromere) the HEp-2 cell should be used as the substrate.

ARA There is a fluorescence of the peritubular fibres (tubular capillaries). Bowman's capsule, as well as the perivascular fibres (arterioles, arteries) of the kidney.

ARIA Anti-ribosomal antibodies (ARIA) are demonstrated by a fluorescence of the gastric main cells, but not of the parietal cells.

ASMA The presence of ASMA is indicated by a fluorescence of the smooth muscle fibres of the blood vessels of kidney and stomach, of muscularis mucosa, tunica muscularis ventriculi, as well as the interglandular fibrillae of stomach mucosa.

APCA Fine granular fluorescence of the parietal cells in the gastric mucous membrane indicates APCA. Since AMA also reacts with parietal cells, anti-mitochondrial antibodies (renal tubules) should be excluded in the ACA assessment.

LKM Two different fluorescent may be seen on the combined tissue sections at the same time:

- A brilliant diffuse cytoplasmatic fluorescence of the hepatocyte cytoplasm, whereby the nucleus region remains dark.
- A comparatively weak fluorescence of the third segment of proximal kidney tubules, whereby the first and second segments display only a weak or no fluorescence.

Note: The distal tubules always appear negative which is where anti-mitochondrial antibodies more often bind.

The appropriate final serum dilution is defined by a positive fluorescence. Weak fluorescence of the cell nuclei with titres between 1:20 and 1:80 or vagueness with respect to the clinical results should be checked by way of monitoring control. In such a case the samples should be collected about every 3 weeks and similarly tested.

Diagnostic Relevance

ABBA

> 1:20 HA generate patterns of fluorescence which may be mistaken for genuine antibodies and though produce false positive results.

AMA

1:20–1:80 A positive reaction is found in several liver diseases.

> 1:160 Indicates biliary cirrhosis. AMA titres remain constant over a long period of time and therapy so that the determination of titre as a measure of therapy control is not useful.

ANA

1:20–1:80 A positive reaction is seen frequently in rheumatoid arthritis (RA) and MCTD. An increase in titre may indicate active SLE. Constant titres are observed in chronic or treated autoimmune diseases.

> 1:160 An acute autoimmune disease is indicated, whereby in 80 % of cases SLE is present.

ARA

> 1:20 ARA are found in Crohn's disease and in children with untreated celiac disease.

ARIA

> 1:20 Occasionally found in SLE patients with renal disorder.

ASMA

1:20–1:80 A positive reaction is found in several liver diseases, viral hepatitis and primary biliary cirrhosis. However, the titres here may fall below the determination border. Low titres may be observed in patients with gallbladder infections, alcoholic cirrhosis, SLE and in 2 % of the normal healthy population.

> 1: 160 Chronically active hepatitis is indicated. In contrast to viral hepatitis the titres fall only slightly and may persist for several years. Patients with infectious mononucleosis may also show high ASMA titres.

APCA

The APCA titre provides no information about the diseases state of the patient. The antibody determination should be evaluated together with the measurement of Intrinsic factor and/or results from histopathology.

LKM

> 1:20 Type II LKM antibodies may be observed in drug-induced hepatitis.

Negative specimen should show no green fluorescence; a green staining without fluorescence may be visible.

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