# MASTAFLUOR™ Combi-3

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## Introduction

Autoimmune diseases are triggered by autoimmune reactions, which are characterised by auto-reactive B and T cells. Autoimmunity is the cause of a loss of self-tolerance of the immune systems. A number of plausible hypotheses are discussed, which also included genetic and/or environmental influences.

**ABBA** Heterophilic antibodies (HA) are antibodies which react with animal tissue, although they are not active against specific tissue antigens. These generate patterns of fluorescence which may be mistaken for genuine antibodies. They are probably active against blood group antigens and thus produce false positive results.

**AMA** Anti-mitochondrial antibodies (AMA) predominantly react with the inner membrane of the mitochondria (rich in phospholipids). AMA mostly appear with diseases such as primary biliary cirrhosis, pseudo-LE syndrome and various forms of chronic aggressive hepatitis. High titre AMA results are mainly found with non-suppurating gallbladder infections or primary biliary cirrhosis (positive results at about 90 %). In these cases antibodies appear before the clinical symptoms and will hardly be influenced by therapy during the course of the disease. Low antibody titres are observed with scleroderma, Sjögren syndrome, rheumatoid arthritis and other autoimmune diseases.

**ANA** The term anti-nuclear antibodies (ANA) comprises all antibodies which react with antigens in the cell nucleus. They are important markers in the diagnosis of the systemic lupus erythematosus (SLE), scleroderma, Sjögren syndrome or various connective tissue diseases, commonly designated as MCTD. Although the ANA detection for these diseases proves highly important, the final diagnosis should be checked against the clinical results. Anti-nuclear antibodies may also occur with non autoimmune based diseases, especially in the case of infections, viral diseases, hepatitis, infectious mononucleosis, leukaemia, lymphoma, melanomas and others. Furthermore, ANA are often detected with chronic hepatitis, primary cirrhosis, thyroiditis and allergic encephalitis. Low ANA titres may occasionally be caused by medication or may be present in healthy people.

**ARA** Anti-reticulon antibodies (ARA) are found in Crohn's disease and in children with untreated celiac disease.

**ARIA** Anti-ribosomal antibodies (ARIA) are occasionally found in SLE patients where renal disorder is involved.

**ASMA** Antibodies against smooth, unstriated muscle occur in various liver diseases, for example acute and

chronic hepatitis, primary biliary cirrhosis and other forms of liver cirrhosis. Furthermore, the detection of ASMA supports the diagnosis of SLE, infectious mononucleosis, breast and ovarian carcinoma and malignant melanomas.

**APCA** Circulating antibodies against the structures of the parietal cell of the gastric mucosa are generally due to pernicious anemia.

They may, however, also be detected with other diseases of the stomach (chronic atrophic gastritis, gastric ulcer), diseases of the thyroid (Hashimoto's thyroiditis, myxedema), and more rarely with hypoferric anemia, diabetes mellitus and in older patients.

**LKM** Liver-kidney-microsomal antibodies (LKM) are predominantly directed at lipoproteins located in the membranes of rough and smooth endoplasmic reticulum.

Two types of LKM (Type I + II) may be differentiated:

- LKM type 1 antibodies are a subgroup of chronic active autoimmune hepatitis which are found mainly in children with being more present in young girls. This type of hepatitis is caused by hepatitis A or other hepatitis virus infection, but never by hepatitis B virus infection.
- LKM type 2 antibodies are found in cases of drug-induced hepatitis.

## Intended Use

MASTFLUOR™ Combi-3 provide an indirect immunofluorescence test system for the detection of specific antibodies on stomach / liver / kidney tissue sections (rat).

## Test Principle

IFA slides are coated with combined sections of kidney stomach tissues. Diluted patient serum or plasma is incubated on the wells of the slide. In the presence of specific antibodies a stable antigen-antibody complex is formed. Unspecific or unbound antibodies are removed in a washing step. Specific immune complexes are then detected by a polyvalent FITC-conjugated anti-human-antibody. Again unspecific or non-bound antibodies are removed in a washing step. Results could be read visually using an immunofluorescence microscope.

Positive reactions show specific fluorescent reactions on stomach, liver and kidney tissue.

## Kit Contents

Content	100 Tests	Specifications
Slides	10 x 10 Wells	Coated with combined tissue section of stomach, liver and kidney sections (rat tissue)
Positive control	0,5 mL	Positive serum, ready to use, contains < 0.1 % NaN <sub>3</sub>
Negative control	0,5 mL	Negative serum, ready to use, contains < 0.1 % NaN <sub>3</sub>
FITC-conjugate	3 mL	Anti-human-globulin polyvalent, FITC-conjugate, ready to use, contains < 0.1 % NaN <sub>3</sub>
Evans Blue	3 mL	Ready to use, contains < 0.1 % NaN <sub>3</sub>
Mounting Medium	3 mL	Ready to use, contains < 0.1 % NaN <sub>3</sub>
PBS	2 Sachets	1 sachet to be diluted in 1 L of distilled water, pH 7.2 ± 0.2

## Materials Required but not Provided

1. Sterile test tubes
2. Micropipettes and tips
3. Staining dish or Coplin jar
4. Moist chamber
5. Volumetric flask for PBS
6. Distilled water or water of higher quality
7. Forceps
8. Transmitted or epi-fluorescence microscope with a 490 nm excitation filter and a 510 nm barrier filter
9. Cover slips
10. Wash bottle

## Stability and Storage

MASTAFLUOR™ Combi-3 is stable until the end of shelf life as indicated on kit label. The kit should be stored at 2–8 °C.

Reconstituted PBS powder should be stored at 2–8 °C for up to 30 days.

Serum and plasma specimens may be stored at 2–8 °C for up to 3 days prior to use or for longer storage at -20 °C. Repeated freezing and thawing of specimens should be avoided.

## Warnings and Precautions

1. The reagents supplied in this kit are for *in vitro* diagnostic use only.
2. Read instructions carefully before conducting the assay. Do not modify the procedure.
3. Do not use kit beyond the expiry date.
4. Wear appropriate protective clothing and use appropriate facilities.
5. Do not mouth pipette.
6. Use disposable plasticware where ever possible. Reusable glassware should be washed thoroughly and rinsed free of detergents before use.
7. Sera provided have been tested for antibodies to HIV and for HBsAg and found to be negative. However, they should be treated as potentially infectious materials and capable of transmitting diseases. No guarantee is given that the sera are free of infection or microbial contamination.
8. Use distilled water or water of higher quality.
9. Do not mix reagents between different lots, as reagents have been calibrated for each batch.
10. Do not cross-contaminate reagents or interchange caps on bottles. Use a separate pipette or pipette tips for each sample and reagent.
11. Do not allow wells to dry out during the assay procedures.
12. Do not expose slides to intense sunlight or similar adverse conditions while incubating.
13. Contaminated plastic ware should be disposed and incinerated. Contaminated glassware should either be autoclaved at 121 °C or decontaminated with a solution of 2.5 % sodium hypochlorite. All liquid waste should be decontaminated appropriately e.g. using sodium hypochlorite.

14. Do not expose the FITC-conjugate at any stage to strong sunlight, UV or fluorescent light. Keep in a dark place whenever possible.
15. Microbial contaminated serum samples should not be used.
16. In the event that haemolysed or lipemic serum must be used, heat inactive at 56 °C before use.
17. Sodium azide is used as a preservative as marked. It may be toxic if ingested. Sodium azide may react with lead and copper plumbing to form highly explosive salts. Always dispose by flushing to drain with plenty of water.

### Test Procedure

1. Allow all materials to reach room temperature (at least 20 °C) prior to use.
2. Reconstitute one sachet of PBS powder with 1 litre of distilled water or water of higher quality.
3. Dilute patient samples 1 : 20 with PBS.  
Example: 10 µL serum + 190 µL of PBS.
4. Controls are supplied ready to use and must not be diluted.
5. Carefully remove the required number of slides from their sachets and mark accordingly. Hold the slide at the edge and do not touch the wells. Place the slides in a moist chamber.
6. Apply 20–25 µL of diluted specimens and ready to use controls to respective wells on the slides according to a prepared work format. Ensure that all the wells are covered and that serum does not escape from the wells.

Direct contact of a pipette tip with the slide surface may result in damage to the antigen substrate and should be avoided.

7. Cover the moist chamber and incubate at room temperature for 30 min.
8. After incubation wash slides carefully with PBS from a wash bottle, taking care not to direct the jet onto the test wells. This may be done by directing the jet of PBS along the centre of the slide, tilting the slide first towards wells 1–5 and then towards 6–10.
9. Immerse the slides in a staining dish or Coplin jar containing PBS and wash for 15 min. A change of the PBS buffer after 5 min increases washing strength.
10. Remove slides one at a time from the staining jar and drain off any excess PBS. Using a blotter, dry the area outside wells. Wipe off any excess PBS from the back of the side. Do not touch the well surface.

11. Immediately transfer the slide to a moist chamber and add 20–25 µL of FITC-conjugate to each well. Ensure that all the wells are covered by the conjugate.
12. Incubate slides at room temperature in a covered moist chamber in the dark for 30 minutes.
13. Switch on microscope 15 minutes prior to use.
14. After incubation wash the slides with PBS as in steps 9 and 10.
15. Optional counter stain: add 3–5 drops of Evans Blue to the PBS wash buffer in the last 5 min washing step. The counter stain increases the contrast of the immunofluorescent pattern.  
  
Briefly wash the slides in PBS to remove all Evans Blue.
16. Apply one small drop of Mounting Medium to each well on a slide. Place a cover slip on the slide and read the results immediately. Take care not to trap air pockets under the cover slip.

Excess Mounting Medium on a slide may result in high background fluorescence due to light scattering. Remove excess Mounting Medium with a paper towel avoiding any direct movement of the cover slip.

17. Examine under a fluorescence microscope at a total magnification of 200x to 400x.

### Interpretation of Results

#### Test Validation

Specimen are assessed by way of comparison with the controls included in the test.

Additional specific controls can be ordered separately.

## Interpretation of Specimen Results

Screening titre: 1:20

The evaluation should always be performed with positive and negative controls

The combined tissue section allows the differentiation of various antibodies within one test area and may thus be applied as a diagnostic test for the following autoimmune antibodies (in case of diverse antibodies it is advisable to look for further diagnostic identification).

**ABBA** The presence of ABBA is indicated by fluorescence of the brush border of the proximal tubules. These are heterophilic antibodies.

**AMA** The presence of anti-mitochondrial antibodies displays a fine granular cytoplasmatic fluorescence of the renal tubules. The distal tubules are richer in mitochondria and therefore display a more intense fluorescence in contrast to the proximal tubules.

**ANA** ANA-positive samples cause a fluorescence of the cell nuclei of the liver tissue. The sample should be assessed as negative if no or a faintly green staining of the nuclei is visible. Highly positive samples should be titrated to determine hidden fluorescence patterns or mixed reactions. For the differentiation of various patterns (e.g. homogenous, nucleolar, speckled, centromere) the HEp-2 cell should be used as the substrate.

**ARA** There is a fluorescence of the peritubular fibres (tubular capillaries). Bowman's capsule, as well as the perivascular fibres (arterioles, arteries) of the kidney.

**ARIA** Anti-ribosomal antibodies (ARIA) are demonstrated by a fluorescence of the gastric main cells, but not of the parietal cells.

**ASMA** The presence of ASMA is indicated by a fluorescence of the smooth muscle fibres of the blood vessels of kidney and stomach, of muscularis mucosa, tunica muscularis ventriculi, as well as the interglandular fibrillae of stomach mucosa.

**APCA** Fine granular fluorescence of the parietal cells in the gastric mucous membrane indicates APCA. Since AMA also reacts with parietal cells, anti-mitochondrial antibodies (renal tubules) should be excluded in the ACA assessment.

**LKM** Two different fluorescent may be seen on the combined tissue sections at the same time:

- A brilliant diffuse cytoplasmatic fluorescence of the hepatocyte cytoplasm, whereby the nucleus region remains dark.

- A comparatively weak fluorescence of the third segment of proximal kidney tubules, whereby the first and second segments display only a weak or no fluorescence.

Note: The distal tubules always appear negative which is where anti-mitochondrial antibodies more often bind.

The appropriate final serum dilution is defined by a positive fluorescence. Weak fluorescence of the cell nuclei with titres between 1:20 and 1:80 or vagueness with respect to the clinical results should be checked by way of monitoring control. In such a case the samples should be collected about every 3 weeks and similarly tested.

## Diagnostic Relevance

### ABBA

> 1:20 HA generate patterns of fluorescence which may be mistaken for genuine antibodies and though produce false positive results.

### AMA

1:20–1:80 A positive reaction is found in several liver diseases.

> 1:160 Indicates biliary cirrhosis. AMA titres remain constant over a long period of time and therapy so that the determination of titre as a measure of therapy control is not useful.

### ANA

1:20–1:80 A positive reaction is seen frequently in rheumatoid arthritis (RA) and MCTD. An increase in titre may indicate active SLE. Constant titres are observed in chronic or treated autoimmune diseases.

> 1:160 An acute autoimmune disease is indicated, whereby in 80 % of cases SLE is present.

### ARA

> 1:20 ARA are found in Crohn's disease and in children with untreated celiac disease.

### ARIA

> 1:20 Occasionally found in SLE patients with renal disorder.

### ASMA

1:20–1:80 A positive reaction is found in several liver diseases, viral hepatitis and primary biliary cirrhosis. However, the titres here may fall below the determination border. Low titres may be observed in patients with gallbladder infections,

alcoholic cirrhosis, SLE and in 2 % of the normal healthy population.

> 1: 160 Chronically active hepatitis is indicated. In contrast to viral hepatitis the titres fall only slightly and may persist for several years. Patients with infectious mononucleosis may also show high ASMA titres.

**APCA** The APCA titre provides no information about the diseases state of the patient. The antibody determination should be evaluated together with the measurement of Intrinsic factor and/or results from histopathology.

#### **LKM**

> 1:20 Type II LKM antibodies may be observed in drug-induced hepatitis.

**Negative specimen** should show no green fluorescence; a green staining without fluorescence may be visible.

## **Limitations / Cross-reactions**

1. No cross-reactivity is known.
2. MASTAFLUOR™ Combi-3 has been demonstrated to be highly sensitive and specific, however, results should be considered along with all serological tests, clinical history and other aspects of patient management to be considered diagnostically significant.

## **Performance**

Inter-assay and intra-assay variances are not detectable in respect to control fluorescence intensity.

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