MASTAFLUOR™ MUMPS

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Literatur

- Bayas JM; Vilella A; Vidal J; Nebot X; Carbo JM; Navarro G. Prat A; Asenjo MA; Salleras L: Susceptibility to measles, rubella and parotitis in young adults. Med Clin (Barc), 106 (15): 561–4 (1996).
- Gut JP; Lablache C; Behr S; Kirn A. Symptomatic mumps virus reinfections. J Med Virol, 45 (1): 17–23 (1995).
- Johnson CE; Kumar ML; Whitwell JK; Staehle BO; Rome LP; Dinakar C; Hurni W; Nalin DR. Antibody persistence after primary measles-mumps-rubella vaccine and response to a second dose given at four to six vs. eleven to thirteen years. Pediatr Infect Dis J, 15 (8): 687–92 (1996).
- King SM; Saunders EF; Petric M; Gold R. Response to measles, mumps and rubella vaccine in paediatric bone marrow transplant recipients. Bone Marrow transplant, 17 (4): 633–6 (1996)
- Matter L; Germann D; Bally F; Schopfer K. Age-stratified seroprevalence of measles, mumps and rubella (MMR) virus infections in Switzerland after the introduction of MMR mass vaccination. Eur J Epidemiol, 13 (1): 61–66, (1997)

Introduction

Mumps (Parotitis) is a common contagious disease with relatively moderate symptoms during childhood, but increasing complications, when adults are infected. The causative agent of mumps is a virus of the Paramyxoviridae family. The virus normally infects children at the age of 4 to 10. The infection is mainly transmitted by the airborne route, but is also spread by various objects contaminated by the patient saliva. The disease shows a seasonal prevalence with the greatest incidence in winter and spring. Both a mumps infection and vaccination lead to a persistent immunity.

The typical symptom associated with a mumps infection is a "parotitis" (swelling of the parotid glands). Additionally pathological involvement of the CNS and various glandular organs (pancreas, thymus, thyroids) is a typical feature of mumps. Mumps induced meningitis is one of the most frequent manifestations of the disease which can also appear without parotitis. The published data about the presence of meningitis vary between 1.4 % and 66 % of the clinically ill patients. In most cases convalescence follows without complications.

The recognition of Mumps in the laboratory is mostly done by detection of the infectious agent itself or by the determination of virus-specific antibodies. The serodiagnosis plays a major part. Besides the classical methods like complement fixation, haemagglutination and neutralisation tests have been introduced a series of modern assays like immunofluorescence, radio immunoassay and enzyme immunoassay. Increasing IgG titers are helpful in the determination of the causative agent because cross reactions with Parainfluenza-Virus type 2 may lead to a wrong interpretation of the patient reports. In the early stage of the disease, significant IgM titers may be monitored.

Intended Use

MASTAFLUOR™ MUMPS provide an indirect immunofluorescence test system for the detection of specific antibodies to Mumps virus antigen.

Test Principle

IFA slides are coated with BMK cells infected with Mumps virus. Diluted patient serum or plasma is incubated on the wells of the slide. In the presence of specific antibodies a stable antigen-antibody complex is formed. Unspecific or un-bound antibodies are removed in a washing step. Specific immune complexes are then detected by either a FITC-conjugated anti human IgG or IgM antibody. Again unspecific or non-bound antibodies are removed in a washing step. Results could be read visually using an immunofluorescence microscope. Positive reactions show specific cytoplasmic inclusions in infected cells.

Kit Contents

Content	100 Tests	50 Tests	Specifications
Slides	10 x 10 Wells	10 x 5 Wells	Coated with Mumps virus infected BMK cells
Positive control	0.5 mL	0.5 mL	IgG or IgM positive serum, ready to use, contains < 0.1 % NaN ₃
Negative control	0.5 mL	0.5 mL	IgG and IgM negative serum, ready to use, contains < 0.1 % NaN ₃
FITC-conjugate	3 mL	2 mL	Anti-human-IgG (γ-chain) FITC-conjugate or Anti-human-IgM (μ-chain) FITC-conjugate; ready to use, contains < 0.1 % NaN ₃
Evans Blue	3 mL	3 mL	ready to use, contains < 0.1 % NaN ₃
Mounting medium	3 mL	3 mL	ready to use, contains < 0.1 % NaN ₃
PBS	2 Sachets	2 Sachets	1 sachet to be diluted in 1 L of distilled water, pH 7.2 ± 0.2

Materials Required but not Provided

1. Sterile test tubes
2. Micropipettes and tips
3. Staining dish or Coplin jar
4. Moist chamber
5. Volumetric flask for PBS
6. Distilled water or water of higher quality
7. Forceps
8. Transmitted or epi-fluorescent microscope with a 490 nm excitation filter and a 510 nm barrier filter
9. Cover slips
10. Wash bottle

Stability and Storage

MASTAFLUOR™ MUMPS is stable until end of shelf life as indicated on kit label. The kit should be stored at 2–8 °C.

Reconstituted PBS powder should be stored at 2–8 °C for up to 30 days.

Serum and plasma specimens may be stored at 2–8 °C for up to 48 hours prior to use or for longer storage at -20 °C. Repeated freezing and thawing of specimens should be avoided.

Warnings and Precautions

1. The reagents supplied in this kit are for *in vitro* diagnostic use only.
2. Read instructions carefully before conducting the assay. Do not modify the procedure.
3. Do not use kit beyond the expiry date.
4. Wear appropriate protective clothing and use appropriate facilities.
5. Do not mouth pipette.
6. Use disposable plasticware where ever possible. Re-usable glassware should be washed thoroughly and rinsed free of detergents before use.
7. Sera provided have been tested for antibodies to HIV and for HBsAg and found to be negative. However, they should be treated as potentially infectious materials and capable of transmitting disease. No guarantee is given that the sera are free of infection or microbial contamination.
8. Use distilled water or water of higher quality.
9. Do not mix reagents between different lots, as reagents have been calibrated for each batch.
10. Do not cross-contaminate reagents or interchange caps on bottles. Use a separate pipette or pipette tips for each sample and reagent.
11. Do not allow wells to dry out during the assay procedures.
12. Do not expose slides to intense sunlight or similar adverse conditions while incubating.
13. Contaminated plastic ware should be disposed of and incinerated. Contaminated glassware should either be autoclaved at 121 °C for 1 hour or decontaminated with a solution of 2.5 % sodium hypochlorite. All liquid waste should be decontaminated appropriately e.g. using sodium hypochlorite.
14. Do not expose the FITC-conjugate at any stage to strong sunlight, UV or fluorescent light. Keep in a dark place whenever possible.
15. Microbial contaminated serum samples should not be used.
16. In the event that haemolysed or lipemic serum must be used, heat inactive at 56 °C before use.
17. Sodium azide is used as a preservative as marked. It may be toxic if ingested. Sodium azide may react with lead and copper plumbing to form highly explosive salts. Always dispose of by flushing to drain with plenty of water.

Test Procedure

1. Allow all materials to reach room temperature (at least 20 °C) prior to use.
2. Reconstitute one sachet of PBS powder with 1 litre of distilled water or water of higher quality.
3. Dilute patient samples 1 : 20 with PBS for IgG and IgM testing: Example: 10 µL serum + 190 µL of PBS.

For IgM determination a preabsorption with Rf absorbant (MASTSORB Order Code: 651003) is recommended to reduce interferences e.g. with rheumatic factors. **Do not absorb the controls.**
4. Controls are supplied ready to use and must not be diluted.
5. Carefully remove the required number of slides from their sachets and mark accordingly. Hold the slide at the edge and do not touch the wells. Place the slides in a moist chamber.
6. Apply 20–25 µL of diluted specimens and ready to use controls to respective wells on the slides according to a prepared work format. Ensure that all the wells are covered and that serum does not escape from the wells.

Direct contact of a pipette tip with the slide surface may result in damage to the antigen substrate and should be avoided.
7. Cover the moist chamber and incubate as follows:
IgG assay: room temperature for 30 min
IgM assay: room temperature for 45 min
8. After incubation wash slides carefully with PBS from a wash bottle, taking care not to direct the jet onto the test wells. This may be done by directing the jet of PBS along the centre of the slide, tilting the slide first towards wells 1–5 and then towards 6–10.
9. Immerse the slides in a staining dish or Coplin jar containing PBS and wash for 15 min with a change of PBS after 5 min.
10. Remove slides one at a time from the staining jar and drain off any excess PBS. Using a blotter, dry the area outside wells. Wipe off any excess PBS from the back of the side. Do not touch the well surface.
11. Immediately transfer the slide to a moist chamber and add 20–25 µL of FITC-conjugate (either IgG or IgM) to each well. Ensure that all the wells are covered by the conjugate.
12. Incubate slides at room temperature in a covered moist chamber in the dark for 30 min (same incubation for IgG and IgM).
13. Switch on microscope 15 min prior to use.
14. After incubation wash the slides with PBS as in steps 9 and 10.
15. Optional counter stain: add 3–5 drops of Evans Blue to the PBS wash buffer in the last 5 min washing step. The counter stain increases the contrast of the immunofluorescent pattern.

Briefly wash the slides in PBS to remove all Evans Blue.

16. Apply one small drop of Mounting Medium to each well on a slide. Place a cover slip on the slide and read the results immediately. Take care not to trap air pockets under the cover slip. Excess Mounting Medium on a slide may result in high background fluorescence due to light scattering. Remove excess Mounting Medium with a paper towel avoiding any direct movement of the cover slip.
17. Examine under a fluorescence microscope at a total magnification of 400x to 800x.

Interpretation of Results

Test Validation

Examine the control samples before reading patient samples. The controls must give the following patterns to be valid.

Sample	Fluorescence Intensity
Positive control IgG	3 + to 4 +
Positive control IgM	1 + to 3 +
Negative control	negative

Specimen sera are assessed by way of comparison with the controls included in the test.

Interpretation of Specimen Results

IgG positive specimen should show a green fluorescence of several coarse or fine speckled inclusions ONLY in the cytoplasm of infected cells. Nuclear fluorescence is not considered as positive.

IgG negative specimen should show no green fluorescence; a green staining without fluorescence may be visible.

IgM positive specimen should show a green fluorescence of several coarse or fine speckled inclusions in the cytoplasm of infected cells. A diffuse fluorescence in the cytoplasm of infected and non-infected cells may be found. The intensity of the fluorescence signal is dependent on the antibody concentration.

IgM negative specimen do not show any green fluorescence.

Limitations / Cross reactions

1. Mumps viruses and viruses of the parainfluenza group do show antigen homologies. Especially with parainfluenza type 2 antibody cross-reactions are found. Positive and borderline results should therefore be re-tested with a fresh serum or plasma sample.
2. MASTAFLUOR™ MUMPS has been demonstrated to be highly sensitive and specific, however, results should be considered along with all serological tests, clinical history and other aspects of patient management to be considered diagnostically significant.

Performance

Data in the table were obtained by comparing MASTAFLUOR™ MUMPS results with data or other assays based on a different method (e.g. ELISA):

	IgG	IgM
Sensitivity	94.6 %	50 %
Specificity	86.2 %	98.9 %

Inter-assay and intra-assay variances are not detectable in respect to control fluorescence intensity.

References

1. A; Vidal J; Nebot X; Carbo JM; Navarro G. Prat A; Asenjo MA; Salleras L: Susceptibility to measles, rubella-Med Clin (Barc), 106 (15): 561-4 (1996).
2. Gut JP; Lablache C; Behr S; Kirn A. Symptomatic mumps virus reinfections. J Med Virol, 45 (1): 17-23 (1995).
3. Johnson CE; Kumar ML; Whitwell JK; Staehle BO; Rome LP; Dinakar C; Hurni W; Nalin DR. Antibody persistence after primary measles-mumps-rubella vaccine and response to a second dose given at four to six vs. eleven to thirteen years. Pediatr Infect Dis J, 15 (8): 687-92 (1996).
4. King SM; Saunders EF; Petric M; Gold R. Response to measles, mumps and rubella vaccine in paediatric bone marrow transplant recipients. Bone Marrow Transplant, 17 (4): 633-6 (1996)
5. Matter L; Germann D; Bally F; Schopfer K. Age-stratified seroprevalence of measles, mumps and rubella (MMR) virus infections in Switzerland after the introduction of MMR mass vaccination. Eur J Epidemiol, 13 (1): 61-66, (1997)
6. Pison Garces FJ; Galbe Sanchez-Ventura J; Arcauz Eguren P; Aguirre y Daban C; Mengual Gil J; Larrad Mur L. Immunity to measles, mumps and rubella in children vaccinated with triple viral vaccine. Aten Primaria, 15 (4): 235-7 (1995)

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