MASTAFLUOR™ VZV

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Introduction

Varicella is an acute viral disease caused by the varicella-zoster virus (VZV) that belongs to the herpes virus family.

The disease begins with fever and malaise and a maculopapular rash that changes within hours to vesicles that remain during 3 to 4 days. Lesions pustulate, crust and heal slowly over a 3 week period. Lesions appear in the scalp, face, neck, trunk, mucous membranes of the mouth and upper respiratory tract. Sometimes, especially in adults, the fever and general clinical status can be severe but the disease is rarely fatal. The most common cause of death in adults is the primary viral pneumonia and in children secondary bacterial infections or central nervous system involvement. Children with acute leukemia or immunocompromised persons are at higher risk of a general dissemination of the disease with fatal outcome. Newborns that are infected after 5 to 10 days and the ones born from mothers that were infected 5 days before or two days after birth can develop severe generalised Varicella with a mortality rate of 30 %. Infection at the beginning of pregnancy rarely gives congenital malformations.

Reactivation of VZV is usually manifested as zoster. Zoster normally begins with unilateral sharp and well localised pain. In this case the titer of the IgG antibodies is increased, whereas the titer of IgM antibodies is raised only in some few cases.

Intended Use

MASTAFUOR™ VZV provide an indirect immunofluorescence test system for the detection of specific antibodies to Mumps virus antigen.

Test Principle

IFA slides are coated with MRC-5 human embryonal lung fibroblasts infected with varicella-zoster virus. Diluted patient serum or plasma is incubated on the wells of the slide. In the presence of specific antibodies a stable antigen-antibody complex is formed. Unspecific or un-bound antibodies are removed in a washing step. Specific immune complexes are then detected by either a FITC-conjugated anti human IgG or IgM antibody. Again unspecific or non-bound antibodies are removed in a washing step. Results could be read visually using an immunofluorescence microscope. Positive reactions show specific cytoplasmic and / or nuclear fluorescence in infected cells.

Kit Contents

Materials Required but not Provided

1. Sterile test tubes
2. Micropipettes and tips
3. Staining dish or Coplin jar

4. Moist chamber
5. Volumetric flask for PBS

Content	100 Tests	50 Tests	Specifications
Slides	10 x 10 Wells	10 x 5 Wells	Coated with VZV infected MRC-5 human embryonal lung fibroblasts
Positive control	0.5 mL	0.5 mL	IgG or IgM positive serum, ready to use, contains < 0.1 % NaN ₃
Negative control	0.5 mL	0.5 mL	IgG and IgM negative serum, ready to use, contains < 0.1 % NaN ₃
FITC-conjugate	3 mL	2 mL	Anti-human-IgG (γ-chain) FITC-conjugate or Anti-human-IgM (μ-chain) FITC-conjugate; ready to use, contains < 0.1 % NaN ₃
Evans Blue	3 mL	3 mL	Ready to use, contains < 0.1 % NaN ₃
Mounting-medium	3 mL	3 mL	Ready to use, contains < 0.1 % NaN ₃
PBS	2 Sachets	2 Sachets	1 sachet to be diluted in 1 L of distilled water, pH 7.2 ± 0.2

6. Distilled water or water of higher quality
7. Forceps
8. Transmitted or epi-fluorescent microscope with a 490 nm excitation filter and a 510 nm barrier filter
9. Cover slips
10. Wash bottle

Stability and Storage

MASTAFUOR™ VZV is stable until end of shelf life as indicated on kit label. The kit should be stored at 2–8 °C .

Reconstituted PBS powder should be stored at 2–8 °C for up to 30 days.

Serum and plasma specimens may be stored at 2–8 °C for up to 48 hours prior to use or for longer storage at -20 °C. Repeated freezing and thawing of specimens should be avoided.

Warnings and Precautions

1. The reagents supplied in this kit are for *in vitro* diagnostic use only.
2. Read instructions carefully before conducting the assay. Do not modify the procedure.
3. Do not use kit beyond the expiry date.
4. Wear appropriate protective clothing and use appropriate facilities.
5. Do not mouth pipette.
6. Use disposable plasticware where ever possible. Re-usable glassware should be washed thoroughly and rinsed free of detergents before use.
7. Sera provided have been tested for antibodies to HIV and for HBsAg and found to be negative. However, they should be treated as potentially infectious materials and capable of transmitting disease. No guarantee is given that the sera are free of infection or microbial contamination.
8. Use distilled water or water of higher quality.
9. Do not mix reagents between different lots, as reagents have been calibrated for each batch.
10. Do not cross-contaminate reagents or interchange caps on bottles. Use a separate pipette or pipette tips for each sample and reagent.
11. Do not allow wells to dry out during the assay procedures.
12. Do not expose slides to intense sunlight or similar adverse conditions while incubating.
13. Contaminated plastic ware should be disposed of and incinerated. Contaminated glassware should either be autoclaved at 121 °C for 1 hour or decontaminated with a solution of 2.5 % sodium hypochlorite. All liquid waste should be decontaminated appropriately e.g. using sodium hypochlorite.
14. Do not expose the FITC-conjugate at any stage to strong sunlight, UV or fluorescent light. Keep in a dark place whenever possible.
15. Microbial contaminated serum samples should not be used.
16. In the event that haemolysed or lipemic serum must be used, heat inactive at 56 °C before use.
17. Sodium azide is used as a preservative as marked. It may be toxic if ingested. Sodium azide may react with lead and copper plumbing to form highly explosive salts. Always dispose of by flushing to drain with plenty of water.

Test Procedure

1. Allow all materials to reach room temperature (at least 20 °C) prior to use.
2. Reconstitute one sachet of PBS powder with 1 litre of distilled water or water of higher quality.
3. Dilute patient samples 1:20 with PBS for IgG and IgM testing: Example: 10 µL serum + 190 µL of PBS.

For IgM determination a preabsorption with Rf absorbent (MASTSORB Order Code: 651003) is recommended to reduce interferences e.g. with rheumatic factors. **Do not absorb the controls.**
4. Controls are supplied ready to use and must not be diluted.
5. Carefully remove the required number of slides from their sachets and mark accordingly. Hold the slide at the edge and do not touch the wells. Place the slides in a moist chamber.
6. Apply 20–25 µL of diluted specimens and ready to use controls to respective wells on the slides according to a prepared work format. Ensure that all the wells are covered and that serum does not escape from the wells.

Direct contact of a pipette tip with the slide surface may result in damage to the antigen substrate and should be avoided.
7. Cover the moist chamber and incubate as follows:
IgG assay: room temperature for 30 min.
IgM assay: room temperature for 45 min.
8. After incubation wash slides carefully with PBS from a wash bottle, taking care not to direct the jet onto the test wells. This may be done by directing the jet of PBS along the centre of the slide, tilting the slide first towards wells 1–5 and then towards 6–10.
9. Immerse the slides in a staining dish or Coplin jar containing PBS and wash for 15 min. A change of the PBS buffer after 10 min increases washing strength.
10. Remove slides one at a time from the staining jar and drain off any excess PBS. Using a blotter, dry the area outside wells. Wipe off any excess PBS from the back of the side. Do not touch the well surface.
11. Immediately transfer the slide to a moist chamber and add 20–25 µL of FITC-conjugate (either IgG or IgM) to each well. Ensure that all the wells are covered by the conjugate.
12. Incubate slides at room temperature in a covered moist chamber in the dark for 30 min (same incubation for IgG and IgM).
13. Switch on microscope 15 min prior to use.
14. After incubation wash the slides with PBS as in steps 9 and 10.
15. Optional counter stain: add 3–5 drops of Evans Blue to the PBS wash buffer in the last 5 min washing step. The counter stain increases the contrast of the immunofluorescent pattern.

Briefly wash the slides in PBS to remove all Evans Blue.

- Apply one small drop of Mounting Medium to each well on a slide. Place a cover slip on the slide and read the results immediately. Take care not to trap air pockets under the cover slip. Excess Mounting Medium on a slide may result in high background fluorescence due to light scattering. Remove excess Mounting Medium with a paper towel avoiding any direct movement of the cover slip.
- Examine under a fluorescence microscope at a total magnification of 400x to 500x.

Interpretation of Results

Test Validation

Examine the control samples before reading patient samples. The controls must give the following patterns to be valid.

Sample	Fluorescence Intensity
Positive control IgG	3 + to 4 +
Positive control IgM	1 + to 3 +
Negative control	negative

Specimen sera are assessed by way of comparison with the controls included in the test.

Interpretation of Specimen Results

IgG positive specimens show a green cytoplasmic and / or nuclear fluorescence in infected cells.

IgG negative specimens show no green fluorescence; a green staining without fluorescence may be visible.

IgM positive specimens show a fluorescence of cytoplasmic inclusions and a more intense fluorescence of the cytoplasm of infected cells. A diffuse fluorescence in the cytoplasm of infected and non-infected cells may be found. The intensity of the fluorescence signal is dependent on the antibody concentration.

IgM negative specimens do not show any green fluorescence.

Limitations / Cross reactions

- Varicella-zoster virus and viruses of the herpes virus group (HSV-1, HSV-2, EBV and CMV) do show antigen homologies. So antibody cross-reactions may be found among those viruses. Positive and borderline results should therefore be re-tested with another test to confirm the very result.
- Borderline results should be re-tested with a fresh collected serum or plasma sample.

MASTAFUOR™ VZV has been demonstrated to be highly sensitive and specific, however, results should be considered along with all serological tests, clinical history and other aspects of patient management to be considered diagnostically significant.

Diagnostic Relevance

Specific IgG and IgM antibodies are found within 3 to 6 days after a primary infection with varicella-zoster virus. IgM antibodies can be detected for about 6 weeks, whereas IgG antibodies may persist at low titers for the rest of the life. So a reinfection with Varicella zoster is found rather rarely.

In case of herpes zoster a rise in IgG and IgM antibody titers is detectable within a few days. IgM antibody titers are found in about 35 % of all case.

Performance

Data in the table were obtained by comparing MASTAFUOR™ VZV results with data of other assays based on a different method (e.g. ELISA):

	IgG	IgM
Sensitivity	97.9 %	75 %
Specificity	33.3 %	100 %

Inter-assay and intra-assay variances are not detectable in respect to control fluorescence intensity.

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