FLOQSwabs[®] hDNA Free

Package insert and How to use guide

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FLOQSwabs[™] is a Copan Italia trademark.

FLOQSwabs[™] are PATENTED. See www.copangroup.com/patents

**hDNA Free FLOQSwabs™ – Buccal FLOQSwab™ with Active Drying System
Information Leaflet & User Guide**

See the glossary of symbols at the end of the information leaflet

INTENDED USE

Buccal FLOQSwab™ with Active Drying System is an hDNA-free buccal swab for taking human DNA samples.

SUMMARY AND PRINCIPLES

One of the routine procedures for genetic investigation is collecting and transporting buccal⁽¹⁾ samples for analysis of the human DNA profile. This can be accomplished using the Buccal FLOQSwab™ with Active Drying System, a product certified for the absence of amplifiable human DNA and detectable DNase and RNase.

The Active Drying System eliminates the need to dry the samples before transporting them; the samples can be stored at ambient temperature (20-25°C) or frozen up to -80°C for maximum 12 months⁽³⁾.


PRODUCT DESCRIPTION

Buccal FLOQSwab™ with Active Drying System is supplied in a tube.

The swab has a standard tip and a breakpoint that allows easily breaking the swab tip inside a test tube after taking the sample. The Buccal FLOQSwab™ shaft is inserted into the tube cap which acts as grip.

The Active Drying System consists of a desiccant located inside the tube cap, which absorbs the moisture of the sample collected aiding its preservation.

Buccal FLOQSwab™ with Active Drying System is subjected to a treatment with ethylene oxide in order to minimize the presence of PCR-amplifiable human DNA.

Product code	Description
50E010D01	Buccal FLOQSwab™ with Active Drying System, standard tip, breakpoint at 20 mm, in tube

REQUIRED MATERIALS BUT NOT INCLUDED

Materials suitable for DNA extraction, amplification and analysis. For these procedures, consult the reference laboratory manuals.

PRODUCT STORAGE

The product is ready to use and does not require additional preparation. It must be stored in its original package at a temperature between 2 and 30°C until the time of use. Do not overheat. Do not incubate or freeze prior to use. If incorrectly stored, its effectiveness will be compromised. Do not use after the expiry date clearly printed on the tube and on each single package.

LIMITATIONS

1. The product is not sterile.
2.  This product is intended for single use. Reuse may lead to a risk of cross-contamination and/or inaccurate results.
3. Buccal FLOQSwab™ is suitable for collecting and transporting buccal samples for human DNA analysis.
4. This device may only be used by adequately trained and qualified personnel.
5. The swab is not intended to be used as buccal sampling system for microbiological investigations.

WARNINGS AND PRECAUTIONS

1. The product is not suitable for applications different from the intended use.
2. Carefully follow the instructions for use.
3. The amount of sample collected for DNA analysis varies considerably from one donor to the next and depends on the collection procedure. For the correct collection procedure, consult the instructions for taking the sample.
4. The product is ready to use: it does not have to be pre-treated or moistened before taking the sample.

5. If the sample has been frozen, bring it to ambient temperature before use.
6. Do not detach the shaft from the grip (tube cap).
7. Do not bend or model the swab before collecting the sample from the patient, as this may cause accidental breakage of the shaft or its detachment from the cap.
8. Adopting handling and analysis procedures different from those described in this information leaflet may compromise product performance and alter the analysis results.
9. All the samples and materials used for the genetic test must be handled in human DNA-free conditions in order to prevent post-collection contamination by personnel with the consequent risk of unreliable results. In order to prevent sample contamination, wear disposable gloves, a mask and a cap, and other protection devices in line with the general precautions.
10. Do not re-pack.
11. The use of this product in combination with extraction and amplification kits must always be validated by the user.
12. The use of this product in combination with rapid diagnostic kits or other instruments must be validated by the user prior to use.
13. Do not use if the swab is visibly damaged (e.g. if the swab tip or shaft is broken).
14. Do not force or excessively press the swab when taking the sample from the patient, as this may result in accidental breakage of the shaft or its detachment from the cap.
15. Unused products, waste and the samples must be disposed of in compliance with the local regulations.

PRODUCT DETERIORATION

Do not use Buccal FLOQSwab™ with Active Drying System if 1. the product shows signs of damage or contamination; 2. the product has passed its expiry date; 3. the swab package is open and/or the protective seal is not intact.

INSTRUCTIONS FOR USE

Sample collection

Correctly collecting the buccal sample from the donor is crucial for the outcome of the DNA analysis. For specific instructions on the sample collection procedures, refer to the procedure described below or to the internal collection procedures.

During sample collection and handling, wear protective clothing, gloves and goggles.

Hold the swab by the grip when taking the sample.


NOTE: The sample must be collected immediately after opening the transparent bag.

NOTE: Do not force or excessively press or bend the swab when taking the sample from the patient, as this may result in accidental breakage of the shaft. The swab shaft has a moulded breakpoint in proximity of the tip, which allows intentionally breaking the swab inside a test tube.

NOTE: The donor must not smoke, eat, drink (except water), brush his or her teeth or chew chewing gum one hour before the sample is taken.

1. Put on gloves, open the bag and remove the plastic tube (Figure 1A).
2. To facilitate opening the tube cap, hold the swab by the plastic tube with one hand and with your thumb press the top of the white cap (Figure 1B).
3. Holding the swab by the grip (white cap), remove the swab from the plastic tube and place the tube on a clean surface (Figure 1C).
4. Holding the grip, insert the swab into the mouth onto the tongue and moisten it with saliva for ten seconds (Figure 1D).
5. Move the swab to one side of the mouth and forcefully rub it between the cheek and the gum for ten seconds. Using the other hand, exercise light pressure on the cheek during collection in order to increase the friction during rubbing (Figure 1E).
6. Repeat step 5 on the opposite side of the mouth (Figure 1F).
7. Once you have taken the sample, put the swab back into the plastic tube and securely close it (Figure 1G).
8. Identify the collected sample. Send the sample to the laboratory for testing as soon as possible.

Figure 1. Specimen Collection

Sample transport and storage

Once the samples have arrived at the laboratory, they can immediately be analysed or stored at ambient temperature (20-25°C) or frozen up to -80°C for maximum 12 months. It is inadvisable to store the samples at a temperature of +4°C after collection, as humidity may accelerate DNA degradation.

Use in the laboratory

Wear protective clothing, gloves and goggles when handling samples for DNA analysis.

When applying molecular methods, take the appropriate precautions in order to prevent cross-contamination. Physical separation of the working areas and a one-way work flow are essential in order to prevent cross-contamination of the amplicon.

Nucleic acid extraction

During sample processing, when handling the applicator, do not touch the area underneath the breakpoint, i.e. the area between the breakpoint and the Nylon flocked swab tip, otherwise you risk contaminating the shaft invalidating the test results.

Acting on the breakpoint of the swab shaft, transfer the swab to an extraction test tube following the standard operating procedures of the laboratory. Add the correct amount of lysis liquid as indicated in the information leaflet of the extraction kit used. Continue according to the extraction and amplification kit procedures.

The device was tested internally using the following extraction methods: ChargeSwitch Forensic DNA Purification Kit di Invitrogen; Chelex 100 e InstaGene[™] matrix di Bio-Rad⁽²⁻⁴⁾; PrepFiler and PrepFiler Express di Thermo Fisher⁽³⁾; Nucleo Spin Tissue di Macherey Nagel; DNA IQ system di Promega; QIAamp DNA Investigator Kit⁽⁵⁾ e QiaAmp mini kit⁽⁶⁾ di Qiagen.

The extracted DNA was tested using the following PCR kits: AmpF[®]STR[®] Yfiler[™] PCR Amplification Kit, Identifier Plus PCR Amplification Kit, NGM PCR amplification kit(2), AmpF[®]STR[®] SGM Plus[™] (6), Quantifiler Trio Human DNA Quantification kit, Quantifiler Duo DNA Quantification kit, Quantifiler Human DNA Quantification kit(3) di Thermo Fisher.

QUALITY CONTROL

Buccal FLOQSwab[™] with Active Drying System is guaranteed for the following characteristics:

- **Analyte release:** the swab is inoculated with a bacterial suspension, transferred to a liquid phase and tested to check complete analyte release. The swab must release $\geq 70\%$ of analyte.
- **Absence of DNase:** Following the standard internal operating procedures, the device is saturated with a solution containing standard DNA ladders. The DNA originating from the solution exposed to the swabs must be intact.
- **Absence of RNase:** Following the standard internal operating procedures, the device is saturated with a solution containing standard RNA ladders. The RNA originating from the solution exposed to the swabs must be intact.
- **Absence of human DNA:** Following the standard internal operating procedures, the device is tested to ascertain the absence of amplifiable human DNA using the Identifier[®] Plus kit by Life Technologies.
- **Effectiveness of the *Active Drying System*:** Following the standard internal operating procedures, the device is inoculated with a known amount of osmotized water. Once enclosed in the relative tube, the swab must dry completely in 24 hours at ambient temperature.

FLOQSwabs[™] is a Copan Italia trademark.

FLOQSwabs[™] is PATENTED. See www.copangroup.com/patents

FLOQSwabs™ hDNA Free – Buccal FLOQSwab™ con Active Drying System Folleto ilustrativo y guía de uso

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BIBLIOGRAPHY

1. **VALUTAZIONE DELL'EFFICIENZA DELLA METODICA DI PRELIEVO SALIVARE MEDIANTE TAMPONI COPAN 4N6FLOQSwabs™ genetics** Corato S., Ponzano E., Giuliodori A., Rodriguez D., Caenazzo L. Presented at XXIV Convegno Nazionale Genetisti Forensi Italiani in Pavia on 20-22 Settembre 2012.
2. **Enhanced semen elution from nylon 4N6FLOQSwabs™ prior DNA analysis** P. MarDñez, MG.Blanco, G.Vallejo National Institute of Toxicology & Forensic Science INTCF José Echegaray 4, 28232 Las Rozas Madrid SPAIN Presented at International Symposium on Human Identification in Phoenix on 29 September-2 October 2014.
3. **Genetics 4N6FLOQSwabs™ are supporting DNA preservation from buccal swabs for long term storage** A. Squassina S. Castriciano Copan Italia S.p.A. Brescia, Italy. Presented at International Society for Forensic Genetics in Krakow on 31 August – 5 September 2015.
4. **FORENSIC FLOCKED COLLECTION DEVICES TO MAXIMIZE CRIME SCENE SAMPLE PROCUREMENT** Squassina, A., Poli, A., Gervasoni, A., Castriciano, S., Paroni, R., Triva, D. Copan Italia S.p.A. Brescia, Italy Presented at International Society for Forensic Genetics in Vienna on 29 August – 3 September 2011.
5. **Improved recovery of DNA with the 4N6FLOQSwab™ system and Nucleic Acid Optimizer (NAO basket™)** A.Ambers, R. Wiley,N. Novroski, B. Frey, A. Machinnis, J. King, B. Budowle
6. **Post-coital vaginal sampling with nylon flocked swabs improves DNA typing** Corina C.G. Benschop, Danielle C. Wiebosch, Ate D. Kloosterman, Titia Sijen.

Tabella dei Simboli / Table of Symbols / Tableau des pictogrammes / Tabelle der verwendeten Symbole / Tabla de símbolos / Tabela de Símbolos

Simbolo / Symbol / Pictogramme Símbolo / Symbol / Símbol	Significato / Meaning / Signification / Bedeutung / Significado
	Fabbricante / Manufacturer / Fabricant / Hersteller / Fabricante / Fabricante
	Marcatura CE / EC marking / Marquage CE / EC –Kennzeichnung / Mercado CE / Marcação CE
	Non riutilizzare / Do no reuse / Ne pas réutiliser / Nicht zur Wiederverwendung / No reutilizar / Não reutilizar
	Numero di catalogo / Catalogue number / Numéro de catalogue / Bestellnummer / Número de catalogo / Número de catálogo
	Limiti di temperatura / Temperature limits / Limites de température / Temperaturbegrenzung / Límites de temperatura / Límites de temperatura
	Utilizzare entro / Use before / À utiliser avant / Verwendbar bis / Usar antes de / Usar antes de
	Consultare le istruzioni per l'uso / Consult the instructions for use / Consulter les instructions d'utilisation / Gebrauchsanweisung beachten / Consultar las instrucciones de uso / Consultar as instruções de utilização
	Codice del lotto (partita) / Batch code (lot) / Code du lot (lot) / Chargenbezeichnung / Código del lote (Lote) / Código do lote (remessa)
	Non utilizzare in caso di confezionamento danneggiato / Do not use if the packaging is damaged / Ne pas utiliser si l'emballage est endommagé / Nicht verwenden, wenn Verpackung beschädigt / No utilizar si el envase está dañado / Não utilizar se a embalagem apresentar danos
	Non sterile / Non-sterile / Non stérile / Nicht steril / No estéril / Não estéril

 Copan Flock Technologies Srl Via F. Perotti, 18 25125 - Brescia, Italy	Copan Flock Technologies Srl Via Perotti 18 25125 Brescia, Italy Tel +39 030 3666100 Fax +39 030 2659932 Email: info@copanflock.com Website: www.copangroup.com	 COPAN Innovating Together™
North American Distributor: Copan Diagnostics Inc. 26055 Jefferson Avenue Murrieta, CA 92562 USA Tel: 951-696-6957 Fax: 951-600-1832 E-mail: customerservice@copanusa.com Website: www.copanusa.com		