**MAST ISOPLEX® DNA / RNA  
Extraction Kit**

MAST ISOPLEX® DNA / RNA Extraction Kit

**REF** 67DNAEX

50 Tests

**Gebrauchsanweisung /  
Instructions for Use /  
Notice d'utilisation**

**Nur zur Verwendung durch Fachpersonal /  
For professional use only /  
Exclusivement pour un usage professionnel**

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## MAST ISOPLEX® DNA / RNA Extraction Kit

Kit for extraction of nucleic acids (NA) for further use in MAST ISOPLEX test systems.

### Intended Use

The kit is suitable for NA extraction of bacteria, viruses and protozoa.

Starting material can be liquid, human swab samples, samples from culture plates as well as stool, urine or blood samples.

### Important Note for Use of these Kit Instructions

Any assay relevant changes to the kit instructions for use (IFU) will lead to a change of the version number at the bottom on the last page. All changes made will be identified on a separate sheet added to the IFU for a period of three months from the date of a change of the version. Please ensure that the latest version of the IFU is used for the assay procedure.

### Principle of the Test

Bead-based methods are simple and reliable methods for the extraction of nucleic acids (NA). First, the sample is lysed. Under optimal conditions, the NA bind selectively to the surface of the magnetic beads, while other components remain in solution. Purified NA can be used directly in various molecular biological applications (sequencing, PCR, etc.). Extraction can be performed manually, semi-automatically or fully automatically.

### Kit Contents

1. **LYSIS** Lysis Solution V  
1 x 15 mL
2. **PROT K** Proteinase K  
for 1 x 1.5 mL, working solution
3. **BINDING** Binding Solution V  
1 x 25 mL
4. **MAG** MAG Suspension F  
1 x 1.1 mL
5. **WASH 1** Washing Solution HS  
1 x 15 mL (final volume 30 mL)
6. **WASH 2** Washing Solution LS  
1 x 12 mL (final volume 60 mL)
7. **WATER** RNase-free Water  
1 x 6 mL
8. **TUBE** RNase free Tubes (1.5 mL)  
50 pieces
9. Instructions for use

### Warning and Precautions

1. Read instructions carefully before conducting the assay. Do not modify the procedure without prior validation.
2. Do not use reagents beyond the expiry date.
3. Comply with the general health and safety guidelines for working with potentially infectious materials. Wear appropriate protective clothing and use appropriate lab facilities.
4. Do not mix reagents between different lots, as reagents have been calibrated for each batch.
5. Do not cross-contaminate reagents or interchange caps on bottles. Use a separate pipette or pipette tips for each sample and reagent.
6. Samples and contaminated disposables should be disposed according to relevant disposal directives and regulations for infectious materials. For the decontamination of surfaces use suitable surface disinfectants. Contaminated glassware should be autoclaved at 121 °C.
7. If buffer bottles are damaged or leaking, wear gloves and protective goggles when discarding the bottles in order to avoid any injuries.
8. Below the European Community risk and safety phrases for components of the MAST ISOPLEX® DNA / RNA Extraction Kit to which they apply are listed.

Binding Solution V: contains 2-propanol; highly flammable, irritant (R11, 36, 67, 7, 16, S24/25/26)

Proteinase K: irritant, sensitizing. Risk and safety phrases: R36/37/38-42; S22-24-26-37/38

Washing Solution HS: contains guanidinium thiocyanate: harmful. Risk and safety phrases: R20/21/22-32; S13-26-36-46

For details please refer to the MSDS.

### Stability and Storage

#### Important

Store lyophilized **Proteinase K** at 4–8 °C!

Divide dissolved **Proteinase K** into aliquots (200 µL) and store at -20 °C. Repeated freezing and thawing will reduce the activity dramatically!

The lyophilized Proteinase K may be used from manufacture until the expiry date indicated on the package, if stored at 4–8 °C.

#### Storage conditions

Store MAG Suspension F at 4–8 °C; avoid freezing.

All other components are stored at room temperature.

### Components not included in the kit

1. RNase free tubes.
2. ddH<sub>2</sub>O (double distilled water) for dissolving **Proteinase K**.
3. 96–99.8 % ethanol, non-denatured or methylated.

## Initial steps

1. Add 15 mL of 96-99.8 % ethanol to the bottle **Washing Solution HS**, mix thoroughly and keep the bottle always firmly closed!
2. Add 48 mL of 96-99.8 % ethanol to the bottle **Washing Solution LS**, mix thoroughly and keep the bottle always firmly closed!
3. Dissolve **Proteinase K** by addition of 1.5 mL of ddH<sub>2</sub>O, mix thoroughly and store as described!

## Recommended steps before starting

1. Ensure that the **Washing Solution HS**, **Washing Solution LS** and, **Proteinase K** have been prepared according to the instruction.
2. Mix **MAG Suspension F** very well by inverting the tube several times.
3. Preheat heating block to 60 °C.

## Sample Preparation

### Liquid / urine / blood samples:

No specific preparation is necessary. Start with step 1 "Isolation of Nucleic Acids".

### Swab samples:

#### • Test systems with liquid media:

The liquid medium can be used directly in step 1 "Isolation of nucleic acids".

#### • Swab samples without liquid medium

Transfer swabs to tubes containing 1 mL (0,9 % NaCl) physiological saline solution or PBS and incubate for 20 minutes, shaking continuously.

Squeeze out and remove swab, then insert 200 µL of liquid sample in step 1 "Isolation of nucleic acids".

### Agar plates:

- Pick a few (3–5) fresh overnight colonies, resuspend in 200 µL PBS and continue with step 1 "Isolation of Nucleic Acids".

### Fecal samples:

- Add approx. 50-100 mg stool sample into a 1.5 mL reaction tube.
- Add 250 µL PBS (not included in kit) and vortex for 10 sec
- Spin at approx. 16000 rpm (approx. 21000 x g) for 3 min
- With 200 µL supernatant, continue as described in step 1 "Isolation of nucleic acids".

### Tissue samples:

- Pipette to 5-10 mg tissue sample 400 µL ddH<sub>2</sub>O (RNase-free) or PBS and homogenize (tissue homogenizer)
- Centrifuge the homogenate for 2 min at 10,000 x g.
- With 200 µL supernatant, continue as described in step 1 "Isolation of nucleic acids".

## Isolation of Nucleic Acids

1. Add **200 µL sample** and **200 µL Lysis Solution V** in a **RNase free tube (1.5 mL)**. Add **20 µL Proteinase K** to the sample and mix well by pipetting up and down.
2. Incubate the sample at room temperature for 15 min.  
**Note:** Mix the sample several times during incubation. No mixing will reduce the lysis efficiency.
3. Add **400 µL Binding Solution V** and **10 µL MAG Suspension F** in the reaction tube.

Mix the sample by pipetting up and down several times. Wait for 30 sec and mix again by pipetting up and down several times.

Place the reaction tube in magnetic rack. Wait for 20 sec and carefully remove the liquid as much as possible.

**Note:** Do not remove the magnetic beads.

4. Take the reaction tube out of magnetic rack and add **500 µL Washing Solution HS**. Mix well by pipetting up and down several times until all magnetic beads are dissolved well.

Place the reaction tube in magnetic rack. Wait for 20 sec and carefully remove the **liquid** as much as possible.

**Note:** Do not remove the magnetic beads.

5. Take the reaction tube out of magnetic rack and add **500 µL Washing Solution LS**. Mix well by pipetting up and down several times until all magnetic beads are dissolved well.

Place the reaction tube in magnetic rack. Wait for 20 sec and carefully remove the liquid as much as possible.

**Note:** Do not remove the magnetic beads.

6. Take the reaction tube out of magnetic rack and add **500 µL Washing Solution LS**. Mix well by pipetting up and down several times until all magnetic beads are dissolved well.

Place the reaction tube in magnetic rack. Wait for 20 sec and carefully remove the liquid as much as possible.

**Note:** Do not remove the magnetic beads.

Remove residual alcohol by pipetting from the reaction tube and also from the lid as much as possible.

7. Take the reaction tube out of magnetic rack and air-dry with open lid at room temperature until alcohol is evaporated completely (approx. 35 to 45 min) or 10–15 min at 60 °C

**Note:** Do not remove the magnetic beads.

8. Resuspend the magnetic beads in **80 µL RNase-free water** by pipetting up and down several times. Wait for 1 min and resuspend again the magnetic bead by pipetting up and down.

Place the tube into the magnetic rack and wait for 20 sec and transfer the liquid sample completely into the new RNase free tube.

**Note:** Avoid carryover of magnetic beads.

**Note:** The extracted nucleic acids are ready to use for downstream application. For short-term store the extracted nucleic acids at +4 °C. For long-term storage at -20 °C is recommended.

## Assay Data

Product specifications:

1. Starting material:  
Up to  $1 \times 10^9$  cells of Gram+ or Gram- bacteria.
2. Time for isolation: Approximately 60 min or approximately 45 min after the lysis step.
3. Binding capacity: > 50 µg NA.
4. The typical yield is sufficient for amplification procedures, up to 35 µg.

## Troubleshooting

Problem / probable cause	Comments and suggestions
<b>Poor lysis of starting material:</b>	
Insufficient disruption or homogenization	After lysis centrifuge lysate to pellet debris and continue with the protocol using the supernatant. Reduce amount of starting material.
<b>Little or no total RNA eluted:</b>	
Insufficient disruption or homogenization	Reduce amount of starting material. Overloading reduces yield!
<b>DNA contamination:</b>	
Too much starting material	Reduce amount of starting material.
Incorrect lysis of starting material	Use the recommended techniques for lysis of cell pellet.
<b>Total RNA degraded:</b>	
RNA source inappropriately handled or stored	Ensure that the starting material is fresh! Ensure that the protocol, especially the first steps, has been performed quickly.
RNase contamination of solutions; Receiver Tubes, etc.	Use sterile, RNase-free filter tips. Before every preparation clean up the pipette, the devices and the working place. Always wear gloves!
<b>Total RNA does not perform well in downstream applications:</b>	
Salt carryover during elution	Ensure that Washing Solution HS and Washing Solution LS are at room temperature. Checkup Washing Solution for salt precipitates. If there are any precipitate dissolves these precipitate by carefully warming.

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**Mast Diagnostica GmbH,**
**Mast Group Ltd.**