eNAT™

**Collection and Preservation System  
Product Insert & How to Use Guide**

ENGLISH

**Copan eNAT™ Collection and Preservation System  
Information Leaflet & User Guide**

See the glossary of symbols at the end of the information leaflet

**INTENDED USE**

Copan eNAT™ Collection and Transport System is intended for the collection, transport and preservation of clinical specimens to be analyzed by nucleic acids amplification techniques.

eNAT™ medium stabilizes and preserves RNA/DNA for prolonged time periods and is compatible with commercial nucleic acid extraction and amplification platforms.

**SUMMARY AND PRINCIPLES**

Clinical specimens stored and transported in eNAT™ medium can be processed, using standard clinical laboratory operating procedures, for the detection of nucleic acids of Viruses, Bacteria, Chlamydia, Protozoa and Mycoplasma with molecular amplification assays.

The primary purpose of nucleic acids amplification techniques is to screen for a wide range of infectious diseases, so nucleic acids integrity of clinical specimens during transport and storage should be preserved. eNAT™ medium contains a detergent and a protein denaturant to prevent microbial proliferation, thus eNAT™ is not intended to be used for culture based techniques.

**PRODUCT DESCRIPTION**

eNAT™ is available in the product configurations indicated in the table below:

CODE	DESCRIPTION
608CS01M	12 mm x 80 mm Screw cap tube containing 1 ml of eNAT™ Transport and Preservation Media Plus One minitip applicator swab with red mark printed break point and flocked nylon fiber tip
608CS01P	12 mm x 80 mm Screw cap tube containing 1 ml of eNAT™ Transport and Preservation Media plus One pernasal applicator swab with red mark printed break point and flocked nylon fiber tip.
608CS01R	12 mm x 80 mm Screw cap tube containing 1 ml of eNAT™ Transport and Preservation Media Plus One regular size applicator swab with red mark printed break point and flocked nylon fiber tip
606CS01L	12 mm x 80 mm Screw cap tube containing 2 ml of eNAT™ Transport and Preservation Media Plus One L shape applicator swab with flocked nylon fiber tip
606CS01M	12 mm x 80 mm Screw cap tube containing 2ml of eNAT™ Transport and Preservation Media Plus One minitip applicator swab with red mark printed break point and flocked nylon fiber tip
606CS01P	12 mm x 80 mm Screw cap tube containing 2 ml of eNAT™ Transport and Preservation Media plus One pernasal applicator swab with red mark printed break point and flocked nylon fiber tip.
606CS01R	12 mm x 80 mm Screw cap tube containing 2 ml of eNAT™ Transport and Preservation Media Plus One regular size applicator swab with red mark printed break point and flocked nylon fiber tip
6E021S	eNAT™ Collection kit comprises: - 12x80 mm screw cap tube filled with 1 ml of eNAT™ transport and preservation medium - one regular applicator FLOQSwab™ - one 2 ml Pasteur pipet

Other product codes may be available. For updates please refer to our website: [www.copangroup.com](http://www.copangroup.com)

Performance testing with eNAT™ was conducted using laboratory strains. Performance testing was not conducted using human specimens.

For eNAT™ codes 608CS01M; 608CS01P; 608CS01R; 606CS01L; 606CS01M; 606CS01P; 606CS01R:

Fifty (50) units are contained in a shelf pack and 10 x 50 units are contained in a box. Each unit consists of a sterile package containing two components: a pre-labeled polypropylene screw-cap tube with conical shaped bottom filled with 1or 2 ml of eNAT™ transport medium and a specimen collection swab which has a tip flocked with soft nylon fiber. Three main types of collection applicators are available; including a tube of medium but each has a different type of swab applicator. One unit type contains a regular size flocked nylon swab applicator intended for the collection of samples from the nose, throat, mouth, vagina, rectum, faeces or wounds, the second unit type contains a minitip size flocked nylon swab applicator intended for the collection of samples from small or less accessible areas such as the eye, ear, nasal passages, urogenital tract and the third type containing a pernasal flocked nylon swab applicator intended for the collection of samples from the nasopharynx. These different types of swab applicators facilitate the collection of specimens from different sites on a patient. Refer to the individual product descriptions for specific information about materials supplied.

For eNAT™ code 6E021S:

Fifty (50) units are contained in a shelf pack and 6 x 50 units are contained in a box. Each unit consists of a sterile peel pouch containing three components: one pre-labelled polypropylene screw-cap tube with conical shaped bottom filled with 1 ml of eNAT™ transport medium, a regular specimen collection FLOQSwab™ and a 2 ml Pasteur pipet.

The Pasteur pipet is designed for the transfer of urine from a standard urine container into the eNAT™ tube. The FLOQSwab™, which has a tip flocked with soft nylon fiber and a red marked molded breakpoint in the shaft of the applicator (Fig.2), is intended for the collection of samples from the nose, throat, vagina, endocervix or wounds. Other sampling sites may be mouth, rectum, faeces.

All collection swab applicators provided with eNAT™ have a molded breakpoint in the shaft of the applicator (Fig 1). Some types of collection swab applicators provided with eNAT may have a red printing in the correspondence of molded breakpoint (Fig.2). The red mark favors the quick and even easier identification of the presence of a break point. Refer to the product configurations described in the table below, for specific information on presence of red mark on shaft break point. After the sample is collected from the patient, the molded breakpoint facilitates easy breakage of the swab applicator into the eNAT™ tube of transport medium eNAT™ tube caps have an internal molded design that is able to capture the swab shaft when it is broken off into the tube and the cap is closed. The action of screwing the cap onto the tube moves the end of the broken swab shaft into a funnel shaped molded docking receptacle in the cap. This molded funnel shape effectively captures the end of the broken applicator shaft and secures it firmly in the dock by friction grip. In the testing laboratory when the swab cap is unscrewed and removed, the swab applicator is attached to the cap. This feature allows the operator to conveniently remove the swab from the transport tube Due to the flexibility of the shaft of minitip, pernasal, urethral and pediatric swabs, the capture cap feature is not applicable, as the broken applicator may not firmly fit into the cap.


Fig. 1


Fig. 2

## REAGENTS

Guanidine thiocyanate  
Tris-EDTA  
HEPES  
Detergent

## REQUIRED MATERIALS THAT ARE NOT INCLUDED

Appropriate materials for molecular testing according to recommended protocols as per laboratory reference manuals. Urine cup for urine collection.

## STORAGE OF THE PRODUCT

This product is ready for use and no further preparation is necessary. The product should be transported and stored in its original container at 5-25°C until use. Do not overheat. Do not incubate or freeze prior to use.

Improper storage will result in a loss of efficacy. Do not use after expiration date, which is clearly printed on the outer box and on each individual specimen transport vial label.

## WARNINGS AND PRECAUTIONS

1. eNAT™ medium is not for external or internal use in humans or animals.
2. This product is for single use only; reuse may cause a risk of infection and/or inaccurate results.
3. Do not use excessive force, pressure or bending when collecting swab samples from patients as this may result in accidental breakage of the swab shaft.
4. Observe approved biohazard precautions and aseptic techniques. To be used only by adequately trained and qualified personnel.
5. Directions should be read and followed carefully.
6. Not suitable for any other application than intended use.
7. This product is intended for professional use or for use in POC
8. The use of eNAT kit for the self- collection of vaginal or cervico/vaginal samples in pregnant women has not been evaluated. In case of pregnancy the self-collection procedure is not recommended.
9. The use of this product in association with diagnostic instrumentation should be previously validated by the user.
10. Do not use if the product is visibly damaged
11. Do not use if the eNAT™ Collection Kit peel pouch is damaged or open.
12. eNAT™ medium does not require any concentration step. Pelleting procedure is not allowed since it may lead to loss of sensitivity.
13. Do not use eNAT™ medium for pre-moistening or pre-wetting the sampling device prior collecting the sample or for rinsing or irrigating the sampling sites.
14. Do not use the same tube for more than one patient. This causes wrong or bad diagnosis.
15. Keep container closed when not in direct use
16. Before transporting, make sure eNAT™ screw cap tube is tightly closed
17. eNAT™ was tested for microbial viability inactivation of Gram positive bacteria, Gram negative bacteria , yeasts and molds like, but not limited to, Candida albicans, Escherichia coli, Staphylococcus aureus, Pseudomonas.aeruginosa, Aspergillus brasiliensis, Shigella sonney,

18. *Salmonella typhimurium*, and *Streptococcus agalactiae*. Starting from an initial strain concentration of  $\geq 10^5$  CFU/ ml, bacteria and yeasts are completely inactivated in  $\leq 30$  minutes. Molds are completely inactivated in  $\leq 1$  hour. However universal precautions for safe handling of biological fluids should be practiced at all times.
19. Dispose of unused reagents, waste and specimens in accordance with local regulations.
20. Avoid contact of eNAT™ medium with skin, mucous membranes. If contact does occur, immediately wash with a large amounts of water.
21. eNAT™ medium contains guanidine thiocyanate. Avoid direct contact between guanidine thiocyanate and sodium hypochlorite (bleach) or other highly reactive reagents such as acids and bases. These mixtures could release noxious gas.
22. Contains Guanidine thiocyanate.

**Warning!**

H302+H332: Harmful if swallowed or if inhaled  
H314: causes severe skin burns and eye damage  
H412: Harmful to aquatic life with long-lasting effects  
EUH032: contact with acids liberates very toxic gas

MSDS available on request from Copan Italia spa via Perotti 10, 25125 Brescia Italy.

**PRODUCT DETERIORATION**

eNAT™ should not be used if (1) there is evidence of damage or contamination to the product, (2) there is evidence of leakage, (3) the expiration date has passed or (4) there are other signs of deterioration.

**PRODUCT STABILITY**


eNAT™ has been tested internally and externally for product stability during its shelf life.

**INSTRUCTIONS FOR USE****Specimen Collection**

Proper specimen collection from the patient is extremely critical for the successful identification of infectious organisms. For specific guidance regarding specimen collection procedures, consult published standard collection manuals. Qualified medical personal only shall collect clinical specimens using proper sampling devices. Vaginal and cervico/vaginal sampling have been validated also for the self-collection procedure performed by non-pregnant patient. The choice of the proper swab for specimen collection must be done by the doctor. Qualification studies for vaginal and cervico/vaginal self-collection were performed by Copan for swabs included in product part number 606CS01R and 608CS01R. In case of self-collection of vaginal or cervico vaginal specimen, patient should follow the specific instructions included in the section: PATIENT SELF-COLLECTION OF VAGINAL OR CERVICO/VAGINAL SWAB SPECIMENS

**Swab specimen collection using eNAT™ system**

1. Open the eNAT™ sample collection pouch, discard the pipet if present and remove the tube. Take the swab out from its pouch.
2. Use the swab to collect the sample from the patient. To avoid risk of contamination, make sure that the swab tip comes into contact only with the sampling site.  
**NOTE:** Swab shafts often exhibit diameter changes to facilitate different sampling requirements. Swab shafts may also have a molded breakpoint point designed for intentional breakage of the swab into a transport tube. In all circumstances when collecting swabs samples from patients, do not use excessive force, pressure or bending of the swab as this may result in accidental breakage of the swab shaft. **Do not use eNAT™ medium for pre-moistening or pre-wetting the sampling device prior collecting the sample or for rinsing or irrigating the sampling site.**
3. Unscrew and remove the cap from eNAT™ tube making sure not to spill the medium.
4. Insert the swab into the tube until the breaking point reaches the level of the opening of the tube.
5. Bend and break the swab at the breaking point holding the tube away from the face.
6. Replace the cap on the tube and close tightly.
7. Write patient information on the tube label or apply patient identification label.
8. Send the sample to the test laboratory.


eNAT™

#### Urine specimen collection

1. Open the eNAT™ sample collection pouch. Remove the eNAT™ tube and discard the pouch containing the FLOQSwab™. Take the Pasteur pipet making sure not to touch the tip. Do not lay the pipet down on a surface.
2. Ask the patient to collect the first 20 to 30 ml of voided urine (the first part of the stream) into a urine cup (not supplied with the kit).
3. Unscrew and remove the cap from eNAT™ tube making sure not to spill the medium.
4. Transfer urine from the cup into the collection tube  
NOTE: In order to avoid excessive medium dilution the volume of liquid specimen to be added to eNAT™ medium should never exceed 1:3 ratio. Maximum fill volume is 6ml.  
Two suggested procedures:
  1. Transfer 2 ml of urine from the cup into the collection tube by using the Pasteur pipet provided with the eNAT™ Collection Kit. The Pasteur pipet has a scale with a mark at every 0.5 ml step. Squeeze the bulb in order to aspirate 2 ml of urine. Make sure to dispose volume exceeding 2 ml back into the cup. Take special care not to introduce contamination in the medium tube.
  2. Transfer 3 ml of urine from the cup into the collection tube by using the Pasteur pipet provided with the eNAT™ Collection Kit. The Pasteur pipet has a scale with a mark at every 0.5 ml step. Squeeze the bulb in order to aspirate 2 ml of urine. Dispense 2 ml urine into collection tube. Make sure to dispose volume exceeding 2 ml back into the cup. Repeat squeezing the bulb in order to aspirate 1 ml of urine. Dispense 1 ml urine into collection tube. Make sure to dispose volume exceeding 1 ml back into the cup. Take special care not to introduce contamination in the medium tube.
5. Replace the cap on the tube and close tightly.  
Mix the urine with the transport medium by vortexing the tube 5 seconds.
6. Write patient information on the tube label or apply patient identification label. Send the sample to the test laboratory

#### Patient self-Collection of Vaginal or Cervico/Vaginal swab specimens

**NOTE:** This procedure is only intended for supervised POC (Point of Care) patient-collection of vaginal or cervico/vaginal specimens. Patient must consult her doctor, nurse or health-care provider and read instructions carefully, before specimen collection . The use of eNAT™ kit for the self- collection of vaginal or cervico/vaginal samples in pregnant women has not been evaluated. In case of pregnancy the self-collection procedure is not recommended.

The choice of the proper swab for specimen collection must be done by the doctor. Qualification studies for vaginal and cervico/vaginal self-collection were performed by Copan for swabs included in product part number 606CS01R and 608CS01R.

*Wash your hands with soap and water before starting, rinse and dry.*

1. Open the eNAT™ sample collection pouch, discard the pipet if present and remove the tube.

Take the swab out from its pouch. Collect the sample with the swab avoiding that the swab tip comes into contact with body areas or surface different from the sampling site.

2. Hold the swab at the red mark molded breakpoint on the shaft

**NOTE:** Swab shaft has a RED MARK which favors the quick and even easier identification of the presence of a break point which was designed for intentional breakage of the swab into the transport tube. Swab shaft exhibits diameter changes to facilitate different sampling requirements. **Do not use eNAT™ medium for pre-moistening or pre-wetting the sampling device prior to collect the sample or for rinsing or irrigating the sampling site.**

3. Get yourself into one of the positions shown on the pictures below:

- A. Stand with the legs apart

- B. Sit on the toilette with your legs wide apart.

Using your free hand, hold back the folds of skin covering the opening of your vagina.

4. Insert the swab into the vaginal opening until the fingers touch the vulva.

5. Gently rotate the swab 3 times against the sides of the vagina. Withdraw the swab carefully.

6. Unscrew and remove the cap from tube making sure not to spill the medium. Insert the swab into the tube until the red mark molded breaking point reaches the level of the opening of the tube.

7. Bend and break the swab at the breaking point holding the tube away from your face.


8. Replace the cap on the tube and close tightly.

9. Return the collection tube to your doctor, nurse or health-care provider.

eNAT™

PI45C Rev.03 Date 2016.07

### **Use in the laboratory**


#### **Processing eNAT™ specimens for molecular testing in the laboratory.**

Specimens received in the laboratory for nucleic acid detection should be processed when received in the laboratory. In case of delay, please refer to the appropriate specimen storage conditions. eNAT™ medium preserves nucleic acids for up to four weeks at room temperature and 4°C and up to 6 months at 20°C. Specimens preserved in eNAT™ medium should be extracted and purified before amplification. eNAT™ medium has been validated with automated platforms like, for example, the easyMAG® (Biomereux), the Abbott m200 platforms (Abbot Diagnostics), QIAasympathy (Qiagen) and the manual method like silica- columns and magnetic beads using the following protocol:

1. Wear gloves and other protection commensurate with universal precautions when handling clinical specimens. Observe other CDC Biosafety Level 2 recommendations
2. When working with NAAT assays, care should be taken to prevent carry over contamination. Spatial separation of working areas and unidirectional workflow are essential to prevent amplicon carry-over.
3. **VORTEXING METHOD:**  
Vortex eNAT™ specimen tube for 10 s  
**NOTE:** if the sample appears too mucousy ,the specimen may remain mainly attached to the swab, extend the vortex time in order to break down mucus clamps and to the sample from the swab.
4. Unscrew the cap and transfer the appropriate amount of sample (e.g. 200ul-400ul or as per the protocol in use for the extraction) directly into the extraction buffer tube.
5. Continue as per extraction and amplification kits procedures.

If unable to use the vortexing method, use the following alternative protocol:

1. **MANUAL SHAKING METHOD:**  
Hold the eNAT™ tube from the cap making sure that it is closed tightly.  
Shake the tube 4 times downward with rapid movements of the wrist (see picture)  
**NOTE:** Inverting the tube up and down, is not recommended. if the sample appears too mucousy, the specimen may remain mainly attached to the swab, it is recommended to extend the shaking times in order to break down mucus clamps and to easier release the sample from the swab
2. Unscrew the cap and transfer the appropriate amount of sample (e.g. 200ul-400ul or as per the protocol in use for the extraction) directly into the extraction buffer tube.
3. Continue as per extraction and amplification kits procedures.

Manual Shaking protocol has been validated using vaginal samples. The use of this protocol with other type of samples, should be previously validated by the user. Other extraction methods may also be applicable prior validation.

For updated list of evaluated kits, please contact Copan Italia.

### **QUALITY CONTROL**

#### **ANTIMICROBIAL ACTIVITY**

eNAT™ medium is routinely tested for its antimicrobial activity against a panel of bacteria (E. coli, S. aureus and C. albicans). A complete microbial viability inactivation is obtained within 30 minutes for these three strains, starting from  $\geq 10^5$  viable CFU inoculated in 1 mL of eNAT™ medium.

### **PERFORMANCE CHARACTERISTICS**

#### **NUCLEIC ACIDS PRESERVATION**

eNAT™ medium preserves nucleic acids for up to 4 weeks at room temperature and 4°C and up to 6 months at -20°C.

Performance testing with Copan eNAT™ was conducted using laboratory strains spiked onto a swab. Performance testing was not conducted using human specimens. Tested strains include among others:

ORGANISM	REF NUMBER	ANALYTE TYPE
INFLUENZA A VIRUS	ATCC VR-822	RNA VIRUS
INFLUENZA B VIRUS	ATCC VR 786	RNA VIRUS
CITOMEGALOVIRUS	ATCC VR-977	DNA VIRUS
HERPES VIRUS TYPE I	ATCC VR-539	DNA VIRUS
CHLAMYDIA TRACHOMATIS	ATCC VR-880	CHLAMYDIA
CHLAMYDIA PNEUMONIAE	ATCC VR-1360	CHLAMYDIA
NEISSERIA GONORRHOEAE	ATCC 43069	GRAM NEGATIVE BACTERIUM
BORDETELLA PERTUSSIS	ATCC 8467	GRAM NEGATIVE BACTERIUM
PSEUDOMONAS AERUGINOSA	ATCC 27853	GRAM NEGATIVE BACTERIUM
MRSA	ATCC 43300	GRAM POSITIVE BACTERIUM
STAPHYLOCOCCUS AERUS	ATCC 6538	GRAM POSITIVE BACTERIUM
MYCOPLASMA HOMINIS	ATCC 23114	MYCOPLASMA
MYCOPLASMA PNEUMONIAE	ATCC 15331	MYCOPLASMA
TRICHOMONAS VAGINALIS	Biomed culture	PROTOZOAN

Results obtained will largely depend on proper and adequate specimen collection, as well as timely transport and processing in the

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Copan Italia SpA
Copan Diagnostics Inc.