

IVD solutions through partnership



# **MASTISOPLEX<sup>®</sup>** *VTEC & E.coli O157*

Individual Kits for the Rapid Detection  
of VTEC and *E.coli* O157

VT1  
VT2

*Escherichia  
coli*

*E.coli*

O157

- Multiplex with novel probe technology
- For the detection of VT1 and VT2
- For the detection of *E.coli* O157
- Inhibition controls included in each test
- Results achieved in 40 minutes

## Introduction

*Escherichia coli* is a Gram-negative, rod-shaped bacterium commonly found in the gastrointestinal tract of humans and warm-blooded animals. While most *E.coli* constitute the commensal flora of the gastrointestinal tract, there are a number of pathogenic *E.coli* strains which can cause a variety of illness in humans and animals.<sup>1</sup>

Verotoxin-producing *Escherichia coli* (VTEC) is one of the major causes of food-borne illness in humans, particularly *E.coli* O157. It is transmitted via the faecal-oral route when food becomes contaminated with the faeces of infected animals, predominantly cattle. The main virulence factors (genes) identified for human pathogenic VTEC are those encoding for verotoxins (VT1, VT2). VTEC causes the onset of bloody diarrhoea and has the potential to cause infections of a more severe nature including haemolytic uraemic syndrome (HUS), which can be problematic for infants often causing renal failure.<sup>1</sup>

Rapid detection and identification of VTEC is essential in order to maintain infection control and carry out vital monitoring. The majority of microbiology laboratories use traditional culture methods to perform a single screen for VTEC-O157 such as Sorbitol MacConkey Agar supplemented with cefixime tellurite (CT-SMAC). However, this can result in other serogroups and Sorbitol-fermenting-O157 strains being overlooked, with an average of 80% of cases going undetected.<sup>2</sup>

### Using novel MAST patented probe technology<sup>3</sup>

**MAST ISOPLEX® VTEC** & **MAST ISOPLEX® E.coli O157** use loop mediated isothermal amplification (LAMP) technology and come in a convenient lyophilised format.


For use with real time thermal cyclers with FAM, TAMRA and CY5 channels

**MAST ISOPLEX® VTEC** kit is optimised for typing of VTEC variants 1 and 2 (VT1 and VT2). VT1 and VT2 genes can be detected simultaneously in a triplex assay format with an inhibition control DNA (IC DNA).

For use with real time thermal cyclers with FAM and TAMRA channels

**MAST ISOPLEX® E.coli O157**. The kit contains reagents for detection of Perosamine synthetase (RfbE) expressed by *E.coli* O157. *E.coli* O157 genes can be detected simultaneously in a duplex assay format with an inhibition control DNA (IC DNA).

## LAMP assay preparation



## Benefits

### Kits contain inhibition controls

- Reducing false-negative results.

### Differentiation between VT1 and VT2

- **MAST ISOPLEX<sup>®</sup> VTEC** detects VT1 and VT2 genes simultaneously

### *E.coli* O157 gene detected

- **MAST ISOPLEX<sup>®</sup> E.coli O157** detects *E.coli* O157 gene

### Highly specific and sensitive

- Highly efficient amplification process can detect femtogram levels of target DNA

### Rapid & easy to use


- Results within 40 minutes, permitting rapid reporting and effective outbreak monitoring

### CE marked kits

- **MAST ISOPLEX<sup>®</sup> VTEC** and **MAST ISOPLEX<sup>®</sup> E.coli O157** can be applied within clinical laboratories.

### Novel probe technology enabling multiplexing

**MAST ISOPLEX<sup>®</sup> LAMP** assays rely on 6 specific primers that recognise 8 locations within a target DNA sequence. **MAST ISOPLEX<sup>®</sup> Probe technology<sup>3</sup>** consists of single gene-specific oligonucleotides labelled with a fluorophore of choice. These oligonucleotides are incorporated into LAMP primer sets against targets of interest and enable accurate, real-time detection and amplification of target DNA.



## Kit Contents

<b>MAST ISOPLEX VTEC</b>	<b>Quantity</b>	<b>MAST ISOPLEX <i>E.coli</i> O157</b>	<b>Quantity</b>
LAMP Pellets (PEL1) <b>Red Cap</b>	X2	LAMP PELLETS (PEL2) <b>Violet Cap</b>	X2
VTEC Primer and Probe Mix with Inhibition Control DNA (PP1) <b>White Cap</b>	X2	<i>E.coli</i> O157 Primer and Probe Mix with Inhibition Control DNA (PP2) <b>Blue Cap</b>	X2
Positive Control DNA <b>Green Cap</b>	X1	Positive Control DNA <b>Brown Cap</b>	X1
Reconstitution Buffer <b>Yellow Cap</b>	X1	Reconstitution Buffer <b>Yellow Cap</b>	X1
Molecular Grade Water <b>Black Cap</b>	X1	Molecular Grade Water <b>Black Cap</b>	X1

## Ordering Information

<b>Order Code</b>	<b>Product</b>	<b>Pack Size</b>
67DNALY3	MAST ISOPLEX® VTEC (DNA/LYO3)	20 Tests
67DNALY4	MAST ISOPLEX® <i>E.coli</i> O157 (DNA/LYO4)	20 Tests

### References

1. Public Health England (www.gov.uk) 2018 - Interim Public Health Operational Guidance for Shiga toxin producing *Escherichia coli* (STEC)
2. Van Duynhoven, Y.T. H. P. et al. 2008. Prevalence, characterisation and clinical profiles of Shiga toxin producing *Escherichia coli* in The Netherlands. *Clin Microbiol Infect* **14**, 437 – 445
3. Suwara MI, Javed S, Gillies EA. Nucleic acid probe with single fluorophore label bound to internal cytosine for use in loop mediated isothermal amplification. World intellectual Property Organization WO2015/063498.

**\* Licensed under International Patent application numbers:**

WO 00/28082, WO 01/34790, WO 01/34838, WO 01/83817, WO 01/77317, WO 02/24902, WO 02/103053 and corresponding patents owned by Eiken Co., Ltd., Japan in other countries.

v 2.0 AH 11/18 (MD2019-02-27)

**United Kingdom**  
**Mast Group Ltd.**  
Mast House, Derby Road, Bootle  
Merseyside L20 1EA

Tel: +44 (0)151 933 7277  
Fax: +44 (0)151 944 1332  
e-mail: sales@mastgrp.com

**Germany**  
**Mast Diagnostica GmbH**  
Feldstraße 20  
DE-23858 Reinfeld

Tel: +49 (0)4533 2007 0  
Fax: +49 (0)4533 2007 68  
e-mail: mast@mast-diagnostica.de

**France**  
**Mast Diagnostic**  
12 rue Jean-Jacques Mention  
CS91106, 80011 Amiens, CEDEX 1

Tél: +33 (0)322 80 80 67  
Fax: +33 (0)322 80 99 22  
e-mail: info@mast-diagnostic.fr