## **MAST ISOPLEX® DNA / RNA Extraction Kit**

MAST ISOPLEX® DNA / RNA Extraction Kit

**REF**

67DNAEXA

960 Tests

### **Gebrauchsanweisung / Instructions for Use / Notice d'utilisation**

**Nur zur Verwendung durch Fachpersonal /  
For professional use only /  
Exclusivement pour un usage professionnel**

	Deutsch	Seiten	02–05
	English	Pages	06–09
	Français	Pages	10–13

## MAST ISOPLEX® DNA / RNA Extraction Kit

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## MAST ISOPLEX® DNA / RNA Extraction Kit

Kit for extraction of nucleic acids (NA) for further use in MAST ISOPLEX test systems or other amplification methods.

### Intended Use

The kit is suitable for NA extraction of bacteria, viruses and protozoa.

Starting material can be liquid, human swab samples, samples from culture plates as well as stool, urine or blood samples. The kit is suitable for automatic processing on KingFisher Flex. On other liquid handling systems like Hamilton systems it can be adapted as well (please contact the corresponding application specialist).

### Important Note for Use of these Kit Instructions

Any assay relevant changes to the kit instructions for use (IFU) will lead to a change of the version number at the bottom on the last page. All changes made will be identified on a separate sheet added to the IFU for a period of three months from the date of a change of the version. Please ensure that the latest version of the IFU is used for the assay procedure.

### Principle of the Test

Bead-based methods are simple and reliable methods for the extraction of nucleic acids (NA). First, the sample is lysed. Under optimal conditions, the NA bind selectively to the surface of the magnetic beads, while other components remain in solution. Purified NA can be used directly in various molecular biological applications (sequencing, PCR, etc.). Extraction can be performed manually, semi-automatically or fully automatically.

### Kit Contents

1. **LYSIS** Lysis Solution V  
1 x 250 mL
2. **PROT K** Proteinase K  
for 16 x 1.5 mL each
3. **BINDING** Binding Solution V  
1 x 650 mL
4. **MAG** MAG Suspension F  
10 x 1.1 mL
5. **WASH 1** Washing Solution HS  
1 x 350 mL (final volume 700 mL)
6. **WASH 2** Washing Solution LS  
2 x 130 mL (final volume 650 mL each)
7. **WATER** RNase-free Water  
4 x 30 mL
8. Instructions for use

### Warning and Precautions

1. Read instructions carefully before conducting the assay. Do not modify the procedure without prior validation.
2. Do not use reagents beyond the expiry date.
3. Comply with the general health and safety guidelines for working with potentially infectious materials. Wear appropriate protective clothing and use appropriate lab facilities.
4. Do not mix reagents between different lots, as reagents have been calibrated for each batch.
5. Do not cross-contaminate reagents or interchange caps on bottles. Use a separate pipette or pipette tips for each sample and reagent.
6. Samples and contaminated disposables should be disposed according to relevant disposal directives and regulations for infectious materials. For the decontamination of surfaces use suitable surface disinfectants. Contaminated glassware should be autoclaved at 121 °C.
7. If buffer bottles are damaged or leaking, wear gloves and protective goggles when discarding the bottles in order to avoid any injuries.
8. Below the European Community risk and safety phrases for components of the MAST ISOPLEX® DNA / RNA Extraction Kit to which they apply are listed.

Binding Solution V: contains 2-propanol; highly flammable, irritant (R11, 36, 67, 7, 16, S24/25/26)

Proteinase K: irritant, sensitizing. Risk and safety phrases: R36/37/38-42; S22-24-26-37/38

Washing Solution HS: contains guanidinium thiocyanate: harmful. Risk and safety phrases: R20/21/22-32; S13-26-36-46

For details please refer to the SDS.

### Stability and Storage

#### Important

Store lyophilized **Proteinase K** at 4–8 °C!

Divide unused, dissolved Proteinase K into aliquots and store at -20 °C. Repeated freezing and thawing will reduce the activity dramatically!

The lyophilized Proteinase K may be used from manufacture until the expiry date indicated on the package, if stored at 4-8 °C.

Store the suspension of Mag F at 4-8 °C, avoid freezing!

All other components are stored at room temperature.

#### Components not included in the kit

1. RNase free tubes.
2. 96–99.8 % ethanol, non-denatured or methylated.

## Initial steps

1. Add 350 mL of 96-99.8 % ethanol to the bottle **Washing Solution HS**, mix thoroughly and keep the bottle always firmly closed!
2. Add 520 mL of 96-99.8 % ethanol to each bottle **Washing Solution LS**, mix thoroughly and keep the bottle always firmly closed!
3. Dissolve **Proteinase K** by addition of 1.5 mL of ddH<sub>2</sub>O, mix thoroughly and store as described!

## Sample Preparation

### Liquid / urine / blood samples:

No specific preparation is necessary. Samples can be used directly for extraction.

### Swab samples:

- **Test systems with liquid media:**

The liquid medium can be used directly for extraction.

- **Swab samples without liquid medium**

Transfer swabs to tubes containing 1 mL (0,9 % NaCl) physiological saline solution or PBS and incubate for 20 minutes, shaking continuously.

Squeeze out and remove swab, then 200 µL of liquid sample can be used directly for extraction.

### Agar plates:

- Pick a few (3-5) fresh overnight colonies, resuspend in 200 µL PBS and used it directly for extraction.

### Fecal samples:

- Add approx. 50-100 mg stool sample into a 1.5 mL reaction tube.
- Add 250 µL PBS (not included in kit) and vortex for 10 sec
- Spin at approx. 16000 rpm (approx. 21000 x g) for 3 min
- 200 µL supernatant can be used directly for extraction.

### Tissue samples:

- Pipette to 5-10 mg tissue sample 400 µL ddH<sub>2</sub>O (RNase-free) or PBS and homogenize (tissue homogenizer)
- Centrifuge the homogenate for 2 min at 10,000 x g.
- 200 µL supernatant can be used directly for extraction.

The prepared samples are extracted on the KingFisher Flex or manually.

## KingFisher Flex: Extraction of Nucleic Acids

### Recommended steps before starting

1. Ensure that the **Washing Solution HS**, **Washing Solution LS** and, **Proteinase K** have been prepared according to the instruction.

2. Mix **MAG Suspension F** very well by vortexing for 1 minute.
3. Avoid freezing and thawing of starting materials.
4. Add Internal control DNA / RNA according the IFU of the amplification test to Lysis Solution V.

### Extraction process

1. Transfer **200 µL Lysis Solution V** (optional with extraction / internal control) into of the Deep Well (DW) Plate labelled with "Lysis Plate".
2. Add **200 µL of the liquid sample** (sample preparation) into the DW Plate containing Lysis Solution V (optional with extraction / internal control).
3. Add **20 µL Proteinase K** to each well used.
4. Labelling and prefilling of Deep Well Plates

Plate	Label	Content
Deep Well	Lysis Plate	Lysed samples (including Lysis Solution V)
Deep Well	Washing Plate 1	500 µL Washing Solution HS
Deep Well	Washing Plate 2	500 µL Washing Solution LS
Deep Well	Washing Plate 3	500 µL Washing Solution LS
96 Plate	Elution Plate	80 µL RNase-free Water
Deep Well	Tip Plate	96 Well Tips

5. Loading Deep Well Plates to KingFisher Flex Turn on and select the protocol "Mast DNA\_RNAex" on the KingFisher Flex instrument and start the run.
6. Follow the instruction and load prefilled Deep Well Plates successive to the sample tray:
  - Tip Comb Plate
  - Elution Plate
  - Washing Plate 3
  - Washing Plate 2
  - Washing Plate 1
  - Lysis Plate (containing sample, Lysis Solution V and Proteinase K)
7. Starting the automated extraction
  - a. The automated extraction process starts with sample lysis. After sample lysis the automated run stops.
  - b. After the device has stopped, take the "Lysis Plate" out of the device and add 10 µL of well mixed MAG Suspension F and 500 µL of Binding Solution V to the lysed samples. **NOTE:** Mix the MAG Solution F well by vortexing for 1 minute.
  - c. After addition of MAG Suspension F and Binding Solution V place the "Lysis Plate" back to the KingFisher Flex and continue the extraction

process by start the KingFisher Flex again (you will find the instruction on the display of the KingFisher Flex).

**IMPORTANT NOTE** After finishing the extraction protocol, the Elution Plate contains the isolated DNA/RNA. Store the DNA/RNA under adequate conditions. We recommend storing the extracted RNA at -80 °C.

If the eluate contains carryover of magnetic particles, place the plate on a magnet or centrifuge the plate at maximum speed for 3 minutes. Pipet the supernatant with DNA/RNA into a new plate.

## Manually: Extraction of Nucleic Acids

### Recommended steps before starting

1. Ensure that the **Washing Solution HS**, **Washing Solution LS** and, **Proteinase K** have been prepared according to the instruction.
2. Mix **MAG Suspension F** very well by vortexing for 1 minute.
3. Avoid freezing and thawing of starting materials.
4. Optional: Heat the heater to 60 °C
5. Add Internal control DNA / RNA according the IFU of the amplification test to Lysis Solution V.

### Extraction process

1. Add **200 µL sample** and **200 µL Lysis Solution V** in a **RNase free tube (1.5 mL)**. Add **20 µL Proteinase K** to the sample and mix well by pipetting up and down.
2. Incubate the sample at room temperature for 15 min.  
**Note:** Mix the sample several times during incubation. No mixing will reduce the lysis efficiency.
3. Add **400 µL Binding Solution V** and **10 µL MAG Suspension F** in the reaction tube.  
Mix the sample by pipetting up and down several times. Wait for 30 sec and mix again by pipetting up and down several times.  
Place the reaction tube in magnetic rack. Wait for 20 sec and carefully remove the liquid as much as possible.  
**Note:** Do not remove the magnetic beads.
4. Take the reaction tube out of magnetic rack and add **500 µL Washing Solution HS**. Mix well by pipetting up and down several times until all magnetic beads are dissolved well.  
Place the reaction tube in magnetic rack. Wait for 20 sec and carefully remove the **liquid** as much as possible.  
**Note:** Do not remove the magnetic beads.
5. Take the reaction tube out of magnetic rack and add **500 µL Washing Solution LS**. Mix well by pipetting up and down several times until all magnetic beads are dissolved well.

Place the reaction tube in magnetic rack. Wait for 20 sec and carefully remove the liquid as much as possible.

**Note:** Do not remove the magnetic beads.

6. Take the reaction tube out of magnetic rack and add **500 µL Washing Solution LS**. Mix well by pipetting up and down several times until all magnetic beads are dissolved well.

Place the reaction tube in magnetic rack. Wait for 20 sec and carefully remove the liquid as much as possible.

**Note:** Do not remove the magnetic beads.

Remove residual alcohol by pipetting from the reaction tube and also from the lid as much as possible.

7. Take the reaction tube out of magnetic rack and air dry with open lid at room temperature until alcohol is evaporated completely (approx. 35 to 45 min) or 10–15 min at 60 °C

**Note:** Do not remove the magnetic beads.

8. Resuspend the magnetic beads in **80 µL RNase-free water** by pipetting up and down several times. Wait for 1 min and resuspend again the magnetic bead by pipetting up and down.

Place the tube into the magnetic rack and wait for 20 sec and transfer the liquid sample completely into the new RNase free tube.

**Note:** Avoid carryover of magnetic beads.

The extracted nucleic acids are ready to use for downstream application. For short-term store the extracted nucleic acids at +4 °C. For long-term storage at -20 °C is recommended.

## Assay Data

Product specifications:

1. Starting material:  
Up to 1 x 10<sup>9</sup> cells of Gram+ or Gram- bacteria.
2. Time for isolation: Approximately 60 min or approximately 45 min after the lysis step.
3. Binding capacity: > 50 µg NA.
4. The typical yield is sufficient for amplification procedures, up to 35 µg.

## Troubleshooting

Problem / probable cause	Comments and suggestions
<b>Poor lysis of starting material:</b>	
Insufficient disruption or homogenization	After lysis centrifuge lysate to pellet debris and continue with the protocol using the supernatant. Reduce amount of starting material.
<b>Little or no total RNA eluted:</b>	
Insufficient disruption or homogenization	Reduce amount of starting material. Overloading reduces yield!
<b>DNA contamination:</b>	
Too much starting material	Reduce amount of starting material.
Incorrect lysis of starting material	Use the recommended techniques for lysis of liquid media.
<b>Total RNA degraded:</b>	
RNA sample inappropriately handled or stored	Ensure that the starting material is fresh! Ensure that the protocol, especially the first steps, has been performed quickly.
RNase contamination of solutions; Receiver Tubes or other materials used.	Use sterile, RNase-free filter tips. Before every preparation clean up the pipette, the devices and the working place. Always wear gloves!
<b>Total RNA does not perform well in downstream applications:</b>	
Salt carryover during elution	Ensure that Washing Solution HS and Washing Solution LS are at room temperature. Checkup Washing Solution for salt precipitates. If there are any precipitate dissolves these precipitate by carefully warming.

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**Mast Diagnostica GmbH,**
**Mast Group Ltd.**