# 4N6

## FLOQSwabs<sup>®</sup>

### Genetics

Package insert and How-to-use guide

## 4N6 FLOQSwabs<sup>®</sup> Genetics

### Product Insert & How to Use Guide

#### INTENDED USE

The 4N6 FLOQSwabs<sup>®</sup> Genetics product line is intended for the collection and transport of swab samples for the detection and the analysis of human DNA for Genetics applications such as human identification testing and reference samples for paternity testing.  
Not intended for diagnostic use.

#### SUMMARY AND PRINCIPLES

One of the routine procedures for genetics applications involves the collection and transportation of swab samples (for example buccal swab (1), vaginal swabs cells, saliva, semen (2), blood) for human DNA profiling. This can be accomplished using **4N6 FLOQSwabs<sup>®</sup> Genetics** which are certified free of amplifiable human DNA, and detectable DNase and RNase.

4N6 FLOQSwabs<sup>®</sup> Genetics can be used to collect and transport specimens for DNA profiling and are available in two main formats:

- Genetics 4N6 FLOQSwabs<sup>®</sup> line, packaged in long tubes or peel pouches.  
Samples collected with the Genetics 4N6 FLOQSwabs<sup>®</sup> product line can be transported directly to the laboratory after proper drying. Proper swab drying prevents contaminant overgrowth, optimizing DNA stabilization. Swab samples should be stored at room temperature (20 – 25°C) or frozen for up to a maximum of 12 months (9)
- 4N6 FLOQSwabs<sup>®</sup> with an Active Drying System  
The Genetics 4N6 FLOQSwabs<sup>®</sup> product line with *Active Drying System* eliminates the need to dry the swab before transport. The swabs can be stored at room temperature (20 – 25°C) or frozen for up to a maximum of 12 months (9).

#### REAGENTS

Some product codes are available with Active Drying System, which consists of a desiccant contained in the tube stopper to absorb the humidity of the collected sample.

#### LIMITATIONS

1. Not for diagnostic use.
2. 4N6 FLOQSwabs<sup>®</sup> Genetics product line is intended for the collection and transport of swab samples for the detection of human DNA, for Genetics applications such as human identification testing, and reference samples for paternity testing.
3. The amount of specimen collected for DNA analysis significantly varies from donor to donor and in relation to collection execution mode. For the proper collection procedure, please refer to Specimen Collection instructions.
4. The product codes supplied in peel pouch format must NOT be reinserted into the peel pouch or envelopes after sample collection. Appropriate transport systems should be adopted.
5. If codes W/O Active Drying System are used for collection, the swab must be air-dried, before reclosing into its own tube or cuvette (it should take approximately 1 hour). Reclosing a wet swab into a tube or a cuvette could result in nucleic acids degradation due to microbial overgrowth.
6. If the sample has been frozen bring it back at room temperature before use
7. The use of handling and testing procedures different from the ones described in this package insert may affect product performance and analysis results.
8. The use of the product in combination with extraction and amplification kits must be always qualified by the user.

#### PRODUCT STORAGE

This product is ready for use and no further preparation is necessary. The product should be stored in its original container at 2 – 30°C until used. Do not overheat. Do not incubate or freeze prior to use. Improper storage will result in a loss of efficacy. Do not use after expiration date, which is clearly printed on the outer box, and on each individual pack.

#### PRODUCT DETERIORATION

4N6 FLOQSwabs<sup>®</sup> Genetics should not be used if 1) there is evidence of damage or contamination to the product, 2) the expiration date is passed, 3) the swab package is open, and/or the security seal is compromised.


#### MATERIALS SUPPLIED

4N6 FLOQSwabs<sup>®</sup> Genetics line is available in different formats shown below.


Regular size tip 4N6 FLOQSwabs<sup>®</sup> individually wrapped in peel pouch (product code 4520CS01)

Regular size tip 4N6 FLOQSwabs<sup>®</sup> individually wrapped in peelpouch with a 2ml cuvette (product code 4520CA)


Regular size tip 4N6 FLOQSwabs<sup>®</sup> individually wrapped in peelpouch with a 2ml cuvette and NAO<sup>®</sup> Basket (product code 4103CS01)

Regular size tip 4N6 FLOQSwabs<sup>®</sup> in a plain dry tube (product code 4508C)

Regular size tip 4N6 FLOQSwabs<sup>®</sup> in a shorter tube format (tube length 109mm) with Active Drying System (product code 4504C), 4N6 FLOQSwabs<sup>®</sup> subungueal shape in a shorter tube format (tube length 109mm) with Active Drying System (product code 40U022D).

All 4N6 FLOQSwabs<sup>®</sup> Genetics are provided with a molded breakpoint which allows easy breakage of the swab tip end into a cuvette and/or NAO<sup>®</sup> Basket.

**MATERIALS REQUIRED BUT NOT SUPPLIED**

Appropriate materials for nucleic acids extractions, amplifications and DNA analysis. Refer to laboratory reference manuals for nucleic acid amplification and analysis.

### DIRECTIONS FOR USE

Copan 4N6 FLOQSwabs<sup>®</sup> Genetics are available in the product configurations indicated in the table below:

Product Code	Description	Packaging	Sample
4520CS01	4N6 FLOQSwabs <sup>®</sup> Genetics, Regular size tip, with breakpoint at 20mm, individually wrapped in peelpouch	100 units/box 10x100 units/carton	Buccal Swabs (Saliva, cells, others)  Vaginal Swabs (Semen, Blood, others)
4520CA	4N6 FLOQSwabs <sup>®</sup> Genetics, Regular size tip, with breakpoint at 20mm, plus 2 ml cuvette individually wrapped in peelpouch	100 units/box 6x100 units/carton	
4103CS01	4N6 FLOQSwabs <sup>®</sup> Genetics, Regular size tip, with breakpoint at 20mm, plus NAO <sup>®</sup> Basket and 2ml cuvette, individually wrapped in peelpouch	100 units/box 6x100 units/carton	
4508C	4N6 FLOQSwabs <sup>®</sup> Regular size tip, with breakpoint at 20mm, in dry tube	100 units/box 10x100 units/carton	
4504C	4N6 FLOQSwabs <sup>®</sup> Regular size tip, with breakpoint at 20mm, in 109mm long tube with Active Drying System	50 units/box 6x50 units/carton	Buccal Swabs (Saliva, cells, others)
40U022D	4N6 FLOQSwabs <sup>®</sup> Genetics with Active Drying System, Subungual shape, breaking point at 20 mm, in 109mm long tube	50 units/box 6x50 units/carton	Subungual swab (blood, cells, others)

Additional product codes may be available. For updates please refer to our website: [www.copangroup.com](http://www.copangroup.com)

The proper taking of a sample from a donor is a crucial aspect in terms of the outcome of the DNA analyses. For specific instructions on sample-taking procedures, please refer to the procedures below or to the internal sample-taking procedures.

Gloves and protective clothing and eyewear should be worn when collecting and handling specimens.

When handling the swab during the sample-taking phase, the operator has to hold the swab from the handle/cap or from the side opposite the tip. Do not touch the tip of the swab.

Identify the collected sample by following the laboratory's internal procedures. Some devices are fitted with a label that makes it possible to enter data, if need be.

#### Taking oral swabs

**NOTE:** Do not use excessive force, pressure or bending when collecting swab samples from humans as this may result in accidental breakage of the swab shaft. Swab shafts often exhibit diameter changes to facilitate different sampling requirements. Swab shaft have a molded breakpoint point close to the tip designed for intentional breakage of the swab into a transport tube.

**NOTE:** The donor must not eat, drink (apart from water), brush his/her teeth, or chew gum one hour before sample collection. The swab must NOT be reinserted into the peel pouch after sample collection.

**NOTE:** for product codes with envelope, the sample MUST NOT be placed back inside the envelope after being taken.

For 4N6 FLOQSwabs<sup>®</sup> Genetics with Active Drying System, codes 4504C:

**NOTE:** for product codes in "4N6 FLOQSwabs<sup>®</sup> format with active drying system", the sample must be collected immediately after opening the aluminium envelope.

- Put on gloves, open the pouch and take out the device
- To facilitate the opening of the device hold the tube with one hand and with the thumb press the portion of the cap above the label
- Holding the device by the cap, remove the cap and place the tube in a vertical position on a clean surface avoiding further manipulation
- Holding the device by the cap, place the swab in the oral cavity placing it on the tongue. Close the mouth and hold the swab in position to moisten it with saliva for 10 seconds
- Position the swab between the cheek and the gum. Using the other hand, exercise light pressure on the cheek and rub the swab vigorously back and forth for 10 times
- Repeat step 5 on the opposite side of the oral cavity
- Insert the cap with the swab inside the test tube and close firmly until you hear a click.

**For Genetics 4N6 FLOQSwabs<sup>®</sup> codes 4520CS01, 4520CA, 4103CS01 (individually wrapped in peel pouch):**

1. Wear gloves, then open the envelope;
2. Remove the swab from the envelope by grabbing it from the side opposite the tip.
3. Place the swab inside the mouth, above the tongue. Close the mouth and keep the swab in position for ten seconds in order to wet the swab with the saliva;
4. Move the swab between cheek and gums. Helping yourself with your other hand, press down lightly on the cheek and rub the swab forcefully back and forth for ten times;
5. Repeat step 4 on the opposite side of the mouth;

Step 3

Step 4

Step 5

6. After taking the sample, let the swab air dry for 1 hour, in an environment protected against contamination.
7. A) for product code 4520CS01, after the swab is dry: by acting on the breaking point on the swab stem, break off the swab inside a test tube (not supplied) in order to take it to the sample analysis laboratory according to laboratory procedures.  
B) for product code 4520CA after the swab is dry: by acting on the breaking point on the swab stem, break off the swab inside the cuvette included in the envelope and close the test tube using the specific cap for transport to the analysis laboratory (See image).

C) for product code 4103CS01, after the swab is dry; by acting on the breaking point near the tip of the swab stem, break off the swab inside the NAO<sup>®</sup> Basket included in the envelope (previously placed in the relative test tube supplied in the package) and close the test tube on the NAO<sup>®</sup> Basket using the special cap for transport to the analysis laboratory.

**For Genetics 4N6 FLOQSwabs<sup>®</sup> code 4508C (in dry tube)**

1. Wear gloves and open the test tube;
2. Remove the swab from the test tube, holding it by its cap.
3. Place the swab inside the mouth, above the tongue. Close the mouth and keep the swab in position for ten seconds in order to wet the swab with the saliva;
4. Move the swab between cheek and gums. Helping yourself with your other hand, press down lightly on the cheek and rub the swab forcefully back and forth for ten times;
5. Repeat step 4 on the opposite side of the mouth;

Step 3


Step 4

Step 5



6. After taking the sample, let the swab air dry for 1 hour, in an environment protected against contamination;
7. Place the swab back into the test tube as shown in the following image and close firmly.

**Vaginal swab specimen collection**

For product codes 4520CS01, 4520CA, 4103CS01, 4508C

1. Put on gloves and open the peel pouch.
2. Remove the 4N6 FLOQSwabs® from the peel pouch. Collect the sample from the patient with the 4N6 FLOQSwabs® as per the following instructions. NOTE: Do not use excessive force, pressure or bending when collecting swab samples from patients as this may result in accidental breakage of the swab shaft. Swab shafts often exhibit configuration to facilitate different sampling requirements. Swab shaft have a molded breakpoint close to the tip designed for intentional breakage of the swab into a transport tube.
3. Insert the swab about 5 cm (2 inches) into the opening of the vagina.
4. Gently rotate the swab 3 times against the sides of the vagina. Withdraw the swab carefully.
5. After collection, allow the swab to air dry for 1 hour in an environment protected from contamination.


For Genetics 4N6 FLOQSwabs® codes 4520CS01:

6. using the breaking point of the swab shaft, break the swab in a cuvette (not provided) to transport to the testing laboratory or process as per laboratory procedures.


For Genetics 4N6 FLOQSwabs® codes 4520CA:

6. using the breaking point close to the tip of the swab shaft, break the swab in the cuvette included in the pouch, close the cap of the cuvette to transport to the testing laboratory. For Product Code 4103CS01: using the breaking point of the swab shaft, break the swab in the NAO® Basket included in the pouch (previously inserted into the relative provided cuvette), close the cap of the tube+NAO® Basket to transport to the testing laboratory.

For Genetics 4N6 FLOQSwabs® code 4508C (in dry tube)


6. Re-insert the swab inside the tube and close firmly

### Subungueal swab specimen collection

Product code 40U022D:

1. Wear gloves and open the aluminum pouch. While holding tightly the cap, open the tube.
2. Place the tip of the swab under fingernail and gently move the swab along the length of the nail. Repeat the collection in the opposite direction as shown in the picture.
3. After collection, replace the swab into the plastic tube and close tightly.

### Processing 4N6 FLOQSwabs<sup>®</sup> Genetics in the laboratory

Once the swab samples arrive at the laboratory, they can be processed immediately or kept at room temperature (20 - 25°C) or frozen for up to a maximum of 12 months. Storing the swabs at +4°C after collection is not recommended, because over time, humidity can give rise to accelerated DNA degradation.

Gloves and protective clothing and eyewear should be worn when collecting and handling forensic specimens.

When working with molecular methods care should be taken to prevent carry over contamination. Spatial separation of working areas and unidirectional workflow are essential to prevent amplicon carry-over.

### Nucleic acids extraction:

During sample processing when handling the 4N6 FLOQSwab<sup>®</sup> applicator, the operator must not touch the area below the breakpoint that is the area from the line to the tip of the nylon flocked swab, as this will lead to contamination of the applicator shaft invalidating the test results.

1. For Product Codes 4520CA: open the cuvette containing the 4N6 FLOQSwabs<sup>®</sup> tip, add directly into the cuvette the appropriate amount of lysing buffer, as per package insert of the extraction kit in use. Continue as per extraction and amplification kits procedures.
2. For Product Code 4103CS01: (In the case of swabs in NAO<sup>®</sup>Basket + cuvette): open the cap of the NAO<sup>®</sup>Basket containing the swabs, and add directly into the NAO<sup>®</sup>Basket the appropriate amount of lysing buffer, as per package insert of the extraction kit in use. The lysis step can proceed, keeping the swabs into NAO<sup>®</sup>Basket, as per package insert of the extraction kit in use. NOTE 1: the maximum volume of lysis buffer that can be added to NAO<sup>®</sup> Basket is 700 µl, the minimum volume recommended is 300 µl, in order to cover most of the swab surface. NOTE 2: at lysis of 70°C (maximum temperature recommended for use), some drops of liquid may come through the valves of the NAO<sup>®</sup> Basket into the cuvette below. This is not a problem as long as the swabs stay immersed in the lysis buffer. At the end of the lysis step, centrifuge the NAO<sup>®</sup> Basket and the cuvette containing the swab for at least 10000 revolutions for at least 1 minute, in order to recover all the liquid containing the lysate sample in the cuvette below. Open the cap of the tube, discard the NAO<sup>®</sup>Basket containing the dry swab and proceed with the cuvette containing the lysate sample, as per package insert of the extraction kit in use.
3. For Product Codes 4508C, 4504C, (4N6 FLOQSwabs<sup>®</sup> in tube) and 4520CS01: Using the breaking point of the swab shaft, transfer the swab to an extraction tube as per laboratory SOP. Add the appropriate amount of lysing buffer, as per package insert of the extraction kit in use. Continue as per extraction and amplification kits procedures. All Genetic swabshave a breaking point compatible with NAO<sup>®</sup>Basket so they can be used in association during the lysing step: in this case, follow the instructions given at previous point n°2.

Genetics 4N6 FLOQSwabs<sup>®</sup> devices have been tested with the following extraction methods: ChargeSwitch Forensic DNA Purification Kit by Invitrogen and Chelex 100 by Biorad (2), PrepFiler and PrepFiler Express by Life Technologies (3, 5); Nucleo Spin Tissue by Macherey Nagel, DNA IQ system by Promega and InstaGene matrix by Bio-Rad (4, 5), MO BIO BiOstic<sup>®</sup> Bacteremia DNA Isolation Kit (MO BIO Laboratories, Inc. Carlsbad, CA) (6).

The extracted DNA has been tested with the following PCR Kits: AmpFtSTR<sup>®</sup> Yfiler<sup>™</sup> PCR Amplification Kit by Life Technologies (2); Identifier Plus by Life Technologies and Quantifiler Trio Human DNA Quantification kit by Life Technologies (3); Human DNA Quantification kit by Life Technologies (4); VeriFiler<sup>™</sup> Express (VFE) PCR amplification Kit by Life Technologies (8).

### QUALITY CONTROL

Genetics 4N6 FLOQSwabs<sup>®</sup> are guaranteed for following features:

- Analyte release: the swab is inoculated with a bacterial suspension, transferred into liquid phase and check for the complete analyte release.  
≥70% of analyte must be released by the swab.
- DNase free: Following internal SOP, 4N6 FLOQSwabs<sup>®</sup> are saturated with a solution containing standard DNA ladder. DNA from solution exposed to the swabs must be intact.
- RNase free: Following internal SOP, 4N6 FLOQSwabs<sup>®</sup> are saturated with a solution containing standard RNA ladder. RNA from solution exposed to the swabs must be intact.
- Human DNA free: Following internal SOP 4N6 FLOQSwabs<sup>®</sup> are tested to be free of any amplifiable human DNA with the Identifier<sup>®</sup> Plus STR kit.
- Active Drying System ability: Following internal SOP, 4N6 FLOQSwabs<sup>®</sup> with Active Drying System are inoculated with a known amount of osmotic water. The swabs, reclosed in their own tubes, must be completely dry after 24 hours at room temperature.

#### WARNING AND PRECAUTIONS

1. Not for diagnostic use.
2. Not suitable for any other application than intended use.
3. To be used only by adequately trained or qualified personnel.
4. Directions for use must be followed carefully.
5. Before transport, where there is no active drying system, make sure that proper drying procedure has been followed. Transport of a not dried sample may result in DNA degradation.
6. Product is ready to use: it is not necessary to pretreat prior sample collection.
7. All specimens and materials used to process them should be treated in human DNA free condition, to avoid post collection contamination by personnel, that could lead to false results: wear disposable gloves, mask, cap and other protection commensurate with universal precautions not to contaminate the samples.
8. This product is for single use only; reuse may cause a risk of infection and/or inaccurate results.
9. Do not re-pack.
10. The use of this product in association with a rapid kit or with instrumentation should be previously validated by the user.
11. Using this product with extraction and amplification kits that are not listed must be validated by the user prior to use.
12. Do not use if the swab is visibly damaged (eg, if the swab tip or rod is broken).
13. Do not use excessive force or pressure when collecting swab samples from patients as this may result in breakage of the swab shaft.
14. Dispose of unused products, waste and specimens in accordance with local regulations.
15. It must be assumed that all specimens may contain infectious micro-organisms; therefore, all specimens must be handled and disposed following appropriate precautions.

4N6 FLOQSwabs<sup>®</sup> and NAO<sup>®</sup> Basket are Copan Italia trademarks.

FLOQSwabs<sup>®</sup> are PATENTED WORLDWIDE No. PCT/EP2004/003392 and in the following countries: European Patent #1608268, US Patent #US8114027B2, Canadian Patent #2515205, Japanese Patent #2007-523663, Australian Patent #2004226798, New Zealand Patent #541560, Chinese Patent #101103931.

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4N6 FLOQSwabs® and NAO® Basket are Copan Italia trademarks.
Copan Italia S.p.A.
Copan Diagnostics Inc.