Copan

HPC235R00 Date 2021.09

eNAT<sup>®</sup>

## Instructions for Use

CE 0123

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## Copan eNAT® Collection and Preservation System Instructions for use

### INTENDED USE

Copan eNAT® Collection and Transport System is intended for the collection, transport and preservation of clinical specimens to be analyzed by nucleic acids amplification techniques.

eNAT® medium stabilizes and preserves RNA/DNA for prolonged time periods and is compatible with commercial nucleic acid extraction and amplification platforms.

### SUMMARY AND PRINCIPLES

Clinical specimens stored and transported in eNAT® medium can be processed, using standard clinical laboratory operating procedures, for the detection of nucleic acids of Viruses, Bacteria, Chlamydia, Protozoa and Mycoplasma with molecular amplification assays.

The primary purpose of nucleic acids amplification techniques is to screen for a wide range of infectious diseases, so nucleic acids integrity of clinical specimens during transport and storage should be preserved <sup>(1)</sup>. eNAT® medium contains a detergent and a protein denaturant to prevent microbial proliferation, thus eNAT® is not intended to be used for culture based techniques.

### PRODUCT DESCRIPTION

eNAT® is available in the product configurations indicated in the table below:

CODE	DESCRIPTION
608CS01M	12 mm x 80 mm Screw cap tube containing 1 ml of eNAT® Transport and Preservation Media Plus One minitip applicator swab with red mark printed break point and flocked nylon fiber tip.
608CS01P	12 mm x 80 mm Screw cap tube containing 1 ml of eNAT® Transport and Preservation Media Plus One flexible applicator swab with red mark printed break point and flocked nylon fiber tip.
608CS01R	12 mm x 80 mm Screw cap tube containing 1 ml of eNAT® Transport and Preservation Media Plus One regular size applicator swab with red mark printed break point and flocked nylon fiber tip.
606CS01L	12 mm x 80 mm Screw cap tube containing 2 ml of eNAT® Transport and Preservation Media Plus One L shape applicator swab with break point and flocked nylon fiber tip.
606CS01M	12 mm x 80 mm Screw cap tube containing 2 ml of eNAT® Transport and Preservation Media Plus One minitip applicator swab with red mark printed break point and flocked nylon fiber tip.
606CS01P	12 mm x 80 mm Screw cap tube containing 2 ml of eNAT® Transport and Preservation Media plus One flexible applicator swab with red mark printed break point and flocked nylon fiber tip.
606CS01R	12 mm x 80 mm Screw cap tube containing 2 ml of eNAT® Transport and Preservation Media Plus One regular size applicator swab with red mark printed break point and flocked nylon fiber tip.

Other product codes may be available. For updates please refer to our website: [www.copangroup.com](http://www.copangroup.com)

Performance testing with eNAT® was conducted using laboratory strains. Performance testing was not conducted using human specimens.

For eNAT® codes 608CS01M; 608CS01P; 608CS01R; 606CS01L; 606CS01M; 606CS01P; 606CS01R:

Fifty (50) units are contained in a Vi-pack shelf pack and 10 x 50 units are contained in a box. Each unit consists of a sterile package containing two components: a pre-labeled polypropylene screw-cap tube with conical shaped bottom filled with 1 or 2 ml of eNAT® transport medium and a specimen collection swab which has a tip flocked with soft nylon fiber. Four main types of collection applicators are available including a tube of medium. One unit type contains a regular size flocked nylon swab applicator intended for the collection of samples from the nose, throat, mouth, vagina, rectum or wounds, the second unit type contains a minitip size flocked nylon swab applicator intended for the collection of samples from small or less accessible areas such as the eye, ear, nasal passages, urogenital tract and the third type containing a pernasal flocked nylon swab applicator intended for the collection of samples from the nasopharynx, the fourth presents an L-shaped conformation for the endocervical / exocervical sampling. These different types of swab applicators facilitate the collection of specimens from different anatomical sites on a patient. Refer to the individual product descriptions for specific information about materials supplied.

All collection swab applicators provided with eNAT® have a molded and red printed break point in the shaft of the applicator (Fig 1) except for product code 606CS01L that has only a molded break point. The red mark favors the quick and even easier identification of the presence of a break point. Refer to the product configurations described in the table for specific information on presence of red mark on shaft break point. After the sample is collected from the patient, the molded break point facilitates easy breakage of the swab applicator into the eNAT® tube of transport medium. eNAT® tube caps have an internal molded design that is able to capture the swab shaft when it is broken off into the tube and the cap is closed. The action of screwing the cap onto the tube moves the end of the broken swab shaft into a funnel shaped molded docking receptacle in the cap. This molded funnel shape effectively captures the end of the broken applicator shaft and secures it firmly in the dock by friction grip. In the testing laboratory when the swab cap is unscrewed and removed, the swab applicator is attached to the cap. This feature allows the operator to conveniently remove the swab from the transport tube due to the flexibility of the shaft of minitip and flexible swabs, the capture cap feature is not applicable, as the broken applicator may not firmly fit into the cap.


Fig.1

**REAGENTS**

Guanidine thiocyanate  
 Tris-EDTA  
 HEPES  
 Detergent

**REQUIRED MATERIALS THAT ARE NOT INCLUDED**

Appropriate materials for molecular testing according to recommended protocols as per laboratory reference manuals.

**STORAGE OF THE PRODUCT**

This product is ready for use and no further preparation is necessary. The product should be transported and stored in its original container at 5-25°C until use. Do not overheat. Do not incubate or freeze prior to use.

Improper storage will result in a loss of efficacy. Do not use after expiration date, which is clearly printed on the outer box and on each individual specimen transport vial label.

**WARNINGS AND PRECAUTIONS**

1. eNAT<sup>®</sup> medium is not for external or internal use in humans or animals.
2. This product is for single use only; reuse may cause a risk of infection and/or inaccurate results.
3. Do not use excessive force, pressure or bending when collecting swab samples from patients as this may result in accidental breakage of the swab shaft.
4. Due to the design of the flexible shaft applicators, the swab will coil when placed in the tube. Therefore it is not recommended to remove the swab from the tube. To process the specimen, collect the liquid using a sterile pipet. If the user must remove the swab, use caution and observe adequate biohazard precaution to protect the operator and the environment in case of splash.
5. Observe approved biohazard precautions and aseptic techniques. To be used only by adequately trained and qualified personnel.
6. Directions should be read and followed carefully. The manufacturer declines any liability arising from use by unauthorized non-qualified persons.
7. Not suitable for any other application than intended use.
8. This product is intended for professional use or for use in POC.
9. The use of eNAT<sup>®</sup> kit for the self- collection of vaginal samples in pregnant women has not been evaluated. In case of pregnancy the self-collection procedure is not recommended.
10. The use of the product with a rapid diagnostic kit, with diagnostic tools or biological sample processing tools must be validated by the user before use.
11. Do not use if the product is visibly damaged.
12. Do not use if the eNAT<sup>®</sup> Collection Kit peel pouch is damaged or open.
13. eNAT<sup>®</sup> medium does not require any concentration step. Pelleting procedure is not allowed since it may lead to loss of sensitivity.
14. Do not use eNAT<sup>®</sup> medium for pre-moistening or pre-wetting the sampling device prior collecting the sample or for rinsing or irrigating the sampling sites.
15. Do not use the same tube for more than one patient. This causes wrong or bad diagnosis.
16. Keep container closed when not in direct use.
17. Before transporting, make sure eNAT<sup>®</sup> screw cap tube is tightly closed.
18. eNAT<sup>®</sup> was tested for microbial viability inactivation of Gram positive bacteria, Gram negative bacteria, yeasts, molds and viruses. However universal precautions for safe handling of biological fluids should always be practiced.
19. Dispose of unused reagents, waste and specimens in accordance with local regulations.
20. Avoid contact of eNAT<sup>®</sup> medium with skin, mucous membranes. If contact does occur, immediately wash with a large amount of water.
21. eNAT<sup>®</sup> medium contains guanidine thiocyanate. Avoid direct contact between guanidine thiocyanate and sodium hypochlorite (bleach) or other highly reactive reagents such as acids and bases. These mixtures could release noxious gas.
- 22.

Danger

Contains Guanidine thiocyanate.  
 H302 Harmful if swallowed  
 H314 Causes severe skin burns and eye damage  
 H412 Harmful to aquatic life with long lasting effects  
 P264 Wash hands thoroughly after handling  
 P273 Avoid release to the environment  
 P280 Wear protective gloves/protective clothing/eye protection/face protection  
 P301+P330+P331 IF SWALLOWED: rinse mouth. Do not induce vomiting

P303+P361+P353+P310 IF ON SKIN (or hair): remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. Immediately call a POISON CENTER or doctor/physician.  
 P305+P351+P338+P310 IF IN EYES: rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do – continue rinsing. Immediately call a POISON CENTER/doctor  
 EUH032 Contact with acids liberates very toxic gas

MSDS available on request from Copan Italia S.p.A. via F. Perotti 10, 25125 Brescia Italy.

#### PRODUCT DETERIORATION

eNAT® should not be used if (1) there is evidence of damage or contamination to the product, (2) there is evidence of leakage, (3) the expiration date has passed or (4) there are other signs of deterioration.

#### INSTRUCTIONS FOR USE

##### Specimen Collection


Proper specimen collection from the patient is extremely critical for the successful identification of infectious organisms. For specific guidance regarding specimen collection procedures, consult published standard collection manuals <sup>(2,3)</sup>. Qualified medical personal only shall collect clinical specimens using proper sampling devices.

The procedures for generic specimen collection with the eNAT system and sampling of endocervical / exocervical samples by qualified personnel (product code 606CS01L) are described below. Vaginal sampling have been validated also for the self-collection procedure performed by non-pregnant patient. The choice of the proper swab for specimen collection must be done by the doctor. Qualification studies for vaginal self-collection were performed by Copan for swabs included in product part number 606CS01R and 608CS01R.

In case of self-collection of vaginal specimen, patient should follow the specific instructions included in the section: PATIENT SELF-COLLECTION OF VAGINAL SWAB SPECIMENS.

##### Swab specimen collection using eNAT® system

1. Open the eNAT® sample collection pouch and remove the tube. Take the swab out from its pouch.
2. Use the swab to collect the sample from the patient. To avoid risk of contamination, make sure that the swab tip comes into contact only with the sampling site.  
**NOTE:** Swab shafts often exhibit diameter changes to facilitate different sampling requirements. Swab shafts may also have a molded break point designed for intentional breakage of the swab into a transport tube. In all circumstances when collecting swabs samples from patients, do not use excessive force, pressure or bending of the swab as this may result in accidental breakage of the swab shaft. **Do not use eNAT® medium for pre-moisting or pre-wetting the sampling device prior collecting the sample or for rinsing or irrigating the sampling site.**
3. Unscrew and remove the cap from eNAT® tube making sure not to spill the medium.
4. Insert the swab into the tube until the breaking point reaches the level of the opening of the tube.
5. Bend and break the swab at the breaking point holding the tube away from the face.
6. Replace the cap on the tube and close tightly.
7. Write patient information on the tube label or apply patient identification label.
8. Send the sample to the test laboratory.

##### Collection of cervical specimens using Copan L-Shape Endocervical/Esocervical product code 606CS01L

1. Prepare patient in position for cervical sample collection.
2. Remove excess mucus if present from the cervical os and surrounding ectocervix using a cotton or polyester swab. Discard this swab.
3. Holding the FLOQSwabs® from the shaft, insert the swab tip (the long section of the L) inside the cervical canal, leaving outside the wider part of the swab (the short section of the L). See Fig. 1
4. Holding the swab in this position, rotate the swab 3 times counterclockwise direction and 3 times clockwise direction, making sure that the long section of the "L" is in contact with the cervical canal and the short section of the "L" rotates around the ectocervix. See Fig. 2  
**NOTE:** Do not use excessive force, pressure or bending when collecting swab samples from patients as this may result in accidental breakage of the swab shaft. Swab shafts often exhibit diameter changes to facilitate different sampling requirements. Swab shafts may also have a moulded break point designed for intentional breakage of the swab into a transport tube.
5. Remove the swab from the cervical canal and place inside the eNAT® tube, rinsing and mixing the swab in the medium to release the sample from the swab.
6. Bend the swab shaft against the side of the tube to break. Discard the broken plastic shaft.
7. Leave the swab inside the eNAT® tube and close tightly.  
**NOTE:** extracting and discarding the swab at the collection site may result in partial sample loss.
8. Label the tube with the patient ID and send the sample to the test laboratory.


Fig. 1


Fig. 2


**Patient self-Collection of Vaginal swab specimens**

**NOTE:** This procedure is only intended for supervised POC (Point of Care) patient-collection of vaginal specimens. Patient must consult her doctor, nurse or health-care provider and read instructions carefully, before specimen collection. The use of eNAT® kit for self- collection of vaginal samples in pregnant women has not been evaluated. In case of pregnancy the self-collection procedure is not recommended.

The choice of the proper swab for specimen collection must be done by the doctor. Studies for qualification of vaginal self-collection were performed by Copan for swabs included in product part number 606CS01R and 608CS01R.

Wash your hands with soap and water before starting, rinse and dry.

1. Open the eNAT® sample collection pouch and remove the tube.  
Take the swab out from its pouch. Collect the sample with the swab avoiding that the swab tip comes into contact with body areas or surface different from the sampling site.
2. Hold the swab at the red mark molded breakpoint on the shaft.  
**NOTE:** Swab shaft has a RED MARK which favors the quick and even easier identification of the presence of a break point which was designed for intentional breakage of the swab into the transport tube. Swab shaft exhibits diameter changes to facilitate different sampling requirements. **Do not use eNAT® medium for pre-moistening or pre-wetting the sampling device prior to collect the sample or for rinsing or irrigating the sampling site.**
3. Get yourself into one of the positions shown on the pictures below:
  - A. Stand with the legs apart
  - B. Sit on the toilette with your legs wide apart.
4. With your free hand, pull back the folds of skin covering the opening of your vagina and insert the swab into the vaginal opening until your fingers touch the vulva.
5. Gently rotate the swab against the walls of the vagina for 10-30 seconds. Withdraw the swab carefully.
6. Unscrew and remove the cap from tube making sure not to spill the medium. Insert the swab into the tube until the red mark molded break point reaches the level of the opening of the tube.
7. Bend and break the swab at the break point holding the tube away from your face.
8. Replace the cap on the tube and close tightly.
9. Return the collection tube to your doctor, nurse or health-care provider.

**Use in the laboratory**

**Processing eNAT® specimens for molecular testing in the laboratory.**

Specimens received in the laboratory for nucleic acid detection should be processed when received in the laboratory. In case of delay, please refer to the appropriate specimen storage conditions. eNAT® medium preserves nucleic acids for up to four weeks at room temperature and at 4°C <sup>(1)</sup> and up to 6 months at -20°C to -80°C.

Specimens preserved in eNAT® medium may need to be extracted and purified before amplification. eNAT® medium has been tested with automated platforms like, for example, NucliSENS® easyMAG® (Biomereux), Abbott m2000 system (Abbot Diagnostics), QIAAsymphony (Qiagen) <sup>(4)</sup>, Microlab NIMBUS (Hamilton) <sup>(5)</sup>, xTAG® Technology (Luminex) <sup>(6)</sup> and other manual extraction methods based on silica columns and magnetic beads. For the complete list of the extraction and purification methods tested with eNAT, contact Copan Italia Customer Care (customercare@copangroup.com).

The following steps must be performed:

1. Wear gloves and other protection commensurate with universal precautions when handling clinical specimens. Observe other CDC Biosafety Level 2 recommendations (7,8,9,10,11).
2. When working with NAAT assays, care should be taken to prevent carry over contamination. Spatial separation of working areas and unidirectional workflow are essential to prevent amplicon carry-over (12).


**VORTEXING METHOD:**

3. Vortex eNAT® specimen tube for 10 s  
**NOTE 1:** If the sample appears too mucousy, the specimen may remain mainly attached to the swab, extend the vortex time in order to break down mucus clumps and to detach the sample from the swab.  
**NOTE 2:** If some foam appears in the tube after vortexing, wait a few seconds before opening it.
4. Unscrew the cap and transfer the appropriate amount of sample (e.g. 200µl - 400µl or as per the protocol in use for the extraction) directly into the extraction buffer tube.  
**NOTE:** If using automated systems, refer to the instructions for use of the instrument.
5. Continue as per extraction and amplification kits procedures.

If unable to use the vortexing method, use the following alternative protocol:

**MANUAL SHAKING METHOD:**

1. Hold the eNAT® tube from the cap making sure that it is closed tightly.
2. Shake the tube 4 times downward with rapid movements of the wrist (see picture)  
**NOTE 1:** Inverting the tube up and down, is not recommended. If the sample appears to be viscous, it is recommended to extend the shaking times in order to break down mucous clumps and to easier release the sample from the swab
3. Unscrew the cap and transfer the appropriate amount of sample (e.g. 200µl - 400µl or as per the protocol in use for the extraction) directly into the extraction buffer tube.  
**NOTE 2:** If using automated systems, refer to the instructions for use of the instrument.
4. Continue as per extraction and amplification kits procedures.

Manual Shaking protocol has been validated using vaginal samples. The use of this protocol with other type of samples, should be previously validated by the user. Other extraction methods may also be applicable prior validation.

For updated list of validated kits, please contact Copan Italia.

**QUALITY CONTROL**

**ANTIMICROBIAL ACTIVITY**

eNAT® medium is routinely tested for its antimicrobial activity against a panel of bacteria (*E. coli*, *S. aureus* and *C. albicans*). A complete microbial viability inactivation is obtained within 30 minutes for these three strains, starting from  $\geq 10^5$  viable CFU inoculated in 1 mL of eNAT® medium.

**PERFORMANCE CHARACTERISTICS**

**NUCLEIC ACIDS PRESERVATION**

eNAT® medium preserves nucleic acids for up to 4 weeks at room temperature and at 4°C and up to 6 months at -20°C to -80°C.

Performance testing with Copan eNAT® was conducted using laboratory strains spiked onto a swab. Performance testing was not conducted using human specimens. Tested strains include among others:

ORGANISM	REF NUMBER	ANALYTE TYPE
INFLUENZA A VIRUS	ATCC® VR-822	RNA VIRUS
INFLUENZA B VIRUS	ATCC® VR 786	RNA VIRUS
CITOMEGALOVIRUS	ATCC® VR-977	DNA VIRUS
HERPES VIRUS TYPE I	ATCC® VR-539	DNA VIRUS
CHLAMYDIA TRACHOMATIS	ATCC® VR-880	CHLAMYDIA
CHLAMYDIA PNEUMONIAE	ATCC® VR-1360	CHLAMYDIA
NEISSERIA GONORRHOEAE	ATCC® 43069	GRAM NEGATIVE BACTERIUM
BORDETELLA PERTUSSIS	ATCC® 8467	GRAM NEGATIVE BACTERIUM
PSEUDOMONAS AERUGINOSA	ATCC® 27853	GRAM NEGATIVE BACTERIUM
MRSA	ATCC® 43300	GRAM POSITIVE BACTERIUM
STAPHYLOCOCCUS AERUS	ATCC® 6538	GRAM POSITIVE BACTERIUM
MYCOPLASMA HOMINIS	ATCC® 23114	MYCOPLASMA
MYCOPLASMA PNEUMONIAE	ATCC® 15331	MYCOPLASMA
TRICHOMONAS VAGINALIS	Biomed culture	PROTOZOAN

**MICROBIAL VIABILITY INACTIVATION**

eNAT® allows to inactivate microbial viability in a short time after inoculum. Starting from an initial strain concentration of  $\geq 10^5$  CFU (or IFU)/ ml, bacteria, yeasts and viruses are completely inactivated in  $\leq 30$  minutes. Molds are completely inactivated in  $\leq 1$  hour. Viability inactivation testing with Copan eNAT® was conducted using highly concentrated laboratory strains. Tested strains include:

Strain	ATCC®	Type	Starting spiked concentration	Growth after 15'	Growth after 30'	Growth after 1 hour
<i>Streptococcus pneumoniae</i>	ATCC® 6305	Gram positive bacteria	$>10^6$ CFU/ml	No growth	No growth	No growth
<i>Streptococcus agalatae</i>	ATCC® 12386		$>10^7$ CFU/ml	No growth	No growth	No growth
<i>Staphylococcus aureus</i>	ATCC® 6538		$>10^7$ CFU/ml	No growth	No growth	No growth
<i>Peptostreptococcus anaerobius</i>	ATCC® 27337		$>10^7$ CFU/ml	No growth	No growth	No growth
<i>Escherichia coli</i>	ATCC® 25922	Gram negative bacteria	$>10^7$ CFU/ml	No growth	No growth	No growth
<i>Pseudomonas aeruginosa</i>	ATCC® 27853		$>10^7$ CFU/ml	No growth	No growth	No growth
<i>Shigella sonnei</i>	ATCC® 9290		$>10^7$ CFU/ml	No growth	No growth	No growth
<i>Salmonella typhimurium</i>	ATCC® 14028		$>10^7$ CFU/ml	No growth	No growth	No growth
<i>Haemophilus influenzae</i>	ATCC® 10211		$>10^6$ CFU/ml	No growth	No growth	No growth
<i>Neisseria gonorrhoeae</i>	ATCC® 43069		$>10^8$ CFU/ml	No growth	No growth	No growth
<i>Candida albicans</i>	ATCC® 10231	fungi	$>10^5$ CFU/ml	No growth	No growth	No growth
<i>Aspergillus brasiliensis</i>	ATCC® 16404		$>10^5$ CFU/ml	140 CFU/100 uL	9 CFU/100 uL	No growth
<i>Chlamydia trachomatis</i>	ATCC® VR-880	Intracellular Gram negative bacteria	$>10^5$ IFU/ml	No IFU detected	No IFU detected	No IFU detected
Influenza B	ATCC® VR-786	RNA virus	$>10^5$ IFU/ml	No IFU detected	No IFU detected	No IFU detected
Influenza A	ATCC® VR-822		$>10^5$ IFU/ml	No IFU detected	No IFU detected	No IFU detected
Herpes simplex virus type 1	ATCC® VR-539	DNA virus	$>10^5$ IFU/ml	No IFU detected	No IFU detected	No IFU detected

CFU= colony forming unit  
IFU= infectious units

Results obtained will largely depend on proper and adequate specimen collection, as well as timely transport and processing in the laboratory.

Please refer to symbol table at the end of the instructions for use

**BIBLIOGRAPHY**

- 1) Mohammad R. Hasan, Rusung Tan, [...], and Peter Tilley. Short-Term Stability of Pathogen-Specific Nucleic Acid Targets in Clinical Samples. *J Clin Microbiol.* 2012 Dec; 50(12): 4147–4150
- 2) Clinical and Laboratory Standards Institute 2005. Collection, transport, preparation, and storage of specimens for molecular methods; approved guideline. MM13-A CLSI, Wayne, PA
- 3) World Health Organization 2019. Guidance on regulations for the transport of infectious substances 2019–2020
- 4) Evaluation of Anatomically Designed Flocked Rectal Swabs for Molecular Detection of Enteric Pathogens in Children Admitted to Hospital with Severe Gastroenteritis in Botswana. David M. Goldfarb, Andrew P. Steenhoff, Jeffrey M. Pernica, Sylvia Chong, Kathy Luinstra, Margaret Mokomane, Loeto Mazhani, Isaac Quaye, Irene Goercke, James Mahony, Marek Smieja. *Journal of Clinical Microbiology.* p.3922-3927, Volume 52, November 2014
- 5) Validation of Copan eNat, a Molecular Transport Medium, For the Collection and Preservation of Urine Specimens For the Detection of STI Infections with the Seegene Anyplex II STI-7 V1.1 Assay, S.Razeti, *BMJ Journals*, 2013.
- 6) COPAN ENAT MEDIUM STABILIZES NUCLEIC ACIDS FOR THE DETECTION OF VIRUSES WITH THE XTAG FAST RESPIRATORY VIRUS PANEL. Alice Squassina, Arnalda Giambra, Santina Castriano, Kathy Luinstra, and Marek Smieja. 26th European Congress of Clinical Microbiology and Infectious Diseases 2016
- 7) J. Michael Miller, Shelley A. Miller, 2017. *A Guide to Specimen Management in Clinical Microbiology*, Third Edition. ASM, Washington DC
- 8) Centers for Disease Control and Prevention (CDC), 2016. *Guide for Shipping Infectious Substances*.
- 9) 42CFR72. Code of Federal Regulations, Title 42, Volume 1, Part 72. October 1, 2007. Interstate Shipment of Etiologic Agents.
- 10) Centers for Disease Control and Prevention (CDC), 2009. *Biosafety in Microbiological and Biomedical Laboratories* 5th Edition. Ex 6
- 11) Hansen DJ. *Healthcare, Laboratories and Biosafety*. Vol 2., 1992. CRC Press. Ex 7
- 12) NHS. UK Standards for Microbiology Investigations- Good practice when performing molecular amplification assays. *Quality Guidance | Q 4 | Issue no: 5 | Issue date: 19.02.18.*

Index of Symbols

Symbol / Simbolo / Símbolo / Symbole / Symbole / Símbolos / Symbol / Symbol / Символ	Meaning / Significato / Signification / Bedeutung / Sens / Significado / Betyder / Betydning / Značenje / Значення
	Manufacturer / Fabricante / Fabricante / Hersteller / Fabricant / Fabricante / Producent / Produsent / Proizvođač / Виробник
	CE marking / Marchio CE / Marca CE / CE-Kennzeichnung / Marquage CE/ Marcação CE / CE-mærke / CE-merke / CE oznaka / Маркування CE
	In vitro diagnostic device / Dispositivo Diagnostico in Vitro / Dispositivo de diagnóstico in vitro / In-vitro-Diagnostikum / Dispositif de diagnostic in vitro/ Dispositivo de diagnóstico in vitro / diagnostisk medicinsk anordning in vitro / In vitro diagnostisk medisinsk utstyr / In vitro dijagnostički uređaj / Медичний виріб для діагностики in vitro
	Identification number of notified body / Numero di identificazione dell'organismo notificato / Identificación del organismo notificado / Identifizierung der benannten Stelle / Identification de l'organisme notifié / Identificação do organismo notificado / Identifikationsnummer for det bemyndigede organ / Identifikasjon for godkjenningsorganet / Ідентифікаційний номер уповноваженого органу
	Sterilized using irradiation / Sterilizzato usando radiazioni ionizzanti / Esterilizado usando radiaciones ionizantes / Sterilisiert mit ionisierenden Strahlungen / Stérilisé à l'aide de radiations ionisantes / Esterilizado por radiação ionizante / Sterilisationsmåde: Bestråling / Steriliseringsmetode med ioniserende stråling / Sterilano upotrebom jonizujućeg zračenja / Стерилізований завдяки опроміненню
	Do not reuse / Non riutilizzare / No reutilizar / Nicht zur Wiederverwendung / Ne pas réutiliser / Não voltar a usar / Brug det ikke igen / Må ikke gjenbrukes / Nemojte ponovo koristiti / Не використовувати повторно
	Catalogue number / Numero di catalogo / Número de catálogo / Bestellnummer / Référence du catalogue / Referência do catálogo / Katalognummer / Katalognr / Kataloški broj / Номер за каталогом
	Temperature limitation / Limiti di temperatura / Límites de temperatura / Temperature Begrenzung / Limites de temperature / Limites de temperatura / Temperaturgrænser / Temperaturgrenser / Opseg temperature / Обмеження температури
	Use by / Utilizzare entro / Fecha de caducidad / Verwendbar bis / Utiliser jusque / Prazo de validade / Anvendes før / Må brukes innen / Uptrebiti do / Використати до
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultar as instruções de utilização / Se vejledningen til brug / Se instruksjoner for bruk / Pročitajte uputstvo za upotrebu / Користуватись інструкціями по використанню
	Batch code (Lot) / Codice del lotto (partita) / Código de lote (Lote) / Chargencode (Chagenbezeichnung) / Code de lot (Lot) / Código do lote (Lote) / Serienummer (parti) / Batch-nummer (parti) / Serijski (lot) broj / Код серії (partii)
	Contains sufficient for <n> tests / Contenido suficiente per <n> test / Contenido suficiente para <n> pruebas / Ausreichend für <n> Tests / Contenu suffisant pour <n> tests / Contémo suficiente para <n> testes / Indhold tilstrækkeligt til <n> prøver / Innhold tilstrekkelig for <n> test / Sadržaj dovoljan za   testova / Придатний для проведення <n> кількості випробувань
	Caution / Attenzione / Precaución / Vorsicht / Attention / Cuidado / Forsigtig / Forsiktig / Pažnja / Увага

**Copan**

**Copan Italia S.p.A.**  
Via F. Perotti, 10  
25125 Brescia, Italy  
Tel +39 030 2687211  
Fax +39 030 2687250

Email: [info@copangroup.com](mailto:info@copangroup.com)  
Website: [www.copangroup.com](http://www.copangroup.com)

**North American Distributor:**

Copan Diagnostics Inc.  
26055 Jefferson Avenue  
Murrieta, CA 92562, USA  
Tel: 951-696-6957  
Fax: 951-600-1832

E-mail: [customerservice@copanusa.net](mailto:customerservice@copanusa.net)  
Website: [www.copanusa.com](http://www.copanusa.com)