# MASTSORB™

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	<b>Content</b>	<b>Page</b>
1.	Intended Use	8
2.	Introduction	8
3.	Principle of the test	8
4.	Kit Contents	8
5.	Warnings and Precautions	9
6.	Storage and Stability	9
7.	Assay Procedure	10
8.	MMASTSORB™ Reagent Performance	11

## 1. Intended Use

The MASTSORB™ reagent has been designed for the elimination of rheumatic factor in serum and plasma samples for the determination of IgM and/or IgA.

This assay is intended for *in vitro* diagnostic use only.

All laboratory test results should be interpreted in conjunction with other clinical data. The clinical judgement and further tests have to be taken into account additionally.

## 2. Introduction

The detection of infectious reagent specific IgM antibodies in patient samples is a common method for the early detection of numerous bacterial and viral infections. The presence of specific IgG and rheumatic factors of the IgM class may interfere with the IGM and/or IgA assay therefore these factors have to be eliminated before IgM/IgA is determined.

Binding of IgM can be reduced if IgG-antibodies of the same specificity are present at high levels in a patient sample. IgG antibodies compete effectively with IgM antibodies for the binding sites (competitive inhibition) and may lead to negative IgM results if IgM is present at low levels only.

Beside the antigen specific IgG also rheumatic factors (RF) of the IgM class may be present. These IgM-RF use antigen specific IgG as a substrate and form immune complexes with the antigen-bound IgG. This complex reacts as a substrate for enzyme- or fluorochrome labelled  $\mu$ -chain-specific IgM conjugate and results false positive for the IgM in the sample.

MASTSORB™ reagent eliminates most of the above mentioned interfering factors.

## 3. Principle of the Test

All IgG subclasses are bound and precipitated by the anti-human IgG antibodies. Human serum IgG in concentrations of up to 15 mg/mL and rheumatoid factors are removed by the absorbent (average serum IgG concentration in adults: 12 mg/mL). While IgG is precipitated recovery of IgM after treatment with MASTSORB™ is almost 100 %

## 4. Kit Contents

The kit contains sufficient reagents for 40 determinations in IFA or 80 determinations in ELISA technique.

1 x 2 mL MASTSORB™ Goat anti-human-IgG, ready to use

## 5. Warnings and Precautions

- For *in vitro* diagnostic use only! Do not ingest or swallow! Laboratory safety precautions should be followed. Do not eat, drink or smoke in the laboratory.
- All sera and plasma or buffers based upon have been tested to HBsAg, HIV and HCV respectively with generally accepted methods and were found negative. Nevertheless, precautions like the use of latex gloves have to be taken.
- Serum and reagent spills have to be wiped off with a disinfecting solution (e.g. sodium hypochlorite 5 %) and have to be disposed of properly.
- All reagents have to be brought to room temperature (18–24 °C) before performing the test.
- Before pipetting all reagents should be mixed thoroughly by gentle tilting or swinging. Vigorous shaking with formation of foam should be avoided.
- It is important to pipet with constant intervals so that all the wells have the same conditions.
- When removing reagents out of the bottles care has to be taken that the stoppers are not contaminated. Further a possible mix-up has to be avoided. The content of the bottles is usually sensitive to oxidation so that they should be opened only for a short time.
- In order to avoid a carry-over or a cross-contamination separate disposable pipet tips have to be used.
- No reagents from different kit lots should be used and they should not be mixed with one another.
- All reagents have to be used within shelf life.
- In accordance with a Good Laboratory Practice (GLP) or following ISO 9001 all laboratory devices employed should be regularly checked regarding the accuracy and precision. This refers e.g. to microliter pipets and washing or reading (ELISA Reader) instrumentation.
- The contact of certain reagents especially the stopping solution and the substrate with skin, eye and mucosa has to be avoided because possible irritations and acid burns could arise and there exists a danger of intoxication.

## 6. Storage and Stability

Store all reagents at 2–8 °C.

The expiry date of the reagent is printed on the label. Do not use the reagents after the expiry date has been exceeded.

## 7. Assay Procedure

Sera containing higher levels of IgG than 15 mg/mL should be treated with a proportional amount of MASTSORB™ reagent.

### 7.1. Preparation of Reagent for IFAs

Allow all reagents and specimens to reach room temperature (RT, 18–24 °C) prior to use and mix well.

- Strictly follow the instructions for reliable test performance. Any changes or modifications are within the responsibility of the user.
- All reagents and samples must be brought to room temperature before use, but should not be left at this temperature for longer than necessary.

#### Assay Steps

Prepare a sufficient amount of tubes for the samples.

Note: Other incubation conditions might be possible. In case of modifications of the recommended test procedure (e.g. incubation temperature 37 °C instead of RT) the user has to validate assay performance. Other volumes can be used if the ration of reagents remains the same.

1. Pipette 30 µL of PBS (pH 7.2) and the ready to use MASTSORB™ into the appropriate tubes, mix thoroughly.
2. Add 20 µL of neat patient sample, mix well
3. Incubate at RT for 30 minutes or over night at 2–8 °C (preferred).
4. Centrifuge the tubes for 5 min at 2000 rpm. The supernatant is used for the determination of IgM or IgA (sample dilution is 1:5) while the precipitate is discarded.

### 7.2. Preparation of Reagent for ELISA

Allow all reagents and specimens to reach room temperature (RT, 18–24 °C) prior to use and mix well.

- Strictly follow the instructions for reliable test performance. Any changes or modifications are within the responsibility of the user.
- All reagents and samples must be brought to room temperature before use, but should not be left at this temperature for longer than necessary.

#### Assay Steps

Prepare a sufficient amount of tubes for the samples.

Note: Other incubation conditions might be possible. In case of modifications of the recommended test procedure (e.g. incubation temperature 37 °C instead of RT) the user has to validate assay performance. Other volumes can be used if the ration of reagents remains the same.

1. Pipette 25 µL of ELISA sample diluent and the ready to use MASTSORB™ into the appropriate tubes, mix thoroughly.
2. Add 10 µL of neat patient sample, mix well
3. Incubate at RT for 30 minutes or over night at 2–8 °C (preferred).
4. Centrifuge the tubes for 5 min at 2000 rpm. The supernatant is used for the determination of IgM or IgA (sample dilution is 1:5) while the precipitate is discarded.
5. Add 25 µL of the supernatant to 400 µL ELISA sample diluent, mix well
6. 100 µL of the reagent mix (final serum dilution 1:102) is used for the ELISA assay.

## 8. MASTSORB™ Reagent Performance

Nephelometric IgG-Determination				EBV-IgG with IFA		
Sample	Neat Serum (mg/mL)	After Ion Exchange-Chromatography (mg/mL)	After MASTSORB™ Reagent (mg/mL)	Neat Serum (Titer)	After Ion Exchange-Chromatography (Titer)	After MASTSORB™ Reagent (Titer)
24	8.80	0.40	0.00	< 10	< 10	< 10
25	19.00	0.66	0.38	1280	10	< 10
26	15.00	0.19	0.18	320	< 10	< 10
27	11.00	0.38	0.00	640	< 10	< 10
28	11.30	0.70	0.38	320	10	< 10
29	11.20	0.87	0.10	320	10	< 10
30	9.20	0.60	0.08	320	>10 but < 60	< 10
31	13.20	0.40	0.15	160		< 10
32	13.80	0.76	0.15	160		< 10
x	12.50	0.55	0.16			
1S	3.17	0.22	0.14			
Recovery		95.6 %	98.7 %			
(+ or -)		1.8 %	1.1 %			

The table above shows the results of a comparison between Ion Exchange Chromatography and MASTSORB™ treatment to reduce IgG in serum

Radial immune diffusion for IgG subclasses showed that treatment with MASTSORB™ eliminates all kinds of IgG subclasses. MASTSORB™ does not show any reactivity with antigens of EBV, HSV, CMV, T. pallidum, B. burgdorferi, T gondii or C. trachomatis.

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