MASTAFLUOR™ Toxoplasma

MASTAFLUOR™ Toxoplasma Screen	REF	631181	UDI-DI: 4250729700088	10 x 5 Tests
MASTAFLUOR™ Toxoplasma Screen	REF	631182	UDI-DI: 4250729700095	10 x 10 Tests
MASTAFLUOR™ Toxoplasma IgG	REF	631183	UDI-DI: 4250729700101	10 x 5 Tests
MASTAFLUOR™ Toxoplasma IgG	REF	631184	UDI-DI: 4250729700118	10 x 10 Tests

**Gebrauchsanweisung /
Instructions for Use /
Notice d'utilisation**

**Nur zur Verwendung durch Fachpersonal /
For professional use only /
Exclusivement pour un usage professionnel**

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MASTAFLUOR™ Toxoplasma

Immunofluorescence assay for the detection of antibodies to *Toxoplasma gondii* in human serum.

Intended Use

Semi-quantitative (titratable) immunofluorescence assay for the detection of IgG antibodies (and additionally IgM antibodies using the Screen assay) directed against *Toxoplasma gondii* tachyzoites in human serum as an aid to diagnosis of immunity against toxoplasmosis (presence of specific antibodies).

The assay is suitable for manual use and is intended for professional in-vitro diagnostic use only. All laboratory test results should be interpreted in conjunction with other clinical data. The clinical judgement and further tests have to be taken into account additionally.

Important Note for Use of these Kit Instructions

Any assay relevant changes to the kit instructions for use (IFU) will lead to a change of the version number at the bottom on the last page. All changes made will be identified on a separate sheet added to the IFU for a period of three months from the date of a change of the version. Please ensure that the latest version of the IFU is used for the assay procedure.

Principle of the Test

Diluted specimens, the negative and positive controls or calibrators are applied to the wells of microscope slides and incubated. All wells are coated with purified antigens of the very pathogen. If specific antibodies are present in the serum they will bind to the fixed pathogen antigens forming a stable antigen-antibody complex. Slides are then washed to remove any unbound material. Complexed antibodies are detected by the addition of a fluorescence labelled anti-human IgG, IgM, or IgA immunoglobulin conjugate. After a further washing step to remove any unbound conjugate, slides are viewed under a fluorescence microscope.

A green fluorescence is observed if pathogen specific antibodies are present in the sample material.

Kit Contents

	100 Tests	50 Tests
1. SLIDE Slides coated with <i>Toxoplasma gondii</i>	10 x 10 Wells	10 x 5 Wells
2. CONTROL+ Positive Control IgG positive serum, human, ready-to-use, contains < 0,1 % NaN ₃	0,5 mL	0,5 mL
3. CONTROL- Negative Control IgG and IgM negative serum, human, ready-to-use, contains < 0,1 % NaN ₃	0,5 mL	0,5 mL
4. CONJ G FITC-conjugate IgG FITC-labelled anti-human IgG (γ-chain) ready-to-use, contains < 0,1 % NaN ₃ or CONJ POLY FITC-conjugate Screen FITC-labelled anti-human IgG/IgM/IgA ready-to-use, contains < 0,1 % NaN ₃	3 mL	2 mL
5. DIL 5x Diluent 5 x concentrated, to be diluted 1:5 before use	2 x 10 mL	1 x 10 mL
6. EVANS BLUE Evans Blue ready-to-use, contains < 0,1 % NaN ₃	3 mL	3 mL
7. MOUNTING MEDIUM Mounting Medium ready-to-use, contains < 0,1 % NaN ₃	3 mL	3 mL
8. BUFFER PBS Dissolve 1 sachet in 1 L distilled water, pH 7,2 ± 0,2	2 Sachets	2 Sachets

9. Instructions for use

Further Abbreviations

- RTU** ready to use

Materials Required but not Provided

- Sterile test tubes.
- Micropipettes and tips.
- Staining dish or Coplin jar.
- Moist chamber.
- Volumetric flask for PBS.
- Distilled or deionised water.
- Forceps.
- Cover slips.
- Wash bottle.
- Fluorescence microscope with a filter combination suitable for FITC (e.g. 490 nm excitation filter and a 510nm barrier filter).

Warning and Precautions

1. The reagents supplied in this kit are for *in vitro* diagnostic use only.
2. Read instructions carefully before conducting the assay. Do not modify the procedure without prior validation.
3. Do not use reagents beyond the expiry date.
4. Comply with the general health and safety guidelines for working with potentially infectious materials. Wear appropriate protective clothing and use appropriate lab facilities.
5. Do not mouth pipette.
6. Use disposable plasticware where ever possible. Re-usable glassware should be washed thoroughly and rinsed free of detergents before use.
7. The control/calibrator materials of human origin provided have been tested for antibodies to HIV, HCV and for HBsAg and found to be negative. However, they should be treated as potentially infectious materials and capable of transmitting diseases. No guarantee is given that the samples are free of infections or microbial contamination.
8. Do not mix reagents between different lots, as reagents have been calibrated for each batch.
9. Do not cross-contaminate reagents or interchange caps on bottles. Use a separate pipette or pipette tips for each sample and reagent.
10. Do not allow wells to dry out during the assay procedures.
11. Do not expose slides to intense sunlight or similar adverse conditions while incubating.
12. Samples and contaminated disposables should be disposed according to relevant disposal directives and regulations for infectious materials. For the decontamination of surfaces use suitable surface disinfectants. Contaminated glassware should be autoclaved at 121 °C.
13. Do not expose the FITC-conjugate at any stage to strong sunlight, UV or fluorescent light. Keep in a dark place whenever possible.
14. Microbial contaminated samples should not be used.
15. Do not use haemolysed, lipaemic or icteric serum.
16. Sodium azide is used as a preservative as marked. It may be toxic if ingested. Sodium azide may react with lead and copper plumbing to form highly explosive salts. Always dispose of by flushing to drain with plenty of water.

Stability and Storage

MASTAFLUOR™ Toxoplasma can be used until the end of expiry date as indicated on kit label. All kits components and reagents should be stored at 2–8 °C and brought to room temperature before use.

After first opening reagents are stable for at least 3 months. An opened slide shall be used the same day.

Reconstituted PBS powder should be stored at 2–8 °C for up to 30 days.

Samples (serum) may be stored according to general recommendation in literature. In general samples can be kept at 2–8 °C for up to 3 days prior to use or for longer storage at -20 °C. Repeated freezing and thawing of specimens should be avoided.

Sample Material

- a) Serum
- b) Samples can be stored at 2–8 °C for up to 3 days. Serum should be aliquoted immediately after sampling and kept at -20°C for longer storage. The samples should not be frozen and thawed repeatedly. After thawing samples should be briefly vortexed carefully before being used in the assay. Lipaemic, haemolytic, icteric or bacterially contaminated samples can cause false positive or false negative results.

Test Procedure

1. Allow all materials to reach room temperature (at least 20 °C) prior to use.
2. Reconstitute one sachet of PBS powder with 1 litre of distilled water or deionised water. It is recommended to use the diluent to prepare the sample dilutions. Dilute the 5x concentrated diluent 1:5 with distilled or deionised water (e.g. 10 mL of diluent + 40 mL of water).
3. Dilute all serum specimens with diluent or PBS buffer as indicated below:
IgG assays: 1:64
e.g. 10 µL serum + 630 µL diluent or PBS
Screen assays: 1:64
e.g. 10 µL serum + 630 µL diluent or PBS
4. All kit controls are ready-to-use and should not be diluted like specimen.
5. Carefully remove the slides from the pouch. Do not touch wells.
6. Apply 20–25 µL of ready-to-use controls and diluted specimens to respective wells on the slides. Ensure that all the wells are covered. Do not scratch the surface of the wells.
7. Add slides to a moist chamber and incubate at room temperature
IgG assays 30 min
Screen assays 30 min
8. After incubation wash slides carefully with PBS from a wash bottle, taking care not to direct the jet onto the test wells. This may be done by directing the jet of PBS along the center of the slide, tilting the slide first towards wells 1–5 and then towards 6–10.
9. Immerse the slides in a staining dish or Coplin jar containing PBS and wash for 15 min at room temperature with a change of PBS after 5 min to increase washing strength.
10. Remove slides from the staining jar and drain off any excess buffer. Using a blotter, dry the area outside wells. Do not touch the well surface.

11. Immediately transfer the slide to a moist chamber and add 20–25 µL of the respective FITC-conjugate to each well. The conjugate shall cover the whole well. Do not allow the wells to dry.
12. Incubate slides at room temperature in a covered moist chamber in the dark for 30 min.
13. After incubation wash the slides with PBS as in steps 9 and 10.

Optionally: Add 3 to 5 drops of Evans Blue for counterstaining to the last wash step. Counterstaining may improve differentiation between negative and positive reactions.

After the Evans Blue counterstain dip the slide in fresh PBS to remove excess Evans Blue.

14. Apply one small drop of Mounting Medium to each well on a slide. Place a cover slip on the slide and read the results immediately. Take care not to trap air pockets under the cover slip.

Note: Remove excess Mounting Medium with a paper towel avoiding any direct movement of the cover slip.

Excess Mounting Medium on a slide may result in high background fluorescence due to light scattering.

15. Examine reactions under a fluorescence microscope at a total magnification of 400x or at 800x.

Interpretation of Results

Controls have to fit the lot specific reactions stated in the QC certificate. The figures in the table below are for orientation only and do not reflect the data of the current lot.

Specimen	Expected Results
Positive control	3 + bis 4 +
Negative Control	negative

For a correct interpretation the results should be compared with positive and negative controls.

If further dilution on the samples or controls are requested on the basis of the screening titer, this can be achieved by serial dilution with PBS on screening titer.

Interpretation of Specimen Results

IgG and Screen Positive Samples: Positive sera show a homogenous green fluorescence of the cytoplasmic membrane of *Toxoplasma gondii* tachyzoites. Cytoplasm may show green staining or even fluorescence together with the cytoplasmic membrane. Cytoplasmic reactions alone must be read as negative.

Negative Samples: The *Toxoplasma* membrane does not show any fluorescence. Fluorescence of cytoplasm only is considered as negative.

Results are defined as positive:

IgG antibody: ≥ 1:64

Screen / polyvalent: ≥ 1:64

For result interpretation recommendations of (inter-) national reference centers should be followed (e.g. www.rki.de).

Limitations / Interferences

1. Autoantibodies (ANA type) may lead to cytoplasmic fluorescence of *Toxoplasma gondii*.
2. The MASTAFLUOR™ *Toxoplasma* test has been demonstrated to be highly sensitive and specific. Like for other lab based assays final diagnosis shall not be made on this assay alone, always interpret the data in connection with the patient's medical history and the clinical picture.

Assay Data

A panel of clinically positive and negatively evaluated samples was selected to determine the performance data.

The MASTAFLUOR™ *Toxoplasma* obtained the following results in this panel:

	IgG	Screen
Sensitivity	90,57 %	95,61 %
Specificity	94,12 %	94,87 %
Positive predictive value	97,96 %	98,20 %
Negative predictive value	76,19 %	88,10 %
Efficiency	91,43 %	95,42 %

No differences in intra-assay and inter-assay variance was found.

Reporting serious incidents

All serious incidents that have occurred in connection with the device must be reported to the manufacturer and the competent authority of the Member State in which the user and / or patient is established.

References

You can find the references at the end of the instruction for use.

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**Mast Diagnostica GmbH,**
Mast Group Ltd.