MASTAZYME™ TG / TPO

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Contents	Page
1. Intended use	10
2. Principle of the Test	10
3. Kit Contents	10
4. Materials Required but not Provided	11
5. Warnings and Precautions	11
6. Storage and Stability	12
7. Specimen Collection and Handling	12
8. Assay Procedure	12
9. Results and Interpretation	13
10. Assay Performance	14
11. References	14

1. Intended use

MASTAZYME™ TG and -TPO ELISA have been designed for the detection and quantification of specific antibodies against TG and TPO respectively in serum and plasma.

This assay is intended for *in-vitro* diagnostic use only.

All laboratory test results should be interpreted in conjunction with other findings. Clinical judgement and additional tests should also be taken into account.

2. Principle of the Test

The ELISA test can be described in four stages:

2.1 Serum incubation

Specific antibodies bind to antigen on the solid phase to form a stable immune complex. After incubation for 30 minutes at room temperature the wells are washed with a prediluted wash buffer to remove all non-reactive serum components.

2.2 Conjugate incubation

The anti-human-IgG horseradish peroxidase conjugate is added to all wells. The conjugate binds to IgG antibodies on the solid phase antigen to form a stable sandwich. After incubating for 30 minutes at room temperature the excess conjugate is removed by washing all wells with wash buffer.

2.3 Substrate reaction and stopping

The TMB substrate is dispensed into each well and the peroxidase enzyme/substrate reaction forms a stable blue chromogen. The reaction and consequent colour development is stopped after a 15 minute incubation at room temperature by adding 0.5 N H₂SO₄ to the wells. The shift in pH also causes the chromogen to change colour from blue to yellow.

2.4 Reading and interpretation

The colour intensity is read in a microtiter plate reader at 450 nm (recommended reference wavelength for bichromatic measurement: 600–690 nm). The intensity of the colour (OD) is directly proportional to the concentration of the specific antibody in the patient sample.

The results can be read from a calibration curve or with an electronic graphing package using a 4 parameter sigmoidal curve.

3 Kit Contents

The kits contains sufficient reagents for 12 x 8 = 96 determinations. Kit components should be stored at 2–8 °C. The expiry date is shown on the individual labels.

12	Microtiter strips	Single strips each with 8 break-apart wells coated with purified TG or TPO respectively.
1 x	Frame holder	
4 x 1.5 mL	Calibrators 1–4	Stabilized human plasma containing antibodies against the above antigens and pre-diluted in buffer. Ready to use (RTU).

	Anti-TG Concentration [IU/mL] as per NIBSC 65/93	Anti-TPO Concentration [IU/mL] as per NIBSC 65/387
Calibrator 1	20	5
Calibrator 2 (Cut-off)	100	60
Calibrator 3	1000	250
Calibrator 4	8000	2000
Positive control	300	500
Negative control	50	20

1 x 1.5 mL	Negative control	Stabilized human plasma, RTU, can be integrated into the calibration curve.
1 x 1.5 mL	Positive control	Stabilized human plasma, RTU, can be integrated into the calibration curve.
2 x 60 mL	Sample diluent	PBS solution with preservative, ready to use.
1 x 12 mL	Enzyme conjugate solution	HRP-labelled goat anti-human-IgG, ready to use.
1 x 12 mL	TMB substrate	3,3',5,5' Tetramethylbenzidine, ready to use.
1 x 12 mL	Stopping solution	0.5 N sulfuric acid, ready to use.
2 x 50 mL	Wash buffer 10 x concentrated	PBS/Tween buffer solution 10x concentrated. To be diluted 1:10 prior to use; the concentrate should be warmed up to 37 °C for 15 min to dissolve any crystals

4. Materials Required but not Provided

- 5 µL-, 100 µL- and 500 µL micropipettes and a multichannel pipette (optional).
- Microtiter plate reader with a 450 nm filter (reference filter 600–690 nm)
- Microtiter plate washer (if washing manually: wash bottle)
- Reagent tubes for serum dilution
- Measuring cylinder
- Distilled water or water of similar or better quality

5. Warning and Precautions

- For *in vitro* diagnostic use only! Do not ingest! Laboratory safety precautions should be followed. Do not eat, drink or smoke in the laboratory.
- All sera, plasma and buffers containing human biological material were found to be negative for Hepatitis B, C and HIV. Nevertheless, precautions such as the use of latex gloves must be taken.
- Serum and reagent spills should be cleaned with a disinfecting solution (e.g. sodium hypochlorite 5 %) and the refuse disposed of appropriately.
- All reagents should be brought to room temperature (18 to 24 °C) before starting the test procedure.
- Before pipetting, all reagents should be mixed thoroughly by gentle agitation. Vigorous shaking leading to formation of foam should be avoided.
- It is important to pipette with constant time intervals, so that all the wells of the microtiter plate have the same reaction conditions.
- When pipetting reagents out of the bottles care has to be taken that the stoppers are not contaminated. Furthermore, disposable pipette tips are strongly recommended to avoid cross-contamination. The contents of the bottles are usually sensitive to oxidation so should be opened only for a short time.
- To avoid cross contamination single use pipet tips should be used
- No reagents from different kit batches should be used and reagents should not be mixed with one another.
- All reagents should be used within the listed shelf life.

- In accordance with Good Laboratory Practice (GLP) or following ISO 9001 all laboratory equipment should be checked regularly regarding accuracy and precision. This refers e.g. to microliter pipets and washing or measurement (ELISA Reader) instrumentation.
- Certain reagents (especially the stop solution and substrate) are irritants- avoid contact with skin, eyes and mucosal membranes. In case of accident, rinse with water and seek medical attention. Clean all equipment after use to avoid secondary contact incidents.

6. Storage and Stability

Store all reagents at 2–8 °C.

The expiry date of each reagent is printed on the individual labels. Do not use any reagents after the expiry date has been exceeded.

The diluted wash buffer is stable for up to 4 weeks when stored at 2–8 °C. However ensure that this is at room temperature before testing.

Kits should be used within three months of opening.

7. Specimen Collection and Handling

Both serum and plasma (EDTA, heparin) can be used for testing. Samples can be stored at 2–8 °C for up to 3 days, although for longer storage periods they should be aliquoted and kept at -20 °C. Repeated freezing and thawing is contra-indicated. Lipaemic, haemolytic, icteric and bacterially contaminated samples can cause false positive or false negative results.

Patient sera should be prediluted 1:101 in sample diluent (e.g. 5 µL serum + 500 µL sample diluent) prior to testing.

8. Assay Procedure

8.1. Preparation of Reagents

Allow all kit components and specimens to reach room temperature (RT, 18–24 °C) prior to use and mix well.

Washing buffer: Dissolve any crystals which may be in the bottle by warming to 37 °C and then mix well.

Dilute the concentrated washing buffer 1:10 with distilled water (e.g. 50 mL buffer concentrate + 450 mL distilled water). Mix thoroughly.

- Strictly follow the instructions for reliable test performance. The user takes sole responsibility for any modifications to the test procedure.
- All reagents and samples must be brought to room temperature before use, but should not be left at this temperature for longer than necessary.
- A standard calibration curve should be established with each assay for quantitative results.
- Replace any unused microtiter strips in the resealable aluminium bag provided and store under dry conditions at 2–8 °C.

8.2. Assay Steps

Prepare a sufficient amount of microtiter wells for the calibrators, controls and samples.

Note: Other incubation conditions are possible. However in case of modifications to the recommended test procedure (e.g. an incubation temperature of 37 °C instead of RT) assay performance should be validated by the user.

1. Pipette 100 µL each of the diluted (1:101) samples and the ready to use calibrators and controls into the appropriate wells.

2. Incubate the test plate at room temperature for 30 minutes.
3. Discard the contents of the microwells and wash 3 times with 300 µL of diluted wash buffer. Afterwards remove residual wash buffer by gentle tapping of the microtiter plate onto a paper towel.
4. Pipette 100 µL of enzyme conjugate solution into each well.
5. Incubate the test plate for 30 minutes at room temperature.
6. Discard the contents of the microwells and wash 3 times with 300 µL of diluted wash buffer. Afterwards remove residual wash buffer by gentle tapping of the microtiter plate onto a paper towel.
7. Pipette 100 µL of TMB substrate into each well.
8. Incubate for 15 minutes in the dark (e.g. drawer) at room temperature.
9. Add 100 µL of stopping solution to each well.
10. After thorough mixing and wiping the bottom of the plate, read the optical density at 450 nm and calculate the results. Blank against air. A bichromatic measurement using a reference wavelength of 600–690 nm is recommended. Concentrations can be plotted with an electronic graphing package or by hand against the calibration curve.

The developed colour is stable for at least 30 minutes. Read optical densities during this time.

9. Results and Interpretation

Qualitative

Calibrator 2 is used for qualitative interpretation of the assays, depending on specificity requirements to include or eliminate borderline results. A quantitative interpretation is not possible when only using calibrator 2 alone.

Quantitative

The results can be calculated using a computer and suitable software program. In this case the results are calculated using a 4 parameter fit algorithm to produce a sigmoidal curve.

Optional: A standard curve is plotted by entering the mean absorbance of each calibrator on the Y-axis and the corresponding concentration on the X-axis using either log/log or linear/log graph paper. The concentration of the controls and specimen can then be read directly from the graph.

Suggested normal values

In a normal range study using serum samples from healthy blood donors and disease-state sera the following normal and elevated ranges have been established. It is recommended that each laboratory should establish its individual normal ranges based on results obtained from the local population.

Anti-TG

negative: < 100 IU/mL
positive: > 100 IU/mL

Anti-TPO

negative: < 60 IU/mL
positive: > 60 IU/mL

Interpretation

Results can be affected by several patient factors. Clinical diagnosis and disease prognosis should not be based on a single test result. Other lab data and clinical observations must be taken into consideration before a conclusive diagnosis can be made.

Interpretation

Autoantibodies to thyroglobulin or thyroid peroxidase (or both) are highly predictive for thyroiditis. Adenocarcinoma may lead to unspecific antibody response which cross-reacts with antigens coated on the solid phase. Thyroid hormone levels should be checked to avoid false-positive results due to hyperthyroidism.

Anti-thyroid antibodies generally vary individually according to age and sex, but only reach cut-off defining levels in case of disease. Cut-off levels should be adjusted according to individual findings if necessary.

There is a strong correlation between purified TPO and TMA (thyroid microsomal antigen) antibody levels. Results obtained with highly purified TPO are more specific than TMA for the diagnosis of thyroiditis (sensitivity 90–100 %). Due to ultrapure TPO antigen no cross-reactivity was observed with anti-myeloperoxidase (so called pANCA) autoantibodies.

Autoimmune based thyroid diseases usually stay at a subclinical level which does not reflect the actual prevalence of autoantibodies directed against thyroid antigens. The numbers of positives is significantly higher than those which can be expected from hospitalised patients.

Summary

For complete thyroid autoimmune diagnostics neither anti-TG nor anti-TPO should be tested alone. Antibodies directed against TPO antigens are more frequently found than those of anti-TG specificity but in few cases, however, positive TG is found together with negative TPO.

Overlapping syndromes have been known within various autoimmune diseases like Sjögren's syndrome, systemic lupus erythematosus (SLE), rheumatoid arthritis or autoimmune haemolytic anaemia, anti-TG antibody concentrations are rather high in this group with a no or very weak or no anti-TPO response. In addition to thyroid disorder anti-TPO or anti-TG are found in other autoimmune diseases, e.g. pernicious anaemia, myasthenia gravis, SLE or rheumatoid arthritis.

Other autoimmune disease expressing organ specific manifestations might give positive results in thyroid autoimmune tests. Overlapping syndromes including the thyroid gland may indicate the need of further testing for other autoantibodies.

10. Assay Performance

The assay characteristics of the MASTAZYME™ TG or TPO ELISAs have been established and assessed according the European IVD directive. Detailed validation data can be provided on special request.

11. References

1. Autoantibodies, Eds. JB Peter, Y Shoenfeld, 1996, Elsevier
2. Furmaniak J., Smith BR, Review: The structure of thyroid autoantigens, Autoimmunity, 7, 63-80 (1990)

Notizen / Notes / Note:

United Kingdom

**Mast Group Ltd.
MAST House Derby Road,
Bootle, Mersey Side L20 1EA**

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**Mast Diagnostica GmbH