# MASTABLOT™ TP

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## Introduction

Syphilis is a highly infectious sexually transmitted disease caused by the spirochaete, *Treponema pallidum*. After infection, a primary lesion or chancre appears on the genitalia at the site of entry of the infecting organism. The disease is systemic from the onset and the course of an untreated case of syphilis can span decades.<sup>1</sup>

Pathogenic treponemas include *T. pallidum*, the cause of syphilis, *T. pertenue*, the cause of yaws and *T. carateum*, the cause of pinta. The only treponeme of importance in Europe is *T. pallidum*. In addition many commensal species of treponemes exist and it is important to differentiate these from *T. pallidum* before a diagnosis of primary syphilis is given.

Commensal treponemas can be cultured in artificial media whereas attempts to culture pathogenic treponemas *in vitro* failed.

The clinical diagnosis of syphilis is confirmed in the laboratory by either demonstrating the presence of *T. pallidum* in the exudates from the lesions, or demonstrating the presence of serum antibodies against the organism. Methods used to measure antibody response to treponemal infection can be divided into two major categories:

1. Tests to measure antibodies against non-specific treponemal antigens i.e. cardiolipin or lipoidal antigen tests. These were formerly called 'Reagin' tests.
2. Tests to measure antibodies against antigens specific for pathogenic treponemas i.e. *T. pallidum* particle agglutination assay (TPPA), the fluorescent treponemal antibody absorbance test (FTA-ABS) or the immuno line / dot blots.

The immuno blot is characterised by the recombinant antigens p15, p17, TmpA and p47, which are highly specific and sensitive for a syphilitic infection. *T. pallidum* blot assays are an accepted reference method as a confirmatory test for diagnosis of syphilis infection.

## Description

MASTABLOT TP is an indirect immunoblot for the detection of IgG and IgM antibodies, respectively, in human serum or plasma. Antibodies bind to purified immobilised antigens which are highly specific for a *T. pallidum* infection.

## Test Principle

Antigens specific for *T. pallidum* are coated onto nitrocellulose strips at defined positions. Each strip contains a cut-off control and one reaction control. Patient serum or plasma specimens are applied together with the strips in the incubation chamber and incubated. If specific antibodies are present in the serum / plasma they will bind to the fixed antigens on the strips forming a stable antigen-antibody complex. Strips are then washed to remove any unbound materials and complexed antibodies are detected by the addition of an alkaline phosphatase labelled anti-human immunoglobulin conjugate. After a further washing step to remove any unbound conjugate, strips are incubated with BCIP/NBT substrate to form a visible dark bluish precipitate. Adding water to this substrate reaction will stop the formation of the precipitating colour.

## Packaging and Ordering Details

MASTABLOT™ TP is presented as complete kits for the performance of 8 tests. Kits are available for the IgG and IgM antibody determination.

### Contents of the kit

- 8 Nitrocellulose strips** coated with recombinant *T. pallidum* antigens (p15, p17, TmpA, p47); depending on the intended use of the kit the strips are specific for IgG or IgM detection.
- 15 mL Sample diluent**  
Green coloured solution, ready to use; contains proclin for conservation.
- 13 mL Conjugate, according to kit specification:**  
Anti-human-IgG H+L chain (goat), labelled with alkaline Phosphatase.  
Anti-human-IgM  $\mu$ -chain (goat) labelled with alkaline Phosphatase.  
Red coloured solution, ready for use.  
The buffer contains 0,02 % methylisothiazolone and 0,02 % bromnitrodioxane.
- 13 mL Chromogenic substrate**  
BCIP/NBT solution, ready to use.
- 15 mL Wash buffer**, 10x concentrated  
The buffer contains proclin for conservation.
- 1 Incubation chamber**
- 1 Interpretation / documentation sheet**
- 1 Instruction leaflet**

### Materials Required but not Provided

1. Sterile test tubes.
2. Micropipettes and tips.
3. Distilled water or water of higher quality.
4. Volumetric flask for dilution of wash buffer concentrate.
5. Forceps.
6. Shaker.
7. Wash bottle.

### Stability and Storage

MASTABLOT™ TP can be used until end of expiry date as indicated on kit label.

All kits components and reagents should be stored at 2–8 °C and brought to room temperature before use.

Diluted wash buffer should be stored at 2–8 °C for up to 30 days.

Serum / plasma specimens may be stored at 2–8 °C for up to 3 days prior to use or for longer storage at -20 °C. Repeated freezing and thawing of specimens should be avoided.

## Warnings and Precautions

1. The reagents supplied in this kit are for *in vitro* diagnostic use only.
2. Read instructions carefully before conducting the assay. Do not modify the procedure.
3. Do not use kit beyond the expiry date.
4. Wear appropriate protective clothing and use appropriate facilities.
5. Do not mouth pipette.
6. Use disposable plastic ware where ever possible. Reusable glassware should be washed thoroughly and rinsed free of detergents before use.
7. Use distilled water or water of higher quality.
8. Do not mix reagents between different lots, as reagents have been calibrated for each batch.
9. Do not cross-contaminate reagents or interchange caps on bottles. Use a separate pipette or pipette tips for each sample and reagent.
10. Do not allow nitrocellulose strips to dry out during the assay procedures.
11. Do not expose slides to intense sunlight or similar adverse conditions while incubating.
12. Contaminated plastic ware should be disposed of and incinerated. Contaminated glassware should either be autoclaved at 121 °C for 1 hour or decontaminated with a solution of 2.5 % sodium hypochlorite. All liquid waste should be decontaminated appropriately e.g. using sodium hypochlorite.
13. Microbial contaminated serum samples should not be used.
14. In the event that haemolysed or lipaemic serum must be used, heat inactive at 56 °C before use.
15. Some reagents contain Proclin as a conservative. Proclin may be toxic if incorporated.

## Test Procedure

1. Optional: Inactivate all sera by heating for 30 min at 56 °C. Do not inactivate plasma.
2. Allow all materials to reach room temperature (at least 20 °C) prior to use.
3. Dilute the wash buffer concentrate 1:10 with distilled water:  
Example:  
25 mL wash buffer + 225 mL of distilled water
4. Important note: Always use the IgG labelled strips for IgG determinations, and IgM-labelled strips for the IgM determinations.  
  
Only touch the strips at the label end. Always use clean disposable gloves when touching the nitrocellulose strips! Alternatively use forceps for touching strips.  
  
Strips must be soaked in the incubation chamber with serum diluent before use. The nitrocellulose side of the strips with the label on top should be upward. For soaking add 1.5 mL of serum diluent to the strip in the incubation chamber. Incubate the strips with the diluent for 10 min ± 2 min.
5. For serum incubation directly add the specimen to the diluent into the incubation chamber.

## Identification of IgG

Sample dilution for **IgG 1:100**

Add 15 µL of serum / plasma to the diluent

- a. **Incubate IgG assays for 60 min** at room temperature on a shaker.
- b. Carefully discard the serum diluent solutions of each incubation chamber. The use of appropriate pipettes or a vacuum pump is helpful.  
  
Immediately add 1.5 mL of 1x wash buffer to each strip. Gently shake for 5 min, discard the wash buffer then repeat the washing step 2–3 times.
- c. After the washing step add 1.5 mL of IgG conjugate to each strip.  
  
The conjugate should cover the strips completely.
- d. Incubate the **assay 60 min** while gently shaking.
- e. After the conjugate reaction wash the strips as described in step b. (Identification of IgG).
- f. Add 1.5 mL of BCIP/NBT substrate solution to each strip.  
  
Gently shake the IgG reactions for 5–8 min at room temperature.  
  
The substrate reaction can be stopped at any time by repeatedly adding tap water or distilled water. It is recommended to stop the reaction immediately after a bluish background colour starts to develop heavily.

## Identification of IgM

Sample dilution for **IgM 1:50**

Add 30 µL of serum / plasma to the diluent

- a. **Incubate IgM assays for 120 min** at room temperature on a shaker.
- b. Carefully discard the serum diluent solutions of each incubation chamber. The use of appropriate pipettes or a vacuum pump is helpful.  
  
Immediately add 1.5 mL of 1x wash buffer to each strip. Gently shake for 5 min, discard the wash buffer then repeat the washing step 2–3 times.
- c. After the washing step add 1.5 mL of IgM conjugate to each strip.  
  
The conjugate should cover the strips completely.
- d. Incubate the **assay 60 min** while gently shaking.
- e. After the conjugate reaction wash the strips as described in step b. (Identification of IgM).
- f. Add 1.5 mL of BCIP/NBT substrate solution to each strip.  
  
Gently shake the IgM reactions for 8-10 min at room temperature.  
  
The substrate reaction can be stopped at any time by repeatedly adding tap water or distilled water. It is recommended to stop the reaction immediately after a bluish background colour starts to develop heavily.

## Interpretation of Results

### 1. Test Validation

A reaction control band as well as a serum control band are coated onto each strip. The reaction control band must show a clear intense signal.

Each strip is coated with a specific cut-off defining control band, too, which is at the lower end of the strip.

The cut-off band is of low reactivity in both IgG and IgM assay versions.

Control band reactivities shall be as follows:

	Expected Result
Serum control band IgG	1+ to 2+ intensity
Serum control band IgM	1+ to 2+ intensity
Reaction control band IgG	2+ to 3+ intensity
Cut-off band IgG	Weak intensity (+)
Reaction control band IgM	2+ to 3+ intensity
Cut-off band IgM	Weak intensity (+)

As long as strips have not dried completely cut-off bands appear to be very weak. With drying the intensity of the cut-off band increases and is best after strips are completely dried through.

A result is only valid if the serum, reaction and cut-off bands are visible. The intensity of each band should be within the specifications listed in the table above or QC document.

### Interpretation of Specimen Results

It is advisable to remove the strips while still wet to fix them onto the lamination. If necessary, apply clear tape over the labeled end to ensure adhesion to the page.

After completely drying the band on the strips can be easily read and interpreted. Drying can be fastened by using a hairdryer or a placing the strips into a drying incubator.

### IgG positive reactivity

**IgG positive** specimens must show at least **1 strong reactive antigen band**, which should belong to one of the antigen specificities: p15, p17, or p47, respectively. The intensity of that strongly reacting band must be far beyond the strength of the cut-off control band.

In most cases further bands are visible, too. These bands must have an intensity equal to or greater than the cut-off band.

The TmpA band may confirm a positive result. However, a TmpA band on its own is not sufficient to be considered a positive result.

Blots with **negative specimens** do not show any reactive antigen bands of the p15, p17, TmpA and p47 specificity. Strips with just 1 single band in the intensity of the cut-off band are considered as positive.

### IgM positive reactivities


**IgM positive** specimens must show at least **1 band** of either the p15, p17 or p47 antigen specificity. The intensity of the band must be equal to or greater than the cut-off control.

The TmpA band may confirm a positive result. However, a TmpA band on its own is not sufficient to be considered a positive result.

Blots with **negative specimens** do not show any reactive antigen bands of the p15, p17, TmpA and p47 specificity. Any antigen band reactivity less than the cut-off band must be read as negative.

### IgG-Blot example

abbreviation: CO = cut-off, SC = Serum control, RC = reactive control

### Correlation of SERODIA TPPA, MASTAFLUOR FTA-ABS and MASTABLOT TP results

1. Weak antigen bands on MASTABLOT TP may occur with samples showing S-TPPA titres of 1:160 and in some case even 1:320. Some of such samples may even give a negative result in MASTABLOT TP. Those samples are negative in an FTA-ABS assay as well.
2. Low reactivity in MASTABLOT IgM does correlate with low reactivity in 19S IgM FTA-ABS (titre ≤ 1:20).

### Limitations

1. The MASTABLOT TP test may occasionally give a false positive result under certain patient conditions or diseases e.g. pregnancy, leprosy and systemic lupus erythematosus.
2. The MASTABLOT test has been demonstrated to be highly sensitive and specific, however results should be considered along with all serological tests, clinical history and other aspects of patient management to be considered diagnostically significant.

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**Mast Group Ltd.
**Mast Diagnostica GmbH