BIOHIT HealthCare

Innovating for Health

ENGLISH

BIOHIT Lactose Intolerance quick test

Quick test for the detection of hypolactasia from biopsy specimen
REF 602 010 (25 tests), REF 602 012 (10 tests)

1. INTENDED USE

The Lactose Intolerance quick test is based on biochemical reactions which qualitatively determine hypolactasia from biopsy specimens.

2. CLINICAL BACKGROUND

Lactose intolerance (hypolactasia) results in the inability to digest lactose, the predominant sugar of milk. This inability results from a decrease in the activity of lactase enzyme, which is produced in the small intestine. Lactase breaks down milk sugar (lactose) into monosaccharides (glucose and galactose), which can then be absorbed into the bloodstream.

Lactase activity begins to decrease after weaning. It is estimated that approximately 15-20% of Western and Northern Europeans, and 90% of Asians, Africans and Native Americans suffer from lactose intolerance. Often the symptoms remain undiagnosed, and, thus lactose intolerance may remain untreated for years. Lactose intolerance may cause the following gastrointestinal symptoms after the intake of lactose: nausea, flatulence, abdominal bloating, cramps and pain, as well as diarrhæa.

3. PRINCIPLE OF THE QUICK TEST

The Lactose Intolerance quick test developed by Biohit for detecting lactose intolerance (hypolactasia of the small intestine) is based on the activity of the lactase enzyme in a biopsy specimen. The biopsy specimen taken from the mucosa of the upper part of the small intestine is examined immediately. The development of color in the test liquid after 20 minutes indicates whether there is functional lactase in the biopsy specimen. Lactase enzyme in the biopsy specimen breaks down the lactose (milk sugar) in the test buffer. As a result blue color develops in the buffer indicating a negative test result. However, if there is no or only a slight color reaction it can be concluded that the patient suffers from lactose intolerance.

4. WARNINGS AND PRECAUTIONS

CAUTION: Handle biopsy specimens as potential biohazardous material. All biopsy specimens should be regarded as potentially contaminated and treated as if they were infectious. Please refer to the U. S. department of Health and Human Services (Bethesda, MD., USA) publication Biosafety in Microbiological and Biomedical Laboratories, 1999, 4th ed. (CDC/NIH) and No. (CDC) 88-8395 on regarding laboratory safety procedures on different diseases or any other local or national regulation. In the United States for research use only. Always use protective gloves when handling patient samples. Read all instructions prior to performing the quick test. Do not use test plates or reagents before the expiry date. Discard the used test plates to biohazardous waste according to local and national regulations. Reagents containing 0.002 % Thiomersal should be disposed of according to local and national regulations.

5. SPECIMEN COLLECTION AND HANDLING

Biopsy specimens are recommended to be taken with e.g. 5 mm forceps from the mucosa of the upper part of the small intestine at any site of the postbulbar mucosa. The biopsy specimen should be collected and stored as follows: If the biopsy is smaller than recommended, there is a risk of a false positive (hypolactasia) result. If the biopsy specimen is too large, there is a risk of a false negative (normolactasia) result. Before performing the test, blood can be removed from the biopsy specimen by placing it briefly on the sterile gauze pad. It is recommended to perform the test immediately.

Do not contaminate the biopsy forceps or biopsy specimen with formalin or other tissue fixatives. It is recommended that the biopsies for LIQT test are collected first, before any biopsies are collected for histology. If the forceps become contaminated with formalin or other tissue fixatives rinse thoroughly in water before collecting the biopsies for LIQT test.

6. KIT CONTENTS AND REAGENT PREPARATION

Each Biohit Lactose Intolerance quick test kit contains reagents for 10 [Cat. no. 602 012] or 25 [Cat. no. 602 010] tests. Store the kit components refrigerated (2-8 °C) until use.

6.1. Test Plates

Contents: 10 (602 012) or 25 (602 010) test plates in a bag, each containing a well for the biopsy specimen. Preparation: Ready for use.

6.2. Substrate Solution, Bottle 1, SUBS
Contents: 1 vial containing 1.6 ml (602 012) or 3.5 ml (602 010) of acetate buffer and lactose with a 0.002% Thiomersal as preservative. Preparation: Ready for use. Stability: Stable until expiration date.

6.3. Chromogen Solution, Bottle 2, CHRO
Contents: 1 vial containing 0.75 ml of chromogen solution in acetic acid with 0.002 % Thiomersal as preservative. Preparation: Ready for use. Stability: Stable until expiration date. **WARNING:** Acetic acid is corrosive and can cause irritation!

6.4. Signal Reaction Solution, Bottle 3, SIGN
Contents: 1 vial containing 1.6 ml (602 012) or 3.5 ml (602 010) of enzyme solution with 0.002 % Thiomersal as preservative. Preparation: Ready for use. Stability: Stable until expiration date.

6.5. Instructions for Use

6.6. Interpretation color chart

7. MATERIALS REQUIRED, BUT NOT PROVIDED
Forceps, Timer, Sterile gauze pads, Gloves

8. STORAGE AND STABILITY

Store the Lactose Intolerance quick test kit refrigerated (2-8 °C). When stored at these temperatures the kit is stable until the expiration date printed on the box label and the label of each individual kit component. Do not freeze or expose the kit to high temperatures or store above 8 °C when not in use. Do not use reagents after the expiration date printed on the label.

9. TEST PROCEDURE

Read the complete assay procedure before starting. Allow all the reagents to reach room temperature (20-25 °C) for at least 15 minutes before use. Mix all the reagents before using by turning the bottles upside down for a few times. Then tap the bottles on the table to ensure the liquid returns to the bottle. When adding the drops into the wells, hold the bottles in vertical upside down position [See figures 1-2].

STEP 1: LACTASE REACTION

1. Open the label covering the well on the plate. Place the biopsy specimen into the well. See figure 1.
2. Add 2 drops (80 µl) of the substrate solution (Bottle 1) into the well. Close the label.
3. Mix the plate by shaking 5-6 times vigorously sideways on the table.
4. Incubate for 15 minutes at room temperature (20-25 °C).

STEP 2: SIGNAL REACTION

1. Open the label again. Add 1 drop (10 µl) of the chromogen solution (Bottle 2) into the well. See figure 2.
2. Immediately add 2 drops (80 µl) of the signal reaction solution (Bottle 3). Close the label.
3. Mix the plate by shaking for 5-6 times vigorously sideways on the table.
4. Incubate for 5 minutes at room temperature (20-25 °C) and interpret the lactase content with the color chart included in the kit package.

10. RESULTS

Presence of functional lactase in the specimen can be interpreted with the color chart given on the label.

See the example below.

- POSITIVE (hypolactasia)
- MILD POSITIVE (mild hypolactasia)
- NEGATIVE (normolactasia)

In case of a negative result hypolactasia can be excluded. For the more detailed interpretation of the results please refer to the interpretation chart included in the kit.

11. LIMITATIONS

The recommended diameter of biopsies is 1.5-2.0 mm in order to avoid imprecision in test results. As with any diagnostic procedure the Biohit Lactose Intolerance quick test results must be interpreted in the light of the patient's clinical presentation and any other information available to the physician.

12. WARRANTY

Biohit shall remedy all defects discovered in any Product (the "Defective Product") that result from unsuitable materials or negligent workmanship and which prevent the mechanical functioning or intended use of the Products including, but not limited to, the functions specified in the Biohit's specifications for the Products. ANY WARRANTY IS DENIED AS VOID IF FAILURE TO HAVE BEEN CAUSED BY MALTREATMENT, MISUSE, ACCIDENTAL DAMAGE, INCORRECT STORAGE OR USE OF THE PRODUCTS FOR OPERATIONS OUTSIDE THEIR SPECIFIED LIMITATIONS OR OUTSIDE THEIR SPECIFICATIONS, CONTRARY TO THE INSTRUCTIONS GIVEN IN THE INSTRUCTION MANUAL. The period of this warranty for the Distributor is defined in the instruction manual of the Products and will commence from the date the relevant Product is shipped by Biohit. This Biohit Diagnostic kit has been manufactured according to ISO 9001 / ISO 13485 quality management protocols and has passed all relevant Quality Assurance procedures related to this product. In case of interpretation disputes the English text applies.

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