# CLASSIQSwabs™

Instructions for Use

CE 0123

## Copan CLASSIQSwabs™ in tube Instructions for use

### INTENDED USE

Copan CLASSIQSwabs™ in tube are sterile and ready-for-use systems intended for clinical samples drawing. The swab is suitable for short contact with the human body. Product should be used as directed.

### SUMMARY AND PRINCIPLES

One of the routine procedures in the diagnosis of infections involves the collection of swab samples. This can be accomplished using the Copan CLASSIQSwabs™.

Copan CLASSIQSwabs™ are available with different size and shape of shaft for the different body collection sites and different types of fiber tips including Viscose, Polyester, Cotton. CLASSIQSwabs™ are provided into polypropylene labeled tube with or without peel pouch.

CLASSIQSwabs™ in tube are available with a single swab or double swab. The double swabs can be removed from the cap if required by the specimen collection protocol or the sample processing.

For ordering details refer to Copan Italia Customer Service, [info@copangroup.com](mailto:info@copangroup.com), or to Copan Italia website, [www.copangroup.com](http://www.copangroup.com). Copan CLASSIQSwabs™ are gamma sterilized.

### TECHNICAL NOTE

A slight yellowing of the tip is a well-known phenomenon. This could be due to many factors: to the type of raw material, to the product sterilization treatment, to the product natural aging. Therefore, product yellowing is not necessarily indicative of product deterioration.

### PRECAUTIONS

- Carefully follow the user's instructions.
- When using a device with plastic or wooden stick, the pressure applied during sampling shall be light since the stick material is breakable.
- The fiber adhesion to the stick is tested for instantaneous sampling: a longer contact between the swab and the sampling area might cause the fiber to come out.
- If the swab is submitted to a chemical or physical sterilizing or micro-biostatic process, its intended functioning could be compromised.
- Product cannot be reused or resterilized.

### STORAGE

Store the Copan Swab devices at 5-25°C. If the product is stored at a different temperature, it would not maintain the stated specifications.

### PRODUCT DETERIORATION

The content of the unopened and undamaged units is guaranteed to be sterile.

CLASSIQSwabs™ should not be used if (1) there is evidence of damage or contamination to the product, (2) the expiration date has passed, (3) the swab package is open or damaged, or (4) there are other signs of deterioration.

### DIRECTIONS FOR USE, SPECIMEN COLLECTION

NOTE: Sterile gloves and protective clothing and eyewear should be worn when collecting and handling microbiology specimens.

- Open the pouch if present
- Open the tube and extract the swab
- Collect the sample. During sampling, the swab end (tip) shall only come in contact with the suspected infection, so as to reduce contamination risks.

NOTE: Do not use excessive force, pressure or bending when collecting swab samples from patients as this may result in accidental breakage of the swab shaft. Swab shafts often exhibit diameter changes to facilitate different sampling requirements. Swab shafts may also have a molded breakpoint point designed for intentional breakage of the swab into a transport tube. In all circumstances when collecting swabs samples from patients, do not use excessive force, pressure or bending of the swab as this may result in accidental breakage of the swab shaft.

- Process the swab according to internal laboratory procedure. Alternatively put the swab into its transport test tube.
- Write the patient's name and data on the tube label.
- Send the sample to the laboratory for immediate testing.


Figure 1 – Directions for use for single swab in tube


Figure 2 – Directions for use for dual swabs in tube

Specimens collected for clinical investigations should be collected and handled following published manuals and guidelines. Specific requirements for the shipment and handling of specimens should be in full compliance with state and federal regulations. Shipping of specimens within medical institutions should comply with internal guidelines of the institution.

#### LIMITATIONS

1. In the laboratory, wear latex gloves and other protection commensurate with universal precautions when handling clinical specimens. Observe Biosafety recommendations when handling or analyzing patient samples.
2. Condition, timing, and volume of specimen collected for clinical investigation are significant variables in obtaining reliable results. Follow recommended guidelines for specimen collection.
3. Copan CLASSIQswabs™ must be used by skilled personnel only.

#### WARNINGS

1. This product is for single use only; reuse may cause a risk of infection and/or inaccurate results.
2. Do not resterilize unused swabs.
3. Do not repack.
4. Not suitable for any other application than intended use.
5. The use of this product in association with a rapid diagnostic kit or with diagnostic instrumentation should be previously validated by the user.
6. Do not use if the swab is visibly damaged (i.e., if the swab tip or swab shaft is broken).
7. Do not use excessive force, pressure or bending when collecting swab samples from patients as this may result in accidental breakage of the swab shaft.
8. Directions for use must be followed carefully. The manufacturer cannot be held responsible for any unauthorized or unqualified use of the product.
9. To be handled by trained personnel only.
10. It must be assumed that all specimens contain infectious micro-organisms; therefore, all specimens must be handled with appropriate precautions. After use, tubes and swabs must be disposed of according to laboratory regulations for infectious waste.

#### WASTE DISPOSAL

Unused swabs may be considered as non-hazardous waste and disposed of accordingly. Dispose of used swabs as well as any other contaminated disposable materials following procedures for infectious or potentially infectious products. It is the responsibility of each laboratory to handle waste and effluents produced according to their nature and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with applicable regulations.

Please refer to symbol table at the end of the instructions for use

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Copan Italia S.p.A.
Copan Diagnostics Inc.