Vitassay qPCR

Campylobacter + Salmonella + Shigella/EIEC

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Termocicladores con bloque de bajo perfil	Termocicladores con bloque de alto perfil
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 ⁽¹⁾
AriaDx Real-Time PCR System	Applied Biosystems
Applied Biosystems	7300 Real-Time PCR System ^{(3) (1)}
7500 Fast Real-Time PCR System ^{(1) (2)}	7500 Real-Time PCR System ⁽¹⁾
7500 Fast Dx Real-Time PCR System ^{(1) (2)}	7900 Real-Time PCR System ⁽³⁾
QuantStudio™ 12K Flex 96-well Fast	ABI PRISM 7000 ⁽³⁾
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7700 ⁽³⁾
QuantStudio™ 7 Flex 96-well Fast	QuantStudio™ 12K Flex 96-well
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
QuantStudio™ 5 Fast Real-Time PCR System	QuantStudio™ 7 Flex 96-well
QuantStudio™ 3 Fast Real-Time PCR System ⁽³⁾	QuantStudio™ 5 Fast Real-Time PCR System
StepOne Plus™ Real-Time PCR System ⁽³⁾	QuantStudio™ 5 Real-Time PCR System
StepOne™ Real-Time PCR System ⁽⁴⁾	QuantStudio™ 3 Real-Time PCR System ⁽³⁾
ViiA™ 7 Fast Real-Time PCR System	ViiA™ 7 Real-Time PCR System
Azure Biosystems	Agilent Technologies
Azure Cielo 3 ⁽⁵⁾	Mx3000P™ Real Time PCR System
Azure Cielo 6	Mx3005P™ Real Time PCR System
BIONEER	Analytik Jena
Exicycler™ 96	qTOWER
Bio-Rad	BIONEER
CFX96™ Real-Time PCR Detection System	Exicycler™ 96
Mini Opticon™ Real-Time PCR Detection System ⁽⁵⁾	BIOER
Roche	QuantGene 9600
LightCycler®480 Real-Time PCR System ^{(1) (6)}	Bio-Rad
LightCycler®96 Real-Time PCR System ⁽¹⁾	CFX96™ Deep Well Real-Time PCR
Cobas z480 Analyzer ^{(1) (6)}	iCycler iQ™ Real-Time PCR
	iCycler iQ™5 Real-Time PCR
	My iQ™ Real-Time PCR Detection System ⁽⁵⁾
	My iQ™ 2 Real-Time PCR Detection System ⁽⁵⁾
Formatos especiales ⁽⁷⁾	DNA-Technology
Bio Molecular Systems	DTlite Real-Time PCR System ⁽⁸⁾
Mic Real Time PCR Cycler	DTprime Real-time Detection Thermal Cycler ⁽⁸⁾
Cepheid	Eppendorf
SmartCycler®	Mastercycler™ ep <i>realplex</i>
Precision System Science Co., Ltd.	Qiagen
geneLEAD VIII System	QIAquant 96
Qiagen	
Rotor-Gene® Q	

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Intended use

Vitassay qPCR *Campylobacter* + *Salmonella* + *Shigella*/EIEC allows the qualitative detection and differentiation of DNA from species of *Campylobacter*, *Salmonella* and/or *Shigella*/ enteroinvasive *Escherichia coli* (EIEC) by real-time PCR in human stool samples. The product is intended for use in the diagnosis of *Campylobacter*, *Salmonella* and/or *Shigella*/EIEC infections alongside the patient's clinical data and other laboratory tests outcomes.

References

Vitassay qPCR *Campylobacter* + *Salmonella* + *Shigella*/EIEC 4x8-well strip, low profile 7041010

Vitassay qPCR *Campylobacter* + *Salmonella* + *Shigella*/EIEC 4x8-well strip, high profile 7042010

Materials/reagents provided

Reference	Reagent/Material	Colour	Amount
7041S010/ 7042S010	<i>Campylobacter</i> + <i>Salmonella</i> + <i>Shigella</i> strips low/high profile	-	4x8-well strip
7C010	<i>Campylobacter</i> + <i>Salmonella</i> + <i>Shigella</i> Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4x8-cap strip

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated on the label.
- The resuspended positive control should be stored at -20°C. To avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

Additional equipment and material required

- Collection and transport system.
- Laboratory freezers (-30°C to -10°C and/or ≤ -70°C).
- DNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)

- Centrifuge for 1.5 mL tubes
- Vortex
- Micropipettes (1-20 μ L, 20-200 μ L)
- Filter tips
- Powder-free disposal gloves

Summary

Campylobacter is 1 of 4 key global causes of diarrheal diseases. It is considered to be the most common bacterial cause of human gastroenteritis in the world. *Campylobacter* is a gram-negative, microaerophilic genus of bacteria of the family *Campylobacteriaceae*. The *Campylobacter* genus consists of a large and diverse group of bacteria currently comprising 26 species, not all of which cause human illness. Approximately 90% of human *Campylobacter* illness is caused by one species, *Campylobacter jejuni*. Less common species, such as *C. coli*, *C. upsaliensis*, *C. fetus*, and *C. lari*, can also infect people.

Anyone can become infected with *Campylobacter*, but infection is more common in males, children younger than 5 years, and people 65 years and older. People can get *Campylobacter* infection by eating raw or undercooked poultry, eating something that touched it, eating other foods (including seafood, meat, and produce), or by contact with animals, and by drinking untreated water. People with *Campylobacter* infection usually have diarrhea (often bloody), fever, and stomach cramps. Nausea and vomiting may accompany the diarrhea. These symptoms usually start 2 to 5 days after the person ingests *Campylobacter* and last about one week. Although people with *Campylobacter* infection usually recover on their own, some need antibiotic treatment. Sometimes *Campylobacter* infections cause complications, such as irritable bowel syndrome, temporary paralysis, and arthritis. In people with weakened immune systems, such as those with a blood disorder, with AIDS, or receiving chemotherapy, *Campylobacter* occasionally spreads to the bloodstream and causes a life-threatening infection.

The burden of foodborne diseases, including *Campylobacteriosis*, is substantial: every year almost 1 in 10 people fall ill. *Campylobacter* causes an estimated 1.5 million illnesses each year in the United States. *Campylobacter* outbreaks are not commonly reported, considering how often people get sick from this bacterium, but the frequency has been increasing. Outbreaks have been associated with unpasteurized dairy products, contaminated water, poultry, and produce.

Campylobacter infection is diagnosed when a laboratory test detects *Campylobacter* bacteria in stool (poop), body tissue, or fluids. The test could be

a culture that isolates the bacteria. As an alternative to culture-based methods, both stool antigen tests and nucleic acid amplification tests (NAATs) are now commercially available. The use of PCR tests for the direct detection of *Campylobacter* in stool is increasing.

Salmonella is 1 of 4 key global causes of diarrhoeal diseases. *Salmonella* is a genus of the family *Enterobacteriaceae*. It is a Gram-negative, non-spore-forming, rod-shaped and facultative anaerobic bacterium. The genus is classified into two species: *Salmonella enterica* and *Salmonella bongori*. Biochemical and genomic analysis of *Salmonella enterica* has led to further classification into subspecies, including *enterica*, *salamae*, *arizonae*, *diarizonae*, *houtenae* and *indica*. The clinically important *Salmonella* species are classified under *Salmonella enterica*, which is further classified into more than 2,579 serovars on the basis of their antigenicity.

The bacterium can be transmitted through faecal–oral routes, where susceptible hosts may acquire *Salmonella* through contaminated foods and water. Infection of humans with *Salmonella* results in three main infectious diseases, namely typhoid fever, paratyphoid fever and NTS infections. Typhoid and paratyphoid fevers are caused by *S. Typhi* and *Salmonella enterica* serovar Paratyphi (*S. Paratyphi*), respectively, and are characterized by gastroenteritis and complications such as septicemia, immunological symptoms, leukopenia and neurological symptoms. NTS infections are restricted to gastroenteritis (nausea, vomiting and diarrhoea) or occasional bacteremia, and are usually non-fatal. Most people who get ill from *Salmonella* have diarrhea, fever, and stomach cramps. Some people may also have nausea, vomiting, or a headache. Symptoms usually begin 6 hours to 6 days after infection and last 4 to 7 days. Most people recover without specific treatment and should not take antibiotics.

Research has estimated that *Salmonella* infection causes 2.8 billion cases of diarrhoea annually worldwide. CDC estimates *Salmonella* bacteria cause about 1.35 million infections, 26,500 hospitalizations, and 420 deaths in the United States every year. Although large *Salmonella* outbreaks usually attract media attention, 60–80% of all salmonellosis cases are not recognized as part of a known outbreak and are classified as sporadic cases.

Diagnosing *Salmonella* infection requires testing a specimen such as stool or blood. The test could be a culture that isolates the bacteria or a culture-independent diagnostic test that detects genetic material of the bacteria such as polymerase chain reaction. The enzyme-linked immunosorbent assay (ELISA) is another useful tool in the diagnosis of *Salmonella* infection.

Shigella bacteria cause an infection called shigellosis. *Shigella* is a member of the Gram-negative *Enterobacteriaceae* family. Four species comprise the genus *Shigella*: *S. dysenteriae*, *S. flexneri*, *S. boydii*, and *S. sonnei*. *Shigella* species are so related to *Escherichia coli*. The differentiation of *Shigella* and *E. coli* is even more complicated with the description of *enteroinvasive E. coli* (EIEC). *Shigella* is very

contagious. People could get sick by touching surfaces contaminated with germs from stool from a sick person, eating food that was prepared by someone who is sick with shigellosis, swallowing recreational water and by exposure to feces through sexual contact.

Most people with *Shigella* infection have diarrhea (sometimes bloody), fever, and stomach cramps. Symptoms usually begin 1–2 days after infection and last 7 days. Most people recover without needing antibiotics. However, people with severe illness and those with underlying conditions that weaken the immune system should be given antibiotics. Some people with shigellosis will not have any symptoms. Possible complications from *Shigella* infections include post-infectious arthritis, bloodstream infections, seizures, and hemolytic-uremic syndrome.

Recent estimates attribute *Shigella* to cause ~125 million diarrhoeal episodes annually, leading to around 160000 deaths, with a third of these associated with young children. *Shigella*, along with enterotoxigenic *Escherichia coli* (*E. coli*), were identified as the predominant bacterial diarrhoeal pathogens in paediatric populations of South Asia and sub-Saharan Africa. CDC estimates about 450,000 cases of shigellosis occur in the United States every year, making it the third most common bacterial enteric disease. Shigellosis does not have a marked seasonality, likely reflecting the importance of person-to-person transmission.

Infection is diagnosed when a laboratory identifies *Shigella* in the stool (poop) of an ill person. The test could be a culture that isolates the bacteria or a rapid diagnostic test that detects genetic material of the bacteria. *Shigella* and EIEC could be separated from other *E. coli* by a PCR targeting the *ipaH*-gene.

Principle of the test

Vitassay qPCR *Campylobacter* + *Salmonella* + *Shigella*/EIEC test is based on the real-time amplification of a conserved region of the *16S rRNA* gene for the identification of *Campylobacter*, *invA* gene for *Salmonella* and *ipaH* gene for *Shigella*/EIEC. After DNA isolation, the presence of *Campylobacter*, *Salmonella* and/or *Shigella*/ enteroinvasive *Escherichia coli* (EIEC) is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

The assay is based on 5' exonuclease activity using two primers and a fluorogenic hydrolysis probe to detect accumulation of the amplified target sequence during the PCR reaction. When the polymerase begins to spread the primers, the probe is hydrolyzed by its exonuclease 5'-3' activity causing the spatial separation of the fluorophore and the quencher. The increase in the resulting fluorescent signal is

proportional to the amount of amplified product in the sample and is detected by means of real-time PCR equipment.

Vitassay qPCR *Campylobacter* + *Salmonella* + *Shigella*/EIEC test is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The amplification of the *Campylobacter* DNA target sequence is detected through the Cy5 channel, *Salmonella* DNA target in FAM channel and *Shigella*/EIEC DNA in ROX channel whereas the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2).

Precautions

- For professional *in vitro* diagnostic use (use by qualified and trained clinical laboratory personnel).
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or is broken or if the foil has been broken or damaged.
- Do not remove desiccant from reagent pouches once is open.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- An appearance of the reaction mixture in stabilized format, normally found at the tube bottom, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or colour different from whitish) does not alter the test functionality.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed. Use separate areas for the patient samples and controls preparation.
- Avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile disposable aerosol resistant pipette tips or RNase/DNase positive displacement pipette tips is recommended.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious and they must be managed according to the national safety legislation and national health waste legislation. Take

necessary precautions during the collection, transport, storage, treatment and disposal of samples.

- Use personal protective equipment (PPE) and biological safety cabinet for handling of potentially infectious samples and reagents according to current recommendations.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

Procedures

Specimen collection, transport, and conservation

Specimen collection, conservation, and transport should be maintained per the conditions validated by the user. Vitassay qPCR *Campylobacter* + *Salmonella* + *Shigella*/EIEC kit has been tested in human stool samples. Other sample types must be validated by the user.

Overall, clinical samples should be collected and labelled appropriately in clean containers with or without transport media (depending on sample type) and be processed as soon as possible to guarantee the test quality. It is recommended to use fresh specimens for the test. Transport must always be carried out in accordance with local and national regulations for the transport of biological samples. The specimens should be transported at room temperature for up to 2h or at 2°C to 8°C for up to 24 hours and for long term transport (more than 24 hours), it is recommended shipping at -20°C or lower. The samples can be stored at 2°C to 8°C for up to 24 hours or frozen at -20°C or ideally at -80°C for conservation during a long time. Repeated freeze-thaw cycles should be avoided to prevent sample and nucleic acids degradation.

The clinical specimens must be collected, transported, and stored according to appropriate laboratory guidelines. For more information, refer to the following guidelines:

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DNA extraction

For nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. For the sample's pretreatment follow the instructions for use of the extraction kit used. The assay has been validated with the following extraction kits:

Maxwell® RSC 16 Viral Total Nucleic Acid Purification Kit, using the Maxwell® RSC 16 instrument (Promega).

Maxwell® RSC Blood DNA Kit, using the Maxwell® 16 instrument (Promega).

MagDEA Dx SV kit, using the magLEAD® 12gC instrument (Precision System Science Co.).

Invisorb® Spin Universal Kit (Invitex).

Positive control preparation

Reconstitute the lyophilized *Campylobacter* + *Salmonella* + *Shigella* Positive Control (red tube) in the 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of extracted DNA, negative control (yellow tube) or positive control (red tube) and add them into the corresponding wells.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Polymerase activation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the Cy5 (*Campylobacter*), FAM (*Salmonella*), ROX (*Shigella*/EIEC) and HEX, JOE or VIC channels (Internal Control (IC)). Depending on the equipment used select the proper detection channels (See Attached II).

Analysis and interpretation of results

The results' analysis is done by the software itself of the used real-time PCR system following manufacturer's instructions. To verify the correct operation of the amplification mix, check the internal control (IC) signal emission.

The threshold value may vary between different instruments due to different signal intensities. Use the positive control amplification curve as a starting point during the reaction validation (prior to interpretation of sample results), to ensure that the threshold falls within the exponential phase of the amplification curves and above any background noise signals. It is recommended to set the threshold values for each channel (target) independently by the end user.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive control used in each run, must show an amplification curve (Ct ≤40) for Cy5 (*Campylobacter*), FAM (*Salmonella*), ROX (*Shigella*/EIEC) and HEX (internal control), which validates the reaction.

Negative control

The negative control included in each run, must show the signal absence (Ct >40 or no signal) in FAM and ROX channels, the signal absence (Ct >37 or no signal) in Cy5, and the signal presence (Ct ≤40) in HEX channel, which validates the reaction.

The experiment seems to be failed if there is amplification signal in negative control and/or signal absence in the positive well for any channel. The assay should be repeated.

The Internal Control (IC) should show an amplification signal (Ct ≤40) in positive and negative control wells.

Note that amplification curves with Ct >37 may appear in channel Cy5 due to ambient *Campylobacter*. Therefore, a cut-off Ct of 37 is set for this target.

The clinical samples test results assessment should be performed once controls' results have been validated. The result interpretation is summarized in the following table:

<i>Campylobacter</i> (Cy5)	<i>Salmonella</i> (FAM)	<i>Shigella</i> /EIEC (ROX)	Internal Control (HEX)	Interpretation
- ¹	+	-	+/-*	<i>Salmonella</i> detected, <i>Shigella</i> /EIEC and <i>Campylobacter</i> not detected
- ¹	-	+	+/-*	<i>Shigella</i> /EIEC detected, <i>Salmonella</i> and <i>Campylobacter</i> not detected
+ ¹	-	-	+/-*	<i>Campylobacter</i> detected, <i>Salmonella</i> and <i>Shigella</i> /EIEC not detected
- ¹	+	+	+/-*	<i>Salmonella</i> and <i>Shigella</i> /EIEC detected, <i>Campylobacter</i> not detected
+ ¹	+	-	+/-*	<i>Salmonella</i> and <i>Campylobacter</i> detected, <i>Shigella</i> /EIEC not detected
+ ¹	-	+	+/-*	<i>Shigella</i> /EIEC and <i>Campylobacter</i> detected, <i>Salmonella</i> not detected
+ ¹	+	+	+/-*	All targets' DNA detected
- ¹	-	-	+ [#]	Targets' DNA not detected
- ¹	-	-	- [#]	Invalid

(+) Positive: Amplification signal (Ct ≤40)

(-) Negative: No amplification signal (Ct ≥40 or no signal)

¹Note that amplification curves with Ct >37 may appear in channel Cy5 due to ambient *Campylobacter*. Therefore, a cut-off Ct of 37 is set for this target.

*The Internal Control (IC) shows or not an amplification signal (Ct ≤40 or no signal). Sometimes, its detection is not necessary because a high copy number of the target can cause preferential amplification of target-specific nucleic acids.

[#]In the case of *Salmonella*, *Campylobacter* and *Shigella*/EIEC target genes negative, IC must show an amplification signal with Ct ≤ 35. If there is an absence of signal or Ct value > 35 of the Internal Control, the result is considered as 'Invalid', and retesting is required. It is recommended to repeat the qPCR diluting the DNA sample 1:10 and/or 1:100, or re-extract and retest to check for possible failure in the extraction procedure and/or inhibition issues.

If the obtained result is confusing or doubtful, it is necessary to check that all the steps have been carried out correctly, to verify the correct performance of each qPCR steps and to review all the parameters, the sigmoid shape of the curve and the fluorescence

intensity. It is also recommended to repeat the assay, preferably in duplicate, depending on the available material (repeat qPCR with the same isolated DNA sample, or re-extract and retest another aliquot of the same specimen or, obtain a new specimen and retest).

The test results must be evaluated by a health professional, together with the medical history, clinical symptoms and/or the results obtained in other diagnostic tests.

Quality Control

To confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

The clinical sensitivity and specificity of Vitassay qPCR *Campylobacter* + *Salmonella* + *Shigella*/EIEC was evaluated using 400 stool samples. The extraction was performed with the Invisorb® Spin Universal Kit (Invitex), and the thermal cycler used was the AriaMx Real-Time PCR System (Agilent Technologies).

The results obtained with Vitassay qPCR *Campylobacter* + *Salmonella* + *Shigella*/EIEC were compared with the values obtained by the routine method (culture + MALDI-TOF). In case of negative results, conventional PCR and sequencing are performed to confirm the negative results and resolve discrepant results. The obtained results are shown in the following table:

Target	TP	TN	FP	FN	SE	SP	PPV	NPV
<i>Campylobacter</i>	53	329	12	0	1 (0.93-1)	0.97(0.94-0.98)	1 (0.99-1)	0.82 (0.70-0.90)
<i>Salmonella</i>	22	377	0	1	0.96 (0.78-0.99)	1 (0.99-1)	1 (0.99-1)	1 (0.85-1)
<i>Shigella</i> /EIEC	2	398	0	0	1 (0.16-1)	1 (0.99-1)	1 (0.16-1)	1 (0.99-1)

TP = True Positive, TN = True Negative, FP = False Positive, FN = False Negative, SE = Sensibility, SP = Specificity, PPV = Positive Predictive Value, NPV = Negative Predictive Value.

These results show a high agreement to detect *Campylobacter*, *Salmonella* and/or *Shigella* / EIEC using Vitassay qPCR *Campylobacter* + *Salmonella* + *Shigella*/EIEC kit.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of *Campylobacter*, *Salmonella* and *Shigella*/EIEC templates ranging from 10^7 to 10^1 copies/reaction. This assay has a detection limit of 10 DNA copies per reaction for *Salmonella* and for *Shigella*/EIEC, and 50 DNA copies per reaction for *Campylobacter*, with a positive rate of $\geq 95\%$.

Analytical specificity

The analytical specificity for *Campylobacter*, *Salmonella* and/or *Shigella*/EIEC detection was tested within the panel of different microorganisms. No cross-reactivity was observed between any of the species:

Cross-reactivity assay		
Sapovirus	<i>Cryptosporidium hominis</i>	<i>Helicobacter pylori</i> J99
<i>Aeromonas caviae</i>	<i>Dientamoeba fragilis</i>	Human Rotavirus A
<i>Aeromonas hydrophila</i> subsp. <i>hydrophila</i> .	Echovirus 30	<i>Klebsiella oxytoca</i>
<i>Arcobacter butzleri</i>	<i>Entamoeba dispar</i>	<i>Listeria monocytogenes</i> serovar 1/2c
Astrovirus genotypes 1-8	<i>Entamoeba histolytica</i> strain DS4-868	Norovirus genotypes I and II.4
<i>Bacteroides fragilis</i>	<i>Enterococcus faecalis</i>	<i>Proteus vulgaris</i>
<i>Blastocystis hominis</i>	Enterovirus 68 and 71	<i>Pseudomonas aeruginosa</i>
<i>Candida albicans</i>	Enteroinvasive <i>Escherichia coli</i>	Adenovirus serotypes 1-5, 8, 15, 31, 40 and 41
<i>Citrobacter freundii</i>	Enteropathogenic <i>Escherichia coli</i>	<i>Serratia liquefaciens</i>
<i>Clostridium difficile</i> O27	<i>Giardia intestinalis</i> strain WB clone C6	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>
<i>Clostridium perfringens</i>	<i>Helicobacter cinaedi</i>	<i>Vibrio parahaemolyticus</i> serotype O1:K1
Coxsackievirus A24, A9 and B3	<i>Helicobacter heilmannii</i>	<i>Yersinia enterocolitica</i> O:3 and O:9
<i>Cryptosporidium parvum</i>	<i>Helicobacter hepaticus</i>	

Analytical reactivity

The Vitassay qPCR Campylobacter + Salmonella + Shigella/EIEC kit's reactivity for *Salmonella* was evaluated against DNA extracted from *Salmonella bongori* serotype Brookfield (66:z41:-); *Salmonella enterica subsp. enterica* serotype 4,5,12:i:1,2 serovar Typhimurium; *Salmonella enterica subsp. enterica* serovars paratyphi A and B; *Salmonella enterica subsp. enterica* serovar enteritidis; *Salmonella enterica subsp. enterica* serovar typhi; *Salmonella enterica subsp. enterica* serovar pullorum; *Salmonella enterica subsp. enterica* serovar gallinarum and *Salmonella enterica* serovar tennessee as templates, showing positive results.

The Vitassay qPCR Campylobacter + Salmonella + Shigella/EIEC kit's reactivity for *Campylobacter* was evaluated against DNA extracted from *Campylobacter jejuni subsp. jejuni*; *Campylobacter upsaliensis*; *Campylobacter lari*; *Campylobacter fetus*; *Campylobacter coli* and *Campylobacter hyointestinalis* as templates, showing positive results.

The Vitassay qPCR Campylobacter + Salmonella + Shigella/EIEC kit's reactivity for *Shigella/EIEC* was evaluated against DNA extracted from *Shigella dysenteriae* serotype 1; *Shigella flexneri* and enteroinvasive *Escherichia coli* as templates, showing positive results.

Compatibles real-time PCR equipment

Vitassay qPCR Campylobacter + Salmonella + Shigella/EIEC has been validated on the following equipment:

- 7500 Fast Real-Time PCR System (Applied Biosystems) ¹
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time Detection Thermal Cycler (DNA-Technology)
- DTPrime Real Time Detection Thermal Cycler (DNA-Technology)
- Linear NEOS-96 Real Time PCR System

¹: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

Limitations

- The test is for professional *in vitro* diagnostic use. This test provides a presumptive diagnosis of *Campylobacter*, *Salmonella* and/or *Shigella/EIEC* infection. All obtained results must be interpreted by a specialist together with other clinical information and laboratory findings available.

- This assay was tried with human faecal samples. The use of other samples has not been established.
- The test quality depends on the sample's quality; proper DNA from clinical specimens must be extracted.
- This test does not provide quantitative values nor indicate the number of organisms present. This is a qualitative test.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by the different bacteria, either by high number of DNA template copies which contains each *Campylobacter* + *Salmonella* + *Shigella* Positive Control vial, samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.
- Detection may be affected by several factors and their combinations which may lead to false negative results, including a) inadequate specimen sampling, shipping, storage, and/or handling; b) procedural errors (including DNA isolation); c) DNA degradation during specimen shipping, storage, and/or preparation; d) potential mutations of the target regions of the pathogens genome covered by this test which may result in DNA being undetectable e) pathogen load below the limit of detection for the assay; f) the presence of qPCR inhibitors or other types of interference (an interference study evaluating the effect of vaccines, antiviral therapeutics, antibiotics, chemotherapeutics or immunosuppressant drugs used to prevent or treat *Salmonella*, *Campylobacter* and/or *Shigella*/EIEC infection was not performed); g) failure to follow the manufacturer's instructions and procedures.
- Detection of bacterial DNA may not indicate the presence of viable and/or infectious bacteria or that these pathogens are the causative agents for clinical symptoms.
- Negative results do not preclude *Salmonella*, *Campylobacter* and/or *Shigella*/EIEC infection and should not be the sole basis of a patient treatment/management decision. Specimen types for the *Salmonella*, *Campylobacter* and/or *Shigella*/EIEC identification and/or stage of infection most suitable for its collection have not been established yet. Consider the collection of multiple specimens from the same patient at different time points, which may increase the probability of detecting the bacteria.
- If the patient clinical data, laboratory tests and epidemiological studies suggest possible *Salmonella*, *Campylobacter* and/or *Shigella*/EIEC infection, and other gastrointestinal illnesses have been discarded, a false negative result might not be discarded, and additional tests should be performed.
- Fluorescence values may vary due to multiple factors such as PCR equipment, extraction system, sample type and its previous treatment, among others.

Attached I. Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 ⁽¹⁾
AriaDx Real-Time PCR System	Applied Biosystems
Applied Biosystems	7300 Real-Time PCR System ^{(3) (1)}
7500 Fast Real-Time PCR System ^{(1) (2)}	7500 Real-Time PCR System ⁽¹⁾
7500 Fast Dx Real-Time PCR System ^{(1) (2)}	7900 Real-Time PCR System ⁽³⁾
QuantStudio™ 12K Flex 96-well Fast	ABI PRISM 7000 ⁽³⁾
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7700 ⁽³⁾
QuantStudio™ 7 Flex 96-well Fast	QuantStudio™ 12K Flex 96-well
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
QuantStudio™ 5 Fast Real-Time PCR System	QuantStudio™ 7 Flex 96-well
QuantStudio™ 3 Fast Real-Time PCR System ⁽³⁾	QuantStudio™ 5 Fast Real-Time PCR System
StepOne Plus™ Real-Time PCR System ⁽³⁾	QuantStudio™ 5 Real-Time PCR System
StepOne™ Real-Time PCR System ⁽⁴⁾	QuantStudio™ 3 Real-Time PCR System ⁽³⁾
ViiA™ 7 Fast Real-Time PCR System	ViiA™ 7 Real-Time PCR System
Azure Biosystems	Agilent Technologies
Azure Cielo 3 ⁽⁵⁾	Mx3000P™ Real Time PCR System
Azure Cielo 6	Mx3005P™ Real Time PCR System
BIONEER	Analytik Jena
Exicycler™ 96	qTOWER
Bio-Rad	BIONEER
CFX96™ Real-Time PCR Detection System	Exicycler™ 96
Mini Opticon™ Real-Time PCR Detection System ⁽⁵⁾	BIOER
Roche	QuantGene 9600
LightCycler®480 Real-Time PCR System ^{(1) (6)}	Bio-Rad
LightCycler®96 Real-Time PCR System ⁽¹⁾	CFX96™ Deep Well Real-Time PCR
Cobas z480 Analyzer ^{(1) (6)}	iCycler iQ™ Real-Time PCR
	iCycler iQ™5 Real-Time PCR
	My iQ™ Real-Time PCR Detection System ⁽⁵⁾
	My iQ™ 2 Real-Time PCR Detection System ⁽⁵⁾
Special Formats ⁽⁷⁾	DNA-Technology
Bio Molecular Systems	DTlite Real-Time PCR System ⁽⁸⁾
Mic Real Time PCR Cycler	DTprime Real-time Detection Thermal Cycler ⁽⁸⁾
Cepheid	Eppendorf
SmartCycler®	Mastercycler™ ep <i>realplex</i>
Precision System Science Co., Ltd.	Qiagen
geneLEAD VIII System	QIAquant 96
Qiagen	
Rotor-Gene® Q	

- (1) A special grid is needed to fit these real-time PCR kits.
- (2) Select Ramp Speed "Standard".
- (3) No ROX caption.
- (4) No Cy5 caption.
- (5) Only FAM and HEX caption.
- (6) Specific compensation color is required.
- (7) The product must be reconstituted following the appropriate procedure (see Test procedure) and transferred to the specific tubes for Mic, SmartCycler®, Rotor-Gene® Q or geneLEAD VIII System.
- (8) See Attached III to configure exposure settings.

Attached II. Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common real-time PCR thermocyclers are specified in the following table:

Thermocycler	Vitassay Channel	Detection Channel	Observations
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Applied Biosystems ABI 7500	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Roche Cobas z 480	FAM	465/510	Colour Compensation required
	HEX	540/580	
	ROX	540/610	
	Cy5	610/670	
Cepheid Smartcycler®	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Agilent Technologies Mx3000P™ Mx 3005P™	FAM	FAM	Passive reference option ROX is none
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Bio Molecular Systems Mic Real Time PCR Cycler	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10).
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Agilent Technologies AriaMx	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Qiagen Rotor-Gene® Q	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise"
	HEX	Yellow	

	ROX	Orange	Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	Cy5	Red	
BIONEER Exicycler™ 96	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Attached III. Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
DTlite Real-Time PCR System (DNA-Technology)	FAM	500
	HEX	500
	ROX	500
	Cy5	500
DTprime Real-time Detection Thermal Cycler (DNA-Technology)	FAM	500*
	HEX	1000
	ROX	1000
	Cy5	1000

*If result in FAM channel is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

Bibliography/Bibliografía

1. Centers for Disease Control and Prevention (CDC). *Campylobacter* (Campylobacteriosis). Available from: <https://www.cdc.gov/campylobacter/> Accessed March 2022
2. World Health Organization (WHO). *Campylobacter*. Available from: <https://www.who.int/news-room/fact-sheets/detail/campylobacter> Accessed March 2022
3. Fitzgerald C. (2015). *Campylobacter*. Clin Lab Med. Jun;35(2):289-98.
4. Gut AM, Vasiljevic T, Yeager T, Donkor ON. (2018). Salmonella infection - prevention and treatment by antibiotics and probiotic yeasts: a review. Microbiology (Reading). Nov;164(11):1327-1344.
5. Centers for Disease Control and Prevention (CDC). *Salmonella*. Available from: <https://www.cdc.gov/salmonella/> Accessed March 2022
6. World Health Organization (WHO). *Salmonella (non-typhoidal)*. Available from: [https://www.who.int/news-room/fact-sheets/detail/salmonella-\(non-typhoidal\)](https://www.who.int/news-room/fact-sheets/detail/salmonella-(non-typhoidal)) Accessed March 2022
7. Centers for Disease Control and Prevention (CDC). *Shigella-Shigellosis*. Available from: <https://www.cdc.gov/shigella/> Accessed March 2022
8. Baker S, The HC. (2018). Recent insights into *Shigella*. Curr Opin Infect Dis. Oct;31(5):449-454.
9. Lampel KA, Formal SB, Maurelli AT. (2018). A Brief History of *Shigella*. EcoSal Plus. Jan;8(1):10.1128/ecosalplus.ESP-0006-2017.
10. Van den Beld MJ, Reubsat FA. (2012). Differentiation between *Shigella*, enteroinvasive *Escherichia coli* (EIEC) and noninvasive *Escherichia coli*. Eur J Clin Microbiol Infect Dis. Jun;31(6):899-904.

Trademarks

All trademarks that may appear in this package insert are the property of their respective owners.

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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