VA Vitassay

Real-Time PCR Kits

**Vitassay qPCR Cryptosporidium+
Entamoeba histolytica + Giardia**

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- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- StepOne™ Real-Time PCR System (Applied Biosystems) ^{II}
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTLite Real-Time PCR System (DNA-Technology)
- DTPRime Real Time Detection Thermal Cyclers (DNA-Technology)
- Rotor-Gene® Q (Qiagen)^I
- SmartCycler® (Cepheid)^I

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Termocicladores con bloque de bajo perfil	Termocicladores con bloque de alto perfil
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 Real-Time System
Applied Biosystems	Applied Biosystems
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR System
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
BIONEER	ViiA™ 7 Real-Time PCR System
Exicycler™ 96	Analytik Jena Biometra
Bio-Rad	TOptical
CFX96™ Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	BIONEER
Cepheid	Exicycler™ 96
SmartCycler®	Bio-Rad
Qiagen	CFX96™ Deep Well Real-Time PCR
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
Roche	iCycler iQ™ 5 Real-Time PCR
LightCycler @480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler @96 Real-Time PCR System	MyiQ™ 2 Real-Time PCR Detection System
Cobas z480 Analyzer	Cepheid
	SmartCycler®
	DNA-Technology
	DTlite Real-Time PCR System *
	DTprime Real-time Detection Thermal Cyclor *
	Eppendorf
	Mastercycler™ ep <i>realplex</i>
	Qiagen
	Rotor-Gene® Q
	Stratagene/ Agilent Technologies
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

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Intended use

Vitassay qPCR Cryptosporidium + Entamoeba histolytica + Giardia allows the detection and differentiation of *Cryptosporidium*, *Entamoeba histolytica* and *Giardia lamblia* by real-time PCR in clinical samples. The product is intended for use in the diagnosis of *Cryptosporidium*, *Entamoeba histolytica* and/or *Giardia lamblia* infections alongside clinical data of the patient and other laboratory tests outcomes.

References

Vitassay qPCR Cryptosporidium + Entamoeba histolytica + Giardia 4x8-well strip, low profile	7041012
Vitassay qPCR Cryptosporidium + Entamoeba histolytica + Giardia 4x8-well strip, high profile	7042012

Reagents provided

In references 7041012 and 7042012:

Reference	Reagent/Material	Colour	Amount
7041S012/ 7042S012	Cryptosporidium + Entamoeba histolytica + Giardia strips low/high profile	-	4 x 8-well strip
7C012	Cryptosporidium + Entamoeba histolytica + Giardia Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4 x 8-cap strip

Material and equipment required, not provided

- DNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Vortex
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. To avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

Summary

Cryptosporidium parvum is a coccidian protozoan parasite responsible for gastrointestinal illnesses in humans. Among all species which have been identified in the *Cryptosporidium* genus, *C. parvum* and *C. hominis* are considered the main species associated with human disease.

The genus *Entamoeba* contains many species, six of which are found in the human intestinal tract: *E. histolytica*, *E. dispar*, *E. moshkovskii*, *E. coli*, *E. hartmanni*, and *E. polecki*. Of these species, only *E. histolytica* is associated with pathological injuries, but the ability of the other species to cause disease is unclear. *Entamoeba dispar* has been traditionally classified as non-pathogenic but some strains have been involved with amebic liver abscess and human colitis.

Giardia lamblia (synonyms are *Giardia intestinalis* and *Giardia duodenalis*) is a common zoonotic protozoan that causes one of the most frequently diagnosed intestinal parasitic disease, giardiasis. A large number of different *G. lamblia* genotypes or assemblages have been detected (A to H); however, humans are mainly infected by assemblages A and B. The largest impact of the disease can be found in developing countries, but the parasite is found worldwide. The main routes of transmission are person-to-person by the fecal-oral route, as well as contaminated food and water.

Test principle

Vitassay qPCR Cryptosporidium + Entamoeba histolytica + Giardia test is based on the real-time amplification of a conserved region of the 18S rRNA gene for the identification of *Cryptosporidium*, *Entamoeba histolytica* and *Giardia lamblia*. After DNA isolation, the presence of *Cryptosporidium*, *Entamoeba histolytica* and/or *Giardia lamblia* is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

Vitassay qPCR Cryptosporidium + Entamoeba histolytica + Giardia test is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The detection channels of the target sequences are described below:

Target	Detection Channel
<i>Cryptosporidium</i>	Cy5
<i>Entamoeba histolytica</i>	ROX
<i>Giardia lamblia</i>	FAM
Internal control	HEX, VIC, or JOE (depending on the equipment used)

Precautions

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or broken or if the foil has been broken or damaged.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles, and mask.
- Do not eat, drink, or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment, and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment, and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

Procedures

Specimen collection, processing, and DNA extraction

For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

Maxwell® RSC Blood DNA Kit, using the Maxwell® 16 instrument (Promega).

QIAamp DNA Mini kit (QIAGEN).

QIAamp DNA Stool Mini Kit (QIAGEN).

RIDA® Xtract (r-Biopharm).

Invisorb® Spin Universal Kit (Invitek).

Positive control preparation

Reconstitute the lyophilized *Cryptosporidium* + *Entamoeba histolytica* + *Giardia* Positive Control (red tube) with the 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of DNA sample, negative and positive controls and add them into the corresponding wells.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the Cy5 (*Cryptosporidium*), ROX (*Entamoeba histolytica*), FAM (*Giardia lamblia*) and HEX, JOE, or VIC channels (Internal Control (IC)). If you use the Applied Biosystems 7500 Fast or the Stratagene Mx3005P™ check that passive reference option ROX is none. (Attached II).

Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive controls used in each run, must show an amplification curve for Cy5 (*Cryptosporidium*), ROX (*Entamoeba histolytica*) and FAM (*Giardia lamblia*), which validates the reaction.

Negative control

The negative controls included in each run, must show the absence of signal in Cy5, ROX and FAM which validates the reaction.

The experiment seems to be failed if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

The result interpretation is summarized in the following table:

Cryptosporidium (Cy5)	Entamoeba histolytica (ROX)	Giardia lamblia (FAM)	Internal Control	Negative Control	Positive Control	Interpretation
+	+	+	+/-	-	+	Cryptosporidium, Entamoeba histolytica and Giardia lamblia Positives
-	-	-	+	-	+	Cryptosporidium, Entamoeba histolytica and Giardia lamblia Negatives
+	-	-	+/-	-	+	Cryptosporidium Positive, Entamoeba histolytica and Giardia lamblia Negatives
+	+	-	+/-	-	+	Cryptosporidium and Entamoeba histolytica Positives, Giardia lamblia Negative
+	-	+	+/-	-	+	Cryptosporidium and Giardia lamblia Positives, Entamoeba histolytica Negative
-	+	-	+/-	-	+	Entamoeba histolytica Positive, Cryptosporidium and Giardia lamblia Negatives
-	+	+	+/-	-	+	Entamoeba histolytica and Giardia lamblia Positives, Cryptosporidium Negative
-	-	+	+/-	-	+	Giardia lamblia Positive, Cryptosporidium and Entamoeba histolytica Negatives
+	+	+	+	+	+	Invalid
-	-	-	-	-	-	Invalid

(+) Positive: Amplification signal

(-) Negative: No amplification signal

If the negative samples do not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction must be repeated due to possible problems caused by PCR inhibitors.

Quality Control

To confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

A total of 172 faecal specimens from symptomatic patients were tested by Real Time PCR using Vitassay qPCR Cryptosporidium + Entamoeba histolytica + Giardia and RIDA®GENE Parasitic Stool Panel (R-biopharm) and FTD Stool parasites (fast-track DIAGNOSTICS) for *Cryptosporidium*, *Entamoeba histolytica* and *Giardia lamblia*.

Vitassay qPCR Cryptosporidium + Entamoeba histolytica + Giardia detected *Cryptosporidium* in 38 positive samples. This test identified even 3 additional low positives that could be confirmed as positive by an additional commercial Real Time PCR Kit (Hadfield et al., 2011).

Vitassay qPCR Cryptosporidium + Entamoeba histolytica + Giardia detected *Entamoeba histolytica* in 8 positive samples (6 clinical specimens and 2 samples which were spiked with culture of *E. histolytica* prior to DNA extraction). This test identified even 1 additional low positive, which had been considered positive for *Entamoeba* by microscopic examination.

Vitassay qPCR Cryptosporidium + Entamoeba histolytica + Giardia detected *Giardia lamblia* in 58 positive samples. Besides, two additional low positive samples were identified by Vitassay test.

The results show a high sensitivity and specificity for detecting *Cryptosporidium*, *Entamoeba histolytica* and *Giardia lamblia* using Vitassay qPCR Cryptosporidium + Entamoeba histolytica + Giardia kit.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of *Cryptosporidium*, *Entamoeba histolytica* and *Giardia lamblia* templates ranging from 10^7 to 10^1 copies/reaction. This assay has a detection limit of ≥ 10 DNA copies per reaction for *G. lamblia* and *E. histolytica* and ≥ 50 DNA copies per reaction for *Cryptosporidium*.

Analytical specificity

The analytical specificity for *Cryptosporidium*, *Entamoeba histolytica* and *Giardia lamblia* was tested within the panel of different microorganisms. No cross-reactivity of *Cryptosporidium*, *Entamoeba histolytica* and *Giardia lamblia* was observed between any of the species:

Cross-reactivity assay		
<i>Shigella flexneri</i>	<i>Proteus vulgaris</i>	<i>Enterococcus faecalis</i>
<i>Shigella dysenteriae</i>	<i>Aeromonas hydrophila</i> subsp. <i>hydrophila</i>	<i>Yersinia enterocolitica</i> O:3
<i>Salmonella typhi</i>	<i>Citrobacter freundii</i>	<i>Yersinia enterocolitica</i> O:9
<i>Salmonella paratyphi A</i>	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	<i>Bacteroides fragilis</i>
<i>Salmonella paratyphi B</i>	<i>Serratia liquefaciens</i>	Adenovirus serotype 40
<i>Salmonella typhimurium</i>	<i>Vibrio parahaemolyticus</i>	Adenovirus serotype 41
<i>Salmonella bongori</i>	<i>Clostridium difficile</i>	Rotavirus A
<i>Salmonella enteritidis</i>	<i>Clostridium perfringens</i>	Norovirus Genotype I and II
<i>Salmonella enterica</i> subsp. <i>enterica</i>	<i>Enterotoxigenic Escherichia coli</i>	Astrovirus Genotype I-VIII
<i>Salmonella pullorum</i>	<i>Enteropathogenic Escherichia coli</i>	<i>Campylobacter lari</i>
<i>Salmonella gallinarum</i>	<i>Klebsiella oxytoca</i>	<i>Campylobacter fetus</i>
<i>Helicobacter pylori</i>	<i>Listeria monocytogenes</i>	<i>Campylobacter coli</i>
<i>Helicobacter hepaticus</i>	<i>Candida albicans</i>	<i>Campylobacter jejuni</i> subsp. <i>jejuni</i>
<i>Helicobacter cinaedi</i>	<i>Arcobacter butzleri</i>	<i>Campylobacter upsaliensis</i>
<i>Helicobacter heilmannii</i>	<i>Pseudomonas aeruginosa</i>	

Analytical reactivity

Analytical reactivity of Vitassay qPCR Cryptosporidium + Entamoeba histolytica + Giardia for *Cryptosporidium* was evaluated against *Cryptosporidium parvum*, showing positive results.

Analytical reactivity of Vitassay qPCR Cryptosporidium + Entamoeba histolytica + Giardia for *Entamoeba histolytica* was evaluated against *Entamoeba histolytica* strain DS4-868, showing positive results.

Analytical reactivity of Vitassay qPCR Cryptosporidium + Entamoeba histolytica + Giardia for *Giardia lamblia* was evaluated against *Giardia lamblia*, showing positive results.

Compatible real-time PCR equipment

Vitassay qPCR Cryptosporidium + Entamoeba histolytica + Giardia has been validated on the following equipment:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- StepOne™ Real-Time PCR System (Applied Biosystems) ^{II}
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTLite Real-Time PCR System (DNA-Technology)
- DTPRime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)^I

- SmartCycler® (Cepheid)^I

^I: For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler tubes.

^{II}: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

Limitations

- This test provides a presumptive diagnosis of *Cryptosporidium*, *Entamoeba histolytica* and/or *Giardia lamblia* infection. All results must be interpreted by a specialist together with other clinical information and laboratory findings available to the physician.
- This assay was tried with human faecal samples. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper DNA from clinical specimens must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by the different parasites, either samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

Attached I. Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 Real-Time System
Applied Biosystems	Applied Biosystems
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR System
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
BIONEER	ViiA™ 7 Real-Time PCR System
Exicycler™ 96	Analytik Jena Biometra
Bio-Rad	TOptical
CFX96™ Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	BIONEER
Cepheid	Exicycler™ 96
SmartCycler®	Bio-Rad
Qiagen	CFX96™ Deep Well Real-Time PCR
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
Roche	iCycler iQ™ 5 Real-Time PCR
LightCycler @480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler @96 Real-Time PCR System	MyiQ™ 2 Real-Time PCR Detection System
Cobas z480 Analyzer	Cepheid
	SmartCycler®
	DNA-Technology
	DTlite Real-Time PCR System *
	DTprime Real-time Detection Thermal Cyclers *
	Eppendorf
	Mastercycler™ ep <i>realplex</i>
	Qiagen
	Rotor-Gene® Q
	Stratagene/ Agilent Technologies
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

* See Attached III to configure exposure settings.

Attached II. Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real-Time PCR thermocyclers are specified in the following table:

Thermocycler	Vitassay Channel	Detection Channel	Observations
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Applied Biosystems ABI 7500	FAM	FAM	Passive reference option ROX is none.
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Color Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Cepheid Smartcycler®	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Agilent Technologies Mx3000P™ Mx3005P™	FAM	FAM	Passive reference option ROX is none.
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Agilent Technologies AriaMx	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Qiagen Rotor-Gene® Q	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range must be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
BIONEER Exicycler™ 96	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Attached III. Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
DTlite Real-Time PCR System (DNA-Technology)	FAM	500
	HEX	500
	ROX	500
	Cy5	500
DTprime Real-time Detection Thermal Cycler (DNA-Technology)	FAM	500
	HEX	1000
	ROX	1000
	Cy5	1000

Bibliography/Bibliografia

1. Fontaine M, Guillot E. Development of a TaqMan quantitative PCR assay specific for *Cryptosporidium parvum*. FEMS Microbiol Lett. 2002; 214(1): 13-17.
2. Limor JR, Lal AA, Xiao L. Detection and differentiation of *Cryptosporidium* parasites that are pathogenic for humans by real-time PCR. J Clin Microbiol. 2002; 40(7): 2335-2338.
3. Verweij JJ, Blangé RA, Templeton K, Schinkel J, Brienen EA, van Rooyen MA, van Lieshout L, Polderman AM. Simultaneous detection of *Entamoeba histolytica*, *Giardia lamblia*, and *Cryptosporidium parvum* in fecal samples by using multiplex real-time PCR. J Clin Microbiol 2004; 42(3): 1220-1223.
4. Qvarnstrom Y, James C, Xayavong M, Holloway BP, Visvesvara GS, Sriram R, da Silva AJ. Comparison of real-time PCR protocols for differential laboratory diagnosis of amebiasis. J Clin Microbiol 2005; 43(11): 5491-5497.
5. Stark D, van Hal S, Fotedar R, Butcher A, Marriott D, Ellis J, Harkness J. Comparison of stool antigen detection kits to PCR for diagnosis of amebiasis. J Clin Microbiol 2008; 46(5): 1678-1681.
6. Jerlström-Hultqvist J, Ankarklev J, Svärd SG. Is human giardiasis caused by two different *Giardia* species? Gut Microbes 2010; 1(6): 379-82.
7. Hadfield SJ1, Robinson G, Elwin K, Chalmers RM. Detection and differentiation of *Cryptosporidium* spp. in human clinical samples by use of real-time PCR. J Clin Microbiol. 2011; 49(3): 918-924
8. Liu J, Gratz J, Amour C, Kibiki G, Becker S, Janaki L, Verweij JJ, Taniuchi M, Sobuz SU, Haque R, Haverstick DM, Houpt ER. A laboratory-developed TaqMan Array Card for simultaneous detection of 19 enteropathogens. J Clin Microbiol 2013; 51(2): 472-80.
9. Lau YL, Anthony C, Fakhurrrazi SA, Ibrahim J, Ithoi I, Mahmud R. Real-time PCR assay in differentiating *Entamoeba histolytica*, *Entamoeba dispar*, and *Entamoeba moshkovskii* infections in Orang Asli settlements in Malaysia. Parasit Vectors. 2013; 6(1):250.

Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™, StepOnePlus™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.


Mx3000P™ and Mx3005P™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid.

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