Vitassay qPCR

Flu A + Flu B

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Intended use

Vitassay qPCR Flu A + Flu B allows the detection and differentiation of Influenza A (Flu A) and Influenza B (Flu B) virus by real-time RT-PCR in respiratory samples. The product is intended for use in Influenza virus infections diagnosis alongside the patient's clinical data and other laboratory tests outcomes.

References

Vitassay qPCR Flu A + Flu B 4x8-well strip, low profile	7041021
Vitassay qPCR Flu A + Flu B 4x8-well strip, high profile	7042021
Vitassay qPCR Flu A + Flu B 96-well plate, low profile	7091021
Vitassay qPCR Flu A + Flu B 96-well plate, high profile	7092021

Materials/reagents provided

Reagents provided in references 7041021 and 7042021:

Reference	Reagent/Material	Colour	Amount
7041S021/ 7042S021	Flu A + Flu B strips low/high profile	-	4 x 8-well strip
7C021	Flu A + Flu B Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4 x 8 cap strip

Reagents provided in references 7091021 and 7092021:

Reference	Reagent/Material	Colour	Amount
7091P021/ 7092P021	Flu A + Flu B Plate low/high profile	-	1 plate
7C021	Flu A + Flu B Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	12 x 8 cap strip

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.

- Keep all reagents in the darkness.

Additional equipment and material required

- RNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Vortex
- Micropipettes (1-20 μ L, 20-200 μ L)
- Filter tips
- Powder-free disposal gloves

Summary

Influenza (flu) is a contagious respiratory illness caused by influenza viruses resulting in significant morbidity, mortality, and economic losses. Influenza viruses are a member of *Orthomyxoviridae* family. There are four types of influenza viruses: A, B, C and D. Influenza A and B viruses are the two main types of influenza viruses that cause seasonal epidemics almost every winter. Worldwide, these annual epidemics are estimated to result in about 3 to 5 million cases of severe illness, and about 290 000 to 650 000 respiratory deaths.

Influenza viruses are single-stranded, negative-sense RNA viruses. The influenza viral genome is composed of eight segments. Influenza A viruses are further classified into subtypes according to the combinations of the hemagglutinin (HA) and the neuraminidase (NA), the proteins on the surface of the virus. Influenza B viruses are further classified into two lineages: B/Yamagata and B/Victoria. Both influenza A and B viruses can be further classified into specific clades and sub-clades.

Flu viruses spread mainly by tiny droplets made when people with flu cough, sneeze or talk. These droplets containing viruses are dispersed into the air and infect persons in close proximity who breathe these droplets in. Spread of the virus is less often by touching contaminated surfaces and then touching the eyes, nose, or mouth. The incubation period is about 2 days but ranges from one to four days. Some people can be infected with the flu virus but have no symptoms. During this time, those people may still spread the virus to others. People with flu are most contagious in the first 3-4 days after their illness begins.

Influenza can cause mild to severe illness, and at times can lead to death. Seasonal influenza is characterized by a sudden onset of some or all these symptoms: fever or chills, cough, sore throat, runny nose, muscle or body aches, headaches, fatigue (tiredness) and some people may have vomiting and diarrhea. People at greater risk of severe disease or complications when infected are pregnant women, children younger than 5 years old, the elderly, individuals with chronic medical conditions and individuals with immunosuppressive conditions (such as HIV/AIDS, receiving chemotherapy).

Complications of flu can include bacterial pneumonia, ear infections, sinus infections and worsening of chronic medical conditions. Other possible serious complications triggered by flu can include inflammation of the heart (myocarditis), brain (encephalitis) or muscle (myositis) tissues, and multi-organ failure.

Laboratory confirmation of influenza virus from throat, nasal and nasopharyngeal secretions or tracheal aspirate or washings is commonly performed using direct antigen detection, virus isolation, or detection of influenza-specific RNA by reverse transcriptase-polymerase chain reaction (RT-PCR).

Principle of the test

Vitassay qPCR Flu A + Flu B test is based on the real-time amplification of specific conserved fragments of the *M1* gene of Influenza A and Influenza B virus. The viral RNA extracted is transcribed into cDNA using a specific primer by reverse transcription step followed immediately in the same well by polymerase chain reaction. The virus presence is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

The assay is based on 5' exonuclease activity that uses two primers and a fluorogenic hydrolysis probe to detect accumulation of the amplified target sequence during the PCR reaction. When the polymerase begins to extend the primers, the probe is hydrolyzed by its 5'-3' exonuclease activity causing the spatial separation of the fluorophore and the quencher. The resulting increase in fluorescent signal is proportional to the amount of amplified product in the sample and is detected by means of a real-time PCR equipment.

Vitassay qPCR Flu A + Flu B test is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The amplification of the Influenza A RNA target sequence is detected through the FAM channel, Influenza B virus RNA target sequence is detected in ROX channel whereas the internal control (IC) is detected in HEX, VIC, or JOE channel (depending on the equipment used select the proper detection channel, see Attached II).

Precautions

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or is broken or if the foil has been broken or damaged.
- Do not remove desiccant from reagent pouches once is open.

- Close protective pouches of reagents promptly with the zip seal after each use.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles, and mask.
- Do not eat, drink, or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment, and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious and they must be managed according to the national safety legislation and national health waste legislation. Take necessary precautions during the collection, storage, treatment, and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

Procedures

Specimen collection, processing, and RNA extraction

Patient samples must be collected, transport and storage according to appropriate laboratory guidelines. For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. Follow the manufacturer's instructions for use. The assay has been validated with the following extraction kits:

Maxwell® RSC 16 Viral Total Nucleic Acid Purification Kit (Promega)

Total Nucleic Acid Isolation (TNAI) Kit (Roche)

Invisorb® Spin Universal Kit (Invitex)

Positive control preparation

Reconstitute the lyophilized Flu A + Flu B Positive Control (red tube) with 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of RNA sample, negative (yellow tube) and positive (red tube) controls and add them to the corresponding wells.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Reverse transcription	45°C	15 min	1
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the FAM (Influenza A), ROX (Influenza B) and HEX, JOE, or VIC channels (Internal Control (IC)). If you use the Applied Biosystems 7500 Fast, the Applied Biosystems StepOne Plus™ or the Stratagene Mx3005P™ Real Time PCR System check that passive reference option ROX is none. (Attached II).

Analysis and interpretation of results

The results' analysis is done by the software itself of the used real-time PCR system following manufacturer's instructions. To verify the correct operation of the amplification mix, check the internal control (IC) signal emission.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive control used in each run must show an amplification curve for FAM (Influenza A) and ROX (Influenza B) channels, which validates the reaction.

Negative control

The negative control included in each run must show the signal' absence in FAM and ROX channels, which validates the reaction.

The experiment seems to be failed if there is amplification signal in negative control or signal' absence in the positive control. The assay should be repeated.

The clinical samples test results assessment should be performed once controls' results have been validated. The result interpretation is summarized in the following table:

Influenza A	Influenza B	Internal control	Negative control	Positive control	Interpretation
+	+	+/-	-	+	Influenza A and Influenza B viruses Positive
-	-	+	-	+	Influenza A and Influenza B viruses Negative
+	-	+/-	-	+	Influenza A Positive, Influenza B negative
-	+	+/-	-	+	Influenza B Positive Influenza A Negative
+	+	+	+	+	Experiment failed
-	-	-	-	-	Experiment failed

(+) Positive: Amplification signal

(-) Negative: No amplification signal

Sometimes, the Internal Control (IC) detection is not necessary since a high copies number of the target can cause a preferential amplification of target-specific nucleic acids.

If the negative samples do not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction must be repeated due to possible problems caused by PCR inhibitors.

If the obtained result is confusing or doubtful, it is necessary to check that all the steps have been carried out correctly, to review all the parameters and the sigmoid shape of the curve. If not resolved, it is recommended to repeat the test preferably in duplicate.

The test results must be evaluated by a health professional, together with the medical history, clinical symptoms and/or the results obtained in other diagnostic tests.

Quality Control

To confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

A total of 186 throat swabs from symptomatic patients were tested by Real Time PCR using Vitassay qPCR Flu A + Flu B and CLART® PneumoVir DNA array assay (Genomica).

Influenza A virus was detected in 107 samples by Vitassay qPCR Flu A + Flu B and in 110 samples by CLART® PneumoVir DNA array assay. For these 3 discordant samples, the low amount of template RNA is below the detection limit of the method used.

For Influenza B virus, a total of 50 clinical samples were positives for Vitassay qPCR Flu A + Flu B whereas CLART® PneumoVir DNA array assay detected 51 positive clinical samples. The low amount of template RNA in this respiratory sample is below the detection limit of the method used.

In addition, Vitassay qPCR Flu A + Flu B was evaluated with the external EQA program QCMD 2017 Influenza virus A and B RNA. That program consisted of a panel composed with 10 samples resuspended in transport medium and they were positive or negative for the studied pathogens. Vitassay qPCR Flu A + Flu B correctly detected all samples from the panel.

These results indicate the high sensitivity and specificity to detect Influenza A and B viruses using the molecular diagnostic test Vitassay qPCR Flu A + Flu B.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of Influenza A and Influenza B virus templates ranging from 10^7 to 10^1 copies/rxn. This assay has a detection limit of 10 viral RNA copies per reaction.

Analytical specificity

The analytical specificity for Influenza A and Influenza B viruses was tested within the panel of following pathogens, where no cross-reactivity was seen between any of the species:

Cross-reactivity assay	
<i>Bordetella pertussis</i>	MERS Coronavirus
<i>Bordetella bronchiseptica</i>	Human parainfluenza 1, 2, 3 and 4 viruses
<i>Bordetella holmesii</i>	Bocavirus
<i>Bordetella parapertussis</i>	<i>Legionella bozemanii</i>
<i>Legionella longbeachae</i>	<i>Legionella micdadei</i>
<i>Legionella pneumophila</i>	<i>Legionella dumoffii</i>
<i>Mycoplasma pneumoniae</i>	<i>Streptococcus pneumoniae</i>
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	<i>Haemophilus influenzae</i> MinnA
Methicillin-resistant <i>Staphylococcus aureus</i>	Influenza A/Mallard/Netherlands/2/2009 (H7N7) virus
<i>Moraxella catarrhalis</i>	Influenza A/South Australia/55/2014 virus
<i>Chlamydia caviae</i>	Influenza A/turkey/Virginia/2002 x PR8-IBCDC-5 (H7N2) virus
<i>Chlamydia psittaci</i> genotypes A and C	Influenza A/chicken/Hong Kong/G9/1997 x PR8-IBCDC-2 (H9N2) virus
<i>Chlamydophila pneumoniae</i>	Influenza A/mallard/Netherlands/12/2000 (H7N7) - IBCDC-1 virus
Human rhinovirus	Influenza A/Anhui/1/2013 (H7N9) virus
Human coronavirus 229E, OC43 and NL63	Influenza A/pheasant/New Jersey/1355/1998 (H5N2)-PR8-IBCDC-4 virus
Human adenovirus types 1-5, 8, 31, 40, 41	Influenza A/Singapore/INFIMH-16-0019/2016 (H3N2) virus
Human metapneumovirus A and B	Influenza A/DE-SH/Reiherente/AR8444/ 2013 (H5N8) virus
Respiratory syncytial virus (RSV)	Influenza A/Uruguay/716/2007 (H3N2)(NYMC-175C) virus
Influenza B/Brisbane/60/2008-like virus	Influenza A/California/7/2009(H1N1) pdm09-like virus
Influenza B/Florida/04/06 virus	Influenza A/Perth/16/2009(H3N2)-like virus
Influenza B/Phuket/3073/2013 virus	Influenza A/New Caledonia/20/99(H1N1) virus
Influenza B/Netherlands/207/06 virus	Influenza A/Netherlands/1250/2016 (H1N1)pdm09 virus (clade 6B.1)
Influenza B/Netherlands/2518/2016 virus (clade 1A)	Influenza A/Michigan/45/2015 (H1N1)pdm09 virus
Influenza B/Netherlands/365/2016 virus (clade 3)	Influenza A/Singapore/GP1908/2015, IVR-180 (H1N1)pdm09 virus
Influenza B/Colorado/6/2017 virus	Influenza A/Thüringen/5/17 (H3N2) virus
Influenza B/Maryland/15/2016 virus	Influenza A/Hong Kong/4801/2014, NYMC X-263B (H3N2) virus
Influenza A/Hong Kong/213/2003 (H5N1) virus	Influenza A/Turkey/Germany R2485+86/2014 (H5N8) virus
Influenza A/Netherlands/2393/2015 (H3N2) virus (clade 3C.2a)	Influenza A/Netherlands/398/2014 (H3N2) virus (clade 3C.3a)
Influenza A/Switzerland/9715293/2013 (H3N2) virus	-

Analytical reactivity

Analytical reactivity of Vitassay qPCR Flu A + Flu B for Influenza A was evaluated against strains: A/New Caledonia/20/99(H1N1) virus, A/California/7/2009(H1N1)pdm09-like virus, A/Netherlands/1250/2016 (H1N1)pdm09 virus (clade 6B.1), A/Michigan/45/2015 (H1N1)pdm09 virus, A/Singapore/GP1908/2015, IVR-180 (H1N1)pdm09 virus, A/Perth/16/2009(H3N2)-like virus, A/Thüringen/5/17 (H3N2) virus, A/Switzerland/9715293/2013 (H3N2) virus, A/Hong Kong/4801/2014, NYMC X-263B (H3N2) virus, A/Netherlands/398/2014 (H3N2) virus (clade 3C.3a), A/Netherlands/2393/2015 (H3N2) virus (clade 3C.2a), A/Hong Kong/213/2003 (H5N1) virus, A/Turkey/Germany R2485+86/2014 (H5N8) virus, A/DE-SH/Reiherente/AR8444/2013 (H5N8) virus, A/Mallard/Netherlands/2/2009 (H7N7) virus, A/Anhui/1/2013 (H7N9) virus, A/turkey/Virginia/2002 x PR8-IBCDC-5 (H7N2), A/chicken/Hong Kong/G9/1997 x PR8-IBCDC-2 (H9N2), A/mallard/Netherlands/12/2000 (H7N7) - IBCDC-1, A/pheasant/New Jersey/1355/1998 (H5N2)-PR8-IBCDC-4, A/Singapore/INFIMH-16-0019/2016 (H3N2), A/South Australia/55/2014 and A/Uruguay/716/2007 (H3N2)(NYMC-175C), showing positive result.

Analytical reactivity of Vitassay qPCR Flu A + Flu B for Influenza B was evaluated against strains: B/Brisbane/60/2008-like virus, B/Netherlands/207/06 virus, B/Florida/04/06 virus, B/Phuket/3073/2013 virus, B/Netherlands/2518/2016 (clade 1A) virus, B/Netherlands/365/2016 (clade 3), B/Colorado/6/2017 and B/Maryland/15/2016, showing positive result.

Compatibles real-time PCR equipment

Vitassay qPCR Flu A + Flu B has been validated on the following equipments:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- StepOne Plus™ Real-Time PCR System (Applied Biosystems) ^{II}
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time PCR System (DNA-Technology)
- DTPrime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)^I

I: For Rotor-Gene® Q thermocycler the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q tubes.

II: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

Limitations

- This test provides a presumptive diagnosis of Influenza A and/or Influenza B virus infections. All results must be interpreted by a specialist together with other clinical information and laboratory findings available to the physician.
- This assay was tried with throat swabs samples. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper RNA from clinical specimens must be extracted.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by the different virus, either samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.

Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 RealTime System
Applied Biosystems	Applied Biosystems
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
BIONEER	ViiA™ 7 Real-Time PCR
Exicycler™ 96	Analytik Jena Biometra
Bio-Rad	TOptical
CFX96™ Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	BIONEER
Cepheid	Exicycler™ 96
SmartCycler®	Bio-Rad
Qiagen	CFX96™ Deep Well Real-Time PCR Detection
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
Roche	iCycler iQ™ 5 Real-Time PCR
LightCycler®480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler®96 Real-Time PCR System	MyiQ™ 2 Real-Time PCR Detection System
Cobas z480 Analyzer	Cepheid
	SmartCycler®
	DNA-Technology
	DTlite Real-Time PCR System
	DTprime Real-time Detection Thermal Cycler
	Eppendorf
	Mastercycler™ ep <i>realplex</i>
	Qiagen
	Rotor-Gene® Q
	Stratagene / Agilent Technologies
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

* See Attached III to configure exposure settings.

Attached II: Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

Thermocycler	Vitassay Channel	Detection Channel	Observations
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Color Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mic Real Time PCR Cycler bms	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10)
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene® Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Attached III: Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
DTIite Real-Time PCR System (DNA-Technology)	FAM	500
	HEX	500
	ROX	500
	Cy5	500
DTprime Real-time Detection Thermal Cycler (DNA-Technology)	FAM	500*
	HEX	1000
	ROX	1000
	Cy5	1000

*If result in FAM channel is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

Bibliography/Bibliografía

1. Centers of Disease Control and Prevention (CDC). Influenza (Flu). Available from: <https://www.cdc.gov/flu/index.htm> Accessed January 2021
2. World Health Organization. Influenza (Seasonal). Available from: [https://www.who.int/news-room/fact-sheets/detail/influenza-\(seasonal\)](https://www.who.int/news-room/fact-sheets/detail/influenza-(seasonal)) Accessed January 2021
3. Wang R, Taubenberger JK. (2010). Methods for molecular surveillance of influenza. Expert review of anti-infective therapy. 8(5):517-527.
4. Lakdawala SS, Lee N, Brooke CB. (2019). Teaching an Old Virus New Tricks: A Review on New Approaches to Study Age-Old Questions in Influenza Biology. *J Mol Biol.* 431(21):4247-4258.
5. Keilman LJ. (2019). Seasonal Influenza (Flu). *Nurs Clin North Am.* 54(2):227-243.

Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™, StepOnePlus™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.


Mx3000P™, Mx3005™ and AriaMx are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid.

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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