# Vitassay qPCR

## RSV A + RSV B

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## Intended use

Vitassay qPCR RSV A + RSV B allows the detection and differentiation of RSV A and RSV B by real-time RT-PCR in clinical samples. The product is intended for use in the diagnosis of RSV infections alongside clinical data of the patient and other laboratory tests outcomes.

## References

Vitassay qPCR RSV A + RSV B 4x8 -well strip, low profile 7041022

Vitassay qPCR RSV A + RSV B 4x8-well strip, high profile 7042022

## Materials/reagents provided

Reference	Reagent/Material	Colour	Amount
7041S022/ 7042S022	RSV A + RSV B Virus strips low/high profile	-	4 x 8-well strip
7C022	RSV A + RSV B Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4 x 8 cap strip

## Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

## Additional equipment and material required

- RNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Vortex
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

## Summary

Human respiratory syncytial viruses (RSV) is the major pathogen responsible for yearly outbreaks of lower respiratory tract disease in young children, but it has also been confirmed as an important cause of respiratory tract illness throughout life. It predictably produces a sizeable outbreak of infections each year.

The respiratory syncytial virus belongs to the *Paramyxoviridae* family, more concretely to the *Pneumovirus* genus. According to differences in the G surface glycoprotein, there are two differentiated groups: A and B. Furthermore, six subgroups can be differentiated in group A and three in group B. Role of strain variation in the epidemiology, severity and clinical impact of an RSV outbreak remains unclear as both can circulate concurrently during seasonal outbreaks (even though the A strains seem to be dominant). RSV is considered to be a ubiquitous virus, which primarily localizes in the upper airway, resulting in cough, rhinorrhea, nasal congestion and fever. In approximately 50% of the cases, several days later, the virus can spread to the lower respiratory tract, causing tachypnea, dyspnea, intercostal retractions, but mostly pneumonia or bronchiolitis (characterized by wheezing and hyperaeration). Transmission of the virus mainly occurs through the upper respiratory tract, primarily the eyes or nose. Being a ubiquitous virus as it is, aspiration or ingestion of contaminated aerosols or respiratory droplets are its main pathways of transmission. Due to this, hospital environments contribute to nosocomial infections.

Real Time PCR is one of the preferred methods when diagnosing RSV, we identify RSV through the targeting of the *N* gene.

## Principle of the test

Vitassay qPCR RSV A + RSV B test is based on the real-time amplification of specific conserved fragments of the *N* gene. The viral RNA extracted is transcribed into cDNA using a specific primer by reverse transcription step followed immediately in the same well by polymerase chain reaction. The presence of the virus is detected by an increment in the fluorescence during the reaction upon hydrolysis of the fluorescent probe.

Vitassay qPCR RSV A + RSV B test is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The amplification of the RSV A RNA target sequence is detected through the Cy5 channel, RSV B virus RNA target in FAM channel whereas the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex II).

## Precautions

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or broken or if the foil has been broken or damaged.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

## Procedures

### Specimen collection, processing and RNA extraction

For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

RIDA® Xtract (r-Biopharm).

Maxwell® 16 Viral Total Nucleic Acid Purification Kit (Promega).

Total Nucleic Acid Isolation (TNAI) Kit (Roche).

Invisorb® Spin Universal kit (Stratec)

### Positive control preparation

Reconstitute the lyophilized RSV A + RSV B Positive Control (red tube) with 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube

thoroughly. After first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

### Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of RNA sample, negative and positive controls and add them into each well.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

### Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Reverse transcription	45°C	15 min	1
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (\*) through the FAM (RSV B), Cy5 (RSV A) and HEX, JOE or VIC channels (Internal Control (IC)). If you use the Applied Biosystems 7500 Fast or the Stratagene Mx3005P™ check that passive reference option ROX is none. (Attached II).

## Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

### Positive control

The positive controls used in each run, must show an amplification curve for FAM (RSV B) and Cy5 (RSV A), which validates the reaction.

### Negative control

The negative controls included in each run, must show the absence of signal in FAM, and Cy5 which validates the reaction.

The experiment seems to be failed if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

The result interpretation is summarized in the following table:

RSV A Virus (Cy5)	RSV B Virus (FAM)	Internal control	Negative control	Positive control	Interpretation
+	+	+/-	-	+	RSV A and RSV B viruses Positive
-	-	+	-	+	RSV A and RSV B viruses Negative
+	-	+/-	-	+	RSV A Positive RSV B negative
-	+	+/-	-	+	RSV B Positive RSV A Negative
+	+	+	+	+	Experiment fail
-	-	-	-	-	Experiment fail

(+) **Positive:** Amplification signal

(-) **Negative:** No amplification signal

If the negative samples do not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

## **Quality Control**

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

## **Performance evaluation**

### **Clinical sensitivity and specificity**

A total of 134 respiratory samples from symptomatic patients were tested by Real Time PCR using Vitassay qPCR RSV A + RSV B, RealStar® RSV RT-PCR Kit (Altona Diagnostics) and CLART® PneumoVir DNA array assay (Genomica).

RSV A was detected in 63 samples by Vitassay qPCR RSV A + RSV B Real Time PCR and RealStar® RSV RT-PCR Kit was positive in 62 clinical samples. The result of the Vitassay test was confirmed by the molecular diagnostic kit CLART® PneumoVir DNA array (GENOMICA).

For RSV B, 61 samples were positives by Vitassay qPCR RSV A + RSV B and RealStar® RSV RT-PCR Kit was able to detect RSV B in 60 samples. The result of the Vitassay test was confirmed by the molecular diagnostic kit CLART® PneumoVir DNA array (GENOMICA).

In addition, Vitassay qPCR RSV A + RSV B was evaluated with QCMD 2017 respiratory syncytial virus RNA EQA programme (RSV17). This consisted of 10 positive or negative samples for the studied pathogens. Vitassay qPCR RSV A + RSV B correctly detected all panel samples.

These results show the high sensitivity and specificity to detect RSV A and B of the Vitassay molecular diagnostic kit qPCR RSV A + RSV B.

### **Analytical sensitivity**

The analytical sensitivity was determined by analysis of 10-fold dilution series of RSV A and RSV B virus templates ranging from  $10^7$  to  $10^1$  copies/rxn. This assay has a detection limit of  $\geq 10$  viral RNA copies per reaction.

## Analytical specificity

The analytical specificity for RSV A and RSV B was tested within the panel of following virus, where no cross-reactivity was seen between any of the species:

Cross-reactivity assay		
<i>Bordetella pertussis</i>	<i>Staphylococcus aureus</i> resistente a la Meticilina	Influenza A/Switzerland/9715293/2013 virus
<i>Legionella bozemanii</i>	<i>Haemophilus influenzae</i>	Influenza A/Turkey/Germany R2485+86/2014 virus
<i>Legionella micdadei</i>	<i>Moraxella catarrhalis</i>	Human parainfluenza 1, 2, 3 and 4 viruses
<i>Legionella dumoffii</i>	Influenza B/Brisbane/60/2008-like virus	Human metapneumovirus A and B
<i>Legionella longbeachae</i>	Influenza B/Florida/04/06 virus	Human coronavirus 229E
<i>Legionella pneumophila</i>	Influenza B/Phuket/3073/2013 virus	Human rhinovirus
<i>Mycoplasma pneumoniae</i>	Influenza A/California/7/2009(H1N1) virus	Human Adenovirus 2 strain Adenoid 6
<i>Streptococcus pneumoniae</i>	Influenza A/Perth/16/2009(H3N2) virus	Human Adenovirus 5
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	Influenza A/New Caledonia/20/99(H1N1) virus	

## Analytical reactivity

Analytical reactivity of Vitassay qPCR RSV A + RSV B Virus Real Time PCR for RSV A and RSV B was evaluated against RSV A and RSV B strains, showing positive result for both.

## Compatibles real-time PCR equipment

Vitassay qPCR RSV A + RSV B has been validated on the following equipments:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) <sup>II</sup>
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time PCR System (DNA-Technology)
- DTPrime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)<sup>I</sup>
- SmartCycler® (Cepheid)<sup>I</sup>

I: For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler® tubes.

II: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

## **Limitations**

- This test provides a presumptive diagnosis of RSV A and/or RSV B virus infection. All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- This assay was tried with nasal and throat swabs, nasopharyngeal aspirates, sputum and bronchoalveolare lavage. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper RNA from clinical specimens must be extracted.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by any of the virus, either samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.

## Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
<b>Agilent Technologies</b>	<b>Abbott</b>
AriaMx/AriaDx Real-Time PCR System	Abbott m2000 RealTime System
<b>Applied Biosystems</b>	<b>Applied Biosystems</b>
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Fast Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Fast Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR System
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
<b>BIONEER</b>	ViiA™ 7 Real-Time PCR System
Exicycler™ 96	<b>Analytik Jena Biometra</b>
<b>Bio-Rad</b>	TOptical
CFX96™ / CFX96™ IVD Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	<b>BIONEER</b>
<b>Cepheid</b>	Exicycler™ 96
SmartCycler®	<b>Bio-Rad</b>
<b>Qiagen</b>	CFX96™ / CFX96™ IVD Real-Time PCR Detection
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
<b>Roche</b>	iCycler iQ™ 5 Real-Time PCR
LightCycler®480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler®96 Real-Time PCR System	MyiQ™ 2Real-Time PCR Detection System
Cobas z480 Analyzer	<b>Cepheid</b>
	SmartCycler®
	<b>DNA-Technology</b>
	DTlite Real-Time PCR System
	DTprime Real-time Detection Thermal Cycler
	<b>Eppendorf</b>
	Mastercycler™ep <i>realplex</i>
	<b>Qiagen</b>
	Rotor-Gene® Q
	<b>Stratagene / Agilent Technologies</b>
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

\* See Attached III to configure exposure settings.

## Attached II: Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

REAL-TIME PCR THERMOCYCLER	Vitassay CHANNEL	DETECTION CHANNEL	OBSERVATIONS
<b>Bio-Rad CFX96™</b>	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>ABI 7500 Applied Biosystems</b>	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>Roche Lightcycler®480II</b>	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
<b>Smartcycler® Cepheid</b>	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
<b>Abbott m2000rt</b>	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>Mx3000P™ Mx 3005P™ Stratagene</b>	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>AriaMx Agilent</b>	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>Rotor-Gene® Q Qiagen</b>	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
<b>Exicycler™ 96 BIONEER</b>	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

### Attached III: Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
<b>DTlite Real-Time PCR System (DNA-Technology)</b>	FAM	500
	HEX	500
	ROX	500
	Cy5	500
<b>DTprime Real-time Detection Thermal Cycler (DNA-Technology)</b>	FAM	500
	HEX	1000
	ROX	1000
	Cy5	1000

## **Bibliography/Bibliografía**

1. Bawage SS, Tiwari PM, Pillai S, Dennis V, Singh SR. Recent advances in diagnosis, prevention, and treatment of human respiratory syncytial virus. *Adv Virol.* 2013;2013:595768.
2. Falsey AR. Respiratory syncytial virus infection in adults. *Semin Respir Crit Care Med.* 2007 Apr;28(2):171-81.
3. Van Woensel JB, Kimpen JL, Brand PL. Respiratory tract infections caused by respiratory syncytial virus in children. Diagnosis and treatment. *Minerva Pediatr.* 2001 Apr;53(2):99-106.
4. Caroline Breese Hall. Respiratory Syncytial Virus. Principles and Practice of Clinical Virology, Fifth Edition. Caroline Breese Hall, et al. Respiratory Syncytial Virus: Concerns and Control. *Pediatrics in review*, 2003; 24 (9): 301-309.

## Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.


Mx3000P™ and Mx3022™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid.

## Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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