Vitassay qPCR

Parainfluenza 1/3 + 2/4

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- Cobas Z480 Analyzer (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time PCR System (DNA-Technology)
- DTPrime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)^I
- SmartCycler® (Cepheid)^I

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Termocicladores con bloque de bajo perfil	Termocicladores con bloque de alto perfil
Agilent Technologies	Abbott
AriaMx/AriaDx Real-Time PCR System	Abbott m2000 RealTime System
Applied Biosystems	Applied Biosystems
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Fast Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Fast Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR System
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
BIONEER	ViiA™ 7 Real-Time PCR System
Exicycler™ 96	Analytik Jena Biometra
Bio-Rad	TOptical
CFX96™ / CFX96™ IVD Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	BIONEER
Cepheid	Exicycler™ 96
SmartCycler®	Bio-Rad
Qiagen	CFX96™ / CFX96™ IVD Real-Time PCR Detection
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
Roche	iCycler iQ™5 Real-Time PCR
LightCycler®480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler®96 Real-Time PCR System	MyiQ™ 2Real-Time PCR Detection System
Cobas z480 Analyzer	Cepheid
	SmartCycler®
	DNA-Technology
	DTlite Real-Time PCR System*
	DTprime Real-time Detection Thermal Cycler*
	Eppendorf
	Mastercycler™ ep <i>realplex</i>
	Qiagen
	Rotor-Gene® Q
	Stratagene / Agilent Technologies
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

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Intended use

Vitassay qPCR Parainfluenza 1/3 + 2/4 allows the detection and differentiation Parainfluenza viruses 1, 2, 3 and 4 by real-time RT-PCR in clinical samples. The product is intended for use in the diagnosis of Parainfluenza 1, Parainfluenza 2, Parainfluenza 3 and Parainfluenza 4 in combination with clinical and epidemiological risk factors and evaluation of infections with Parainfluenza subtype infections alongside clinical data of the patient and other laboratory tests outcomes.

References

Vitassay qPCR Parainfluenza 1/3 + 2/4 4x8 -well strip, low profile 7041025

Vitassay qPCR Parainfluenza 1/3 + 2/4 4x8-well strip, high profile 7042025

Materials/reagents provided

Reference	Reagent/Material	Colour	Amount
7041S025A/ 7042S025A	Parainfluenza 1+3 strips low/high profile	-	4 x 8-well strip
7041S025B/ 7042S025B	Parainfluenza 2+4 strips low/high profile	-	4 x 8-well strip
7C025	Parainfluenza viruses Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	8 x 8 cap strip

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

Additional equipment and material required

- RNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Vortex
- Micropipettes (1-20 μ L, 20-200 μ L)
- Filter tips
- Powder-free disposal gloves

Summary

Parainfluenza viruses (PIV or HPIV in humans) are from the Paramyxoviridae family and are divided into 4 types, Parainfluenza 1, Parainfluenza 2, Parainfluenza 3 and Parainfluenza 4. Parainfluenza viruses are medium-sized enveloped viruses whose genomes consist of a single-stranded RNA of negative sense, which encodes 6 structural proteins. In the viral envelope these viruses have two glycoproteins: the fusion protein (F) and the viral binding protein Hemagglutinin-Neuraminidase (HN), which also causes hemadsorption and haemagglutination.

HPIV is usually transmitted through the air due to coughing and sneezing, personal contact and / or touching objects / surfaces contaminated with the virus to later touch the nose, mouth or eyes. After infection, the virus is incubated for 2-7 days until the symptoms develop.

Type 1 and type 2 lead to respiratory diseases in the upper and lower tracts (cold and flu) being the most frequent type 1 in children. Type 3 is associated with diseases of the lower tract such as bronchiolitis, bronchitis and pneumonia, and type 4 is less frequent and is associated with mild and severe respiratory diseases.

The traditional method for the diagnosis of this virus is with immunofluorescence, but it has the disadvantage that it is a very slow method. That with the addition that antigenic methods are less sensitive and specific than other tools such as real-time PCR, so the latter is considered one of the best methods of choice.

Principle of the test

Vitassay qPCR Parainfluenza 1/3 + 2/4 test is based on the real-time amplification of specific conserved fragments of the *hemagglutinin-neuraminidase* gene for HPIV 1, 2, 3 and 4. The viral RNA extracted is transcribed into cDNA using a specific primer by reverse transcription step followed immediately in the same well by polymerase chain

reaction. The presence of the virus is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

Vitassay qPCR Parainfluenza 1/3 + 2/4, it is a ready-to-use test that contains in each well all the necessary reagents in stabilized format to carry out the PCR in real time. Each kit includes two types of strips and each of them corresponds to a different test. The first strip contains the multiplex reaction mixture for the detection of HPIV 1 and HPIV 3. After the amplification reaction, HPIV 1 is detected in the Cy5 channel and HPIV 3 is detected in the FAM channel. The second strip contains the multiplex reaction mixture for subtyping HPIV 2 and HPIV 4. After the amplification reaction, HPIV 2 is detected in the Cy5 channel and HPIV 4 is detected in the FAM channel. Both reaction strips contain an internal control (IC) which is detected in the HEX, VIC or JOE channel (depending on the equipment used).

Precautions

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or broken or if the foil has been broken or damaged.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

Procedures

Specimen collection, processing and RNA extraction

For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

Maxwell® 16 Viral Total Nucleic Acid Purification Kit (Promega).

Total Nucleic Acid Isolation (TNAI) Kit (ROCHE).

Invisorb® Spin Universal Kit (Stratec).

Positive control preparation

Reconstitute the lyophilized Parainfluenza 1/3 + 2/4 Positive Control (red tube) in the 200 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of RNA sample, negative and positive controls and add them into each well.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Reverse transcription	45°C	15 min	1
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the FAM (Parainfluenza 3 and Parainfluenza 4), Cy5 (Parainfluenza 1 and Parainfluenza 2) and HEX, JOE or VIC channels (Internal Control (IC)). If you use the Applied Biosystems 7500 Fast or the Stratagene Mx3005P™ check that passive reference option ROX is none. (Attached II)

Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive controls used in each run, must show an amplification curve for FAM (Parainfluenza 3 and Parainfluenza 4), Cy5 (Parainfluenza 1 and Parainfluenza 2) and HEX, JOE or VIC channels (Internal Control (IC)), which validates the reaction.

Negative control

The negative controls included in each run, must show the absence of signal in FAM, Cy5 and HEX, JOE or VIC channels which validates the reaction.

The experiment seems to be failed if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

The result interpretation is summarized following:

Interpretation of results Vitassay qPCR Parainfluenza 1/3 + 2/4 strips:

First multiplex reaction mix for the detection of Parainfluenza 1 and/or Parainfluenza 3					
Parainfluenza 1 (Cy5)	Parainfluenza 3 (FAM)	Internal control	Negative control	Positive control	Interpretation
+	+	+/-	-	+	Parainfluenza 1 and 3 Positive
-	-	+	-	+	Parainfluenza 1 and 3 Negative
+	-	+/-	-	+	Parainfluenza 1 Positive and Parainfluenza 3 Negative
-	+	+/-	-	+	Parainfluenza 3 Positive and Parainfluenza 1 Negative
+	+	+	+	+	Experiment fail
-	-	-	-	-	Experiment fail
Second multiplex reaction mix for the detection of Parainfluenza 2 and/or Parainfluenza 4					
Parainfluenza 2 (Cy5)	Parainfluenza 4 (FAM)	Internal control	Negative control	Positive control	Interpretation
+	+	+/-	-	+	Parainfluenza 2 and 4 Positive
-	-	+	-	+	Parainfluenza 2 and 4 Negative
+	-	+/-	-	+	Parainfluenza 2 Positive and Parainfluenza 4 Negative
-	+	+/-	-	+	Parainfluenza 4 Positive and Parainfluenza 2 Negative
+	+	+	+	+	Experiment fail
-	-	-	-	-	Experiment fail

(+) Positive: Amplification signal
(-) Negative: No amplification signal

If the negative samples do not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

Quality Control

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

To analyze the sensitivity and clinical specificity, 78 samples of throat swabs were used for the detection of HPIV-1, HPIV-3, HPIV-2 and HPIV-4. It was carried out with the Vitassay qPCR Parainfluenza 1/3 + 2/4 molecular diagnostic tests and the results obtained were compared with those obtained by a CLART® PneumoVir molecular detection method (Genomics)

HPIV-1 was not detected in any sample by the CLART® PneumoVir molecular detection method (Genomics) while Vitassay qPCR Parainfluenza 1/3 + 2/4 detected 1 positive sample (the low amount of template RNA detected in this sample is below the detection limit of the method used). The result of the Vitassay test was confirmed by the commercial RIDA®GENE Parainfluenza kit (R-Biopharm).

HPIV-3 was detected in 40 positive samples by both kits. In addition, Vitassay qPCR Parainfluenza 1/3 + 2/4 detected one more positive sample (the low amount of template RNA detected in this sample is below the detection limit of the method used). The result of the Vitassay test was confirmed by the commercial RIDA®GENE Parainfluenza kit (R-Biopharm).

HPIV-2 was not detected in any sample using both kits.

HPIV-4 was detected in 10 positive samples by both kits. Apart from this, Vitassay qPCR Parainfluenza 1/3 + 2/4 detected 1 more positive sample (the low amount of template RNA detected in this sample is below the detection limit of the method used). Furthermore, this sample was evaluated with a commercial kit (RIDA®GENE Parainfluenza (R-Biopharm)), confirming the results.

The results give a high specificity and sensitivity to detect Parainfluenza using Vitassay and qPCR Parainfluenza 1/3 + 2/4.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of Parainfluenza 1, Parainfluenza 2, Parainfluenza 3 and Parainfluenza 4 templates ranging from 10⁷ to 10¹ copies/rxn. This assay has a detection limit of ≥10 viral RNA copies per reaction.

Analytical specificity

The analytical specificity for Parainfluenza 1/3 + 2/4 was tested within the panel of following virus, where no cross-reactivity was seen between any of the species:

Cross-reactivity testing		
<i>Bordetella pertussis</i>	<i>Moraxella catarrhalis</i>	Influenza A/Anhui/1/2013 (H7N9) virus
<i>Bordetella bronchiseptica</i>	<i>Haemophilus influenzae MinnA</i>	Influenza B/Brisbane/60/2008 - like virus
<i>Bordetella holmesii</i>	<i>Chlamydomphila pneumoniae</i>	Influenza B/Florida/04/06 virus
<i>Bordetella parapertussis</i>	<i>Chlamydia caviae</i>	Influenza B/Phuket/3073/2013 virus
<i>Legionella bozemanii</i>	<i>Chlamydia psittaci</i> genotypes A and C	MERS Coronavirus
<i>Legionella micdadei</i>	Influenza A/New Caledonia/20/99(H1N1) virus	Human metapneumovirus A and B
<i>Legionella dumoffii</i>	Influenza A/California/7/2009 (H1N1) virus	Human coronavirus 229E, OC43 and NL63
<i>Legionella longbeachae</i>	Influenza A/Michigan/45/2015 (H1N1)pdm09 virus	Influenza A/Singapore/GP1908/2015 virus
<i>Legionella pneumophila</i>	Influenza A/Perth/16/2009 (H3N2) virus	Human rhinovirus
<i>Mycoplasma pneumoniae</i>	Influenza A/Thüringen/5/17 (H3N2) virus	Human Adenovirus
<i>Streptococcus pneumoniae</i>	Influenza A/Switzerland/9715293/2013 (H3N2) virus	Respiratory Syncytial virus (RSV)
<i>Staphylococcus aureus subsp. aureus</i>	Influenza A/Turkey/Germany R2485+86/ 2014 (H5N8) virus	Human Bocavirus
Influenza A/Hong Kong/4801/2014 (H3N2) virus	Influenza A/DE-SH/Reiherente/AR8444/ 2016 (H5N8) virus	Human parainfluenza 1, 2, 3 and 4 viruses

Analytical reactivity

Vitassay qPCR Parainfluenza 1/3 + 2/4 has been evaluated against Parainfluenza 1, Parainfluenza 2, Parainfluenza 3 and Parainfluenza 4, with a positive result.

Compatibles real-time PCR equipment

Vitassay qPCR Parainfluenza 1/3 + 2/4 has been validated on the following equipments:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time PCR System (DNA-Technology)
- DTPrime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)^I
- SmartCycler® (Cepheid)^I

I: For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler tubes.

II: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

Limitations

- All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- This assay was tried with throat swab samples dissolved in transport medium. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper RNA from clinical specimens must be extracted.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by the different virus, either samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.

Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers
Agilent Technologies
AriaMx/AriaDx Real-Time PCR System
Applied Biosystems
7500 Fast Real-Time PCR System
7500 Fast Dx Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast
QuantStudio™ 6 Flex 96-well Fast
QuantStudio™ 7 Flex 96-well Fast
QuantStudio™ 5 Fast Real-Time PCR System
QuantStudio™ 3 Fast Real-Time PCR System
StepOne Plus™ Real-Time PCR System
StepOne™ Real-Time PCR System
ViiA™ 7 Fast Real-Time PCR System
BIONEER
Exicycler™ 96
Bio-Rad
CFX96™/ CFX96™ IVD Real-Time PCR Detection System
Mini Opticon™ Real-Time PCR Detection System
Cepheid
SmartCycler®
Qiagen
Rotor-Gene® Q
Roche
LightCycler®480 Real-Time PCR System
LightCycler®96 Real-Time PCR System
Cobas z480 Analyzer

High profile Block Thermocyclers
Abbott
Abbott m2000 RealTime System
Applied Biosystems
7300 Real-Time PCR System
7500 Real-Time PCR System
7900 HT Real-Time PCR System
ABI PRISM 7000
ABI PRISM 7700
QuantStudio™ 12K Flex 96-well
QuantStudio™ 6 Flex 96-well
QuantStudio™ 7 Flex 96-well
QuantStudio™ 5 Real-Time PCR System
QuantStudio™ 3 Real-Time PCR System
ViiA™ 7 Real-Time PCR System
Analytik Jena Biometra
TOptical
qTOWER 2.0
BIONEER
Exicycler™ 96
Bio-Rad
CFX96™ / CFX96™ IVD Real-Time PCR Detection
iCycler iQ™ Real-Time PCR
iCycler iQ™ 5 Real-Time PCR
MyiQ™ Real-Time PCR Detection System
MyiQ™ 2Real-Time PCR Detection System
Cepheid
SmartCycler®
DNA-Technology
DTlite Real-Time PCR System
DTprime Real-time Detection Thermal Cycler
Eppendorf
Mastercycler™ep <i>realplex</i>
Qiagen
Rotor-Gene® Q
Stratagene / Agilent Technologies
Mx3000P™ Real Time PCR System
Mx3005P™ Real Time PCR System

* See Attached III to configure exposure settings.

Attached II: Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

REAL-TIME PCR THERMOCYCLER	Vitassay CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Attached III: Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
DTlite Real-Time PCR System (DNA-Technology)	FAM	500
	HEX	500
	ROX	500
	Cy5	500
DTprime Real-time Detection Thermal Cycler (DNA-Technology)	FAM	500
	HEX	1000
	ROX	1000
	Cy5	1000

Bibliography/Bibliografía

1. R. Pechhini *et al.* Parainfluenza virus as a cause of acute respiratory infection in hospitalized children. *The Brazilian Journal of Infectious Diseases* 2015; 19(4): 358-362
2. L. Magnusson. Analysis of the Human Parainfluenza Virus Replication Promoter. *UW-L Journal of Undergraduate Research X* 2007.
3. R. Vainionpaa *et al.* Biology of Parainfluenza Viruses. *Clinical Microbiology Reviews*, Apr. 1994, p. 265-275.
4. .J. Henrickson *et al.* Parainfluenza viruses. *Clinical Microbiology Reviews* 2003; 16(2): 242–264.
5. M. E. Terlizzi *et al.* Quantitative RT Real Time PCR and indirect immunofluorescence for the detection of human parainfluenza virus 1, 2, 3. *Journal of Virological Methods* 160 (2009) 172–177.
6. K.E. Templeton *et al.* Rapid and Sensitive Rapid and sensitive method using multiplex real-time PCR for diagnosis of infections by influenza A and influenza B viruses, respiratory syncytial virus, and parainfluenza viruses 1, 2, 3, and 4. *Journal of Clinical Microbiology* 2004; 42(4): 1564–1569.

Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.


Mx3000P™ and Mx3005™ are registered trademarks of Agilent Technologies.

Mastecycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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