# Vitassay qPCR

**AMB**

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Adenovirus	Metapneumovirus	Bocavirus
<i>Bordetella pertussis</i>	<i>Bordetella pertussis</i>	<i>Bordetella pertussis</i>
<i>Legionella bozemanii</i>	<i>Legionella bozemanii</i>	<i>Legionella bozemanii</i>
<i>Legionella micdadei</i>	<i>Legionella micdadei</i>	<i>Legionella micdadei</i>
<i>Legionella dumoffii</i>	<i>Legionella dumoffii</i>	<i>Legionella dumoffii</i>
<i>Legionella longbeachae</i>	<i>Legionella longbeachae</i>	<i>Legionella longbeachae</i>
<i>Legionella pneumophila</i>	<i>Legionella pneumophila</i>	<i>Legionella pneumophila</i>
<i>Mycoplasma pneumoniae</i>	<i>Mycoplasma pneumoniae</i>	<i>Mycoplasma pneumoniae</i>
<i>Streptococcus pneumoniae</i>	<i>Streptococcus pneumoniae</i>	<i>Streptococcus pneumoniae</i>
<i>Staphylococcus aureus</i> subsp. <i>Aureus</i>	<i>Staphylococcus aureus</i> subsp. <i>Aureus</i>	<i>Staphylococcus aureus</i> subsp. <i>Aureus</i>
<i>Staphylococcus aureus</i> resistente a la meticilina	<i>Staphylococcus aureus</i> resistente a la meticilina	<i>Staphylococcus aureus</i> resistente a la meticilina
<i>Haemophilus influenzae</i>	<i>Haemophilus influenzae</i>	<i>Haemophilus influenzae</i>
<i>Moraxella catarrhalis</i>	<i>Moraxella catarrhalis</i>	<i>Moraxella catarrhalis</i>
Influenza B/Brisbane/60/2008	Influenza B/Brisbane/60/2008	Influenza B/Brisbane/60/2008
Influenza B/Florida/04/06	Influenza B/Florida/04/06	Influenza B/Florida/04/06
Influenza B/Phuket/3073/2013	Influenza B/Phuket/3073/2013	Influenza B/Phuket/3073/2013
Influenza A/Switzerland/9715293/2013	Influenza A/Switzerland/9715293/2013	Influenza A/Switzerland/9715293/2013
Influenza A/Turkey/Germany R2485+86/2014	Influenza A/Turkey/Germany R2485+86/2014	Influenza A/Turkey/Germany R2485+86/2014
Influenza A/California/7/2009(H1N1) virus	Influenza A/California/7/2009(H1N1) virus	Influenza A/California/7/2009(H1N1) virus
Influenza A/Perth/16/2009(H3N2) virus	Influenza A/Perth/16/2009(H3N2) virus	Influenza A/Perth/16/2009(H3N2) virus
Influenza A/New Caledonia/20/99(H1N1) virus	Influenza A/New Caledonia/20/99(H1N1) virus	Influenza A/New Caledonia/20/99(H1N1) virus
Rinovirus humano	Rinovirus humano	Rinovirus humano
Coronavirus 229E, OC43 y NL63 humano	Coronavirus 229E, OC43 y NL63 humano	Coronavirus 229E, OC43 y NL63 humano
Human parainfluenza 1, 2, 3, 4	Human parainfluenza 1, 2, 3, 4	Human parainfluenza 1, 2, 3, 4
MERS Coronavirus	MERS Coronavirus	MERS Coronavirus
Respiratory syncytial virus (RSV)	Respiratory syncytial virus (RSV)	Respiratory syncytial virus (RSV)
Metapneumovirus A y B humano	Adenovirus	Adenovirus
Bocavirus	Bocavirus	Metapneumovirus A y B humano
Enterovirus D68 and A71	Enterovirus D68 and A71	Enterovirus D68 and A71
Echovirus 11 and 30	Echovirus 11 and 30	Echovirus 11 and 30
Coxsackievirus A9, A21, A24, B3 and B4	Coxsackievirus A9, A21, A24, B3 and B4	Coxsackievirus A9, A21, A24, B3 and B4

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- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) <sup>II</sup>
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time PCR System (DNA-Technology)
- DTPRime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)<sup>I</sup>
- SmartCycler® (Cepheid)<sup>I</sup>

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## Intended use

Vitassay qPCR AMB allows the detection and differentiation of human Adenovirus, Metapneumovirus and Bocavirus by real-time RT-PCR in clinical samples. The product is intended for use in the diagnosis of Adenovirus, Metapneumovirus and Bocavirus infections alongside clinical data of the patient and other laboratory tests outcomes.

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## Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

## Additional equipment and material required

- RNA/DNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 ml tubes
- Vortexer
- Micropipettes (1-20 µl, 20-200 µl)
- Filter tips
- Powder-free disposal gloves

## Summary

Acute viral respiratory tract infections are the leading cause of hospitalization for infants and young children in developed countries, whilst a major cause of death in developing countries. As coinfections are very normal when diagnosing this type of infections, Adenovirus and two recently discovered viruses: Metapneumovirus and Bocavirus, are some of the most common pathogens found.

Bocavirus belongs to the *Parvoviridae* family and the *Bocavirus* genus, having four different genotypes. These being: HBoV 1, HBoV 2, HBoV 3 and HBoV 4. HBoV1 is mainly found in respiratory samples and is related with acute wheezing illness, whilst HBoV 2, 3 and 4 are found in fecal samples and related with gastroenteritis. Depending in which way it is transmitted, through respiration or direct ingestion; it will probably develop into one type of disease or the other. Also, no seasonal prevalence has been observed for HBoV, being found throughout the year.

Metapneumovirus mainly affects children, but in adults appears with flu-like symptoms, and is a major cause of morbidity in elderly and immunocompromised patients. It belongs to the *Paramyxoviridae* family, more concretely to the *Metapneumovirus* genus, possessing two genotypes: A and B. Then, each genotype having two subtypes, these being: A1, A2 (which coincidentally is the most common among hosts), B1 and B2. Lastly, the A2 subtype is also divided into A2a and A2b sub-genotypes. The virus is characterized by a late winter-early spring epidemiologic pattern, even though summer infections cannot be excluded.

Human adenoviruses (HAdVs) are a diverse group of double-stranded DNA viruses responsible for a wide variety of human ailments. Belonging to the *Adenoviridae* family and *Mastadenovirus* genus, they are categorized by species (AdA, AdB1, AdB2, AdC, AdD, AdE, and AdF) and each species can be then categorized by serotype (Ad1 to Ad51). As of 2007, 51 serotypes are recognized. AdB1, AdC and AdE can cause widespread outbreaks with severe clinical presentations. Ad3 and Ad7 from AdB1 are often responsible for the most severe respiratory symptomatology. AdC (1, 2, 5, and 6) are associated with febrile respiratory illness in children and illnesses in immunocompromised patients. With this given, serotype directly correlates with the severity and symptomology of disease.

Real Time PCR has proven to be able to cope with the complications that traditional techniques cannot and is capable of revealing coinfections of multiple Ad serotypes, differentiates Bocavirus, Metapneumovirus and Adenovirus from RSV or influenzas. And it is all thanks to the targeting of the *Hexon* gene for Adenovirus, *N* gene for Metapneumovirus and the *NS1* gene for the Bocavirus.

## Principle of the test

Vitassay qPCR AMB test is based on the real-time amplification of specific conserved fragments of the *hexon* gene (Adenovirus) and the *NS1* gene (Bocavirus). The viral RNA extracted is transcribed into cDNA using a specific primer by reverse transcription step followed immediately in the same well by polymerase chain reaction. The presence of the virus is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe. After RNA isolation, the detection of Metapneumovirus is done in one step real time RT format where the reverse transcription and the subsequent amplification of specific target sequence occur in the same reaction well. The isolated RNA target is transcribed generating complementary DNA by reverse transcriptase which is followed by the amplification of a conserved region of *N* gene of Metapneumovirus.

Vitassay qPCR AMB test is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The amplification of the Adenovirus DNA target sequence is detected through the FAM channel, Metapneumovirus RNA target in ROX channel and Bocavirus DNA in Cy5 channel whereas the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Attached II).

## Precautions

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or broken or if the foil has been broken or damaged.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.

- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

## **Procedures**

### **Specimen collection, processing and RNA/DNA extraction**

For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

Maxwell® 16 Viral Total Nucleic Acid Purification Kit (Promega).

RIDA® Xtract (r-Biopharm).

Total Nucleic Acid Isolation (TNAI) Kit (Roche).

Invisorb® Spin Universal Kit (Stratec).

### **Positive control preparation**

Reconstitute the lyophilized AMB Positive Control (red tube) in the 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

### **Reaction setup**

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of RNA/DNA sample, negative and positive controls and add them into each well.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

## Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Reverse transcription	45°C	15 min	1
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (\*) through the FAM (Adenovirus), ROX (Metapneumovirus), Cy5 (Bocavirus) and HEX, JOE or VIC channels (Internal Control (IC)). If you use the Applied Biosystems 7500 Fast or the Stratagene Mx3005P™ check that passive reference option ROX is none. (Attached II)

### Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

#### Positive control

The positive controls used in each run, must show an amplification curve for FAM (Adenovirus), ROX (Metapneumovirus) and Cy5 (Bocavirus), which validates the reaction.

#### Negative control

The negative controls included in each run, must show the absence of signal in FAM, ROX and Cy5 which validates the reaction.

The experiment seems to be failed if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

The result interpretation is summarized in the following table:

Adenovirus	Metapneumovirus	Bocavirus	Internal Control	Negative Control	Positive Control	Interpretation
+	+	+	+/-	-	+	Adenovirus, Metapneumovirus and Bocavirus Positives
-	-	-	+	-	+	Adenovirus, Metapneumovirus and Bocavirus Negatives
+	-	-	+/-	-	+	Adenovirus Positive, Metapneumovirus and Bocavirus Negatives
+	+	-	+/-	-	+	Adenovirus and Metapneumovirus Positives, Bocavirus Negative
+	-	+	+/-	-	+	Adenovirus and Bocavirus Positives, Metapneumovirus Negative
-	+	-	+/-	-	+	Metapneumovirus Positive, Adenovirus and Bocavirus Negatives
-	+	+	+/-	-	+	Metapneumovirus and Bocavirus Positives, Adenovirus Negative
-	-	+	+/-	-	+	Bocavirus Positive, Adenovirus and Metapneumovirus Negatives
+	+	+	+	+	+	Invalid
-	-	-	-	-	-	Invalid

**(+) Positive:** Amplification signal

**(-) Negative:** No amplification signal

If the negative samples do not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

## **Quality Control**

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

## **Performance evaluation**

### **Clinical sensitivity and specificity**

A total of 142 throat swabs from symptomatic patients were tested by Real Time PCR using Vitassay qPCR AMB and the results were compared with those obtained using FTD Adenovirus/Metapneumovirus/Bocavirus (Fast-Track). Adenovirus was detected in 68 samples with Vitassay qPCR AMB and in 69 samples using FTD Adenovirus/Metapneumovirus/Bocavirus (Fast-Track). This sample of difference could not be detected with Vitassay qPCR AMB because the amount of template DNA is below the detection limit of the method used.

Metapneumovirus was detected in 39 samples with Vitassay qPCR AMB and in 37 samples using FTD Adenovirus/Metapneumovirus/Bocavirus (Fast-Track). The amount of RNA template of this non-detected sample was below the detection limit of the method used.

Bocavirus was detected in 34 samples with Vitassay qPCR AMB and in 37 samples using FTD Adenovirus/Metapneumovirus/Bocavirus (Fast-Track). The amount of DNA template of these non-detected samples was below the detection limit of the method used.

These results indicate the high sensitivity and specificity to detect Adenovirus, Metapneumovirus and Bocavirus using the molecular diagnostic test Vitassay qPCR AMB.

### **Analytical sensitivity**

The analytical sensitivity was determined by analysis of 10-fold dilution series of Adenovirus, Metapneumovirus and Bocavirus templates ranging from  $10^7$  to  $10^1$  copies/rxn. This assay has a detection limit of  $\geq 10$  viral RNA/DNA copies per reaction.

### **Analytical specificity**

The analytical specificity for Adenovirus, Metapneumovirus and Bocavirus was tested within the panel of following virus, where no cross-reactivity was seen between any of the species:

Cross-reactivity assay		
Adenovirus	Metapneumovirus	Bocavirus
<i>Bordetella pertussis</i>	<i>Bordetella pertussis</i>	<i>Bordetella pertussis</i>
<i>Legionella bozemanii</i>	<i>Legionella bozemanii</i>	<i>Legionella bozemanii</i>
<i>Legionella micdadei</i>	<i>Legionella micdadei</i>	<i>Legionella micdadei</i>
<i>Legionella dumoffii</i>	<i>Legionella dumoffii</i>	<i>Legionella dumoffii</i>
<i>Legionella longbeachae</i>	<i>Legionella longbeachae</i>	<i>Legionella longbeachae</i>
<i>Legionella pneumophila</i>	<i>Legionella pneumophila</i>	<i>Legionella pneumophila</i>
<i>Mycoplasma pneumoniae</i>	<i>Mycoplasma pneumoniae</i>	<i>Mycoplasma pneumoniae</i>
<i>Streptococcus pneumoniae</i>	<i>Streptococcus pneumoniae</i>	<i>Streptococcus pneumoniae</i>
<i>Staphylococcus aureus</i> subsp. <i>Aureus</i>	<i>Staphylococcus aureus</i> subsp. <i>Aureus</i>	<i>Staphylococcus aureus</i> subsp. <i>Aureus</i>
Meticilin resistant <i>Stahylococcus aureus</i>	Meticilin resistant <i>Stahylococcus aureus</i>	Meticilin resistant <i>Stahylococcus aureus</i>
<i>Haemophilus influenzae</i>	<i>Haemophilus influenzae</i>	<i>Haemophilus influenzae</i>
<i>Moraxella catarrhalis</i>	<i>Moraxella catarrhalis</i>	<i>Moraxella catarrhalis</i>
Influenza B/Brisbane/60/2008	Influenza B/Brisbane/60/2008	Influenza B/Brisbane/60/2008
Influenza B/Florida/04/06	Influenza B/Florida/04/06	Influenza B/Florida/04/06
Influenza B/Phuket/3073/2013	Influenza B/Phuket/3073/2013	Influenza B/Phuket/3073/2013
Influenza A/Switzerland/9715293/2013	Influenza A/Switzerland/9715293/2013	Influenza A/Switzerland/9715293/2013
Influenza A/Turkey/Germany R2485+86/2014	Influenza A/Turkey/Germany R2485+86/2014	Influenza A/Turkey/Germany R2485+86/2014
Influenza A/California/7/2009(H1N1) virus	Influenza A/California/7/2009(H1N1) virus	Influenza A/California/7/2009(H1N1) virus
Influenza A/Perth/16/2009(H3N2) virus	Influenza A/Perth/16/2009(H3N2) virus	Influenza A/Perth/16/2009(H3N2) virus
Influenza A/New Caledonia/20/99(H1N1) virus	Influenza A/New Caledonia/20/99(H1N1) virus	Influenza A/New Caledonia/20/99(H1N1) virus
Human Rhinovirus	Human Rhinovirus	Human Rhinovirus
Human coronavirus 229E, OC43 and NL63	Human coronavirus 229E, OC43 and NL63	Human coronavirus 229E, OC43 and NL63
Human parainfluenza 1, 2, 3, 4	Human parainfluenza 1, 2, 3, 4	Human parainfluenza 1, 2, 3, 4
MERS Coronavirus	MERS Coronavirus	MERS Coronavirus
Respiratory syncytial virus (RSV)	Respiratory syncytial virus (RSV)	Respiratory syncytial virus (RSV)
Human Metapneumovirus A y B	Adenovirus	Adenovirus
Bocavirus	Bocavirus	Human Metapneumovirus A y B
Enterovirus D68 and A71	Enterovirus D68 and A71	Enterovirus D68 and A71
Echovirus 11 and 30	Echovirus 11 and 30	Echovirus 11 and 30
Coxsackievirus A9, A21, A24, B3 and B4	Coxsackievirus A9, A21, A24, B3 and B4	Coxsackievirus A9, A21, A24, B3 and B4

## Analytical reactivity

Analytical reactivity of Vitassay qPCR AMB for Adenovirus was evaluated against Human Adenovirus type 1, Human Adenovirus 2 strain Adenoid 6, Human Adenovirus Type 3, strain GB, Human Adenovirus Type 4, strain RI-67, Human Adenovirus Type 5, Human Adenovirus Type 6, Human Adenovirus Type 8, Human Adenovirus Type 15 strain CH. 38, Human Adenovirus Type 31 strain 1315/63, Human Adenovirus 40 strain Dugan and Human Adenovirus 41 strain Tak, showing positive result.

Analytical reactivity of Vitassay qPCR AMB for Metapneumovirus was evaluated against Human Metapneumovirus genotype A and Human Metapneumovirus genotype B, showing positive result.

Analytical reactivity of Vitassay qPCR AMB for Bocavirus was evaluated against Human Bocavirus showing positive result.

## Compatibles real-time PCR equipment

Vitassay qPCR AMB has been validated on the following equipments:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) <sup>II</sup>
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time PCR System (DNA-Technology)
- DTPrime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)<sup>I</sup>
- SmartCycler® (Cepheid)<sup>I</sup>

I: For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler® tubes.

II: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

## Limitations

- This test provides a presumptive diagnosis of Adenovirus, Metapneumovirus and Bocavirus infection. All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- This assay was tried with throat swab and CSF samples. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper RNA/DNA from clinical specimens must be extracted.

- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by the different virus, either samples containing high concentrations of target RNA/DNA or contamination due to PCR products from previous reactions.

## Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
<b>Agilent Technologies</b>	<b>Abbott</b>
AriaMx/AriaDx Real-Time PCR System	Abbott m2000 RealTime System
<b>Applied Biosystems</b>	<b>Applied Biosystems</b>
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Fast Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Fast Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR System
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
<b>BIONEER</b>	ViiA™ 7 Real-Time PCR System
Exicycler™ 96	<b>Analytik Jena Biometra</b>
<b>Bio-Rad</b>	TOptical
CFX96™/ CFX96™ IVD Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	<b>BIONEER</b>
<b>Cepheid</b>	Exicycler™ 96
SmartCycler®	<b>Bio-Rad</b>
<b>Qiagen</b>	CFX96™ / CFX96™ IVD Real-Time PCR Detection
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
<b>Roche</b>	iCycler iQ™ 5 Real-Time PCR
LightCycler®480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler®96 Real-Time PCR System	MyiQ™ 2Real-Time PCR Detection System
Cobas z480 Analyzer	<b>Cepheid</b>
	SmartCycler®
	<b>DNA-Technology</b>
	DTlite Real-Time PCR System
	DTprime Real-time Detection Thermal Cycler
	<b>Eppendorf</b>
	Mastercycler™ep <i>realplex</i>
	<b>Qiagen</b>
	Rotor-Gene® Q
	<b>Stratagene / Agilent Technologies</b>
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

\* See Attached III to configure exposure settings.

## Attached II: Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

REAL-TIME PCR THERMOCYCLER	Vitassay CHANNEL	DETECTION CHANNEL	OBSERVATIONS
<b>Bio-Rad CFX96™</b>	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>ABI 7500 Applied Biosystems</b>	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>Roche Lightcycler®480II</b>	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
<b>Smartcycler® Cepheid</b>	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
<b>Abbott m2000rt</b>	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>Mx3000P™ Mx 3005P™ Stratagene</b>	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>AriaMx Agilent</b>	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>Rotor-Gene® Q Qiagen</b>	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
<b>Exicycler™ 96 BIONEER</b>	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

### Attached III: Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
<b>DTlite Real-Time PCR System (DNA-Technology)</b>	FAM	500
	HEX	500
	ROX	500
	Cy5	500
<b>DTprime Real-time Detection Thermal Cycler (DNA-Technology)</b>	FAM	500
	HEX	1000
	ROX	1000
	Cy5	1000

## Bibliography/Bibliografía

1. David Metzgar, *et al.* PCR Analysis of Egyptian Respiratory Adenovirus Isolates, Including Identification of Species, Serotypes, and Coinfections. *Journal of Clinical Microbiology*, 2005; 5743-5752.
2. Gregory C. Gray, *et al.* Genotype Prevalence and Risk Factors for Severe Clinical Adenovirus Infection, United States 2004-2006. *Clinical Infectious Diseases*, 2007; 45(9): 1120-1131.
3. Mbayame Ndiaye Niang, *et al.* Respiratory viruses in patients with influenza-like illness in Senegal: Focus on human respiratory adenoviruses. *PLOS ONE*, 2017; 12 (3): e0174287.
4. Gregory C. Gray, *et al.* Adenovirus Transmission—Worthy of Our Attention. *The Journal of Infectious Diseases*, 2006; 194 (7): 871-873.
5. Magdalena Kendall Scott, *et al.* Human Adenovirus Associated with Severe Respiratory Infection, Oregon, USA, 2013-2014. *Emerging Infectious Diseases*, 2016 Jun; 102-108.
6. John V. Williams, *et al.* Human Metapneumovirus and Lower Respiratory Tract Disease in Otherwise Healthy Infants and Children. *The New England Journal of Medicine*, 2004; 350: 443-450.
7. Guy Boivin, *et al.* Virological Features and Clinical Manifestations Associated with Human Metapneumovirus: A New Paramyxovirus Responsible for Acute Respiratory-Tract Infections in All Age Groups. *The Journal of Infectious Diseases*, 2002; 186 (9): 1330-1334.
8. Wenhua Kong, *et al.* Circulation of Human Metapneumovirus Among Children With Influenza-Like Illness in Wuhan, China. *Journal of Medical Virology*, 2016; 88: 774-781.
9. Ann R. Falsey, *et al.* Human Metapneumovirus Infections in Young and Elderly Adults. *The Journal of Infectious Diseases*, 2003; 187 (5): 785-790.
10. Bernadette G. van den Hoogen, *et al.* Prevalence and Clinical Symptoms of Human Metapneumovirus Infection in Hospitalized Patients. *The Journal of Infectious Diseases*, 2003; 188 (10): 1571-1577.
11. Sai-Zhen Zeng, *et al.* Clinical Features of Human Metapneumovirus Genotypes in Children With Acute Lower Respiratory Tract Infection in Changsha, China. *Journal of Medical Virology*, 2015; 87: 1839-1845.
12. Amber K. Haynes, *et al.* Human Metapneumovirus Circulation in the United States, 2008 to 2014. *PEDIATRICS*, 2016; Volume 137: number 5.
13. Tobias Allander, *et al.* Human Bocavirus and Acute Wheezing in Children. *Clinical Infectious Diseases*, 2007; 44 (7): 904-910.
14. Nathalie Bastien, *et al.* Human Bocavirus Infection, Canada. *Emerging Infectious Diseases*, 2006; 12 (5): 848-850.

15. Xiaoming Ma, et al. Detection of Human Bocavirus in Japanese Children with Lower Respiratory Tract Infections. *Journal of Clinical Microbiology*, 2006; 44 (3): 1132-1134.
16. Diego Vicente, et al. Human Bocavirus, a Respiratory and Enteric Virus. *Emerging Infectious Diseases*, 2007; 13 (4): 636-637.
17. Xiaochun Wang, et al. Complete Genomes of Three Human Bocavirus Strains from Children with Gastroenteritis and Respiratory Tract Illnesses in Jiangsu, China. *Journal of Virology*, 2012; 86 (24): 13826-13827.
18. Marcello Guido, et al. Human bocavirus: Current knowledge and future challenges. *World Journal of Gastroenterology*, 2016; 22 (39): 8684-8697.

## Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.


Mx3000P™ and Mx3029™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid.

## Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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