Vitassay qPCR

Blastocystis hominis + Dientamoeba fragilis

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- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- StepOne™ Real-Time PCR System (Applied Biosystems) ^{II}
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)

- DTlite Real-Time PCR System (DNA-Technology)
- DTPrime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)¹
- SmartCycler® (Cepheid)¹

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| Termocicladores con bloque de bajo perfil | Termocicladores con bloque de alto perfil |
|--|---|
| Agilent Technologies | Abbott |
| AriaMx Real-Time PCR System | Abbott m2000 RealTime System |
| Applied Biosystems | Applied Biosystems |
| 7500 Fast Real-Time PCR System | 7300 Real-Time PCR System |
| 7500 Fast Dx Real-Time PCR System | 7500 Real-Time PCR System |
| QuantStudio™ 12K Flex 96-well Fast | 7900 HT Real-Time PCR System |
| QuantStudio™ 6 Flex 96-well Fast | ABI PRISM 7000 |
| QuantStudio™ 7 Flex 96-well Fast | ABI PRISM 7700 |
| QuantStudio™ 5 Real-Time PCR System | QuantStudio™ 12K Flex 96-well |
| QuantStudio™ 3 Real-Time PCR System | QuantStudio™ 6 Flex 96-well |
| StepOne Plus™ Real-Time PCR System | QuantStudio™ 7 Flex 96-well |
| StepOne™ Real-Time PCR System | QuantStudio™ 5 Real-Time PCR |
| ViiA™ 7 Fast Real-Time PCR System | QuantStudio™ 3 Real-Time PCR System |
| BIONEER | ViiA™ 7 Real-Time PCR |
| Exicycler™ 96 | Analytik Jena Biometra |
| Bio-Rad | TOptical |
| CFX96™ Real-Time PCR Detection System | qTOWER.2.0 |
| Mini Opticon™ Real-Time PCR Detection System | BIONEER |
| Cepheid | Exicycler™ 96 |
| SmartCycler® | Bio-Rad |
| Qiagen | CFX96™ Deep Well Real-Time PCR Detection |
| Rotor-Gene® Q | iCycler iQ™ Real-Time PCR |
| Roche | iCycler iQ™ 5 Real-Time PCR |
| LightCycler® 480 Real-Time PCR System | MyiQ™ Real-Time PCR Detection System |
| LightCycler® 96 Real-Time PCR System | MyiQ™ 2Real-Time PCR Detection System |
| Cobas z480 Analyzer | Cepheid |
| | SmartCycler® |
| | DNA-Technology |
| | DTlite Real-Time PCR System* |
| | DTprime Real-time Detection Thermal Cycler* |
| | Eppendorf |
| | Mastercycler™ ep <i>realplex</i> |
| | Qiagen |
| | Rotor-Gene® Q |
| | Stratagene / Agilent Technologies |
| | Mx3000P™ Real Time PCR System |
| | Mx3005P™ Real Time PCR System |

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Intended use

Vitassay qPCR *Blastocystis hominis* + *Dientamoeba fragilis* allows the detection and differentiation of *Blastocystis hominis* and *Dientamoeba fragilis* by real-time PCR in human stool samples. The product is intended for use in the diagnosis of *Blastocystis hominis* and/or *Dientamoeba fragilis* infections alongside clinical data of the patient and other laboratory tests outcomes.

References

Vitassay qPCR *Blastocystis hominis* + *Dientamoeba fragilis* 4x8-well strip, low profile
7041030

Vitassay qPCR *Blastocystis hominis* + *Dientamoeba fragilis* 4x8-well strip, high profile
7042030

Materials/reagents provided

| Reference | Reagent/Material | Colour | Amount |
|-----------------------|---|--------|-----------------|
| 7041S030/ 7042S030 | <i>Blastocystis hominis</i> + <i>Dientamoeba fragilis</i> strips low/high profile | - | 4x8-well strip |
| 7C030 | <i>Blastocystis hominis</i> + <i>Dientamoeba fragilis</i> Positive Control | red | 1 vial |
| 7001A | PCR grade water | white | 1 vial x 1 mL |
| 7002B | Resuspension Buffer | green | 1 vial x 1,8 mL |
| 7003N | Negative control | yellow | 1 vial x 1 mL |
| 7004O | Optical caps | - | 4x8 cap strip |

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

Additional equipment and material required

- DNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Vortexer
- Micropipettes (1-20 µL, 20-200 µL)

- Filter tips
- Powder-free disposal gloves

Summary

Blastocystis hominis and *Dientamoeba fragilis* are frequently encountered together. They are both associated with symptoms such as: abdominal pain, bloating, cramps, discomfort, nausea or alteration of bowel habits, resembling irritable bowel syndrome. Although *Dientamoeba fragilis* can be symptomatic or asymptomatic, when symptomatic might be ignored due to the presence of other suspect pathogens. For *Blastocystis hominis* it is truly unknown whether it is a pathogen, a commensal or a circumstantial pathogen. Recently, various groups of *Blastocystis* have been isolated from humans, each proceeding from a different host, these are: *B. galli* (from chickens), *B. anatis* (from ducks), *B. anseri* (from geese) and *B. lapemi* (from sea snakes). This has raised the possibility for more than one species of *Blastocystis* infecting humans. The principal problem when determining whether *Blastocystis hominis* is a pathogen or not, is the unknown pathogenicity of other concomitant organisms. It is believed that both *Blastocystis* and *Dientamoeba* are transmitted through the oral-fecal route. *Dientamoeba fragilis* is characterized for having a binucleate nucleus and two SSU rDNA genotypes, genotype 1 and genotype 2. Given the fragile nature of *Dientamoeba fragilis*, the binucleate structure cannot be appreciated correctly through saline preparations. And once its trophozoites degenerate they become difficult to recognize. Real Time PCR has demonstrated excellent sensitivity and specificity when diagnosing *Dientamoeba fragilis* and *Blastocystis hominis*, through the targeting of the 5,8S rRNA gene and 18S RNA gene, respectively.

Principle of the test

Vitassay qPCR Blastocystis hominis + Dientamoeba fragilis test is based on the real-time amplification of a conserved region of the 5,8S rRNA gene for the identification of *Blastocystis hominis* and 18S rRNA gene for the identification of *Dientamoeba fragilis*. After DNA isolation, the presence of *Blastocystis hominis* and *Dientamoeba fragilis* is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

Vitassay qPCR Blastocystis hominis + Dientamoeba fragilis test is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The amplification of the *Dientamoeba fragilis* DNA target sequence is detected through the FAM channel and *Blastocystis hominis* DNA target in ROX channel whereas the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Attached II).

Precautions

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or broken or if the foil has been broken or damaged.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

Procedures

Specimen collection, processing and DNA extraction

For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

Maxwell® RSC Blood DNA Kit, Maxwell® 16 instrument (Promega).

QIAamp DNA Mini kit (QIAGEN).

QIAamp DNA Stool Mini Kit (QIAGEN).

RIDA® Xtract (r-Biopharm).

Invisorb® Spin Universal Kit (Stratec).

UltraClean® Tissue & Cells DNA Isolation Kit (Mobio).

Positive control preparation

Reconstitute the lyophilized *Blastocystis hominis* + *Dientamoeba fragilis* Positive Control (red tube) with 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of DNA sample, negative and positive controls and add them into each well.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

| Step | Temperature | Time | Cycles |
|---|-------------|--------|--------|
| Initial denaturation | 95°C | 2 min | 1 |
| Denaturation | 95°C | 10 sec | 45 |
| Annealing/Extension (Data collection*) | 60°C | 50 sec | |

Set the fluorescence data collection during the extension step (*) through the FAM (*Dientamoeba fragilis*), ROX (*Blastocystis hominis*) and HEX, JOE or VIC channels (Internal Control (IC)). If you use the Applied Biosystems 7500 Fast or the Stratagene Mx3005P™ check that passive reference option ROX is none (Attached II).

Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive controls used in each run, must show an amplification curve for FAM (*Dientamoeba fragilis*) and ROX (*Blastocystis hominis*), which validates the reaction.

Negative control

The negative controls included in each run, must show the absence of signal in ROX and FAM which validates the reaction.

The experiment seems to be failed if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

The result interpretation is summarized in the following table:

| <i>Blastocystis hominis</i> (ROX) | <i>Dientamoeba fragilis</i> (FAM) | Internal Control | Negative Control | Positive Control | Interpretación |
|-----------------------------------|-----------------------------------|------------------|------------------|------------------|--|
| + | + | +/- | - | + | <i>Blastocystis hominis</i> and <i>Dientamoeba fragilis</i> Positives |
| - | - | + | - | + | <i>Blastocystis hominis</i> and <i>Dientamoeba fragilis</i> Negatives |
| + | - | +/- | - | + | <i>Blastocystis hominis</i> Positive, <i>Dientamoeba fragilis</i> Negative |
| - | + | +/- | - | + | <i>Dientamoeba fragilis</i> Positive, <i>Blastocystis hominis</i> Negative |
| + | + | + | + | + | Experiment fail |
| - | - | - | - | - | Experiment fail |

(+) **Positive:** Amplification signal

(-) **Negative:** No amplification signal

If the negative samples not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

Quality Control

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

A total of 299 faecal specimens from symptomatic patients were tested by Real Time PCR using Vitassay qPCR Blastocystis hominis + Dientamoeba fragilis and RIDA®GENE Dientamoeba fragilis (R-biopharm) and these results were compared with those obtained by routine parallel diagnosis of parasites using microscopy for *Blastocystis hominis* and *Dientamoeba fragilis*.

Vitassay qPCR Blastocystis hominis + Dientamoeba fragilis detected *Blastocystis hominis* in 47 positive samples, 23 of which were evaluated and confirmed as positive for *Blastocystis hominis* by microscopic examination.

Vitassay qPCR Blastocystis hominis + Dientamoeba fragilis detected *Dientamoeba fragilis* in 79 positive samples, 72 of them were evaluated using PCR-in house, with 69 samples showing positive result for *Dientamoeba fragilis*. In addition, 45 samples were evaluated and confirmed as positive by RIDA®GENE Dientamoeba fragilis (R-biopharm)

Vitassay qPCR Blastocystis hominis + Dientamoeba fragilis shows a high sensitivity and specificity to detect *Blastocystis hominis* and *Dientamoeba fragilis*.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of *Blastocystis hominis* and *Dientamoeba fragilis* templates ranging from 10^7 to 10^1 copies/rxn. This assay has a detection limit of ≥ 10 DNA copies per reaction.

Analytical specificity

The analytical specificity for *Blastocystis hominis* and *Dientamoeba fragilis* was tested within the panel of different microorganisms.

No cross-reactivity of *Blastocystis hominis* and *Dientamoeba fragilis* was seen between any of the species:

| Cross-reactivity assay | | |
|---|--|--|
| <i>Shigella flexneri</i> | <i>Proteus vulgaris</i> | <i>Enterococcus faecalis</i> |
| <i>Shigella dysenteriae</i> | <i>Aeromonas hydrophila</i> subsp. <i>hydrophila</i> | <i>Yersinia enterocolitica</i> O:3 |
| <i>Salmonella typhi</i> | <i>Citrobacter freundii</i> | <i>Yersinia enterocolitica</i> O:9 |
| <i>Salmonella paratyphi A</i> | <i>Staphylococcus aureus</i> subsp. <i>aureus</i> | <i>Bacteroides fragilis</i> |
| <i>Salmonella paratyphi B</i> | <i>Serratia liquefaciens</i> | Adenovirus serotype 40 |
| <i>Salmonella typhimurium</i> | <i>Vibrio parahaemolyticus</i> | Adenovirus serotype 41 |
| <i>Salmonella bongori</i> | <i>Clostridium difficile</i> | Rotavirus A |
| <i>Salmonella enteritidis</i> | <i>Clostridium perfringens</i> | Norovirus Genotype I y II |
| <i>Salmonella enterica</i> subsp. <i>entérica</i> | <i>Enterotoxigenic Escherichia coli</i> | Astrovirus Genotype I-VIII |
| <i>Salmonella pullorum</i> | <i>Enteropathogenic Escherichia coli</i> | <i>Campylobacter lari</i> |
| <i>Salmonella gallinarum</i> | <i>Klebsiella oxytoca</i> | <i>Campylobacter fetus</i> |
| <i>Helicobacter pylori</i> | <i>Listeria monocytogenes</i> | <i>Campylobacter coli</i> |
| <i>Helicobacter hepaticus</i> | <i>Candida albicans</i> | <i>Campylobacter jejuni</i> subsp. <i>jejuni</i> |
| <i>Helicobacter cinaedi</i> | <i>Arcobacter butzleri</i> | <i>Campylobacter upsaliensis</i> |
| <i>Helicobacter heilmannii</i> | <i>Giardia intestinalis</i> | |
| <i>Pseudomonas aeruginosa</i> | <i>Entamoeba histolytica</i> | |

Analytical reactivity

Analytical reactivity of Vitassay qPCR *Blastocystis hominis* + *Dientamoeba fragilis* for *Blastocystis hominis* was evaluated against *Blastocystis hominis*, showing positive results.

Analytical reactivity of Vitassay qPCR *Blastocystis hominis* + *Dientamoeba fragilis* for *Dientamoeba fragilis* was evaluated against *Dientamoeba fragilis*, showing positive results.

Compatibles real-time PCR equipment

Vitassay qPCR *Blastocystis hominis* + *Dientamoeba fragilis* has been validated on the following equipments:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- StepOne™ Real-Time PCR System (Applied Biosystems) ^{II}
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time PCR System (DNA-Technology)
- DTPRime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)^I
- SmartCycler® (Cepheid)^I

I: For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler tubes.

II: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

Limitations

- This test provides a presumptive diagnosis of *Blastocystis hominis* and *Dientamoeba fragilis* infection. All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- This assay was tried with human faecal samples. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper DNA from clinical specimens must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by the different parasites, either samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

| Low profile Block Thermocyclers |
|--|
| Agilent Technologies |
| AriaMx Real-Time PCR System |
| Applied Biosystems |
| 7500 Fast Real-Time PCR System |
| 7500 Fast Dx Real-Time PCR System |
| QuantStudio™ 12K Flex 96-well Fast |
| QuantStudio™ 6 Flex 96-well Fast |
| QuantStudio™ 7 Flex 96-well Fast |
| QuantStudio™ 5 Real-Time PCR System |
| QuantStudio™ 3 Real-Time PCR System |
| StepOne Plus™ Real-Time PCR System |
| StepOne™ Real-Time PCR System |
| ViiA™ 7 Fast Real-Time PCR System |
| BIONEER |
| Exicycler™ 96 |
| Bio-Rad |
| CFX96™ Real-Time PCR Detection System |
| Mini Opticon™ Real-Time PCR Detection System |
| Cepheid |
| SmartCycler® |
| Qiagen |
| Rotor-Gene® Q |
| Roche |
| LightCycler®480 Real-Time PCR System |
| LightCycler®96 Real-Time PCR System |
| Cobas z480 Analyzer |

| High profile Block Thermocyclers |
|--|
| Abbott |
| Abbott m2000 RealTime System |
| Applied Biosystems |
| 7300 Real-Time PCR System |
| 7500 Real-Time PCR System |
| 7900 HT Real-Time PCR System |
| ABI PRISM 7000 |
| ABI PRISM 7700 |
| QuantStudio™ 12K Flex 96-well |
| QuantStudio™ 6 Flex 96-well |
| QuantStudio™ 7 Flex 96-well |
| QuantStudio™ 5 Real-Time PCR |
| QuantStudio™ 3 Real-Time PCR System |
| ViiA™ 7 Real-Time PCR |
| Analytik Jena Biometra |
| TOptical |
| qTOWER 2.0 |
| BIONEER |
| Exicycler™ 96 |
| Bio-Rad |
| CFX96™ Deep Well Real-Time PCR Detection |
| iCycler iQ™ Real-Time PCR |
| iCycler iQ™ 5 Real-Time PCR |
| MyiQ™ Real-Time PCR Detection System |
| MyiQ™ 2Real-Time PCR Detection System |
| Cepheid |
| SmartCycler® |
| DNA-Technology |
| DTlite Real-Time PCR System |
| DTprime Real-time Detection Thermal Cycler |
| Eppendorf |
| Mastercycler™ ep <i>realplex</i> |
| Qiagen |
| Rotor-Gene® Q |
| Stratagene / Agilent Technologies |
| Mx3000P™ Real Time PCR System |
| Mx3005P™ Real Time PCR System |

* See Attached III to configure exposure settings.

Attached II: Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

| REAL-TIME PCR THERMOCYCLER | Vitassay CHANNEL | DETECTION CHANNEL | OBSERVATIONS |
|--------------------------------------|------------------|-------------------|---|
| Bio-Rad CFX96™ | FAM | FAM | |
| | HEX | HEX | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |
| ABI 7500 Applied Biosystems | FAM | FAM | Passive reference option ROX is none |
| | HEX | VIC | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |
| Roche Lightcycler®480II | FAM | 465/510 | Colour Compensation required |
| | HEX | 533/580 | |
| | ROX | 533/610 | |
| | Cy5 | 618/660 | |
| Smartcycler® Cepheid | FAM | Channel 1 | |
| | HEX | Channel 2 | |
| | ROX | Channel 3 | |
| | Cy5 | Channel 4 | |
| Abbott m2000rt | FAM | FAM | |
| | HEX | VIC | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |
| Mx3000P™ Mx 3005P™ Stratagene | FAM | FAM | Passive reference option ROX is none |
| | HEX | VIC | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |
| AriaMx Agilent | FAM | FAM | |
| | HEX | HEX | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |
| Rotor-Gene®Q Qiagen | FAM | Green | In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition". |
| | HEX | Yellow | |
| | ROX | Orange | |
| | Cy5 | Red | |
| Exicycler™ 96 BIONEER | FAM | FAM | |
| | HEX | JOE | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |

Attached III: Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

| Thermocycler | Vitassay channel | Exposure values |
|--|------------------|-----------------|
| DTlite Real-Time PCR System (DNA-Technology) | FAM | 500 |
| | HEX | 500 |
| | ROX | 500 |
| | Cy5 | 500 |
| DTprime Real-time Detection Thermal Cycler (DNA-Technology) | FAM | 500 |
| | HEX | 1000 |
| | ROX | 1000 |
| | Cy5 | 1000 |

Bibliography/Bibliografía

1. Prospective Study of the Prevalence, Genotyping, and Clinical Relevance of *Dientamoeba fragilis* Infections in an Australian Population. D. Stark, et al. *Journal of Clinical Microbiology*, 2005; 2718-2723.
2. A review of *Dientamoeba fragilis* carriage in humans: Several reasons why this organism should be considered in the diagnosis of gastrointestinal illness. Joel L. N. Barrat, et al. *Gut Microbes*, 2011; 2:1: 3-12.
3. Emerging from Obscurity: Biological, Clinical, and Diagnostic Aspects of *Dientamoeba fragilis*. Eugene H. Johnson, et al. *Clinical Microbiology Reviews*, 2004: 553-570.
4. *Blastocystis hominis* and *Dientamoeba fragilis* in patients fulfillin irritable bowel syndrome criteria. Javed Yakoob, et al. *Parasitology Research*, 2010; 107: 679-684.
5. Epidemiology and Pathogenicity of *Blastocystis hominis*. Patrick W. Doyle, et al. *Journal of Clinical Microbiology*, 1990; 116-121.
6. Clinical Significance of *Blastocystis hominis*. S. M. Hussain Qadri, et al. *Journal of Clinical Microbiology*, 1989; 2407-2409.
7. *Blastocystis hominis* Revisited. D. J. Stenzel, et al. *Clinical Microbiology Reviews*, 1996; 563-584.

Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™, StepOnePlus™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.


Mx3000P™ and Mx3029™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid.

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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