# Vitassay qPCR

## Panel I / Respiratory viruses

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Termocicladores con bloque de bajo perfil	Termocicladores con bloque de alto perfil
<b>Agilent Technologies</b>	<b>Abbott</b>
AriaMx Real-Time PCR System	Abbott m2000 RealTime System
<b>Applied Biosystems</b>	<b>Applied Biosystems</b>
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
<b>BIONEER</b>	ViiA™ 7 Real-Time PCR
Exicycler™ 96	<b>Analytik Jena Biometra</b>
<b>Bio-Rad</b>	TOptical
CFX96™ Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	<b>BIONEER</b>
<b>Cepheid</b>	Exicycler™ 96
SmartCycler®	<b>Bio-Rad</b>
<b>Qiagen</b>	CFX96™ Deep Well Real-Time PCR Detection
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
<b>Roche</b>	iCycler iQ™ 5 Real-Time PCR
LightCycler® 480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler® 96 Real-Time PCR System	MyiQ™ 2 Real-Time PCR Detection System
Cobas z480 Analyzer	<b>Cepheid</b>
	SmartCycler®
	<b>DNA-Technology</b>
	DTlite Real-Time PCR System*
	DTprime Real-time Detection Thermal Cycler*
	<b>Eppendorf</b>
	Mastercycler™ ep <i>realplex</i>
	<b>Qiagen</b>
	Rotor-Gene® Q
	<b>Stratagene / Agilent Technologies</b>
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

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## Intended use

Vitassay qPCR Panel I / Respiratory viruses allows the detection and differentiation of Flu A, Flu B, and Respiratory Syncytial (RSV) viruses and subtyping of Flu A (H1N1)pdm09, H3N2, H5N1, and H7N9 by real-time RT-PCR in clinical samples. The product is intended for use in the diagnosis of Flu A, Flu B and/or RSV in combination with clinical and epidemiological risk factors and evaluation of infections with Influenza A subtypes: (H1N1)pdm09, H3N2, H5N1, and H7N9 infections alongside clinical data of the patient and other laboratory tests outcomes.

## References

Vitassay qPCR Panel I / Respiratory viruses 8x8 -well strip, low profile 7041033

Vitassay qPCR Panel I / Respiratory viruses 8x8-well strip, high profile 7042033

## Materials/reagents provided

Reference	Reagent/Material	Colour	Amount
7041S027/ 7042S027	Flu A + Flu B + RSV strips low/high profile	-	4 x 8-well strip
7041S032/ 7042S032	Type II / Flu strips low/high profile	-	4 x 8-well strip
7011S032C/ 7012S032C	Neuroamidinase confirmation strips low/high profile	-	1 x 8-well strip
7C033	Panel I / Respiratory viruses Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	9 x 8 cap strip

## Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

## Additional equipment and material required

- RNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Vortex
- Micropipettes (1-20  $\mu$ L, 20-200  $\mu$ L)
- Filter tips
- Powder-free disposal gloves

## Summary

Influenza A and B viruses, which belong to the *Orthomyxoviridae* family, cause the majority of viral lower respiratory tract infections; elderly and compromised individuals are especially at risk of developing severe illness and complications. Influenza viruses are a significant cause of morbidity and mortality worldwide.

Seasonal influenza epidemics impose a heavy burden on society, with 3–5 million cases and 250-500 000 deaths worldwide every year. The resulting economic impact is large and includes both direct and indirect costs. Also, influenza A can be divided into numerous subgroups, but the most important ones are the ones that can cause illness in humans, them being, H1N1, H1N2, H2N2, H3N2 and H5N1.

Influenza can be transmitted either by maintaining direct contact with infected individuals, contaminated objects or by inhalation of contaminated aerosols. After an incubation period of one to two days, the illness has an abrupt onset, with chills, fever, myalgia, headache, and anorexia. Pneumonia is the most important complication of influenza. There may be primary viral pneumonia, secondary bacteria pneumonia, or a combination of both.

The human respiratory syncytial virus (RSV) is an enveloped, single stranded linear RNA genome virus. RSV, which belongs to the genus *Orthopneumovirus*, is divided into two major groups, A and B, based on antigenic and genomic differences.

RSV infection most commonly causes a cold-like illness. But it can also cause bronchitis, croup, and lower respiratory infections like bronchiolitis and pneumonia. Of every 100 infants and young children with RSV infection, 25 to 40 (25% to 40%) will show signs of pneumonia or bronchiolitis.

RSV is transmitted via large nasopharyngeal secretion droplets coming from infected individuals. These droplets enter via the mucus membranes of the eyes, nose and mouth following close contact, or self-inoculation after touching contaminated surfaces.

Several different types of laboratory tests are available for diagnosis of Influenza or RSV infection. Laboratory diagnosis of influenza has been accomplished using cell culture, rapid antigen testing, and, more recently, real-time PCR, concretely qPCR that targets the *M1* gene. Whereas for RSV, antigen detection tests and cultures are generally reliable when diagnosing young children, but they are less useful in older children and adults. Use of highly sensitive RT-PCR assays, such as the ones that target the N-gene, should be considered, particularly, when testing older children and adults, as they may have low viral loads in the obtained respiratory specimens.

Of the total positive samples found for Influenza A, 72% were confirmed as subtype H1N1) pdm09 and 28% as subtype H3N2. The number of positive samples for the other two subtypes (H5N1 and H7N9) is decreasing with the years, the last reported case of H5N1 was in 2017 and that of H7N9 in 2018. Therefore, the number of times it will be necessary to perform the Confirmation PCR with the *Neuroaminidase* confirmation strip will be less than 10% of Influenza cases detected with this test.

### **Principle of the test**

Vitassay qPCR Panel I / Respiratory viruses test is based on the real-time amplification of specific conserved fragments of the *M1* gene for Flu A and Flu B and the *N* gene for RSV (Flu A, Flu B & RSV, strip 1), as well as, the *hemagglutinin* gene for Influenza A subtyping (H1N1)pdm09, H3N2, H5N1 and H7N9 (Type II / Flu, strip 2) and of specific conserved fragments of the *neuroaminidase* gene for Influenza A subtyping H5N1 and H7N9 (*Neuroaminidase* confirmation, strip 3). The viral RNA extracted is transcribed into cDNA using a specific primer by reverse transcription step followed immediately in by polymerase chain reaction. The presence of the RNA viral is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

Vitassay qPCR Panel I / Respiratory viruses test is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. Each kit includes three kind of strips and each one corresponds to one different assay. The first strip contains the multiplex reaction mix for the detection of Flu A, Flu B and RSV. Influenza A RNA targets are amplified and detected in channel FAM, Influenza B RNA in channel ROX, and RSV RNA in channel Cy5. The second strip contains the multiplex reaction mix for the detection of subtypes (H1N1)pdm09 and H3N2 and the screening of H5N1 and H7N9 (Type II / Flu, strip 2). The amplification of the Influenza A (H1N1)pdm09 RNA target sequence is detected through the FAM channel, H3N2 RNA target in ROX channel, H7N9 RNA target in Cy5 channel and H5N1 RNA target in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Attached 2). The third strip contains the multiplex reaction mix for the confirmation of the H5N1 and H7N9 subtypes (*Neuroaminidase* confirmation, strip 3). This reaction mix should be used when amplification for H5N1 and H7N9 in the second reaction mix is observed. The amplification of the H5N1 RNA target sequence is

detected through the FAM channel, the H7N9 RNA target in ROX channel and the internal control (IC) is detected in HEX, VIC or JOE channel (depending on the equipment used).

## **Precautions**

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or broken or if the foil has been broken or damaged.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

## **Procedures**

### **Specimen collection, processing and RNA extraction**

For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

RIDA® Xtract (r-Biopharm).

Maxwell® 16 Viral Total Nucleic Acid Purification Kit (Promega).

Total Nucleic Acid Isolation (TNAI) Kit (Roche).

Invisorb® Spin Universal Kit (Stratec).

## Positive control preparation

Reconstitute the lyophilized Panel I / Respiratory viruses Positive Control (red tube) in the 200 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

## Reaction setup

The third strip (*Neuroaminidase* confirmation, strip 3) should only be used when amplification of H5N1 and / or H7N9 is observed in the second strip. If amplification of these subtypes is not observed in the second strip, do not use the third strip. If you observe amplification of these subtypes proceed as indicated below. In addition to this, do not repeat the analysis performed with the first reaction mix

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run. (the wells of the third strip (*Neuroaminidase* confirmation, strip 3) only if necessary, as explained above).
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of RNA sample, negative and positive controls and add them into each well.
- Cover the wells with the caps provided. Spin down briefly(optional).
- Place the strips in the Real-time PCR instrument.

## Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Reverse transcription	45°C	15 min	1
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (\*) in the first strip through the FAM (Influenza A), HEX, JOE or VIC (Internal Control (IC)), ROX (Influenza B), and Cy5 channels (RSV); in the second strip through the FAM (Influenza A H1N1), HEX, JOE or VIC (H5N1), ROX (H3N2), and Cy5 channels (H7N9) and in the third reaction mix through the FAM (H5N1), ROX (H7N9) and HEX, JOE or VIC (Internal Control (IC)). If you use the Applied Biosystems 7500 Fast or the Stratagene Mx3005P™ check that passive reference option ROX is none. (Attached II)

### **Analysis and interpretation of results**

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

#### **Positive control**

The positive controls used in each run, must show an amplification curve in the channels of the respiratory virus panel.

#### **Negative control**

The negative controls included in each run, must show the absence of signal in the well of the negative control.

The experiment seems to be failed if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

#### **Interpretation of results for Flu A + Flu B + RSV 8-well strips:**

- A sample is considered positive for Influenza A if there is an amplification signal in FAM channel, the Ct value obtained is less than 40 and the internal control shows or not an amplification signal. Sometimes, the detection of internal control is not necessary because a high copy number of targets can cause preferential amplification of target-specific nucleic acids.
- A sample is considered positive for Influenza B if there is an amplification signal in ROX channel, the Ct value obtained is less than 40 and the internal control shows or not an amplification signal. Sometimes, the detection of internal control is not necessary because a high copy number of targets can cause preferential amplification of target-specific nucleic acids.
- A sample is considered positive for RSV if there is an amplification signal in CY5 channel, the Ct value obtained is less than 40 and the internal control shows or not an amplification signal. Sometimes, the detection of internal

control is not necessary because a high copy number of targets can cause preferential amplification of target-specific nucleic acids.

- A sample is considered negative, if the sample shows no amplification signal in the detection system but the internal control is positive (HEX channel). An inhibition of the PCR reaction can be excluded by the amplification of internal control.
- The experiment is considered failed if there is an amplification signal in the Negative Control well and/or there is not amplification signal in the Positive Control well. We recommend to repeat the assay again.
- In case of absence of internal control signal in sample wells we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of inhibition.

**Interpretation of results for Type II / Flu 8-well strips:**

(H1N1) pdm09	H5N1	H3N2	H7N9	Negative Control	Positive Control	Interpretation
+	+	+	+	-	+	Influenza A(H1N1)pdm09 and H3N2 subtypes positives; H5N1 and H7N9 subtypes possible positive*
-	-	-	-	-	+	Influenza A (H1N1)pdm09, H3N2, H5N1, and H7N9 subtypes negatives
+	-	-	-	-	+	Influenza A (H1N1)pdm09 positive;H3N2, H5N1 and H7N9 negatives
+	+	-	-	-	+	Influenza A (H1N1)pdm09 positive and H5N1 possible positive; H3N2 and H7N9 negatives*
+	-	+	-	-	+	Influenza A (H1N1)pdm09 and H3N2 positives;H5N1 and H7N9 negatives
+	-	-	+	-	+	Influenza A (H1N1)pdm09 positive and H7N9 possible positive; H3N2 and H5N1 negatives*
+	+	+	-	-	+	Influenza A (H1N1)pdm09 and H3N2 positives and H5N1 possible positive; H7N9 negative*
+	+	-	+	-	+	Influenza A (H1N1)pdm09 positive, H5N1 and H7N9 possible positive; H3N2 negative*
+	-	+	+	-	+	Influenza A (H1N1)pdm09 and H3N2 positives, H7N9 possible positive; H5N1 negative*
-	+	-	-	-	+	Influenza A H5N1 possible positive; Influenza A (H1N1)pdm09, H3N2 and H7N9 negatives*
-	+	+	-	-	+	Influenza A H5N1 possible positive and H3N2 positive; Influenza (H1N1)pdm09 and H7N9 negatives*
-	+	-	+	-	+	Influenza A H5N1 and H7N9 possible positive;(H1N1)pdm09 and H3N2 negatives*
-	+	+	+	-	+	Influenza A H3N2 Positive, H5N1 and H7N9 possible positive; and A (H1N1)pdm09 negative*
-	-	+	-	-	+	Influenza A H3N2 positive; (H1N1)pdm09, H5N1, and H7N9 negatives
-	-	+	+	-	+	Influenza A H3N2 positive and H7N9 possible positive; (H1N1)pdm09 and H5N1 negatives*
-	-	-	+	-	+	Influenza A H7N9 possible positive; (H1N1)pdm09, H3N2, and H5N1 negatives*
-	-	-	-	-	-	Experiment fail
+	+	+	+	+	+	Experiment fail

**(+) Positive:** Amplification signal

**(-) Negative:** No amplification signal

\* Note that if you observe amplification in the HEX and/or Cy5 channels you should use the *Neuroaminidase* Confirmation strip, to confirm the H5N1 and / or H7N9 subtypes

## The result interpretation for *Neuroaminidase* confirmation 8-well strips:

H5N1	H7N9	Internal Control	Negative Control	Positive Control	Interpretation
+	+	+/-	-	+	H5N1 and H7N9 Positives
-	-	+	-	+	H5N1 and H7N9 Negatives
+	-	+/-	-	+	H5N1 Positive and H7N9 Negative
-	+	+/-	-	+	H7N9 Positive and H5N1 Negative
-	-	-	-	-	Experiment fail
+	+	+	+	+	Experiment fail

(+) **Positive:** Amplification signal

(-) **Negative:** No amplification signal

If the negative samples do not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

### Quality Control

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in strips 1 and 3. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

### Performance evaluation

#### Clinical sensitivity and specificity

To evaluate the clinical sensitivity and the specificity, two clinical assays were carried out. In the first one, 256 throat swabs samples were analyzed with Vitassay qPCR Panel I / Respiratory viruses and CLART® PneumoVir DNA array (Genomica).

Vitassay qPCR Panel I / Respiratory viruses was able to detect Influenza A virus in 105 while CLART® PneumoVir DNA array was able to detect influenza A 107 samples. For this two discordant samples, the amount of RNA was below the detection limit of the method used.

Influenza B virus was detected in 51 samples with Vitassay qPCR Panel I / Respiratory viruses and in 52 samples by CLART® PneumoVir DNA array. For this discordant sample, the amount of RNA template was below the detection limit of the method used.

For RSV, a total of 41 clinical samples were positives for both assays.

For the second analysis, 127 throat swabs and samples in transport medium from symptomatic patients were tested by Real Time PCR using Vitassay qPCR Panel I / Respiratory viruses, CLART® PneumoVir DNA array (Genomica) and Cobas® Influenza A/B & RSV (Roche).

Influenza A (H1N1)pdm09 was detected in 82 samples by Vitassay qPCR Panel I / Respiratory viruses and in 85 samples by CLART® PneumoVir DNA + Cobas® Influenza A/B & RSV. For these 3 discordant samples, the low amount of template RNA is below the detection limit of the method used.

Influenza A (H3N2) was detected in 33 samples by Vitassay qPCR Panel I / Respiratory viruses and in 34 samples by CLART® PneumoVir DNA + Cobas® Influenza A/B & RSV. For this discordant sample, the low amount of template RNA is below the detection limit of the method used.

Vitassay qPCR Panel I / Respiratory viruses (Neuroaminidase confirmation) was evaluated with 59 samples dissolved in transport medium from EQAs programs. The results were compared with the final reports. Four positive samples of H5N1 and two positive samples of H7N9 were detected.

The specificity of Vitassay qPCR Panel I / Respiratory viruses shows a high sensitivity and specificity to detect Influenza A, Influenza B, and Respiratory Syncytial (RSV) viruses and subtyping of Influenza A (H1N1)pdm09, H3N2, H5N1 and H7N9.

### **Analytical sensitivity**

The analytical sensitivity was determined by analysis of 10-fold dilution series of Influenza A, Influenza B, RSV, (H1N1)pdm09, H3N2, H5N1 and H7N9 templates ranging from  $10^7$  to  $10^1$  copies/rxn. This assay has a detection limit of  $\geq 10$  viral RNA copies per reaction.

### **Analytical specificity**

The analytical specificity for Panel I / Respiratory viruses was tested within the panel of following microorganisms, where no cross-reactivity was seen between any of the species:

Cross-reactivity testing		
<i>Bordetella pertussis</i>	Human coronavirus 229E, OC43 and NL63	Influenza A/DE-SH/Reiherente/AR8444/ 2013 (H5N8) virus
<i>Bordetella bronchiseptica</i>	Human rhinovirus	Influenza A/Anhui/1/2013 (H7N9) virus
<i>Bordetella holmesii</i>	Bocavirus	Influenza A/Mallard/Netherlands/2/2009 (H7N7) virus
<i>Bordetella parapertussis</i>	Human Adenovirus types 1-5, 8, 31, 40 and 41	Influenza A/turkey/Virginia/2002 x PR8-IBCDC-5 (H7N2)
<i>Legionella bozemanii</i>	MERS Coronavirus	Influenza A/chicken/Hong Kong/G9/1997 x PR8-IBCDC-2 (H9N2)
<i>Legionella micdadei</i>	Influenza A/New Caledonia/20/99(H1N1) virus	Influenza A/mallard/Netherlands/12/2000 (H7N7) - IBCDC-1
<i>Legionella dumoffii</i>	Influenza A/California/7/2009(H1N1)pdm09-like virus	Influenza A/pheasant/New Jersey/1355/1998 (H5N2)-PR8-IBCDC-4
<i>Legionella longbeachae</i>	Influenza A/Netherlands/1250/2016 (H1N1)pdm09 virus (clade 6B.1)	Influenza A/Singapore/INFIMH-16-0019/2016 (H3N2)
<i>Legionella pneumophila</i>	Influenza A/Michigan/45/2015 (H1N1)pdm09 virus	Influenza A/South Australia/55/2014 (H3N2)
<i>Mycoplasma pneumoniae</i>	Influenza A/Singapore/GP1908/2015, IVR-180 (H1N1)pdm09 virus	Influenza A/Uruguay/716/2007 (H3N2) (NYMC-175C)
<i>Streptococcus pneumoniae</i>	Influenza A/Perth/16/2009(H3N2)-like virus	Influenza B/Brisbane/60/2008-like virus
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	Influenza A/Thüringen/5/17 (H3N2) virus	Influenza B/Florida/04/06 virus
<i>Haemophilus influenzae</i> MinnA	Influenza A/Switzerland/9715293/2013 (H3N2) virus	Influenza B/Phuket/3073/2013 virus
<i>Moraxella catarrhalis</i>	Influenza A/Hong Kong/4801/2014, NYMC X-263B (H3N2) virus	Influenza B/Netherlands/2518/2016 (clade 1A) virus
<i>Chlamydia caviae</i>	Influenza A/Netherlands/398/2014 (H3N2) virus (clade 3C.3a)	Influenza B/Netherlands/365/2016 (clade 3)
<i>Chlamydia psittaci</i> genotype A and C	Influenza A/Netherlands/2393/2015 (H3N2) virus (clade 3C.2a)	Influenza B/Colorado/6/2017
<i>Chlamydophila pneumoniae</i>	Influenza A/Hong Kong/213/2003 (H5N1) virus	Influenza B/Maryland/15/2016
Human parainfluenza 1, 2, 3 and 4 viruses	Influenza A/Turkey/Germany R2485+86/2014 (H5N8) virus	Respiratory syncytial virus (RSV)
Human metapneumovirus A and B		

## Analytical reactivity

Analytical reactivity of Vitassay qPCR Panel I / Respiratory viruses for Flu A was evaluated against strains: A/New Caledonia/20/99(H1N1) virus, A/California/7/2009(H1N1)pdm09-like virus, A/Netherlands/1250/2016 (H1N1)pdm09 virus (clade 6B.1), A/Michigan/45/2015 (H1N1)pdm09 virus, A/Singapore/GP1908/2015, IVR-180 (H1N1)pdm09 virus, A/Perth/16/2009(H3N2)-like virus, A/Thüringen/5/17 (H3N2) virus, A/Switzerland/9715293/2013 (H3N2) virus, A/Hong Kong/4801/2014, NYMC X-263B (H3N2) virus, A/Netherlands/398/2014 (H3N2) virus (clade 3C.3a), A/Netherlands/2393/2015 (H3N2) virus (clade 3C.2a), A/Hong Kong/213/2003 (H5N1) virus, A/Turkey/Germany R2485+86/2014 (H5N8) virus, A/DE-SH/Reiherente/AR8444/2013 (H5N8) virus, A/Mallard/Netherlands/2/2009 (H7N7) virus, A/Anhui/1/2013 (H7N9) virus, A/turkey/Virginia/2002 x PR8-IBCDC-5 (H7N2), A/chicken/Hong Kong/G9/1997 x PR8-IBCDC-2 (H9N2), A/mallard/Netherlands/12/2000 (H7N7) - IBCDC-1, A/pheasant/New Jersey/1355/1998 (H5N2)-PR8-IBCDC-4, A/Singapore/INFIMH-16-0019/2016 (H3N2), A/South Australia/55/2014 and A/Uruguay/716/2007 (H3N2)(NYMC-175C) virus showing positive results.

Analytical reactivity of Vitassay qPCR Panel I / Respiratory viruses for Flu B was evaluated against strains: B/Brisbane/60/2008-like virus, B/Florida/04/06 virus, B/Phuket/3073/2013 virus, B/Netherlands/2518/2016 (clade 1A) virus, B/Netherlands/365/2016 (clade 3), B/Colorado/6/2017 and B/Maryland/15/2016, showing positive results.

Analytical reactivity of Vitassay qPCR Panel I / Respiratory viruses for RSV was evaluated against strains: Human Respiratory Syncytial Virus (RSV) A and B, showing positive results.

Analytical reactivity of Vitassay qPCR Panel I / Respiratory viruses for Influenza A (H1N1)pdm09 was evaluated against the strains: A/Netherlands/1250/2016 (H1N1)pdm09 virus (clade 6B.1), A/California/7/2009(H1N1)pdm09-like virus, A/Michigan/45/2015 (H1N1)pdm09 virus and A/Singapore/GP1908/2015 IVR-180 (H1N1)pdm09 virus, showing positive result.

Analytical reactivity of Vitassay qPCR Panel I / Respiratory viruses for Influenza A (H3N2) was evaluated against strains: A/Perth/16/2009(H3N2)-like virus, A/Thüringen/5/17 (H3N2) virus, A/Switzerland/9715293/2013 (H3N2) virus, A/Hong Kong/4801/2014 NYMC X-263B (H3N2) virus, A/Netherlands/398/2014 (H3N2) virus (clade 3C.3a), A/Netherlands/2393/2015 (H3N2) virus (clade 3C.2a), A/Singapore/INFIMH-16-0019/2016 (H3N2), Influenza A/South Australia/55/2014 (H3N2) and A/Uruguay/716/2007 (H3N2)(NYMC-175C)) viruses, showing positive result.

Analytical reactivity of Vitassay qPCR Panel I / Respiratory viruses for Influenza A (H5N1) was evaluated against strains A/Hong Kong/213/2003 (H5N1) virus, showing positive result.

Analytical reactivity of Vitassay qPCR Panel I / Respiratory viruses for Influenza A (H7N9) was evaluated against strains: A/Anhui/1/2013(H7N9) virus, showing positive result.

### **Compatibles real-time PCR equipment**

Vitassay qPCR Panel I / Respiratory viruses has been validated on the following equipments:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems)<sup>II</sup>
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTPRime Real Time Detection Thermal Cyclers (DNA-Technology)
- DTlite Real-Time PCR System (DNA-Technology)
- Rotor-Gene® (Qiagen)<sup>I</sup>
- SmartCycler® (Cepheid)<sup>I</sup>

I: For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler® tubes.

II: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

### **Limitations**

- All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- This assay was tried with throat swab samples dissolved in transport medium. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper RNA from clinical specimens must be extracted.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by the different virus, either samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.

## Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
<b>Agilent Technologies</b>	<b>Abbott</b>
AriaMx Real-Time PCR System	Abbott m2000 RealTime System
<b>Applied Biosystems</b>	<b>Applied Biosystems</b>
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
<b>BIONEER</b>	ViiA™ 7 Real-Time PCR
Exicycler™ 96	<b>Analytik Jena Biometra</b>
<b>Bio-Rad</b>	TOptical
CFX96™ Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	<b>BIONEER</b>
<b>Cepheid</b>	Exicycler™ 96
SmartCycler®	<b>Bio-Rad</b>
<b>Qiagen</b>	CFX96™ Deep Well Real-Time PCR Detection
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
<b>Roche</b>	iCycler iQ™ 5 Real-Time PCR
LightCycler®480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler®96 Real-Time PCR System	MyiQ™ 2Real-Time PCR Detection System
Cobas z480 Analyzer	<b>Cepheid</b>
	SmartCycler®
	<b>DNA-Technology</b>
	DTlite Real-Time PCR System
	DTprime Real-time Detection Thermal Cycler
	<b>Eppendorf</b>
	Mastercycler™ep <i>realplex</i>
	<b>Qiagen</b>
	Rotor-Gene® Q
	<b>Stratagene / Agilent Technologies</b>
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

\* See Attached III to configure exposure settings.

## Attached II: Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

REAL-TIME PCR THERMOCYCLER	Vitassay CHANNEL	DETECTION CHANNEL	OBSERVATIONS
<b>Bio-Rad CFX96™</b>	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>ABI 7500 Applied Biosystems</b>	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>Roche Lightcycler®480II</b>	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
<b>Smartcycler® Cepheid</b>	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
<b>Abbott m2000rt</b>	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>Mx3000P™ Mx 3005P™ Stratagene</b>	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>AriaMx Agilent</b>	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>Rotor-Gene® Q Qiagen</b>	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
<b>Exicycler™ 96 BIONEER</b>	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

**Attached III: Optical measurement exposure setting**

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

<b>Thermocycler</b>	<b>Vitassay channel</b>	<b>Exposure values</b>
<b>DTlite Real-Time PCR System (DNA-Technology)</b>	FAM	500
	HEX	500
	ROX	500
	Cy5	500
<b>DTprime Real-time Detection Thermal Cycler (DNA-Technology)</b>	FAM	500
	HEX	1000
	ROX	1000
	Cy5	1000

## Bibliography/Bibliografía

1. Methods for molecular surveillance of influenza. Wang R, Taubenberger JK. Expert review of anti-infective therapy. 2010;8(5):517-527.
2. WHO information for molecular diagnosis of influenza virus—update. World Health Organization. Available: [http://www.who.int/influenza/gisrs\\_laboratory/molecular\\_diagnosis/en/](http://www.who.int/influenza/gisrs_laboratory/molecular_diagnosis/en/). Accessed 2015 Dec 30.
3. Kuo RL, Yang SL, Liu YC, Chen LT, Mok CK, Kuo SM, Shih SR, Tsao KC. Influenza A/B virus detection and influenza A virus subtyping with emphasis on the novel H7N9 virus by using multiplex real-time RT-PCR. J Virol Methods. 2014 Nov; 208:41-6.
4. Bawage SS, Tiwari PM, Pillai S, Dennis V, Singh SR. Recent advances in diagnosis, prevention, and treatment of human respiratory syncytial virus. Adv Virol. 2013;2013:595768.
5. Falsey AR. Respiratory syncytial virus infection in adults. Semin Respir Crit Care Med. 2007 Apr;28(2):171-81.
6. Van Woensel JB, Kimpen JL, Brand PL. Respiratory tract infections caused by respiratory syncytial virus in children. Diagnosis and treatment. Minerva Pediatr. 2001 Apr;53(2):99-106.
7. I.H. Brown, et al. Multiple genetic reassortment of avian and human influenza A viruses in European pigs, resulting in the emergence of an H1N2 virus of novel genotype. Journal of General Virology, 1998; 79: 2947-2955.

## Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.


Mx3000P™ and Mx3029™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid.

## Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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