Vitassay qPCR

Genital ulcer

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- LightCycler Z480, 480II (Roche)
- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- StepOne™ Real-Time PCR System (Applied Biosystems) ^{II}
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time PCR System (DNA-Technology)
- DTPPrime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)^I
- SmartCycler® (Cepheid)^I

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Termocicladores con bloque de bajo perfil	Termocicladores con bloque de alto perfil
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 RealTime System
Applied Biosystems	Applied Biosystems
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
BIONEER	ViiA™ 7 Real-Time PCR
Exicycler™ 96	Analytik Jena Biometra
Bio-Rad	TOptical
CFX96™ Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	BIONEER
Cepheid	Exicycler™ 96
SmartCycler®	Bio-Rad
Qiagen	CFX96™ Deep Well Real-Time PCR Detection
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
Roche	iCycler iQ™ 5 Real-Time PCR
LightCycler® 480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler® 96 Real-Time PCR System	MyiQ™ 2Real-Time PCR Detection System
Cobas z480 Analyzer	Cepheid
	SmartCycler®
	DNA-Technology
	DTlite Real-Time PCR System*
	DTprime Real-time Detection Thermal Cycler*
	Eppendorf
	Mastercycler™ ep <i>realplex</i>
	Qiagen
	Rotor-Gene® Q
	Stratagene / Agilent Technologies
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

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Intended use

Vitassay qPCR Genital ulcer allows the detection and differentiation of Herpes simplex virus 1, Herpes simplex virus 2 and *Treponema pallidum* by real-time PCR in clinical samples. The product is intended for use in the diagnosis of Herpes simplex virus 1, Herpes simplex virus 2 and *Treponema pallidum* infections alongside clinical data of the patient and other laboratory tests outcomes.

References

Vitassay qPCR Genital ulcer 4x8-well strip, low profile	7041036
Vitassay qPCR Genital ulcer 4x8-well strip, high profile	7042036

Materials/reagents provided

Reference	Reagent/Material	Colour	Amount
7041S036/ 7042S036	Genital ulcer strips low/high profile	-	4x8-well strip
7C036	Genital ulcer Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4x8 cap strip

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

Additional equipment and material required

- DNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Vortexer
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

Summary

Genital, anal and oropharyngeal ulcers are caused by herpes simplex viruses type 1 and 2 and *Treponema pallidum*. Their detection is important since these ulcers are a risk factor for the transmission and acquisition of HIV. There are different types of tests for the detection of these pathogens, such as PCR and serology, although real-time PCR has proven to be the most specific and sensitive tool for detecting them.

Treponema pallidum is a bacterium which causes syphilis (venereal and endemic), yaws and pinta. Although only venereal syphilis is sexually transmitted, because it requires close personal contact, due to the poor tolerance of these to desiccation and elevated temperature.

Herpes simplex viruses type 1 and 2 cause genital herpes and persistent infections. They cannot be distinguished clinically; the same treatment is used for infections caused by both viruses. The dose and frequency will be determined by the location of the lesions and the chronicity of the infection. Infection with herpes simplex virus type 1 and 2 occurs through the inoculation of viral particles in susceptible mucosae. Thanks to real-time PCR, they can be detected and differentiated.

Principle of the test

Vitassay qPCR Genital ulcer test is based on the real-time amplification of specific conserved fragments of the *US4* gene for Herpes simplex virus 1, *US6* gene for Herpes simplex virus 2 and 16S rRNA gene for *Treponema pallidum*. After DNA extraction, the presence of the Herpes simplex virus 1, Herpes simplex virus 2 and *Treponema pallidum* is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

Vitassay qPCR Genital ulcer test is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The amplification of the Herpes simplex virus 1 DNA target sequence is detected through the FAM channel, Herpes simplex virus 2 virus DNA target in ROX channel and *Treponema pallidum* DNA in Cy5 channel whereas the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Attached II).

Precautions

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.

- Do not use the kit if the desiccant from the different reagents is not present or broken or if the foil has been broken or damaged.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

Procedures

Specimen collection, processing and DNA extraction

For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

ZP02012 MagPurix Viral/Pathogen Nucleic Acids Extraction Kit B, using the MagPurix 12A instrument (Zinexts Life Science Corp.).

Invisorb® Spin Universal Kit (Stratec).

Positive control preparation

Reconstitute the lyophilized Genital ulcer Positive Control (red tube) with the 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of RNA sample, negative and positive controls and add them into each well.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the FAM (Herpes simplex virus 1), ROX (Herpes simplex virus 2), Cy5 (*Treponema pallidum*) and HEX, JOE or VIC channels (Internal Control (IC)). If you use the Applied Biosystems 7500 Fast or the Stratagene Mx3029P™ check that passive reference option ROX is none. (Attached II)

Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive controls used in each run, must show an amplification curve for FAM (Herpes simplex virus 1), ROX (Herpes simplex virus 2) and Cy5 (*Treponema pallidum*), which validates the reaction.

Negative control

The negative controls included in each run, must show the absence of signal in FAM, ROX and Cy5 which validates the reaction.

The experiment seems to be failed if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

The result interpretation is summarized in the following table:

Herpes simplex virus 1	Herpes simplex virus 2	<i>Treponema pallidum</i>	Internal Control	Negative Control	Positive Control	Interpretation
+	+	+	+/-	-	+	Herpes simplex virus 1, Herpes simplex virus 2 and <i>Treponema pallidum</i> Positives
-	-	-	+	-	+	Herpes simplex virus 1, Herpes simplex virus 2 and <i>Treponema pallidum</i> Negatives
+	-	-	+/-	-	+	Herpes simplex virus 1 Positive, Herpes simplex virus 2 and <i>Treponema pallidum</i> Negatives
+	+	-	+/-	-	+	Herpes simplex virus 1 and Herpes simplex virus 2 Positives, <i>Treponema pallidum</i> Negative
+	-	+	+/-	-	+	Herpes simplex virus 1 and <i>Treponema pallidum</i> Positives, Herpes simplex virus 2 Negative
-	+	-	+/-	-	+	Herpes simplex virus 2 Positive, Herpes simplex virus 1 and <i>Treponema pallidum</i> Negatives
-	+	+	+/-	-	+	Herpes simplex virus 2 and <i>Treponema pallidum</i> Positives, Herpes simplex virus 1 Negative
-	-	+	+/-	-	+	<i>Treponema pallidum</i> Positive, Herpes simplex virus 1 and Herpes simplex virus 2 Negatives
+	+	+	+	+	+	Experiment fail
-	-	-	-	-	-	Experiment fail

(+) **Positive:** Amplification signal

(-) **Negative:** No amplification signal

If the negative samples not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

Quality Control

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

A total of 33 samples in transport medium (belonging to EQA programs) of genital, anal or oropharyngeal pathogens from symptomatic patients were analyzed by Vitassay qPCR Genital ulcer and FTD Genital ulcer (Fast Track Diagnostics). Vitassay qPCR Genital ulcer detected Herpes Simplex Virus 1 in 8 samples, Herpes Simplex Virus 2 in 9 samples and *Treponema pallidum* in 8 samples.

In addition, Vitassay qPCR Genital ulcer was evaluated with the following programmes: QCMD 2017 Central nervous system I EQA pilot study, QCMD 2017 Syphilis EQA programme (SYPH17) and QCMD 2017 enterovirus RNA EQA programme (EVRNA17). The programmes consisted of panels of samples (10 samples, 8 samples and 10 samples, respectively) dissolved in transport medium. These samples were positive or negative to the analyzed pathogens. Vitassay qPCR Genital ulcer detected correctly all the samples from the different panels.

The results show a high sensitivity and specificity to detect Herpes Simplex Virus 1, Herpes Simplex Virus 2 and *Treponema pallidum* using the molecular diagnostic test Vitassay qPCR Genital ulcer.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of Herpes simplex virus 1, Herpes simplex virus 2 and *Treponema pallidum* templates ranging from 10^7 to 10^1 copies/rxn. This assay has a detection limit of ≥ 10 viral DNA copies per reaction.

Analytical specificity

The analytical specificity for Herpes simplex virus 1, Herpes simplex virus 2 and *Treponema pallidum* was tested within the panel of following pathogens, where no cross-reactivity was seen between any of the species:

Cross-reactivity assay		
<i>Acinetobacter baumannii</i>	<i>Enterococcus faecium</i>	<i>Proteus mirabilis</i>
<i>Aspergillus fumigatus</i>	<i>Escherichia coli</i> 0.1285; O18:H7:K1	<i>Pseudomonas aeruginosa</i>
<i>Bacteroides fragilis</i>	<i>Gardnerella vaginalis</i>	<i>Serratia marcescens</i>
<i>Candida albicans</i>	<i>Haemophilus influenzae</i>	<i>Staphylococcus aureus</i>
<i>Candida parapsilosis</i>	<i>Haemophilus ducrey</i>	<i>Stenotrophomonas maltophilia</i>
<i>Candida tropicalis</i>	<i>Klebsiella oxytoca</i>	<i>Streptococcus pneumoniae</i>
<i>Candida glabrata</i>	<i>Klebsiella pneumoniae</i>	<i>Streptococcus agalactiae</i>
<i>Candida krusei</i>	<i>Listeria ivanovii</i>	<i>Trichomonas vaginalis</i>
<i>Candida dubliniensis</i>	<i>Listeria monocytogenes</i>	<i>Ureaplasma parvum</i>
<i>Chlamydia trachomatis</i> (LGV)	<i>Listeria Innocua</i>	<i>Ureaplasma urealyticum</i>
<i>Chlamydia trachomatis</i> (SW)	<i>Mycoplasma genitalium</i>	Cytomegalovirus
<i>Enterobacter cloacae</i>	<i>Mycoplasma hominis</i>	Hepatitis A
<i>Enterobacter aerogenes</i>	<i>Neisseria gonorrhoeae</i>	Human Papillomavirus genotypes HP16 and HPV18
<i>Enterococcus faecalis</i>	<i>Neisseria meningitidis</i>	

Analytical reactivity

Analytical reactivity of Vitassay qPCR Genital ulcer for Herpes simplex virus 1 was evaluated against HSV 1 MacIntyre strain, showing positive result.

Analytical reactivity of Vitassay qPCR Genital ulcer for Herpes simplex virus 2 was evaluated against HSV 2 MS strain, showing positive result.

Analytical reactivity of Vitassay qPCR Genital ulcer for *Treponema pallidum* was evaluated against *Treponema pallidum* showing positive result.

Compatibles real-time PCR equipment

Vitassay qPCR Genital ulcer has been validated on the following equipments:

- LightCycler Z480, 480II (Roche)
- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- StepOne™ Real-Time PCR System (Applied Biosystems) ^{II}
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time PCR System (DNA-Technology)
- DTPrime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)^I
- SmartCycler® (Cepheid)^I

I: For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler tubes.

II: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

Limitations

- This test provides a presumptive diagnosis of Herpes simplex virus 1, Herpes simplex virus 2 and *Treponema pallidum* infection. All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- This assay was tried with samples in transport medium, endocervical swabs, urine, plasma, cell lysates and synthetic CSF. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper DNA from clinical specimens must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by the different pathogens, either samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers
Agilent Technologies
AriaMx Real-Time PCR System
Applied Biosystems
7500 Fast Real-Time PCR System
7500 Fast Dx Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast
QuantStudio™ 6 Flex 96-well Fast
QuantStudio™ 7 Flex 96-well Fast
QuantStudio™ 5 Real-Time PCR System
QuantStudio™ 3 Real-Time PCR System
StepOne Plus™ Real-Time PCR System
StepOne™ Real-Time PCR System
ViiA™ 7 Fast Real-Time PCR System
BIONEER
Exicycler™ 96
Bio-Rad
CFX96™ Real-Time PCR Detection System
Mini Opticon™ Real-Time PCR Detection System
Cepheid
SmartCycler®
Qiagen
Rotor-Gene® Q
Roche
LightCycler®480 Real-Time PCR System
LightCycler®96 Real-Time PCR System
Cobas z480 Analyzer

High profile Block Thermocyclers
Abbott
Abbott m2000 RealTime System
Applied Biosystems
7300 Real-Time PCR System
7500 Real-Time PCR System
7900 HT Real-Time PCR System
ABI PRISM 7000
ABI PRISM 7700
QuantStudio™ 12K Flex 96-well
QuantStudio™ 6 Flex 96-well
QuantStudio™ 7 Flex 96-well
QuantStudio™ 5 Real-Time PCR
QuantStudio™ 3 Real-Time PCR System
ViiA™ 7 Real-Time PCR
Analytik Jena Biometra
TOptical
qTOWER 2.0
BIONEER
Exicycler™ 96
Bio-Rad
CFX96™ Deep Well Real-Time PCR Detection
iCycler iQ™ Real-Time PCR
iCycler iQ™ 5 Real-Time PCR
MyiQ™ Real-Time PCR Detection System
MyiQ™ 2Real-Time PCR Detection System
Cepheid
SmartCycler®
DNA-Technology
DTlite Real-Time PCR System
DTprime Real-time Detection Thermal Cycler
Eppendorf
Mastercycler™ ep <i>realplex</i>
Qiagen
Rotor-Gene® Q
Stratagene / Agilent Technologies
Mx3000P™ Real Time PCR System
Mx3005P™ Real Time PCR System

* See Attached III to configure exposure settings.

Attached II: Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

REAL-TIME PCR THERMOCYCLER	Vitassay CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Attached III: Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
DTlite Real-Time PCR System (DNA-Technology)	FAM	500
	HEX	500
	ROX	500
	Cy5	500
DTprime Real-time Detection Thermal Cycler (DNA-Technology)	FAM	500
	HEX	1000
	ROX	1000
	Cy5	1000

Bibliography/Bibliografía

1. J. D. Radolf et al. *Treponema pallidum*, the syphilis spirochete: making a living as a stealth pathogen. *Nature Reviews Microbiology* 2016; 14(12): 744–759.
2. M. Lieveld et al. A high resolution melting (HRM) technology-based assay for cost-efficient clinical detection and genotyping of herpes simplex virus (HSV)-1 and HSV-2. *Journal of Virological Methods* 2017; Volume 248, 181-186.
3. M. Costa-Silva et al. Cross-sectional study of *Treponema pallidum* PCR in diagnosis of primary and secondary syphilis. *International Journal of Dermatology* 2017; doi: 10.1111/ijd.13823.
4. M.A. Minaya et al. Molecular evolution of herpes simplex virus 2 complete genomes: Comparison between primary and recurrent infections. *American Society for Microbiology* 2017; doi:10.1128/JVI.00942-17.
5. C.O. Onyango et al. Evaluation of a TaqMan Array Card for Detection of Central Nervous System Infections. *Journal of Clinical Microbiology* 2017 Jul;55(7): 2035-2044.
6. Blanco DR, Miller JN, Lovett MA. Surface Antigens of the Syphilis Spirochete and Their Potential as Virulence Determinants. *Emerging Infectious Diseases*. 1997;3(1):11-20.
7. Genital Herpes Simplex Virus (HSV) Module. Genital and Perirectal Herpes Simplex Virus Infection Faculty Notes. CDC. STD Curriculum for Clinical Educators, 2014.

Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.


Mx3000P™ and Mx3029™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid.

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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