VA Vitassay

Real-Time PCR Kits

Vitassay qPCR HPV

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- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- CFX96 TM Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTLite Real-Time PCR System (DNA-Technology)
- DTPRime Real Time Detection Thermal Cyclers (DNA-Technology)
- Rotor-Gene® Q (Qiagen)^I
- SmartCycler® (Cepheid)^I

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Termocicladores con bloque de bajo perfil	Termocicladores con bloque de alto perfil
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 Real-Time System
Applied Biosystems	Applied Biosystems
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR System
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
BIONEER	ViiA™ 7 Real-Time PCR System
Exicycler™ 96	Analytik Jena Biometra
Bio-Rad	TOptical
CFX96™ Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	BIONEER
Cepheid	Exicycler™ 96
SmartCycler®	Bio-Rad
Qiagen	CFX96™ Deep Well Real-Time PCR
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
Roche	iCycler iQ™ 5 Real-Time PCR
LightCycler @480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler @96 Real-Time PCR System	MyiQ™ 2 Real-Time PCR Detection System
Cobas z480 Analyzer	Cepheid
	SmartCycler®
	DNA-Technology
	DTlite Real-Time PCR System *
	DTprime Real-time Detection Thermal Cyclor *
	Eppendorf
	Mastercycler™ ep <i>realplex</i>
	Qiagen
	Rotor-Gene® Q
	Stratagene/ Agilent Technologies
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

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Intended use

Vitassay qPCR HPV allows the detection and differentiation of human papilloma virus 16 and 18 by real-time PCR in clinical samples. The product is intended for use in the diagnosis of human papillomavirus 16 and 18 infections alongside clinical data of the patient and other laboratory tests outcomes.

References

Vitassay qPCR HPV 4x8-well strip, low profile	7041037
Vitassay qPCR HPV 4x8-well strip, high profile	7042037

Reagents provided

In references 7041037 and 7042037:

Reference	Reagent/Material	Color	Amount
7041S037/ 7042S037	HPV strips low/high profile	-	4 x 8-well strip
7C037	HPV Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4 x 8-cap strip

Material and equipment required, not provided

- DNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 ml tubes
- Vortex
- Micropipettes (1-20 µl, 20-200 µl)
- Filter tips
- Powder-free disposal gloves

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. To avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

Summary

Human papilloma virus (HPV) is responsible for the most frequent sexually transmitted infection (STI) of the reproductive tract. This virus is so common that almost all sexually active women and men contract it at some point in their lives.

More than 200 types of HPV are known and have been classified predominantly in three genera: Alpha-papillomavirus, Beta-papillomavirus and Gamma-papillomavirus. Among the 65 HPV types belonging to Alpha-HPV, there are 15 genotypes most frequently associated with cervical cancer cases, called high-risk HPV (HPV-16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, 68 and 82), type 16 is the responsible for more than 50% of cases of squamous cell carcinoma.

The transmission is produced by sexual contacts and the organs most susceptible to the infection are the cervix and the pectinate line of the anal canal. HPV infections are frequently in sheet, in which cases the viral DNA can be recovered from the cervix, vulva, vagina, anal canal, penis and scrotum. Most HPV infections do not cause symptoms or diseases and resolve spontaneously. But when HPV does not go away, it can cause health problems such as genital warts, respiratory papillomatosis and cervical, anogenital and oropharyngeal cancers.

Tumors of female genital tract account for a fifth of women's tumors according to global estimates. The most frequent tumor is that of cervical cancer (9.9%), followed by cancer of the ovary (1%), endometrium (4%) and cancers of the vagina and vulva (1%). These data are of particular importance in developing countries, due to the lack of adequate screening programs for cervical cancer. In particular, HPV 16 and 18 are the ones that most often can lead to precancerous lesions. In addition, both men and women can develop mouth / throat and anus / rectum cancer caused by HPV infections. Men can also develop penile cancer. Currently, there are vaccines that could prevent infection with the types of HPV that most often cause cancer.

A variety of diagnostic methods of different sensitivity and specificity have been developed to detect HPV in skin, oral and anogenital samples (mostly, cervical scrapings and biopsy material). Since HPV cannot be cultivated efficiently and the clinical performance of serological assays is poor, the diagnosis of HPV infection is based on molecular tools. Currently, real-time PCR offers high sensitivity and can be performed in different types of samples.

Test principle

Vitassay qPCR HPV is based on the real-time amplification of specific conserved fragments of the *L2* gene for human papilloma virus 16 and *L1* gene for human papilloma virus 18. After DNA extraction, the presence of HPV16 and HPV18 viruses is detected by an increase in observed fluorescence during the reaction, after hydrolysis of the fluorescent probe.

Vitassay qPCR HPV is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The detection channels of the target sequences are described below:

Target	Detection Channel
Human papilloma virus 16	FAM
Human papilloma virus 18	HEX, VIC, or JOE (depending on the equipment used)
Internal Control (IC)	Cy5

Precautions

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or is broken or damaged.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles, and mask.
- Do not eat, drink, or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment, and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment, and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

Procedures

Specimen collection, transport, and DNA extraction

For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

- Maxwell® 16 Viral Total Nucleic Acid Purification Kit (Promega)
- Invisorb® Spin Universal Kit (Invitek)
- QIAamp Viral RNA Mini Kit, utilizando QIAcube instrument (Qiagen)
- NucleoMag® Pathogen (Macherey Nagel)
- NucleoSpin RNA Virus (Macherey Nagel)
- MagDEA Dx SV kit con magLEAD® 6gC instrument (Precision System Science Co.)
- NX-48 Urine/Swab DNA Kit utilizando Nextractor® NX-48 system, (Genolution)

- ZP02004 MagPurix Tissue DNA Extraction Kit, utilizando MagPurix 12A instrument (Zinexts Life Science Corp)

Positive control preparation

Reconstitute the lyophilized HPV Positive Control (red tube) with 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template, and it is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

Reaction setup

- Separate the number of reactions required including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of DNA sample, negative control (yellow tube) or positive control reconstituted (red tube) and add them into the corresponding wells.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the channels FAM (Human papilloma virus 16), HEX, JOE or VIC (Human papilloma virus 18), and Cy5 (Internal Control (IC)). If you use the Applied Biosystems 7500 Fast or the Stratagene Mx3005P™ check that passive reference option ROX is none. (Attached II).

Analysis and interpretation of results

The results analysis is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive control used in each run must show an amplification curve in FAM (Human papilloma virus 16) and HEX, VIC, or JOE (Human papilloma virus 18) channels, which validates the reaction.

Negative control

The negative control included in each run must show signal' absence in FAM, and HEX, VIC or JOE channels, which validates the reaction.

The experiment seems to be failed if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

The result interpretation is summarized in the following table:

HPV 16 (FAM)	HPV 18 (HEX)	Internal Control	Negative Control	Positive Control	Interpretación
+	+	+/-	-	+	Human papillomavirus 16 and Human papilloma virus 18 Positive
-	-	+	-	+	Human papillomavirus 16 and Human papilloma virus 18 Negatives
+	-	+/-	-	+	Human papillomavirus 16 Positive and Human papilloma virus 18 Negative
-	+	+/-	-	+	Human papillomavirus 16 Negative and Human papilloma virus 18 Positive
+	+	+	+	+	Experiment fail
-	-	-	-	-	Experiment fail

(+) **Positive:** Amplification signal

(-) **Negative:** No amplification signal

If the negative samples do not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction must be repeated due to possible problems caused by PCR inhibitors.

Quality Control

To confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

A total of 94 clinical samples (32 liquid cytology and 62 endocervical / vaginal swabs) from symptomatic patients were analyzed by Vitassay qPCR HPV and CLART®HPV4" Genotyping of human Papillomavirus via genomic identification for *in vitro* diagnosis (Genomica) or "Anyplex™ II HPV28 Detection" using DPO™ technology and melting curve analysis method of TOCE™ technology" (Seegene)).

Both kits detected the presence of human papilloma virus 16 in 43 samples and human papillomavirus 18 in 5 samples. In addition, it was observed that 36% of the samples were coinfections of HPV16 / HPV18 (20/34).

In addition to this, Vitassay qPCR HPV test was evaluated with QCMD 2017 Human Papillomavirus DNA EQA programme. The panel consisted on 12 vials containing cell suspensions with various concentrations and subtypes of humanpapillomavirus or samples negative for HPV. All samples were detected correctly confirming the sensitivity and specificity of the evaluated molecular diagnostic kit.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of Human papillomavirus 16 and Human papillomavirus 18 standards ranging from 10^7 to 10^1 copies/reaction.

Vitassay qPCR HPV has a detection limit (LoD) of ≥ 10 DNA copies per reaction.

Analytical specificity

The analytical specificity of the Human papillomavirus 16 and 18 assay was tested within the panel of following pathogens, where no cross-reactivity was seen between any of the species.

Cross-reactivity assay		
<i>Candida parapsilosis</i>	<i>Listeria monocytogenes</i>	Human papillomavirus 18
<i>Candida tropicalis</i>	<i>Listeria innocua</i>	Human papillomavirus 6
<i>Candida glabrata</i>	<i>Mycoplasma genitalium</i>	Human papillomavirus 11
<i>Candida krusei</i>	<i>Mycoplasma hominis</i>	Human papillomavirus 31
<i>Candida dubliniensis</i>	<i>Neisseria gonorrhoeae</i>	Human papillomavirus 33
<i>Candida albicans</i>	<i>Neisseria meningitidis</i>	Human papillomavirus 39
<i>Aspergillus fumigatus</i>	<i>Proteus mirabilis</i>	Human papillomavirus 40
<i>Bacteroides fragilis</i>	<i>Pseudomonas aeruginosa</i>	Human papillomavirus 42
<i>Acinetobacter baumannii</i>	<i>Serratia marcescens</i>	Human papillomavirus 44
<i>Chlamydia trachomatis</i> (LGV)	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	Human papillomavirus 45
<i>Chlamydia trachomatis</i> (SW)	<i>Stenotrophomonas maltophilia</i>	Human papillomavirus 51
<i>Chlamydia trachomatis</i> genovar F	<i>Streptococcus agalactiae</i>	Human papillomavirus 52
<i>Enterobacter cloacae</i>	<i>Streptococcus pneumoniae</i>	Human papillomavirus 53
<i>Enterobacter aerogenes</i>	<i>Treponema pallidum</i>	Human papillomavirus 54
<i>Enterococcus faecalis</i>	<i>Trichomonas vaginalis</i>	Human papillomavirus 56
<i>Enterococcus faecium</i>	<i>Ureaplasma parvum</i>	Human papillomavirus 58
<i>E. coli</i> 0.1285;O18:H7:K1	<i>Ureaplasma urealyticum</i>	Human papillomavirus 59
<i>Gardnerella vaginalis</i>	Cytomegalovirus AD-169	Human papillomavirus 61
<i>Haemophilus influenzae</i> MinnA	Hepatitis A	Human papillomavirus 66
<i>Haemophilus ducreyi</i> class 1	Herpes simplex virus 1	Human papillomavirus 68
<i>Klebsiella oxytoca</i>	Herpes simplex virus 2	Human papillomavirus 70
<i>Klebsiella pneumoniae</i>	Human papillomavirus 16	Human papillomavirus 73
<i>Listeria ivanovii</i>		

Analytical reactivity

The reactivity of Vitassay qPCR HPV for Human papillomavirus 16 was evaluated against Human papillomavirus 16, showing positive result.

The reactivity of Vitassay qPCR HPV for Human papillomavirus 18 was evaluated against Human papillomavirus 18, showing positive result.

Compatibles real-time PCR equipment

Vitassay qPCR HPV has been validated on the following equipment:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTLite Real-Time PCR System (DNA-Technology)
- DTPPrime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)^I
- SmartCycler® (Cepheid)^I

I: For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® and/or SmartCycler tubes.

I: When using the Applied Biosystems 7500 Fast with strips, it is recommended to place a plate holder for the tubes (Ref. PN 4388506).

Limitations

- This test provides a presumptive diagnosis of HPV 16 and HPV 18 infection. All results must be interpreted by a specialist together with other clinical information and laboratory findings available to the physician.
- This assay was tried with liquid cytology and endocervical/vaginal swabs samples. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper DNA from clinical specimens must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by the different pathogens, either samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

Attached I. Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 Real-Time System
Applied Biosystems	Applied Biosystems
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR System
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
BIONEER	ViiA™ 7 Real-Time PCR System
Exicycler™ 96	Analytik Jena Biometra
Bio-Rad	TOptical
CFX96™ Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	BIONEER
Cepheid	Exicycler™ 96
SmartCycler®	Bio-Rad
Qiagen	CFX96™ Deep Well Real-Time PCR
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
Roche	iCycler iQ™ 5 Real-Time PCR
LightCycler @480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler @96 Real-Time PCR System	MyiQ™ 2 Real-Time PCR Detection System
Cobas z480 Analyzer	Cepheid
	SmartCycler®
	DNA-Technology
	DTlite Real-Time PCR System *
	DTprime Real-time Detection Thermal Cyclor *
	Eppendorf
	Mastercycler™ ep <i>realplex</i>
	Qiagen
	Rotor-Gene® Q
	Stratagene/ Agilent Technologies
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

* See Attached III to configure exposure settings.

Attached II. Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common real-time PCR thermocyclers are specified in the following table:

Thermocycler	Vitassay Channel	Detection Channel	Observations
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Applied Biosystems ABI 7500	FAM	FAM	Passive reference option ROX is none.
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Color Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Cepheid Smartcycler®	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Agilent Technologies Mx3000P™ Mx3005P™	FAM	FAM	Passive reference option ROX is none.
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Agilent Technologies AriaMx	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Qiagen Rotor-Gene® Q	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range must be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
BIONEER Exicycler™ 96	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Attached III. Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
DTlite Real-Time PCR System (DNA-Technology)	FAM	500
	HEX	500
	ROX	500
	Cy5	500
DTprime Real-time Detection Thermal Cycler (DNA-Technology)	FAM	500
	HEX	1000
	ROX	1000
	Cy5	1000

Bibliography/Bibliografia

1. Y. Zhao *et al.* A novel multiplex real-time PCR assay for the detection and quantification of HPV16/18 and HSV1/2 in cervical cancer screening. *Molecular and Cellular Probes* 2012; 66-72.
2. B. Biesaga *et al.* Comparison of the sensitivity and specificity of real-time PCR and in situ hybridization in HPV16 and 18 detection in archival cervical cancer specimens. *Folia Histochemica et Cytobiologica* 2012; 239-247.
3. S. Szostek *et al.* Physical state of human papillomavirus type 16 in cervical intraepithelial lesions and cancers determined by two different quantitative real-time PCR methods. *Acta Biochimica Polonica* 2015; 923-928.
4. Z. Chen *et al.* Classification and evolution of human papillomavirus genome variants: Alpha-5 (HPV26, 51, 69, 82), Alpha-6 (HPV30, 53, 56, 66), Alpha-11 (HPV34, 73), Alpha-13 (HPV54) and Alpha-3 (HPV61). *Virology* 2018; 86-101.
5. L. Zhang *et al.* Comparison of the performance in detection of HPV infections between the High-risk HPV Genotyping Real Time PCR and the PCR-Reverse Dot Blot assays. *Journal of Medical Virology*. DOI 10.1002/jmv.24931.
6. World Health Organization. Human papillomavirus (HPV) and cervical cancer. <http://www.who.int/mediacentre/factsheets/fs380/en/>
7. Centers for Disease Control and Prevention. HPV. <https://www.cdc.gov/hpv/parents/whatishpv.html>.

Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.


Mx3000P™ and Mx3005™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid.

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