Vitassay qPCR

Viral Meningitis

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<i>E. coli</i> 0.1285;O18:H7:K1	<i>Staphylococcus aureus</i> subsp <i>aureus</i>	Coxsackievirus type B3
<i>Listeria innocua</i> Serotype 6a/strain CCUG 15531	<i>Streptococcus agalactiae</i> Z019	Parechovirus Type 3
<i>Listeria ivanovii</i> Serovar 5/strain CCUG 15528	<i>Streptococcus pneumoniae</i> Z022	Parvovirus B19
<i>Listeria monocytogenes</i> Serotype 1/2b	<i>Plasmodium falciparum</i> 3D7	St Louis Encephalitis Virus strain 17D
<i>Listeria monocytogenes</i> Serovar 4b/Strain CIP 59.53	<i>Cryptococcus gatti</i> Z156	Tick-Borne encephalitis virus strain Neudorfl
<i>Bacteroides fragilis</i>	<i>Haemophilus influenzae</i> MinnA	Chikungunya virus (S27 Petersfield)
<i>Pseudomonas aeruginosa</i>	<i>Neisseria meningitidis</i> serogroup A	Dengue 1 (Hawaii A strain)
<i>Yersinia enterocolitica</i> O:3	<i>Treponema pallidum</i>	Dengue 2 (New Guinea C strain)
<i>Klebsiella oxytoca</i>	<i>Toxoplasma gondii</i> Type II	Dengue 3 (H87)
<i>Enterobacter cloacae</i> Serotype Cloaca A	Enterovirus type 68	Dengue 4 (H241)
<i>Enterobacter aerogenes</i> Serotype Cloaca B	Enterovirus type 71	West Nile Virus (Heja)
<i>Enterococcus durans</i>	Epstein-Barr virus	West Nile Virus (NY99)
<i>Enterococcus faecalis</i>	BK Virus Type Ib2	West Nile Virus (Ug37)
<i>Enterococcus faecium</i> serotype 11	BK Virus Type IV	Yellow Fever Virus (17D strain)
<i>Acinetobacter baumannii</i>	JC Virus Type 2B	Echovirus Type 30
<i>Candida albicans</i>	JC Virus Type 1A	Echovirus Type 11
<i>Aspergillus fumigatus</i>	Citomegalovirus strain AD-169	HHV6 strain Z29
<i>Pneumocystis jirovecii</i> Type A1	Coxsackievirus type A24	HHV6 Type A
<i>Pneumocystis jirovecii</i> g885652	Coxsackievirus type A9	HHV6 Type B

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- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time PCR System (DNA-Technology)
- DTPrime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)^I
- SmartCycler® (Cepheid)^I

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Intended use

Vitassay qPCR Viral Meningitis allows the detection and differentiation of herpes simplex virus 1 (HSV-1), herpes simplex virus 2 (HSV-2) and/or varicella zoster virus (VZV) by real-time PCR in clinical samples. The product is intended for use in the diagnosis of herpes simplex virus 1, herpes simplex virus 2 and/or varicella zoster virus infections alongside clinical data of the patient and other laboratory tests outcomes.

References

Vitassay qPCR Viral Meningitis 4x8-well strip, low profile	7041045
Vitassay qPCR Viral Meningitis 4x8-well strip, high profile	7042045

Materials/reagents provided

Code	Reagent/Material	Colour	Amount
7041S045/ 7042S045	Viral Meningitis strips low/high profile	-	4x8-well strip
7C045	Viral Meningitis Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yelow	1 vial x 1 mL
7004O	Optical caps	-	4x8 cap strip

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

Additional equipment and material required

- DNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Vortexer
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

Summary

HSV-1, HSV-2 and VZV are highly contagious infections which are common and endemic throughout the world. These pathogens are double-stranded DNA viruses belonging to *Herpesviridae* family. Most HSV-1 infections are acquired during childhood, and infection is lifelong. The vast majority of HSV-1 infections are oral herpes but a proportion of HSV-1 infections are genital herpes. Symptoms of oral herpes include painful blisters or open sores in or around the mouth. HSV-1 is mainly transmitted by oral-to-oral contact to cause oral herpes infection, via contact with the HSV-1 virus in sores, saliva, and surfaces in or around the mouth. However, HSV-1 can also be transmitted to the genital area through oral-genital contact to cause genital herpes. These can be asymptomatic or can have mild symptoms that go unrecognized. When symptoms do occur, genital herpes is characterised by one or more genital or anal blisters or ulcers.

People who already have HSV-1 infection are not at risk of getting it again, but they are still at risk of acquiring herpes simplex virus type 2 (HSV-2) genital infection. This virus is almost exclusively sexually transmitted, causing genital herpes. The infection is also lifelong and incurable. More women are infected with HSV-2 than men; this is because sexual transmission of this pathogen is more efficient from men to women than from women to men. Genital herpes infections often have no symptoms, or mild symptoms that go unrecognized. Typically, about 10-20% of people with HSV-2 infection report a prior diagnosis of genital herpes. When symptoms do occur, genital herpes is characterised by one or more genital or anal blisters or open sores. In addition to genital ulcers, symptoms of new genital herpes infections often include fever, body aches, and swollen lymph nodes. HSV-2 is mainly transmitted during sex, through contact with genital surfaces, skin, sores or fluids of someone infected with the virus.

In rare circumstances, HSV-1 or HSV-2 infections can be transmitted from a mother to her infant during delivery.

Varicella-zoster virus causes both varicella (chickenpox) and herpes zoster (HZ). VZV produces a generalized vesicular rash on the dermis in normal children, usually before 10 years of age. After primary infection, the virus persists in latent form and may emerge as herpes zoster (shingles). Varicella is highly contagious. The virus can be spread from person to person by direct contact, inhalation of aerosols from vesicular fluid of skin lesions of acute varicella or zoster. A person with varicella is contagious from 1-2 days before rash onset until the lesions have crusted. Varicella is often diagnosed clinically on the basis of a generalized maculopapulovesicular rash. For laboratory confirmation, skin lesions are the preferred specimen. Vesicular swabs or scrapings and scabs from crusted lesions can be used to identify varicella-zoster virus DNA by PCR.

Real-time PCR assays have proven to be a tool for the detection of herpes simplex virus types 1 and 2 and Varicella-zoster virus to shorten detection time and improve sensitivity.

Principle of the test

Vitassay qPCR Viral Meningitis test is based on the real-time amplification of a conserved region of the *US4* gene for herpes simplex virus type 1, *US6* gene for herpes simplex virus type 2 and *ORF29* gene for varicella zoster virus. After DNA isolation, the presence of herpes simplex virus 1, herpes simplex virus 2 and/or varicella zoster virus is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

Vitassay qPCR Viral Meningitis test is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The amplification of the herpes simplex virus type 1 DNA target sequence is detected through the FAM channel, varicella zoster virus DNA target in Cy5 channel, herpes simplex virus type 2 DNA target in ROX channel whereas the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex II).

Precautions

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or broken or if the foil has been broken or damaged.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.

- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

Procedures

Specimen collection, processing and DNA extraction

For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

ZP02012 MagPurix Viral/Pathogen Nucleic Acids Extraction Kit B, using the MagPurix 12A instrument (Zinexts Life Science Corp.).

NucleoMag® Pathogen (Macherey Nagel).

NucleoSpin RNA Virus (Macherey Nagel).

Invisorb® Spin Universal Kit (Strattec).

MagDEA Dx SV kit, using the magLEAD® 6gC instrument (Precision System Science Co.)

NX-48 Urine/Swab DNA Kit using the Nextractor® NX-48 system, (Genolution).

Positive control preparation

Reconstitute the lyophilized Viral Meningitis Positive Control (red tube) with the 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of DNA sample, negative and positive controls and add them into each well.
- Cover the wells with the caps provided. Spin down briefly (optional).

- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the FAM (herpes simplex virus 1), ROX (herpes simplex virus 2), Cy5 (varicella zoster virus) and HEX, JOE or VIC channels (Internal Control (IC)). If you use the Applied Biosystems 7500 Fast or the Stratagene Mx3005P™ check that passive reference option ROX is none. (Attached II).

Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive controls used in each run, must show an amplification curve for FAM (herpes simplex virus 1), ROX (herpes simplex virus 2) and Cy5 (varicella zoster virus), which validates the reaction.

Negative control

The negative controls included in each run, must show the absence of signal in Cy5, ROX and FAM which validates the reaction.

The experiment seems to be failed if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

The result interpretation is summarized in the following table:

Herpes simplex virus 1 (FAM)	Herpes simplex virus 2 (ROX)	Varicella zoster virus (Cy5)	Internal control	Negative control	Positive control	Interpretation
+	+	+	+/-	-	+	Herpes virus 1, Herpes virus 2 and Varicella zoster virus Positive
-	-	-	+	-	+	Herpes virus 1, Herpes virus 2 and Varicella zoster virus Negative
+	-	-	+/-	-	+	Herpes virus 1 Positive, Herpes virus 2 and Varicella zoster virus Negative
+	+	-	+/-	-	+	Herpes virus 1 and Herpes virus 2 Positive, and Varicella zoster virus Negative.
+	-	+	+/-	-	+	Herpes virus 1 and Varicella zoster virus Positive, and Herpes virus 2 Negative.
-	+	-	+/-	-	+	Herpes virus 2 Positive, Herpes virus 1 and Varicella zoster virus Negative.
-	+	+	+/-	-	+	Herpes virus 2 and Varicella zoster virus Positive, Herpes virus 1 Negative.
-	-	+	+/-	-	+	Varicella zoster virus Positive, Herpes virus 1 and Herpes virus 2 Negative.
+	+	+	+	+	+	Invalid
-	-	-	-	-	-	Invalid

(+) Positive: Amplification signal

(-) Negative: No amplification signal

If the negative samples do not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

Quality Control

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

Vitassay qPCR Viral Meningitis was evaluated using different QCMD panels and INSTAND panels. These panels consist of 91 clinical specimens (samples in transport medium, swabs, urine and plasma). The results were compared with the QCMD and INSTAND final reports. 17 of 91 samples from the panels were positive for herpes simplex virus type 1, 16 of 91 samples were positives for herpes simplex virus type 2 and 25 of 91 samples were positives for varicella zoster virus Vitassay qPCR Viral Meningitis.

The results show a high sensitivity and specificity to detect herpes simplex virus type 1, herpes simplex virus type 2 and varicella zoster virus using Vitassay qPCR Viral Meningitis.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of herpes simplex virus type 1, herpes virus simplex type 2 and varicella zoster virus templates ranging from 10^7 to 10^1 copies/rxn. This assay has a detection limit of ≥ 10 DNA copies per reaction for herpes simplex virus type 2 and varicella zoster virus and a detection limit of ≥ 50 DNA copies per reaction for herpes simplex virus type 1.

Analytical specificity

The analytical specificity for herpes simplex virus type 1, herpes simplex virus type 2 and varicella zoster virus was tested within the panel of different microorganisms.

No cross-reactivity of herpes simplex virus type 1, herpes simplex virus type 2 and varicella zoster virus was seen between any of the species:

Cross-reactivity assay		
<i>E. coli</i> 0.1285;O18:H7:K1	<i>Staphylococcus aureus</i> subsp <i>aureus</i>	Coxsackievirus type B3
<i>Listeria innocua</i> Serotype 6a/strain CCUG 15531	<i>Streptococcus agalactiae</i> Z019	Parechovirus Type 3
<i>Listeria ivanovii</i> Serovar 5/strain CCUG 15528	<i>Streptococcus pneumoniae</i> Z022	Parvovirus B19
<i>Listeria monocytogenes</i> Serotype 1/2b	<i>Plasmodium falciparum</i> 3D7	St Louis Encephalitis Virus strain 17D
<i>Listeria monocytogenes</i> Serovar 4b/Strain CIP 59.53	<i>Cryptococcus gatti</i> Z156	Tick-Borne encephalitis virus strain Neudorfl
<i>Bacteroides fragilis</i>	<i>Haemophilus influenzae</i> MinnA	Chikungunya virus (S27 Petersfield)
<i>Pseudomonas aeruginosa</i>	<i>Neisseria meningitidis</i> serogroup A	Dengue 1 (Hawaii A strain)
<i>Yersinia enterocolitica</i> O:3	<i>Treponema pallidum</i>	Dengue 2 (New Guinea C strain)
<i>Klebsiella oxytoca</i>	<i>Toxoplasma gondii</i> Type II	Dengue 3 (H87)
<i>Enterobacter cloacae</i> Serotype Cloaca A	Enterovirus type 68	Dengue 4 (H241)
<i>Enterobacter aerogenes</i> Serotype Cloaca B	Enterovirus type 71	West Nile Virus (Heja)
<i>Enterococcus durans</i>	Epstein-Barr virus	West Nile Virus (NY99)
<i>Enterococcus faecalis</i>	BK Virus Type Ib2	West Nile Virus (Ug37)
<i>Enterococcus faecium</i> serotype 11	BK Virus Type IV	Yellow Fever Virus (17D strain)
<i>Acinetobacter baumannii</i>	JC Virus Type 2B	Echovirus Type 30
<i>Candida albicans</i>	JC Virus Type 1A	Echovirus Type 11
<i>Aspergillus fumigatus</i>	Citomegalovirus strain AD-169	HHV6 strain Z29
<i>Pneumocystis jirovecii</i> Type A1	Coxsackievirus type A24	HHV6 Type A
<i>Pneumocystis jirovecii</i> g885652	Coxsackievirus type A9	HHV6 Type B

Analytical reactivity

Analytical reactivity of Vitassay qPCR Viral Meningitis for herpes simplex virus type 1 was evaluated against HSV 1 MacIntyre strain and HSV-1 (95/1906) strain showing positive results.

Analytical reactivity of Vitassay qPCR Viral Meningitis for herpes simplex virus type 2 was evaluated against HSV 2 MS strain, showing positive results.

Analytical reactivity of Vitassay qPCR Viral Meningitis for varicella zoster virus was evaluated against Varicella-Zoster Virus Ellen, Varicella Zoster Virus OKA and Varicella Zoster virus (9/84), showing positive results.

Compatibles real-time PCR equipment

Vitassay qPCR Viral Meningitis has been validated on the following equipments:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time PCR System (DNA-Technology)
- DTPRime Real Time Detection Thermal Cyclers (DNA-Technology)
- Rotor-Gene® Q (Qiagen)^I
- SmartCycler® (Cepheid)^I

I: For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler tubes.

II: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

Limitations

- This test provides a presumptive diagnosis of herpes simplex virus type 1, herpes simplex virus type 2 and varicella zoster virus infection. All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- This assay was tried with samples in transport medium, urogenital swabs, urine, plasma, CFS and epithelial cells of chickenpox lesions or crusts. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper DNA from clinical specimens must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by the different pathogens, either samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 RealTime System
Applied Biosystems	Applied Biosystems
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
BIONEER	ViiA™ 7 Real-Time PCR
Exicycler™ 96	Analytik Jena Biometra
Bio-Rad	TOptical
CFX96™ Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	BIONEER
Cepheid	Exicycler™ 96
SmartCycler®	Bio-Rad
Qiagen	CFX96™ Deep Well Real-Time PCR Detection
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
Roche	iCycler iQ™ 5 Real-Time PCR
LightCycler®480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler®96 Real-Time PCR System	MyiQ™ 2Real-Time PCR Detection System
Cobas z480 Analyzer	Cepheid
	SmartCycler®
	DNA-Technology
	DTlite Real-Time PCR System
	DTprime Real-time Detection Thermal Cyclers
	Eppendorf
	Mastercycler™ ep <i>realplex</i>
	Qiagen
	Rotor-Gene® Q
	Stratagene / Agilent Technologies
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

* See Attached III to configure exposure settings.

Attached II: Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

REAL-TIME PCR THERMOCYCLER	Vitassay CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Attached III: Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
DTlite Real-Time PCR System (DNA-Technology)	FAM	500
	HEX	500
	ROX	500
	Cy5	500
DTprime Real-time Detection Thermal Cycler (DNA-Technology)	FAM	500
	HEX	1000
	ROX	1000
	Cy5	1000

Bibliography/Bibliografía

1. <https://www.cdc.gov/chickenpox/hcp/clinical-overview.html>
2. <http://www.who.int/news-room/fact-sheets/detail/herpes-simplex-virus>
3. G. Freer et al. Varicella-zoster virus infection: natural history, clinical manifestations, immunity and current and future vaccination strategies. *New microbiológica*. 2018. 41, 2, 95-105.
4. <https://wwwnc.cdc.gov/travel/yellowbook/2018/infectious-diseases-related-to-travel/varicella-chickenpox>.

Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.


Mx3000P™ and Mx3029™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid.

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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