Vitassay qPCR

SARS-CoV-2

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Intended use

Vitassay qPCR SARS-CoV-2 allows the qualitative detection of SARS-CoV-2 by real-time RT-PCR in respiratory samples (nasopharyngeal and oropharyngeal swabs) and saliva samples. The product is intended for use in SARS-CoV-2 infections diagnosis alongside the patient's clinical data and other laboratory tests outcomes.

References

Vitassay qPCR SARS-CoV-2 4 x 8-well strip, low profile	7041046
Vitassay qPCR SARS-CoV-2 4 x 8-well strip, high profile	7042046
Vitassay qPCR SARS-CoV-2 8 x 8-well strip, low profile	7081046
Vitassay qPCR SARS-CoV-2 8 x 8-well strip, high profile	7082046
Vitassay qPCR SARS-CoV-2 96-well plate, low profile	7091046
Vitassay qPCR SARS-CoV-2 96-well plate, high profile	7092046

Materials/reagents provided

Reagents provided in references 7041046, 7042046, 7081046 and 7082046:

Reference	Reagent/Material	Colour	Amount
7041S046/ 7042S046	SARS-CoV-2 strips	-	4/8 x 8-well strip
7C046	SARS-CoV-2 Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4/8 x 8 cap strip

Reagents provided in references 7091046 and 7092046:

Reference	Reagent/Material	Colour	Amount
7091P046/ 7092P046	SARS-CoV-2 Plate	-	1 plate
7C046	SARS-CoV-2 Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	12 x 8 cap strip

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. To avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

Additional equipment and material required

- Real-time PCR instrument (thermocycler) (Attached I)
- RNA extraction kit
- Centrifuge for 1.5 mL tubes
- Vortex
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

Summary

Coronaviruses are unsegmented single-stranded RNA viruses belonging to family *Coronaviridae* of the order *Nidovirales*. Six kinds of human CoVs have been previously identified: *NL63*, *229E*, *OC43* and *HKU1* which cause upper respiratory tract disease and common cold symptoms and the severe acute respiratory syndrome coronavirus (SARS-CoV), and the Middle East respiratory syndrome coronavirus (MERS-CoV) which are highly pathogenic in humans, with high rates of severe pneumonia and fatal outcomes.

In December 2019, a group of pneumonia cases was reported at a wholesale seafood market in Wuhan, Hubei province, which was found to be caused by previously unknown Coronaviruses. Airway epithelial cells from infected patients were used to isolate a novel coronavirus, temporarily named 2019-nCoV. Later, it was found that the new coronavirus is related to the SARS-CoV but is sufficiently divergent from SARS-CoV to be considered a new human-infecting betacoronavirus. Therefore, the International Committee for the classification of viruses designated the name of this coronavirus as the severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2). The World Health Organization has named the disease caused by the SARS-CoV-2 as coronavirus disease 2019 (COVID-19).

People with COVID-19 have had a wide range of symptoms reported, ranging from mild symptoms to severe illness. Symptoms include fever or chills, cough, shortness of breath or difficulty breathing, fatigue, muscle or body aches, headache, loss of taste or smell, sore throat, congestion or runny nose, nausea or vomiting and diarrhea. Symptoms may appear 2-14 days after exposure to the virus. Older people, and those with underlying

medical problems like cardiovascular disease, diabetes, chronic respiratory disease, and cancer are more likely to develop serious illness and infection may progress to pneumonia, acute respiratory distress syndrome and multi-organ failure. The fatality rate ranges from 3% to 4%.

The COVID-19 may be transmitted from person to person through several different routes. The SARS-CoV-2 spreads primarily through droplets of saliva or discharge from the nose when an infected person coughs or sneezes but also by direct contact with an infected subject or indirect contact (through hand-mediated transfer of the virus from contaminated fomites to the mouth, nose, or eyes). Transmission of this virus is occurring from person to person, even during the asymptomatic incubation period.

WHO recommend, for all suspect cases, collection of upper respiratory tract (URT) specimens (nasopharyngeal and oropharyngeal) for testing by RT-PCR and, where clinical suspicion remains and URT specimens are negative, to collect specimens from the lower respiratory tract (LRT) when readily available (expectorated sputum, or endotracheal aspirate/bronchoalveolar lavage in ventilated patient). Additional clinical specimens may be collected as COVID-19 virus has been detected in blood, stool, urine, and saliva.

Diagnostic testing for the SARS-CoV-2 is currently undertaken using chest scan, whole genome sequencing and real-time reverse transcriptase-polymerase chain reaction (rRT-PCR).

Principle of the test

Vitassay qPCR SARS-CoV-2 is based on the real-time amplification of specific conserved fragments of the *ORF1ab* and *N* genes encoded by the SARS-CoV-2. The viral RNA extracted is transcribed into cDNA using a specific primer by reverse transcription step followed immediately in the same well by polymerase chain reaction. The SARS-CoV-2 presence is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

The assay is based on 5' exonuclease activity that uses two primers and a fluorogenic hydrolysis probe to detect accumulation of the amplified target sequence during the PCR reaction. When the polymerase begins to extend the primers, the probe is hydrolyzed by its 5'-3' exonuclease activity causing the spatial separation of the fluorophore and the quencher. The resulting increase in fluorescent signal is proportional to the amount of amplified product in the sample and is detected by means of a real-time PCR equipment.

Vitassay qPCR SARS-CoV-2 is a ready-to used test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The amplification of the target sequence *ORF1ab* gene is detected through the FAM channel, the amplification of the

target sequence *N* gene is detected through the ROX channel whereas the internal control (IC) in HEX, VIC, or JOE channel (depending on the equipment used).

Precautions

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or is broken or if the foil has been broken or damaged.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles, and mask.
- Do not eat, drink, or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment, and reagents to the area in which the previous step was performed. Use separate areas for the patient samples and controls preparation.
- Avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious and they must be managed according to the national safety legislation and national health waste legislation. Take necessary precautions during the collection, transport, storage, treatment, and disposal of samples.
- Use personal protective equipment (PPE) and biological safety cabinet for handling of potentially infectious samples and reagents according to current recommendations.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

Procedures

Specimen collection, processing, and RNA extraction

Patient samples must be collected, transport and storage according to appropriate laboratory guidelines. For pre-treatment and nucleic acid isolation, it is recommended to

use your optimized manual or automatic system compatible with real-time PCR technology. For respiratory samples, the assay has been validated with the following extraction kits:

Maxwell® RSC 16 Viral Total Nucleic Acid Purification Kit (Promega).

Total Nucleic Acid Isolation (TNAI) Kit (Roche).

MagDEA Dx SV kit, using magLEAD® 12gC instrument (Precision System Science Co.).

MagCore® Viral Nucleic Acid Extraction kit, using MagCore® HF16 automated Nucleic Acid Extractor System.

mSample Preparation Systems RNA, using Abbott m2000 RealTime System (Abbott Molecular).

QIAamp Viral RNA Mini Kit (Qiagen).

EMAG®Nucleic Acid Extraction System (bioMérieux, Inc).

For saliva samples, the assay has been validated with the following extraction kits:

QIAamp Viral RNA Mini Kit (Qiagen).

Maxwell® RSC 16 Viral Total Nucleic Acid Purification Kit, (Promega).

Positive control preparation

Reconstitute the lyophilized SARS-CoV-2 Positive Control (red tube) with the 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.

- Pipette 5 µL of RNA sample, negative (yellow tube) and positive (red tube) controls and add them to the corresponding wells.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Reverse transcription	45°C	15 min	1
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the FAM (*ORF1ab* gene), ROX (*N* gene) and HEX, JOE, or VIC channels (Internal Control (IC)). If you use the Applied Biosystems 7500 Fast, the Applied Biosystems StepOne™ or the Stratagene Mx3005P™ check that passive reference option ROX is none (attached II).

Analysis and interpretation of results

The results' analysis is done by the software itself of the used real-time PCR system following manufacturer's instructions. To verify the correct operation of the amplification mix, check the internal control (IC) signal emission.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive control used in each run must show an amplification curve for SARS-CoV-2 (Ct ≤40) in FAM and ROX channels, which validates the reaction.

Negative control

The negative control included in each run must show signal' absence for SARS-CoV-2 (Ct ≥40 or no signal) in FAM and ROX channels, which validates the reaction.

The Internal Control (IC) should show an amplification signal (Ct ≤40 in Positive and Negative Controls wells).

The experiment seems to be failed if there is amplification signal in negative control or signal' absence in the positive control for any channel. The assay should be repeated.

The clinical samples test results assessment should be performed once controls' results have been validated. The result interpretation is summarized in the following table:

ORF1ab gene (FAM)	N gene (ROX)	Internal Control (HEX)	Interpretation	
+	+	+/-*	Valid	SARS-CoV-2 RNA Detected
+	-	+/-*	Inconclusive	If only one SARS-CoV-2 target gene amplifies, repeat test depending on the available material: a) obtain a new specimen, re-extract, and retest (ideally) or, b) re-extract and retest another aliquot of the same specimen or, c) repeat RT-qPCR with the same isolated RNA sample. After retesting one time, if at least one target gene is positive, the sample should be considered SARS-CoV-2 positive.
-	+	+/-*	Inconclusive	
-	-	+ [#]	Valid	SARS-CoV-2 RNA Not detected [#]
-	-	- [#]	Invalid	Test failed [#]

(+) Positive: Amplification signal

(-) Negative: No amplification signal

* Sometimes, the Internal Control (IC) detection is not necessary since a high copies number of the target can cause a preferential amplification of target-specific nucleic acids.

[#] In the case of negative SARS-CoV-2 target genes detection, IC must show an amplification signal with Ct less than 40. If there is a signal' absence or Ct value ≥ 40 of Internal Control, the result is considered 'Invalid', and retesting is required. It is recommended to repeat the RT-qPCR by diluting the RNA sample 1:10 and/or 1:100, or to re-extract and retest to check for possible failure in the extraction procedure and/or inhibition problems.

If the obtained result is confusing or doubtful, it is necessary to check that all the steps have been carried out correctly, to verify the correct performance of each RT-qPCR steps and to review all the parameters, the sigmoid shape of the curve and the fluorescence intensity.

The test results must be evaluated by a health professional, together with the medical history, clinical symptoms and/or the results obtained in other diagnostic tests.

Quality Control

To confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

A total of 100 respiratory samples (nasopharyngeal swab) from symptomatic patients with suspicion of COVID-19, were tested by Real Time PCR using Vitassay qPCR SARS-CoV-2 and the results were compared with those obtained using a molecular detection method by the Spanish National Reference Center (ISCIII). SARS-CoV-2 was detected in the 2 positive samples by Vitassay qPCR SARS-CoV-2.

In another Vitassay qPCR SARS-CoV-2 clinical evaluation, 409 respiratory samples (nasopharyngeal swabs) from patients with suspected COVID-19 disease were analyzed by real-time PCR using Vitassay qPCR SARS-CoV-2 and the results were compared with those previously obtained with the molecular detection method Abbott Real Time SARS-CoV-2 assay (Abbott Molecular). SARS-CoV-2 was detected in 145 positive samples by Vitassay qPCR SARS-CoV-2. Eight discordant samples were obtained but the low amount of template RNA in these respiratory samples is below the detection limit of the method used. These samples were also negative by a test developed according to the 2019-Novel Coronavirus (2019-nCoV) Real-time rRT-PCR panel primers and probes (US CDC protocol).

Besides, the Vitassay qPCR SARS-CoV-2 clinical performance was conducted by testing in parallel 28 anonymized saliva specimens and oropharyngeal swabs from patients with signs and symptoms of COVID-19 infection. The values obtained with the oropharyngeal samples were taken as the reference values. SARS-CoV-2 was detected in 18 saliva samples of the 19 positive samples detected in oropharyngeal smears by Vitassay qPCR SARS-CoV-2. It was demonstrated that the results obtained when analyzing the saliva samples in parallel, were as expected.

These results show high agreement to detect SARS-CoV-2 using the molecular diagnostic test Vitassay qPCR SARS-CoV-2.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of SARS-CoV-2 template ranging from 10^7 to 10^1 copies/rxn. This assay has a detection limit of 20 genome copies per reaction for *ORF1ab* and *N* genes with a positive rate of $\geq 95\%$.

Analytical specificity

The analytical specificity for SARS-CoV-2 was tested within the panel of following microorganisms, where no cross-reactivity was seen between any of the species:

SARS-CoV-2		
<i>Bordetella bronchiseptica</i>	Influenza A/California/7/2009(H1N1)pdm09	<i>Legionella longbeachae</i>
<i>Bordetella holmesii</i>	Influenza A/Michigan/45/2015 (H1N1)pdm09 virus	<i>Legionella micdadei</i>
<i>Bordetella parapertussis</i>	Influenza A/Singapore/GP1908/2015, IVR-180 (H1N1)pdm09 virus	<i>Legionella pneumophila</i>
<i>Bordetella pertussis</i>	Influenza A/Virginia/ATCC2/2009 (H1N1)pdm09	Human metapneumovirus A, B and 8 Peru6-2003
<i>Candida albicans</i> 3147	Influenza A/Victoria/210/2009 (H3N2)	<i>Moraxella catarrhalis</i>
<i>Chlamydia caviae</i>	Influenza A/Thüringen/5/17 (H3N2) virus	<i>Mycoplasma pneumoniae</i>
<i>Chlamydia psittaci</i> genotype A and C	Influenza A/Switzerland/9715293/2013 (H3N2) virus	<i>Mycobacterium tuberculosis</i> not rifampin resistant
<i>Chlamydophila pneumoniae</i> CM-1 and TWAR strain 2023	Influenza A/Hong Kong/4801/2014, NYMC X-263B (H3N2) virus	Human parainfluenza 1, 2, 3 and 4 viruses
Human coronavirus 229E, OC43, NL63 and HKU1	Influenza A/South Australia/55/2014, IVR-175 (H3N2) virus	<i>Pneumocystis jirovecii</i>
MERS Coronavirus	Influenza A/Wisconsin/15/2009 (H3N2)	<i>Pseudomonas aeruginosa</i> PAO1LAC
SARS Coronavirus Strain Frankfurt 1	Influenza A/DE-SH/Reiherente/AR8444/ 2016 (H5N8) virus	Human rhinovirus type C and 17
Enterovirus 68, 71 and D68 US/MO/1418947	Influenza A/Anhui/1/2013 (H7N9) virus	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>
Enterovirus Echovirus 11 and 30	Influenza B/Brisbane/60/2008	<i>Staphylococcus epidermidis</i> strain RP62A
Enterovirus Coxsackievirus A24, A9 and B3	Influenza B/Florida/04/06 virus	<i>Streptococcus pneumoniae</i> Z022 and strain TIGR4
<i>Haemophilus influenzae</i> MinnA and Rd	Influenza B/Phuket/3073/2013 virus	<i>Streptococcus pyogenes</i> typing strain T1
Human Adenovirus types 1-5, 8, 14, 15, 31, 40 and 41	<i>Legionella bozemanii</i>	<i>Streptococcus salivarius</i> subsp. <i>salivarius</i> strain DSM 13084
Human Bocavirus	<i>Legionella dumoffii</i>	Respiratory syncytial virus (RSV) A and B

Analytical reactivity

The reactivity of Vitassay qPCR SARS-CoV-2 was confirmed by real-time amplification using RNA from Human 2019-nCoV strain BetaCoV/Germany/BavPat1/2020 p.1, Human 2019-nCoV strain 2019-nCoV/Italy-INMI1, SARS-CoV-2 strain 2019nCoV/USA-WA1/2020, SARS-CoV-2 strain hCoV-19/Netherlands/Diemen_1363454/2020 and synthetic RNA controls for two SARS-CoV-2 virus variants: MT007544.1 (SARS-CoV2 isolate Australia/VIC01/2020) and MN908947.3 (SARS-CoV-2 isolate Wuhan-Hu-1), showing positive result.

Compatibles real-time PCR equipment

Vitassay qPCR SARS-CoV-2 has been validated on the following equipments:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- CFX96 TM Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time PCR System (DNA-Technology)
- DTPrime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)^I
- Mic Real Time PCR Cycler bms^I

I: For Rotor-Gene® Q and Mic Real Time PCR Cycler bms. thermocyclers the product should be reconstituted following the procedure and transferred into specific tubes of each equipment.

II: When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder for the tubes (Ref. PN 4388506).

Limitations

- This test provides a presumptive diagnosis of SARS-CoV-2 infection. All results must be interpreted by a specialist together with other clinical information and laboratory findings available to the physician.
- This assay was tried with respiratory (nasopharyngeal and oropharyngeal swabs) and saliva samples. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper RNA from clinical specimens must be extracted.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.

- There is a possibility of false positive results due to cross-contamination by SARS-CoV-2, either high number of cDNA template copies which contains each SARS-CoV-2 Positive Control vial, samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.
- Vitassay qPCR SARS-CoV-2, designed for the specific SARS-CoV-2 detection, mostly showed no significant combined homologies with the human genome, other coronaviruses, or human microflora that would predict potential false positive RT-qPCR results (except for some N and/or ORF1ab sequences from SARS-CoV, and other coronaviruses identified in bats and pangolin).
- Detection of SARS-CoV-2 may be affected by several factors and their combinations which may lead to false negative results, including a) inadequate specimen sampling, shipping, storage, handling; b) procedural errors (including RNA isolation); c) RNA degradation during specimen shipping, storage, and/or preparation; d) potential mutations of the target regions of the SARS-CoV-2 genome covered by this test which may result in RNA being undetectable e) pathogen load below the limit of detection for the assay; f) the presence of retrotranscription and/or Real Time amplification inhibitors or other types of interference (an interference study evaluating the effect of vaccines, antiviral therapeutics, antibiotics, chemotherapeutics or immunosuppressant drugs related to COVID-19 was not performed); g) failure to follow the manufacturer's instructions and procedures.
- The single-gene amplification or even random (irregular) positive results is suggestive of slightly different target regions' amplification yield. A negative result for *ORF1ab* gene and the only amplification of the *N* gene may occur in specimens with low pathogen load, due to subtle differences in clinical sensitivity. In case of an ambiguous result, it is recommended to seek confirmatory testing by a reference laboratory if clinically indicated.
- Detection of viral RNA (*ORF1ab* and/or *N* gene sequences) may not indicate the presence of viable and/or infectious virus or that SARS-CoV-2 is the causative agent for clinical symptoms.
- Negative results or single target amplification do not preclude SARS-CoV-2 virus infection and should not be the sole basis of a patient treatment/management decision. Specimens types for the SARS-CoV-2 identification and/or stage of infection most suitable for its collection have not been established yet. Consider the collection of multiple specimens from the same patient at different time points, which may increase the probability of detecting the virus.
- If the patient clinical data, laboratory tests and epidemiological studies suggest possible SARS-CoV-2 infection, and other respiratory illnesses have been discarded, a false negative result might not be discarded, and additional tests should be performed.

Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 RealTime System
Applied Biosystems	Applied Biosystems
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
BIONEER	ViiA™ 7 Real-Time PCR
Exicycler™ 96	Analytik Jena Biometra
Bio-Rad	TOptical
CFX96™ Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	BIONEER
Bio Molecular Systems	Exicycler™ 96
Mic Real Time PCR Cyclor	Bio-Rad
Cepheid	CFX96™ Deep Well Real-Time PCR Detection
SmartCycler®	iCycler iQ™ Real-Time PCR
Precision System Science Co., Ltd.	iCycler iQ™ 5 Real-Time PCR
geneLEAD VIII System	MyiQ™ Real-Time PCR Detection System
Qiagen	MyiQ™ 2Real-Time PCR Detection System
Rotor-Gene® Q	Bio Molecular Systems
Roche	Mic Real Time PCR Cyclor
LightCycler®480 Real-Time PCR System	Cepheid
LightCycler®96 Real-Time PCR System	SmartCycler®
Cobas z480 Analyzer	DNA-Technology
	DTlite Real-Time PCR System*
	DTprime Real-time Detection Thermal Cyclor*
	Eppendorf
	Mastercycler™ ep <i>realplex</i>
	Qiagen
	Rotor-Gene® Q
	Precision System Science Co., Ltd.
	geneLEAD VIII System
	Stratagene / Agilent Technologies
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

* See Attached III to configure exposure settings.

Attached II: Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

Thermocycler	Vitassay Channel	Detection Channel	Observations
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Color Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Roche Cobas z 480	FAM	465/510	Color Compensation required
	HEX	540/580	
	ROX	540/610	
	Cy5	610/670	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mic Real Time PCR Cycler bms	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10)
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene® Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Attached III: Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
DTIite Real-Time PCR System (DNA-Technology)	FAM	500
	HEX	500
	ROX	500
	Cy5	500
DTprime Real-time Detection Thermal Cycler (DNA-Technology)	FAM	500*
	HEX	1000
	ROX	1000
	Cy5	1000

*If result in FAM channel is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

Bibliography/Bibliografía

1. Wang H, Li X, Li T, et al. (2020). The genetic sequence, origin, and diagnosis of SARS-CoV-2. *Eur J Clin Microbiol Infect Dis*. 39(9):1629-1635.
2. Ogimi C, Kim YJ, Martin ET, Huh HJ, Chiu CH, Englund JA. (2020). What's New With the Old Coronaviruses? *J Pediatric Infect Dis Soc*. 9(2):210-217.
3. Centers of Disease Control and Prevention (CDC). 2019 Novel Coronavirus, Symptoms. Available from <https://www.cdc.gov/coronavirus/2019-ncov/symptoms-testing/symptoms.html> Accessed December 2020.
4. World Health Organization. Coronavirus. https://www.who.int/health-topics/coronavirus#tab=tab_1 Accessed September 2020.
5. Lu R, Zhao X, Li J, et al. (2020). Genomic characterisation and epidemiology of 2019 novel coronavirus: implications for virus origins and receptor binding. *Lancet*. 395(10224):565-574.
6. Rothe, C., Schunk, M., Sothmann, P., Bretzel, G., Froeschl, G., Wallrauch, C., ... & Hoelscher, M. (2020). Transmission of 2019-nCoV Infection from an Asymptomatic Contact in Germany. *New England Journal of Medicine*.
7. Liu Z, Chu R, Gong L, Su B, Wu J. (2020). The assessment of transmission efficiency and latent infection period on asymptomatic carriers of SARS-CoV-2 infection [published online ahead of print, 2020 Jun 13]. *Int J Infect Dis*. S1201-9712(20)30471-9.
8. Kumar M, Taki K, Gahlot R, Sharma A, Dhangar K. (2020). A chronicle of SARS-CoV-2: Part-I - Epidemiology, diagnosis, prognosis, transmission and treatment. *Sci Total Environ*. 734:139278.
9. Chu DKW, Pan Y, Cheng SMS, et al. (2020). Molecular Diagnosis of a Novel Coronavirus (2019-nCoV) Causing an Outbreak of Pneumonia. *Clin Chem*. 66(4):549-555.
10. Lv DF, Ying QM, Weng YS, et al. (2020). Dynamic change process of target genes by RT-PCR testing of SARS-Cov-2 during the course of a Coronavirus Disease 2019 patient. *Clin Chim Acta*. 506:172-175.
11. Centers of Disease Control and Prevention (CDC). 2019 Novel Coronavirus, Real-time rRT-PCR Panel Primers and Probes. Available from <https://www.cdc.gov/coronavirus/2019-ncov/downloads/rt-pcr-panel-primer-probes.pdf> Accessed December 2020
12. European Centre for Disease Prevention and Control. Transmission of COVID-19. Available from: <https://www.ecdc.europa.eu/en/covid-19/latest-evidence/transmission> Accessed December 2020
13. Sharma A, Tiwari S, Deb MK, Marty JL. (2020). Severe acute respiratory syndrome coronavirus-2 (SARS-CoV-2): a global pandemic and treatment strategies. *Int J Antimicrob Agents*. 56(2):106054.

14. World Health Organization. Clinical management of COVID-19. Interim guidance 27 May 2020. Available from: <https://www.who.int/publications/i/item/clinical-management-of-covid-19> Accessed December 2020
15. World Health Organization. Laboratory testing for coronavirus disease 2019 (COVID-19) in suspected human cases. Interim guidance 19 March 2020. Available from <https://www.who.int/publications/i/item/laboratory-testing-for-2019-novel-coronavirus-in-suspected-human-cases-20200117> Accessed December 2020
16. To KK, Tsang OT, Yip CC, Chan KH, Wu TC, Chan JM, et al. (2020) Consistent Detection of 2019 Novel Coronavirus in Saliva. *Clin Infect Dis.* 71(15):841-843.
17. Peng L, Liu J, Xu W, Luo Q, Chen D, Lei Z, Huang Z, Li X, Deng K, Lin B, Gao Z. (2020) SARS-CoV-2 can be detected in urine, blood, anal swabs, and oropharyngeal swabs specimens. *J Med Virol.* 92(9):1676-1680.

Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™, StepOnePlus™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.

Mx3000P™, Mx3005™ and AriaMx are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.


Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid.

MagDEA® Dx is PSS patented technology (Japanese Patent: No.3115501; US Patent: No.5,702,950; European Patent: No.687501)

EMAG®Nucleic Acid Extraction System is a registered trademark of bioMérieux, Inc.

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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