Vitassay qPCR

Vaginosis

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Termocicladores con bloque de bajo perfil	Termocicladores con bloque de alto perfil
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 RealTime System
Applied Biosystems	Applied Biosystems
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
BIONEER	ViiA™ 7 Real-Time PCR
Exicycler™ 96	Analytik Jena Biometra
Bio-Rad	TOptical
CFX96™ Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	BIONEER
Cepheid	Exicycler™ 96
SmartCycler®	Bio-Rad
Qiagen	CFX96™ Deep Well Real-Time PCR Detection
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
Roche	iCycler iQ™ 5 Real-Time PCR
LightCycler® 480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler® 96 Real-Time PCR System	MyiQ™ 2 Real-Time PCR Detection System
Cobas z480 Analyzer	Cepheid
	SmartCycler®
	DNA-Technology
	DTlite Real-Time PCR System*
	DTprime Real-time Detection Thermal Cycler*
	Eppendorf
	Mastercycler™ ep <i>realplex</i>
	Qiagen
	Rotor-Gene® Q
	Stratagene / Agilent Technologies
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

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Intended use

Vitassay qPCR Vaginosis allows the detection and differentiation of *Candida albicans*, *Gardnerella vaginalis* and *Trichomonas vaginalis* by real-time PCR in clinical samples. The product is intended for use in the diagnosis of *Candida albicans*, *Gardnerella vaginalis* and *Trichomonas vaginalis* infections alongside clinical data of the patient and other laboratory tests outcomes.

References

Vitassay qPCR Vaginosis 4x8-well strip, low profile 7041051

Vitassay qPCR Vaginosis 4x8-well strip, high profile 7042051

Materials/reagents provided

Reference	Reagent/Material	Colour	Amount
7041S051/ 7042S051	Vaginosis strips low/high profile	-	4x8-well strip
7C051	Vaginosis Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4x8 cap strip

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

Additional equipment and material required

- DNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Vortexer
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

Summary

Candida albicans is the most common human opportunistic fungal pathogen. *Candida albicans* can be isolated from the gastrointestinal tract, and oral and vaginal mucosa of many healthy individuals. Although, a number of fungal species belonging to the genus *Candida* can cause acute vulvovaginal candidiasis (VVC), *Candida albicans* is the most prevalent etiological agent, particularly for the most severe chronic condition known as recurrent vulvovaginal candidiasis (RVVC). Although most cases of VVC are still caused by *Candida albicans*, non-*albicans Candida* (NAC) species are increasingly being implicated in the etiology of VVC with the most common including *C. glabrata*, *C. krusei*, *C. parapsilosis*, *C. tropicalis*, *C. dubliniensis*, *C. africana*, *C. duobushaemulonii*, and *C. auris*.

Vulvovaginal candidiasis is a commonly diagnosed vaginitis. Symptoms are nonspecific and include vaginal discharge, soreness, irritation, burning, and dyspareunia. In addition to the morbidity associated with VVC, *C. albicans* colonization has been associated with increased HIV acquisition in longitudinal studies. Vulvovaginal candidiasis is relatively common with 75% of women experiencing at least one episode of VVC in their lifetime, and 5% developing RVVC. Although *Candida* species are not regarded as classical sexually transmitted pathogens, there is some evidence that this organism can be transmitted by both oral and vaginal intercourse.

Gardnerella vaginalis is implicated as one of the causative agents of bacterial vaginosis (BV), but it can also be isolated from the vagina of healthy women. One explanation is that certain strains of *G. vaginalis* are more pathogenic than others. Recent advances in BV pathogenesis research have determined the existence of 13 different species within the genus *Gardnerella*. Bacterial vaginosis is a vaginal dysbiosis and is linked with a variety of heterogeneous pathogenic anaerobic bacteria which results in the decrease in number of healthy *Lactobacillus* species.

Epidemiological data strongly suggest that BV is a sexually transmitted infection (STI) with an incubation period of around 4 days. The infection is most common in women of age from 15-44 years. Approximately 50% to 75% of women with BV are asymptomatic. When symptoms are present, the most common presenting symptoms include an off-white, thin, homogenous vaginal discharge and/or an unpleasant vaginal odor or "fishy smell" that is particularly noticeable after sexual intercourse and around the time of menstruation. The prevalence rate of BV is about 21.2 million (29.2%) worldwide and BV prevalence rate is higher in African countries in parallel to Asia and Europe population.

BV disease is linked with a compelling amount of obstetric and gynecological complications which includes spontaneous abortion, preterm delivery, subclinical pelvic inflammatory disease, post-Caesarean delivery wound infections, chorioamnionitis, premature rupture of membranes, postpartum endometritis, and postsurgical infections.

Moreover, bacterial vaginosis is associated with an increased risk of acquisition of sexually transmitted infections including human immunodeficiency virus (HIV).

Trichomonas vaginalis is a flagellated protozoan parasite of the human genital tract which causes Trichomoniasis, the most prevalent non-viral sexually transmitted infection worldwide. The protist genus *Trichomonas* contains this globally distributed endemic human pathogen called *Trichomonas vaginalis*. The parasite resides in the female lower genital tract and the male urethra and prostate. Transmission occurs almost exclusively via sexual contact.

The World Health Organization (WHO) estimated that 248 million people were infected with *T. vaginalis* in 2005, and by 2008 this number had increased by 11.8% to 276.4 million. The prevalence and incidence of vaginal *T. vaginalis* infection are higher in African regions and in the Americas than in other parts of the world. The prevalence of *T. vaginalis* infection for men it was calculated as one tenth of the prevalence in women.

In women, the disease may range from asymptomatic (up to 50% of women) to severe with serious sequelae. Women infected with *T. vaginalis* are at high risk of vaginitis, endometritis, atypical pelvic inflammatory disease, and adverse pregnancy outcomes. This parasite can also result in serious consequences, such as infertility, premature rupture of placental membranes, premature delivery, low-birth-weight infants, and neonatal death. Classical symptoms include a malodorous and purulent discharge which results in local pain and irritation. Men are often asymptomatic carriers of *T. vaginalis*, although dysuria and discharge have also been reported and *T. vaginalis* increases the risk of infertility and prostate cancer in men. Furthermore, *T. vaginalis* infection has been linked with increased risk and transmission of HIV for both men and women.

Principle of the test

Vitassay qPCR Vaginosis test is based on the real-time amplification of *T. vaginalis*-specific 2-kb repeated sequence for *Trichomonas vaginalis*, a conserved region of 5.8S rRNA gene for *Candida albicans* and 16S rRNA gene for *Gardnerella vaginalis*. After DNA extraction, the presence of *Candida albicans*, *Gardnerella vaginalis* and *Trichomonas vaginalis* is detected by an increase in observed fluorescence during the reaction upon hydrolysis of the fluorescent probe.

Vitassay qPCR Vaginosis test is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The amplification of *Trichomonas vaginalis* DNA target sequence is detected through the FAM channel, *Candida albicans* DNA target in ROX channel and *Gardnerella vaginalis* DNA in Cy5 channel whereas the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Attached II).

Precautions

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or broken or if the foil has been broken or damaged.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

Procedures

Specimen collection, processing and DNA extraction

For pre-treatment and nucleic acid isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

ZP02014 Mag Purix Plant DNA Extraction Kit, using the MagPurix 12A instrument (Zinexts Life Science Corp.).

Invisorb® Spin Universal Kit (Invitex).

Positive control preparation

Reconstitute the lyophilized Vaginosis Positive Control (red tube) with the 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of DNA sample, negative and positive controls and add them to the corresponding wells.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the FAM (*Trichomonas vaginalis*), ROX (*Candida albicans*), Cy5 (*Gardnerella vaginalis*) and HEX, JOE or VIC channels (Internal Control (IC)). If you use the Applied Biosystems 7500 Fast or the Stratagene Mx3005P™ Real Time PCR System check that passive reference option ROX is none. (See Attached II)

Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive controls used in each run, must show an amplification curve for FAM (*Trichomonas vaginalis*), ROX (*Candida albicans*) and Cy5 (*Gardnerella vaginalis*), which validates the reaction.

Negative control

The negative controls included in each run, must show the absence of signal in FAM, ROX and Cy5 which validates the reaction.

The experiment seems to be failed if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

The result interpretation is summarized in the following table:

<i>Trichomonas vaginalis</i>	<i>Candida albicans</i>	<i>Gardnerella vaginalis</i>	Internal Control	Negative Control	Positive Control	Interpretation
+	+	+	+/-	-	+	<i>Trichomonas vaginalis</i> , <i>Candida albicans</i> and <i>Gardnerella vaginalis</i> Positive
-	-	-	+	-	+	<i>Trichomonas vaginalis</i> , <i>Candida albicans</i> and <i>Gardnerella vaginalis</i> Negative
+	-	-	+/-	-	+	<i>Trichomonas vaginalis</i> Positive, <i>Candida albicans</i> and <i>Gardnerella vaginalis</i> Negative
+	+	-	+/-	-	+	<i>Trichomonas vaginalis</i> and <i>Candida albicans</i> Positive, and <i>Gardnerella vaginalis</i> Negative
+	-	+	+/-	-	+	<i>Trichomonas vaginalis</i> and <i>Gardnerella vaginalis</i> Positive, and <i>Candida</i> <i>albicans</i> Negative
-	+	-	+/-	-	+	<i>Candida albicans</i> Positive, <i>Trichomonas vaginalis</i> and <i>Gardnerella vaginalis</i> Negative
-	+	+	+/-	-	+	<i>Candida albicans</i> and <i>Gardnerella vaginalis</i> Positive, <i>Trichomonas</i> <i>vaginalis</i> Negative
-	-	+	+/-	-	+	<i>Gardnerella vaginalis</i> Positive, <i>Trichomonas</i> <i>vaginalis</i> and <i>Candida</i> <i>albicans</i> Negative
-	-	-	-	-	+	Experiment fail
+	+	+	+	+	-	Experiment fail

(+) Positive: Amplification signal

(-) Negative: No amplification signal

If the negative samples do not show a positive result for the internal control, they should be retested from the diluted original sample 1:10 or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

A positive test result for *C. albicans* and/or *G. vaginalis* does not necessarily indicate vulvovaginal candidiasis or bacterial vaginosis (respectively), since those microorganisms can be part of the vaginal microbiota in healthy women. The results should be evaluated by a professional alongside clinical data of the patient and other laboratory tests outcomes.

Quality Control

In order to confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

Vitassay qPCR Vaginosis was evaluated with the following EQA programmes: QCMD 2017 *Candida* spp. EQA Programme and QCMD 2017 Immunocompromised EQA Pilot Study for *Candida albicans* and QCMD 2017 Sexually Transmitted Infections I EQA Pilot Study and QCMD 2018 *Trichomonas vaginalis* EQA Pilot Study for *Gardnerella vaginalis* and *Trichomonas vaginalis*.

The samples were positive or negative to the analysed pathogens. Vitassay qPCR Vaginosis detected correctly all the samples (7 *Candida albicans* positive samples, 2 *Gardnerella vaginalis* positive samples and 9 *Trichomonas vaginalis* positive samples). Negative samples were reported as negative. The results were compared with the corresponding reports.

The results show a high sensitivity and specificity to detect *Candida albicans*, *Gardnerella vaginalis* and *Trichomonas vaginalis* using the molecular diagnostic test Vitassay qPCR Vaginosis.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of *Candida albicans*, *Gardnerella vaginalis* and *Trichomonas vaginalis* templates ranging from 10^7 to 10^1 copies/rxn. This assay has a detection limit of ≥ 10 viral DNA copies per reaction.

Analytical specificity

The analytical specificity for *Candida albicans*, *Gardnerella vaginalis* and *Trichomonas vaginalis* was tested within the panel of following pathogens, where no cross-reactivity was seen between any of the species:

Cross-reactivity assay		
<i>Candida albicans</i>	<i>Escherichia coli</i> 0.1285; O18:H7:K1	Human papillomavirus genotypes HPV16 and HPV18
<i>Candida parapsilosis</i>	<i>Gardnerella vaginalis</i>	HSV-1 and HSV-2
<i>Candida tropicalis</i>	<i>Haemophilus influenzae</i>	<i>Serratia marcescens</i>
<i>Candida glabrata</i>	<i>Haemophilus ducrey</i>	<i>Staphylococcus aureus</i>
<i>Candida krusei</i>	<i>Klebsiella oxytoca</i>	<i>Stenotrophomonas maltophilia</i>
<i>Candida dubliniensis</i>	<i>Klebsiella pneumoniae</i>	<i>Streptococcus pneumoniae</i>
<i>Aspergillus fumigatus</i>	<i>Listeria ivanovii</i>	<i>Streptococcus agalactiae</i>
<i>Bacteroides fragilis</i>	<i>Listeria monocytogenes</i>	<i>Trichomonas vaginalis</i>
<i>Acinetobacter baumannii</i>	<i>Listeria Innocua</i>	<i>Treponema pallidum</i>
<i>Chlamydia trachomatis</i> (LGV)	<i>Mycoplasma genitalium</i>	<i>Ureaplasma parvum</i>
<i>Chlamydia trachomatis</i> (SW)	<i>Mycoplasma hominis</i>	<i>Ureaplasma urealyticum</i>
<i>Enterobacter cloacae</i>	<i>Neisseria gonorrhoeae</i>	Hepatitis A
<i>Enterobacter aerogenes</i>	<i>Neisseria meningitidis</i>	Cytomegalovirus
<i>Enterococcus faecalis</i>	<i>Proteus mirabilis</i>	<i>Pseudomonas aeruginosa</i>
<i>Enterococcus faecium</i>		

Analytical reactivity

Analytical reactivity of Vitassay qPCR Vaginosis for *Candida albicans* was evaluated against *Candida albicans*, showing positive result.

Analytical reactivity of Vitassay qPCR Vaginosis for *Gardnerella vaginalis* was evaluated against *Gardnerella vaginalis*, showing positive result.

Analytical reactivity of Vitassay qPCR Vaginosis for *Trichomonas vaginalis* was evaluated against *Trichomonas vaginalis* showing positive result.

Compatibles real-time PCR equipment

Vitassay qPCR Vaginosis has been validated on the following equipments:

- Cobas z480 Analyzer (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTprime Real-time Detection Thermal Cycler (DNA-Technology)
- DTlite Real-Time PCR System (DNA-Technology)
- Rotor-Gene® Q (Qiagen) ^I
- SmartCycler® (Cepheid) ^I

I: For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler tubes.

II: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

Limitations

- All results must be interpreted together with other clinical information and laboratory findings available to the physician.
- This assay was tried with DNA extracted from endocervical, vaginal and vagino-rectal smears, urine, plasma, bronchoalveolar lavage (BAL), sputum, saline, and samples in transport medium. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper DNA from clinical specimens must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target template DNA below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by the different pathogens, either samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers
Agilent Technologies
AriaMx Real-Time PCR System
Applied Biosystems
7500 Fast Real-Time PCR System
7500 Fast Dx Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast
QuantStudio™ 6 Flex 96-well Fast
QuantStudio™ 7 Flex 96-well Fast
QuantStudio™ 5 Real-Time PCR System
QuantStudio™ 3 Real-Time PCR System
StepOne Plus™ Real-Time PCR System
StepOne™ Real-Time PCR System
ViiA™ 7 Fast Real-Time PCR System
BIONEER
Exicycler™ 96
Bio-Rad
CFX96™ Real-Time PCR Detection System
Mini Opticon™ Real-Time PCR Detection System
Cepheid
SmartCycler®
Qiagen
Rotor-Gene® Q
Roche
LightCycler®480 Real-Time PCR System
LightCycler®96 Real-Time PCR System
Cobas z480 Analyzer

High profile Block Thermocyclers
Abbott
Abbott m2000 RealTime System
Applied Biosystems
7300 Real-Time PCR System
7500 Real-Time PCR System
7900 HT Real-Time PCR System
ABI PRISM 7000
ABI PRISM 7700
QuantStudio™ 12K Flex 96-well
QuantStudio™ 6 Flex 96-well
QuantStudio™ 7 Flex 96-well
QuantStudio™ 5 Real-Time PCR
QuantStudio™ 3 Real-Time PCR System
ViiA™ 7 Real-Time PCR
Analytik Jena Biometra
TOptical
qTOWER 2.0
BIONEER
Exicycler™ 96
Bio-Rad
CFX96™ Deep Well Real-Time PCR Detection
iCycler iQ™ Real-Time PCR
iCycler iQ™ 5 Real-Time PCR
MyiQ™ Real-Time PCR Detection System
MyiQ™ 2Real-Time PCR Detection System
Cepheid
SmartCycler®
DNA-Technology
DTlite Real-Time PCR System
DTprime Real-time Detection Thermal Cycler
Eppendorf
Mastercycler™ ep <i>realplex</i>
Qiagen
Rotor-Gene® Q
Stratagene / Agilent Technologies
Mx3000P™ Real Time PCR System
Mx3005P™ Real Time PCR System

* See Attached III to configure exposure settings.

Attached II: Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

REAL-TIME PCR THERMOCYCLER	Vitassay CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene® Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Attached III: Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
DTlite Real-Time PCR System (DNA-Technology)	FAM	500
	HEX	500
	ROX	500
	Cy5	500
DTprime Real-time Detection Thermal Cycler (DNA-Technology)	FAM	500
	HEX	1000
	ROX	1000
	Cy5	1000

Bibliography/Bibliografía

1. Reiter S, Kellogg Spadt S. (2019). Bacterial vaginosis: a primer for clinicians. *Postgrad Med.* 131(1):8-18.
2. Javed A, Parvaiz F, Manzoor S. (2019). Bacterial vaginosis: An insight into the prevalence, alternative treatments regimen and it's associated resistance patterns. *Microb Pathog.* 127:21-30.
3. Muzny CA, Taylor CM, Swords WE, et al. (2019). An Updated Conceptual Model on the Pathogenesis of Bacterial Vaginosis. *J Infect Dis.* 220(9):1399-1405.
4. Potter RF, Burnham CD, Dantas G. (2019). In Silico Analysis of Gardnerella Genomespecies Detected in the Setting of Bacterial Vaginosis. *Clin Chem.* 65(11):1375-1387.
5. Peters A, Das S, Raidal SR. (2020). Diverse Trichomonas lineages in Australasian pigeons and doves support a columbid origin for the genus Trichomonas. *Mol Phylogenet Evol.* 143:106674.
6. Noh CS, Kim SS, Park SY, Moon HS, Hong Y, Ryu JS. (2019). Comparison of Two PCR Assays for Trichomonas vaginalis. *Korean J Parasitol.* 57(1):27-31.
7. Bouchemal K, Bories C, Loiseau PM. (2017). Strategies for Prevention and Treatment of Trichomonas vaginalis Infections. *Clin Microbiol Rev.* 30(3):811-825.
8. Edwards T, Burke P, Smalley H, Hobbs G. (2016). Trichomonas vaginalis: Clinical relevance, pathogenicity and diagnosis. *Crit Rev Microbiol.* 42(3):406-417.
9. Mielczarek E, Blaszkowska J. (2016). Trichomonas vaginalis: pathogenicity and potential role in human reproductive failure. *Infection.* 44(4):447-458.
10. Šoba B, Skvarč M, Matičič M. (2015). Trichomoniasis: a brief review of diagnostic methods and our experience with real-time PCR for detecting infection. *Acta Dermatovenerol Alp Pannonica Adriat.* 24(1):7-10.
11. Brown SE, Schwartz JA, Robinson CK, et al. (2019). The Vaginal Microbiota and Behavioral Factors Associated with Genital Candida albicans Detection in Reproductive-Age Women. *Sex Transm Dis.* 46(11):753-758.

12. Cassone A. (2015). Vulvovaginal *Candida albicans* infections: pathogenesis, immunity and vaccine prospects. *BJOG*. 122(6):785-794.
13. Kim J, Sudbery P. (2011). *Candida albicans*, a major human fungal pathogen. *J Microbiol*. 49(2):171-177.
14. Makanjuola O, Bongomin F, Fayemiwo SA. (2018). An Update on the Roles of Non-*albicans* *Candida* Species in Vulvovaginitis. *J Fungi (Basel)*. 4(4):121. Published 2018 Oct 31.

Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.


Mx3000P™ and Mx3029™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid.

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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