Vitassay qPCR

SARS-CoV-2 (ORF1ab, E & N)

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Intended use

Vitassay qPCR SARS-CoV-2 (ORF1ab, E & N) allows the qualitative detection of SARS-CoV-2 RNA by real-time RT-PCR in respiratory samples. This product is intended to aid in the diagnosis of infections caused by SARS-CoV-2 virus, alongside the patient's clinical data, epidemiological risk factors, and other laboratory tests outcomes.

References

Vitassay qPCR SARS-CoV-2 (ORF1ab, E & N) 4x8-well strip, low profile	7041060
Vitassay qPCR SARS-CoV-2 (ORF1ab, E & N) 4x8-well strip, high profile	7042060
Vitassay qPCR SARS-CoV-2 (ORF1ab, E & N) 96-well plate, low profile	7091060
Vitassay qPCR SARS-CoV-2 (ORF1ab, E & N) 96-well plate, high profile	7092060

Materials/reagents provided

Reagents provided in references 7041060 and 7042060:

Reference	Reagent/Material	Colour	Amount
7041S060/ 7042S060	SARS-CoV-2 (ORF1ab, E & N) strips low/high profile	-	4 x 8-well strip
7C060	SARS-CoV-2 (ORF1ab, E & N) Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4 x 8-cap strip

Reagents provided in references 7091060 and 7092060:

Reference	Reagent/Material	Colour	Amount
7091P060/ 7092P060	SARS-CoV-2 (ORF1ab, E & N) Plate low/high profile	-	1 plate
7C060	SARS-CoV-2 (ORF1ab, E & N) Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	12 x 8-cap strip

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. To avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

Additional equipment and material required

- Collection and transport system.
- Laboratory freezers (- 30°C to - 10°C and/or ≤ -70°C).
- RNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Vortex
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

Summary

Coronaviruses are unsegmented single-stranded RNA viruses belonging to family Coronaviridae of the order Nidovirales. Six kinds of human CoVs have been previously identified: NL63, 229E, OC43 and HKU1 which cause upper respiratory tract disease and common cold symptoms and the severe acute respiratory syndrome coronavirus (SARS-CoV), and the Middle East respiratory syndrome coronavirus (MERS-CoV) which are highly pathogenic in humans, with high rates of severe pneumonia and fatal outcomes.

In December 2019, a group of pneumonia cases was reported at a wholesale seafood market in Wuhan, Hubei province, which was found to be caused by previously unknown Coronaviruses. Airway epithelial cells from infected patients were used to isolate a novel coronavirus, temporarily named 2019-nCoV. Later, it was found that the new coronavirus is related to the SARS-CoV but is sufficiently divergent from SARS-CoV to be considered a new human-infecting betacoronavirus. Therefore, the International Committee for the classification of viruses designated the name of this coronavirus as the severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2). The World Health Organization has named the disease caused by the SARS-CoV-2 as coronavirus disease 2019 (COVID-19).

People with COVID-19 have had a wide range of symptoms reported, ranging from mild symptoms to severe illness. Symptoms include fever or chills, cough, shortness of breath

or difficulty breathing, fatigue, muscle or body aches, headache, loss of taste or smell, sore throat, congestion or runny nose, nausea or vomiting and diarrhea. Symptoms may appear 2-14 days after exposure to the virus. Older people, and those with underlying medical problems like cardiovascular disease, diabetes, chronic respiratory disease, and cancer are more likely to develop serious illness and infection may progress to pneumonia, acute respiratory distress syndrome and multi-organ failure. The fatality rate ranges from 3% to 4%.

The COVID-19 may be transmitted from person to person through several different routes. The SARS-CoV-2 spreads primarily through droplets of saliva or discharge from the nose when an infected person coughs or sneezes but also by direct contact with an infected subject or indirect contact (through hand-mediated transfer of the virus from contaminated fomites to the mouth, nose, or eyes). Transmission of this virus is occurring from person to person, even during the asymptomatic incubation period.

WHO recommend, for all suspect cases, collection of upper respiratory tract (URT) specimens (nasopharyngeal and oropharyngeal) for testing by RT-PCR and, where clinical suspicion remains and URT specimens are negative, to collect specimens from the lower respiratory tract (LRT) when readily available (expectorated sputum, or endotracheal aspirate/bronchoalveolar lavage in ventilated patient). Additional clinical specimens may be collected as COVID-19 virus has been detected in blood, stool, urine, and saliva.

Diagnostic testing for the SARS-CoV-2 is undertaken using chest scan, whole genome sequencing and real-time reverse transcriptase-polymerase chain reaction (RT-PCR).

Principle of the test

Vitassay qPCR SARS-CoV-2 (ORF1ab, E & N) is based on the real-time amplification of a conserved region of SARS-CoV-2 *ORF1ab*, *E* and *N* genes. After RNA extraction, the SARS-CoV-2 presence is detected by an increase in observed fluorescence during the reaction, after hydrolysis of the fluorescent probe.

The assay is based on 5' exonuclease activity using two primers and a fluorogenic hydrolysis probe to detect accumulation of the amplified target sequence during the PCR reaction. When the polymerase begins to spread the primers, the probe is hydrolysed by its exonuclease 5'-3' activity causing the spatial separation of the fluorophore and the quencher. The increase in the resulting fluorescent signal is proportional to the amount of amplified product in the sample and is detected by means of real-time PCR equipment.

Vitassay qPCR SARS-CoV-2 (ORF1ab, E & N) is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an endogenous internal control allows the extraction process control and the detection of a possible inhibition reaction. The assay uses a human housekeeping gene

as an endogenous internal control (IC) (*RNase P* gene present in human DNA) that is expected to be present in all nucleated human cells. The amplification of the target sequence *ORF1ab* gene is detected through the FAM channel, the amplification of the target sequence *E* gene is detected through the HEX, VIC, or JOE channel (depending on the equipment used), the amplification of the target sequence *N* gene is detected through the Cy5 channel whereas the endogenous internal control (IC) in ROX channel.

Precautions

- For professional *in vitro* diagnostic use (use by qualified and trained clinical laboratory personnel).
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or is broken or if the foil has been broken or damaged.
- Do not remove desiccant from reagent pouches once is open.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles, and mask.
- Do not eat, drink, or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment, and reagents to the area in which the previous step was performed. Use separate areas for the patient samples and controls preparation.
- Avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile disposable aerosol resistant pipette tips or RNase/DNase positive displacement pipette tips is recommended.
- An appearance of the reaction mixture in stabilized format, normally found at the tube bottom, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or color different from whitish) does not alter the test functionality.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious and they must be managed according to the national safety legislation and national health waste legislation. Take necessary precautions during the collection, transport, storage, treatment, and disposal of samples.
- Use personal protective equipment (PPE) and biological safety cabinet for handling of potentially infectious samples and reagents according to current recommendations.

- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

Procedures

Specimen collection, processing, and RNA extraction

For specimen collection, conservation and transport, user-validated conditions must be followed. The Vitassay qPCR SARS-CoV-2 (ORF1ab, E & N) has been tested on nasopharyngeal samples collected with synthetic fiber swabs with plastic and deposited immediately into a sterile transport tube in different transport media as Viral Transport Media (Vircell®), BD™ universal viral transport medium, Universal Transport Medium and preservation medium (Biocomma®). Other sample types should be validated by the user.

Overall, respiratory samples should be collected and labelled appropriately in clean containers with or without transport media (depending on sample type) and be processed as soon as possible to guarantee the test quality. It is recommended to use fresh specimens for the test. Transport must always be carried out in accordance with local and national regulations for the transport of biological samples. The specimens should be transported at 2°C to 8°C for up to 72 hours and for long term transport (more than 72 hours), it is recommended shipping at -20°C or lower. The samples can be stored at 2°C to 8°C for up to 72 hours or frozen at -20°C or ideally at -70°C for conservation during a long time. Repeated freeze-thaw cycles should be avoided to prevent sample and nucleic acids degradation.

The clinical specimens must be collected, transported, and stored according to appropriate laboratory guidelines. For more information, refer to the following guidelines:

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For nucleic acid isolation from clinical specimens, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. For the sample's pretreatment follow the instructions for use of the extraction kit used. The following extraction kits have been validated:

MagMAX™ Viral/Pathogen II (MVP II) Nucleic Acid Isolation Kit using the KingFisher Flex System (ThermoFisher).

MagDEA Dx SV kit, using the magLEAD® 12gC instrument (Precision System Science Co).

Maxwell® RSC 16 Viral Total Nucleic Acid (TNA) Purification Kit (Promega).

Positive control preparation

Reconstitute the lyophilized SARS-CoV-2 (ORF1ab, E & N) Positive Control (red tube) with 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template, and it is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of RNA sample, negative (yellow tube) or positive (red tube) controls and add them into the corresponding wells.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Reverse transcription	45°C	15 min	1
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the channels FAM (*ORF1ab* gen), HEX, JOE, or VIC (*E* gen), Cy5 (*N* gen) and ROX (Endogenous Internal Control (IC)). Depending on the equipment used select the proper detection channel (Attached II).

Analysis and interpretation of results

The results analysis is done by the software itself of the used real-time PCR system following manufacturer's instructions. To verify the correct operation of the amplification mix, and the extraction procedure, check the endogenous internal control (IC) signal emission.

Use the positive control amplification curve as a starting point during reaction validation (prior to interpretation of sample results), to ensure that the threshold falls within the exponential phase of the amplification curves and above any background noise signals. The threshold value may vary between different instruments due to different signal intensities. It is recommended to set the threshold values for each channel (target) independently by the end user.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive control used in each run must show an amplification curve ($C_t \leq 40$) in FAM, ROX, Cy5 and HEX, VIC, or JOE channels, which validates the reaction.

Negative control

The negative control included in each run must show signal' absence ($C_t > 40$ or no signal) in FAM, Cy5, ROX and HEX, VIC, or JOE channels, which validates the reaction.

The positive control includes the human *housekeeping RNase P* gene target; therefore, amplification signals are observed in all target channels, including the Endogenous Internal Control.

The experiment seems to be failed if there is amplification signal in negative control or signal absence in the positive control for any channel. The assay should be repeated.

The clinical samples test results assessment should be performed once controls' results have been validated. The result interpretation is summarized in the following table:

ORF1ab gen (FAM)	E gen (HEX)	N gen (Cy5)	Endogenous Internal Control (ROX)	Interpretation	
+	+	+	+/-*	Valid	SARS-CoV-2 RNA Detected. If one of the genes is not amplified, this may be due to: Low concentration of the target gene below the detection limit, genetic mutations in the specific target region or other factors.
+	+	-	+/-*	Valid	
+	-	+	+/-*	Valid	
-	+	+	+/-*	Valid	
+	-	-	+/-*	Inconclusive	If only one SARS-CoV-2 target gene amplifies, repeat test depending on the available material: a) obtain a new specimen, re-extract and retest (ideally) or, b) re-extract and retest another aliquot of the same specimen or, c) repeat RT-qPCR with the same isolated RNA sample. After retesting one time, if at least one target gene is positive, the sample should be considered SARS-CoV-2 positive.
-	-	+	+/-*	Inconclusive	
-	+	-	+/-*	Inconclusive	
-	-	-	+ [#]	Valid	Sarbecovirus RNA is detected but no specific SARS-CoV-2 RNA is detected. Repeat the test. In case of obtaining the same result, additional confirmation is required to differentiate between SARS-CoV-1, SARS-CoV-2 and other Sarbecoviruses whose infection in humans is unknown.
-	-	-	- [#]	Invalid	Failed test [#]

(+) Positive: Amplification signal (Ct ≤40)

(-) Negative: No amplification signal (Ct >40 or no signal)

* Sometimes, the Endogenous Internal Control (IC) detection is not necessary since a high copies number of the target can cause a preferential amplification of target-specific nucleic acids. The Endogenous Internal Control (IC) shows or not an amplification signal (Ct ≤40 or no signal).

In the case of negative SARS-CoV-2 target genes detection, IC must show an amplification signal with Ct ≤35. The Ct value could be very variable due to the Endogenous Internal Control is a human *housekeeping* gene that should be present in all human nucleated cells in the original sample. If there is a signal' absence or Ct value > 35 of Endogenous Internal Control, the result is considered 'invalid', and retesting is required. It is recommended to repeat the RT-qPCR by diluting the RNA sample 1:10 and/or 1:100, or to re-extract and retest to check for possible failure in the extraction procedure and/or inhibition problems.

If the obtained result is confusing or doubtful, it is necessary to check that all the steps have been carried out correctly, to verify the correct performance of each RT-qPCR steps and to review all the parameters, the sigmoid shape of the curve and the fluorescence intensity. It is also recommended to repeat the assay, preferably in duplicate, depending on the available material (repeat RT-qPCR with the same isolated RNA sample, or re-extract and retest another aliquot of the same specimen or, obtain a new specimen and retest).

The test results must be evaluated by a health professional, together with the medical history, clinical symptoms and/or the results obtained in other diagnostic tests.

Quality Control

To confirm the appropriate performance of the molecular diagnostic technique, an Endogenous Internal Control (IC) is included in each reaction. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

For the Vitassay qPCR SARS-CoV-2 (ORF1ab, E & N) kit's clinical performance evaluation, a total of 1003 remnants of nasopharyngeal samples from patients with clinical manifestations and/or suspected SARS-CoV-2 infection, previously characterised by molecular methods and/or sequencing in the case of those positive for SARS-CoV-2, were used. After swab collection, the sample was inoculated in transport media and analysed by routine method with different molecular assays (Panther Fusion SARS-CoV-2 assay (Hologic®), Simplexa™ COVID-19 Direct Kit (DiaSorin Molecular), Allplex™ SARS-CoV-2 Assay (Seegene) and Cobas® SARS-CoV-2 Real Time RT-PCR test (Roche)) and sequencing for initial diagnosis. The remnant was frozen at -20°C until further analysis with the Vitassay kit. Extraction was performed with the MagMAX™ Viral/Pathogenic Nucleic Acid Isolation Kit using the KingFisher Flex System

(ThermoFisher) and the thermocycler used was the CFX96™ Real-Time PCR Detection System (Bio-Rad).

Initial characterisation showed 556 negative and 447 positive samples for SARS-CoV-2, among which there were different variants, including 14 samples of the initial Wuhan strain and different lineages of variants later emerged.

Vitassay qPCR SARS-CoV-2 (ORF1ab, E & N) detected 465 positive and 538 negative samples, of which 24 were discrepant results: 21 false positives and 3 false negatives. After this comparison, sensitivity, and specificity values with 95% confidence interval were calculated, where sensitivity was 0.99 (0.98-0.99), and specificity 0.96 (0.94-0.97).

In conclusion, Vitassay qPCR SARS-CoV-2 (ORF1ab, E & N) has a high sensitivity and specificity for the SARS-CoV-2 detection, when compared to sequenced samples.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of standards from the pathogens ranging from 10^7 to 10^1 copies/reaction. Vitassay qPCR SARS-CoV-2 (ORF1ab, E & N) has a detection limit of ≥ 0.625 genome copies per reaction for *ORF1ab* gene, ≥ 1.25 genome copies per reaction for *E* gene and ≥ 1.25 genome copies per reaction for *N* gene with a positive rate of 95%.

Analytical specificity

The analytical specificity for SARS-CoV-2 detection was tested within the panel of following microorganisms, where no cross-reactivity was observed between any of the species:

Cross-reactivity testing					
<i>Chlamydia caviae</i>	-	Influenza A/Hong Kong/4801/2014, NYMC X-263B (H3N2) virus	-	<i>Mycobacterium tuberculosis</i> not rifampin resistant	-
<i>Chlamydia psittaci</i> genotypes A and C	-	Influenza A/South Australia/55/2014, IVR-175 (H3N2) virus	-	Human Parainfluenza virus 1, 2, 3 and 4	-
<i>Chlamydophila pneumoniae</i> CM-1	-	Influenza A/DE-SH/Reiherente/AR8444/2016 (H5N8) virus	-	<i>Pneumocystis jirovecii</i> Type A1 and g885652	-
Human Coronavirus 229E, OC43, NL63 and HKU1	-	Influenza A/Anhui/1/2013 (H7N9) virus	-	Human Rhinovirus type C and 17	-
MERS-CoV strain Florida/USA-2_Saudi Arabia_2014	-	Influenza A/PR/8/34 (H1N1) virus	-	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-
Human Adenovirus types 1-5, 8, 15, 31, 40 and 41	-	Influenza A/California/7/2009(H1N1)pdm09 virus	-	<i>Legionella micdadei</i>	-
Bocavirus	-	Influenza A/Michigan/45/2015 (H1N1)pdm09 virus	-	<i>Legionella pneumophila</i>	-
<i>Bordetella bronchiseptica</i>	-	Influenza A/Singapore/GP1908/2015, IVR-180 (H1N1)pdm09 virus	-	Human Metapneumovirus A and B	-
<i>Bordetella holmesii</i>	-	Influenza A/Victoria/210/2009 (H3N2) virus	-	<i>Moraxella catarrhalis</i>	-
<i>Bordetella parapertussis</i>	-	Influenza A/Thüringen/5/2017 (H3N2) virus (Clade 3C2a.1)	-	<i>Mycoplasma pneumoniae</i>	-
<i>Bordetella pertussis</i>	-	Influenza A/Switzerland/9715293/2013 (H3N2) virus	-	<i>Mycobacterium tuberculosis</i> TMC 331	-
SARS Coronavirus strain Frankfurt 1	-	Influenza B/Brisbane/60/2008 virus	-	<i>Staphylococcus epidermidis</i>	-
Enterovirus 68 and 71	-	Influenza B/Florida/04/06 virus	-	<i>Streptococcus pneumoniae</i> Z022	-
Enterovirus Echovirus 30	-	Influenza B/Phuket/3073/2013 virus	-	<i>Streptococcus pyogenes</i>	-
Enterovirus Coxsackievirus A24, A9 and B3	-	<i>Legionella bozemanii</i>	-	<i>Streptococcus salivarius</i>	-
<i>Haemophilus influenzae</i> MinnA	-	<i>Legionella dumoffii</i>	-	Respiratory syncytial virus (RSV) A and B	-
Influenza A/New Caledonia/20/99(H1N1) virus	-	<i>Legionella longbeachae</i>	-	Respiratory syncytial virus (RSV) strain Long	-

Analytical reactivity

The Vitassay qPCR SARS-CoV-2 (ORF1ab, E & N) kit's reactivity for SARS-CoV-2 was evaluated against the following strains: Human 2019-nCoV strain BetaCoV/Germany/BavPat1/2020 p.1; Human 2019-nCoV strain 2019-nCoV/Italy-INMI1; Human 2019-nCoV strain 2019nCoV/USA-WA1/2020; SARS-CoV-2 BetaCoV/Berlin/ChVir1670/2020_IsolatBER, SARS-CoV-2 BetaCoV/Munich/ChVir984/2020; SARS-CoV-2 BetaCoV/Baden-Wuerttemberg/1/ChVir1577/2020_IsolatBER and synthetic RNA controls for seven SARS-CoV-2 virus variants: MT007544.1 (SARS-CoV-2/human/AUS/VIC01/2020); MN908947.3 (SARS-CoV-2 isolate Wuhan-Hu-1); B.1.1.7_710528 (Twist Synthetic SARS-CoV-2 RNA Control 14); B.1.1.7_601443 (Twist Synthetic SARS-CoV-2 RNA Control 15); B.1.351 (Twist Synthetic SARS-CoV-2 RNA Control 16); P.1 (Twist Synthetic SARS-CoV-2 RNA Control 17) and B.1.617.1 (Twist Synthetic SARS-CoV-2 RNA Control 18), showing positive results.

Compatibles real-time PCR equipment

Vitassay qPCR SARS-CoV-2 (ORF1ab, E & N) has been validated on the following equipment:

- Cobas z480 Analyzer (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^I
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTprime Real-time Detection Thermal Cycler (DNA-Technology)
- DTlite Real-Time PCR System (DNA-Technology)
- Rotor-Gene® Q (Qiagen) ^{II}
- CFX Opus 96 Touch™ Real-Time PCR Detection System (Bio-Rad)
- LightCycler 480 Instrument II (Roche)
- NEOS-96 Real Time PCR System (Linear)

^I: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

^{II}: For Rotor-Gene® Q thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q tubes.

Limitations

- All obtained results should be interpreted by a specialist in conjunction with available clinical information and laboratory findings.
- This assay has been validated with RNA extracted from nasopharyngeal samples. The use of other samples has not been established.
- The test quality depends on the sample's quality; RNA must be properly extracted from clinical samples.
- This test does not provide quantitative values nor indicate the number of organisms present. This test is a qualitative assay.
- A low copy number of the target template may be detected below the limit of detection, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination, either by SARS-CoV-2 (ORF1ab, E & N) positive control during reconstitution, by samples containing high concentrations of target template, or by PCR products of previous reactions.
- Detection may be affected by several factors and their combinations which may lead to false negative results, including a) inadequate specimen sampling, shipping, storage, handling; b) procedural errors (including RNA isolation); c) RNA degradation during specimen shipping, storage, and/or preparation; d) pathogen load below the limit of detection for the assay; e) the presence of retrotranscription and/or Real-Time amplification inhibitors or other types of interference (an interference study evaluating the effect of vaccines, antiviral therapeutics, antibiotics, chemotherapeutics or immunosuppressant drugs used to prevent the infection or used during the treatment of the infection was not performed); f) Mutations or polymorphisms in primer or probes binding regions that may affect the detection of new or unknown SARS-CoV-2 variants; g) failure to follow the manufacturer's instructions and procedures. If in doubt, refer to section Analysis and interpretation of results to check the correct interpretation of the results.
- Detection of viral RNA may not indicate the presence of viable and/or infectious virus or that SARS-CoV-2 is the causative agent for clinical symptoms.
- Negative results do not preclude SARS-CoV-2 virus infection and should not be the sole basis of a patient treatment/management decision. Optimal specimen types and/or stage of infection most suitable for its collection have not been established yet. Consider the collection of multiple specimens from the same patient at different time points, which may increase the probability of detecting the virus.
- If the patient clinical data, laboratory tests and epidemiological studies suggest possible SARS-CoV-2 infection, and other respiratory illnesses have been

discarded, a false negative result might not be discarded, and additional tests should be performed.

- Some samples may fail to exhibit *RNase P* amplification curves due to low human cell numbers in the original clinical sample. A negative endogenous IC signal does not preclude the SARS-CoV-2 RNA presence in a clinical specimen.
- Fluorescence values may vary due to multiple factors such as PCR equipment, extraction system, sample type and its previous treatment, among others.

Attached I. Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 ⁽¹⁾
AriaDx Real-Time PCR System	Applied Biosystems
Applied Biosystems	7300 Real-Time PCR System ^{(3) (1)}
7500 Fast Real-Time PCR System ^{(1) (2)}	7500 Real-Time PCR System ⁽¹⁾
7500 Fast Dx Real-Time PCR System ^{(1) (2)}	7900 Real-Time PCR System ⁽³⁾
QuantStudio™ 12K Flex 96-well Fast	ABI PRISM 7000 ⁽³⁾
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7700 ⁽³⁾
QuantStudio™ 7 Flex 96-well Fast	QuantStudio™ 12K Flex 96-well
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
QuantStudio™ 5 Fast Real-Time PCR System	QuantStudio™ 7 Flex 96-well
QuantStudio™ 3 Fast Real-Time PCR System ⁽³⁾	QuantStudio™ 5 Fast Real-Time PCR System
StepOne Plus™ Real-Time PCR System ⁽³⁾	QuantStudio™ 5 Real-Time PCR System
StepOne™ Real-Time PCR System ⁽⁴⁾	QuantStudio™ 3 Real-Time PCR System ⁽³⁾
ViiA™ 7 Fast Real-Time PCR System	ViiA™ 7 Real-Time PCR System
Azure Biosystems	Agilent Technologies
Azure Cielo 3 ⁽⁵⁾	Mx3000P™ Real Time PCR System
Azure Cielo 6	Mx3005P™ Real Time PCR System
BIONEER	Analytik Jena
Exicycler™ 96	qTOWER
Bio-Rad	BIONEER
CFX96™ Real-Time PCR Detection System	Exicycler™ 96
Mini Opticon™ Real-Time PCR Detection System ⁽⁵⁾	BIOER
Roche	QuantGene 9600
LightCycler @480 Real-Time PCR System ^{(1) (6)}	Bio-Rad
LightCycler @96 Real-Time PCR System ⁽¹⁾	CFX96™ Deep Well Real-Time PCR
Cobas z480 Analyzer ^{(1) (6)}	iCycler iQ™ Real-Time PCR
	iCycler iQ™ 5 Real-Time PCR
	My iQ™ Real-Time PCR Detection System ⁽⁵⁾
	My iQ™ 2 Real-Time PCR Detection System ⁽⁵⁾
Special Formats ⁽⁷⁾	DNA-Technology
Bio Molecular Systems	DTlite Real-Time PCR System ⁽⁸⁾
Mic Real Time PCR Cyclor	DTprime Real-time Detection Thermal Cyclor ⁽⁸⁾
Cepheid	Eppendorf
SmartCycler®	Mastercycler™ ep <i>realplex</i>
Precision System Science Co., Ltd.	Qiagen
geneLEAD VIII System	QIAquant 96
Qiagen	
Rotor-Gene® Q	

- (1) A special grid is needed to fit these real-time PCR kits.
- (2) Select Ramp Speed "Standard".
- (3) No ROX caption.
- (4) No Cy5 caption.
- (5) Only FAM and HEX caption.
- (6) Specific compensation color is required.
- (7) The product must be reconstituted following the appropriate procedure (see Test procedure) and transferred to the specific tubes for Mic, SmartCycler®, Rotor-Gene® Q or geneLEAD VIII System.
- (8) See Attached III to configure exposure settings.

Attached II. Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR thermocyclers are specified in the following table:

Thermocycler	Vitassay Channel	Detection Channel	Observations
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Applied Biosystems ABI 7500	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Roche Cobas z 480	FAM	465/510	Colour Compensation required
	HEX	540/580	
	ROX	540/610	
	Cy5	610/670	
Cepheid Smartcycler®	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Agilent Technologies Mx3000P™ Mx 3005P™	FAM	FAM	Passive reference option ROX is none
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Bio Molecular Systems Mic Real Time PCR Cycler	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10).
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Agilent Technologies AriaMx	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Qiagen Rotor-Gene® Q	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquaring". The fluorescence Target Sample Range must be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
BIONEER Exicycler™ 96	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Attached III. Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
DTIite Real-Time PCR System (DNA-Technology)	FAM	500
	HEX	500
	ROX	500
	Cy5	500
DTprime Real-time Detection Thermal Cycler (DNA-Technology)	FAM	500*
	HEX	1000
	ROX	1000
	Cy5	1000

*If result in FAM channel is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

Bibliography/Bibliography


1. Wang H, Li X, Li T, et al. (2020). The genetic sequence, origin, and diagnosis of SARS-CoV-2. *Eur J Clin Microbiol Infect Dis*. 39(9):1629-1635.
2. Ogimi C, Kim YJ, Martin ET, Huh HJ, Chiu CH, Englund JA. (2020). What's New With the Old Coronaviruses? *J Pediatric Infect Dis Soc*. 9(2):210-217.
3. Centers of Disease Control and Prevention (CDC). 2019 Novel Coronavirus, Symptoms. Available from <https://www.cdc.gov/coronavirus/2019-ncov/symptoms-testing/symptoms.html> Accessed December 2020.
4. World Health Organization. Coronavirus. https://www.who.int/health-topics/coronavirus#tab=tab_1 Accessed September 2020.
5. Lu R, Zhao X, Li J, et al. (2020). Genomic characterisation and epidemiology of 2019 novel coronavirus: implications for virus origins and receptor binding. *Lancet*. 395(10224):565-574.
6. Rothe, C., Schunk, M., Sothmann, P., Bretzel, G., Froeschl, G., Wallrauch, C., ... & Hoelscher, M. (2020). Transmission of 2019-nCoV Infection from an Asymptomatic Contact in Germany. *New England Journal of Medicine*.
7. Liu Z, Chu R, Gong L, Su B, Wu J. (2020). The assessment of transmission efficiency and latent infection period on asymptomatic carriers of SARS-CoV-2 infection [published online ahead of print, 2020 Jun 13]. *Int J Infect Dis*. S1201-9712(20)30471-9.
8. Kumar M, Taki K, Gahlot R, Sharma A, Dhangar K. (2020). A chronicle of SARS-CoV-2: Part-I - Epidemiology, diagnosis, prognosis, transmission and treatment. *Sci Total Environ*. 734:139278.
9. Chu DKW, Pan Y, Cheng SMS, et al. (2020). Molecular Diagnosis of a Novel Coronavirus (2019-nCoV) Causing an Outbreak of Pneumonia. *Clin Chem*. 66(4):549-555.
10. Lv DF, Ying QM, Weng YS, et al. (2020). Dynamic change process of target genes by RT-PCR testing of SARS-Cov-2 during the course of a Coronavirus Disease 2019 patient. *Clin Chim Acta*. 506:172-175.
11. Centers of Disease Control and Prevention (CDC). 2019 Novel Coronavirus, Real-time rRT-PCR Panel Primers and Probes. Available from <https://www.cdc.gov/coronavirus/2019-ncov/downloads/rt-pcr-panel-primer-probes.pdf> Accessed December 2020

12. European Centre for Disease Prevention and Control. Transmission of COVID-19. Available from: <https://www.ecdc.europa.eu/en/covid-19/latest-evidence/transmission> Accessed December 2020
13. Sharma A, Tiwari S, Deb MK, Marty JL. (2020). Severe acute respiratory syndrome coronavirus-2 (SARS-CoV-2): a global pandemic and treatment strategies. *Int J Antimicrob Agents*. 56(2):106054.
14. World Health Organization. Clinical management of COVID-19. Interim guidance 27 May 2020. Available from: <https://www.who.int/publications/i/item/clinical-management-of-covid-19> Accessed December 2020
15. World Health Organization. Laboratory testing for coronavirus disease 2019 (COVID-19) in suspected human cases. Interim guidance 19 March 2020. Available from <https://www.who.int/publications/i/item/laboratory-testing-for-2019-novel-coronavirus-in-suspected-human-cases-20200117> Accessed December 2020
16. To KK, Tsang OT, Yip CC, Chan KH, Wu TC, Chan JM, et al. (2020) Consistent Detection of 2019 Novel Coronavirus in Saliva. *Clin Infect Dis*. 71(15):841-843.
17. Peng L, Liu J, Xu W, Luo Q, Chen D, Lei Z, Huang Z, Li X, Deng K, Lin B, Gao Z. (2020) SARS-CoV-2 can be detected in urine, blood, anal swabs, and oropharyngeal swabs specimens. *J Med Virol*. 92(9):1676-1680.

Trademarks

All trademarks that may appear in this package insert are the property of their respective owners.

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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