Vitassay qPCR

NG CiproRes

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- 7500 Fast Real-Time PCR System (Applied Biosystems) ^I
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- CFX Opus 96 Touch™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTPrime Real Time Detection Thermal Cycler (DNA-Technology)
- DTlite Real-Time PCR System (DNA-Technology)
- Rotor-Gene® Q (Qiagen) ^{II}
- Cobas z480 Analyzer (Roche)
- LightCycler® 480 Instrument II (Roche)
- MIC qPCR Cycler (Bio Molecular Systems) ^{II}
- NEOS-96 Real Time PCR System (Linear)

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Intended use

Vitassay qPCR NG CiproRes allows the identification and qualitative detection of DNA from *N. gonorrhoeae* and/or a specific point mutation in *N. gonorrhoeae* (in the *gyrA* gene, changing the serine in position 91 of the wild type genotype to phenylalanine) that generates resistance to ciprofloxacin (CIP), from urethral urine, urethral swabs, endocervical swabs and rectal swabs samples already characterized as *N. gonorrhoeae* positive (using molecular assays) of patients suspected of gonorrhea disease. This product is intended to aid in the diagnosis of infections caused by *N. gonorrhoeae* and its possible resistance or sensitivity to CIP, alongside the patient's clinical data and other laboratory tests outcomes.

References

Vitassay qPCR NG CiproRes 4x8-well strip, low profile	7041068
Vitassay qPCR NG CiproRes 4x8-well strip, high profile	7042068

Materials/reagents provided

Reference	Reagent/Material	Colour	Amount
7041S068/ 7042S068	NG CiproRes strips low/high profile	-	4 x 8-well strip
7C068	NG CiproRes Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4 x 8-cap strip

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated on the label.
- The resuspended positive control should be stored at -20°C. To avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

Additional equipment and material required

- Collection and transport system.
- Laboratory freezers (- 30°C to - 10°C and/or ≤ -70°C).
- DNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Vortex
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

Summary

Sexually Transmitted Infections (STIs) embrace many clinical syndromes and infections caused by pathogens that can be acquired and transmitted through sexual activity. There is an increase in all STIs worldwide, and antibiotic resistance is also on the rise.

Neisseria gonorrhoeae is a Gram-negative diplococcus, an obligate human pathogen which is the etiological agent of gonorrhoea. *N. gonorrhoeae* belongs to the genus *Neisseria*, in which is also included *Neisseria meningitidis*, which causes bacterial meningitis. *N. gonorrhoeae* mainly colonizes the genital mucosa, but it can also colonize the ocular, nasopharyngeal and anal mucosa. There are at least eight non-pathogenic commensal *Neisseria* spp. located in the human nasal and oropharyngeal flora. As both the commensal and pathogenic *Neisseria* spp. are in the same niches, it is often difficult to differentiate colonization factors from the virulence factors that are necessary to elicit host damage.

The gonorrhoea is a sexually transmitted infection, which remains a major global public health concern. This disease has a high morbidity burden, with more than 106 million new cases being diagnosed every year worldwide. The morbidity is increasing exponentially due to *Neisseria gonorrhoeae* facility to acquire resistance to all antimicrobials introduced for its treatment.

Complications from untreated, ascending, genital-tract infections in women can include pelvic inflammatory disease, infertility, and ectopic pregnancy. The pathogen cannot survive outside the host, so its transmission requires a sexual contact. If the infection is disseminated, can cause arthritis and endocarditis. Also, maternal transmission can lead to neonatal blindness. Female genital infections are mostly asymptomatic and male genital infection are mostly symptomatic, which is based on the fact that in men is easier to observe the symptoms such as purulent exudate from the penis and painful urination, whereas in women the clinical manifestation is less painful and specific and could be mistaken for bacterial vaginosis.

Seeing as drug resistance appears to be unstoppable, new treatment options are necessary to control the disease. There are three approaches for the development of new therapies against drug-resistant *N. gonorrhoeae*: novel combinations of already existing antibiotics, development of new antibiotics and development of alternative therapies.

Ciprofloxacin is a second-generation fluoroquinolone antibiotic that is commonly used to treat mild to moderate urinary and respiratory tract infections. It is also used against infectious diarrhea, typhoid fever, uncomplicated gonorrhea, treatment of *Neisseria meningitidis* nasal carriage and prophylaxis against anthrax. The fluoroquinolones work by inhibiting the bacterial enzymes DNA gyrase and topoisomerase IV, which are needed for mRNA production (transcription) and DNA replication. The advantage is that human cells lack DNA gyrases, and human topoisomerases are not affected by fluoroquinolone suppression.

Ciprofloxacin resistance in *N. gonorrhoeae* is usually caused by mutations in the *gyrA* and *parC* genes, which code for topoisomerases II and IV, respectively. Resistance is normally determined by simultaneous combinations of alterations in both proteins, but it depends on geographical region. S91F or T, A92P, D95N or A or G, 197M and Q102H are the most commonly documented substitutions in the *N. gonorrhoeae* GyrA protein.

Principle of the test

Vitassay qPCR NG CiproRes is based on the real-time amplification of a conserved region of *porA* gene (*N. gonorrhoeae*) and a specific point mutation in the *gyrA* gene (*N. gonorrhoeae*) that generates resistance to ciprofloxacin (CIP). After nucleic acids isolation, the pathogens presence is detected by an increase in observed fluorescence during the reaction, after hydrolysis of the fluorescent probe.

The assay is based on 5' exonuclease activity using two primers and a fluorogenic hydrolysis probe to detect accumulation of the amplified target sequence during the PCR reaction. When the polymerase begins to spread the primers, the probe is hydrolyzed by its exonuclease 5'-3' activity causing the spatial separation of the fluorophore and the quencher. The increase in the resulting fluorescent signal is proportional to the amount of amplified product in the sample and is detected by means of real-time PCR equipment.

Vitassay qPCR NG CiproRes is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible inhibition reaction. After the amplification reaction, *N. gonorrhoeae ciprofloxacin resistant* is detected in FAM channel, *N. gonorrhoeae* is detected in ROX channel, the Internal Control (IC) is detected in Cy5 channel, and the *N. gonorrhoeae* ciprofloxacin sensitive is detected in HEX, VIC, or JOE channels (depending on the equipment used).

Precautions

- For professional *in vitro* diagnostic use (use by qualified and trained clinical laboratory personnel).
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or is broken or if the foil has been broken or damaged.
- Do not remove desiccant from reagent pouches once is open.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles, and mask.
- Do not eat, drink, or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment, and reagents to the area in which the previous step was performed. Use separate areas for the patient samples and controls preparation.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious and they must be managed according to the national safety legislation and national health waste legislation. Take necessary precautions during the collection, transport, storage, treatment, and disposal of samples.
- An appearance of the reaction mixture in stabilized format, normally found at the tube bottom, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or color different from whitish) does not alter the test functionality.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile disposable aerosol resistant pipette tips or RNase/DNase positive displacement pipette tips is recommended.
- Use personal protective equipment (PPE) and biological safety cabinet for handling of potentially infectious samples and reagents according to current recommendations.

Procedures

Specimen collection, transport, and conservation

For specimen collection, conservation, and transport, user-validated conditions must be followed. The Vitassay qPCR NG CiproRes has been tested in urethral urine, endocervical swabs, urethral swabs, and rectal swabs already characterized as positive for *N. gonorrhoeae* and collected in sterile tubes (APTIMA® Urine Specimen Collection Kit – HOLOGIC). Other sample types should be validated by the user.

Overall, clinical samples should be collected appropriately in clean containers, correctly labelled and promptly processed to ensure the test's quality. It is recommended to use fresh specimens for the test. Transport must always be carried out in accordance with local and national regulations for the transport of biological samples. The specimens should be transported at room temperature for up to 1 hour. For long term transport (more than 1 hour), we recommend shipping at -20°C or lower. The samples can be stored frozen at -20°C or lower (at -80°C ideally). Repeated freeze-thaw cycles should be avoided to prevent sample and nucleic acids degradation.

The clinical specimens must be collected, transported, and stored according to appropriate laboratory guidelines. For more information, refer to the following guidelines:

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DNA extraction

For nucleic acid isolation from clinical samples, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. For the sample's pretreatment follow the instructions for use of the extraction kit used. The assay has been validated with the following extraction kits:

- MagDEA Dx SV Kit, using the magLEAD® 12gC instrument (Precision System Science Co.)
- Invisorb® Spin Universal Kit (Invitex)
- MagCore® Viral Nucleic Acid Extraction Kit, using the MagCore® HF16 Nucleic Acid extractor (RBC Biotech)

Positive control preparation

Reconstitute the lyophilized NG CiproRes Positive Control (red tube) with 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run for each assay.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of DNA sample, negative (yellow tube) or positive (red tube) controls and add them into the corresponding wells.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Polymerase activation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	63°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the FAM (*N. gonorrhoeae* ciprofloxacin resistant), ROX (*N. gonorrhoeae*), HEX, JOE, or VIC (*N. gonorrhoeae* ciprofloxacin sensitive), and Cy5 (Internal Control) channels. Depending on the equipment used select the proper detection channel (Attached II).

Analysis and interpretation of results

The results analysis is done by the software itself of the used real-time PCR system following manufacturer's instructions. To verify the correct operation of the amplification mix, check the internal control (IC) signal emission.

Use the positive control amplification curve as a starting point during reaction validation (prior to interpretation of sample results), to ensure that the threshold falls within the exponential phase of the amplification curves and above any background noise signals. The threshold value may vary between different instruments due to different signal intensities. It is recommended to set the threshold values for each channel (target) independently by the end user.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive control used in each run must show an amplification curve ($Ct \leq 40$) in FAM, ROX, and HEX, VIC, or JOE channels, which validates the reaction.

Negative control

The negative control included in each run must show signal' absence ($Ct > 40$ or no signal) in FAM, ROX, and HEX, VIC, or JOE, which validates the reaction.

The experiment seems to be failed if there is amplification signal in negative control or signal absence in the positive control for any channel. The assay should be repeated.

The Internal Control (IC) should show an amplification signal ($Ct \leq 40$) in positive control and negative control wells.

The clinical samples test results assessment should be performed once controls' results have been validated. The result interpretation is summarized in the following table:

NG CiproRes				Interpretation
CIP resistance (FAM)	CIP sensitivity (HEX)	<i>N. gonorrhoeae</i> (ROX)	Internal Control (Cy5)	
-	-	+	+/- ¹	<i>N. gonorrhoeae</i> DNA detected ³
+	-	+	+/- ¹	<i>N. gonorrhoeae</i> CIP resistant DNA detected
-	+	+	+/- ¹	<i>N. gonorrhoeae</i> CIP sensitive DNA detected
+	+	+	+/- ¹	<i>N. gonorrhoeae</i> CIP resistant and CIP sensitive DNA detected
-	-	-	+ ²	Targets not detected ²
-	-	-	- ²	Test failure ²

(+) Positive: Amplification signal (Ct ≤40)

(-) Negative: No amplification signal (Ct >40 or no signal)

¹ Sometimes, the Internal Control (IC) detection is not necessary since a high copies number of the target can cause a preferential amplification of target-specific nucleic acids. The Internal Control (IC) shows or not an amplification signal (Ct ≤40 or no signal).

² In the case of negative pathogens target genes detection, IC must show an amplification signal with Ct ≤ 35. If there is a signal' absence or Ct value >35 of Internal Control, the result is considered 'Invalid', and retesting is required. It is recommended to repeat the qPCR by diluting the DNA sample 1:10 and/or 1:100, or to re-extract and retest to check for possible failure in the extraction procedure and/or inhibition problems.

³ Low concentration of *N. gonorrhoeae* (Ct>35) may prevent detection of CIP resistance/sensitivity gene.

If the obtained result is confusing or doubtful, it is necessary to check that all the steps have been carried out correctly, to verify the correct performance of each qPCR steps, to review all the parameters, the sigmoid shape of the curve and the fluorescence intensity. It is also recommended to repeat the assay, preferably in duplicate, depending on the available material (repeat qPCR with the same isolated DNA sample, or re-extract and retest another aliquot of the same specimen or, obtain a new specimen and retest).

The test results must be evaluated by a health professional, together with the medical history, clinical symptoms and/or the results obtained in other diagnostic tests.

Quality Control

To confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Moreover, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

The Vitassay qPCR NG CiproRes clinical performance was assessed in a study with 124 leftovers of swabs and first urine (urethral urine) from patients of both sexes with an initial diagnosis of *N. gonorrhoeae* infection using laboratory routine methods. All samples were initially identified as *N. gonorrhoeae* and they were analysed with Vitassay qPCR kit and another molecular kit and/or sequencing to assess the ciprofloxacin resistance. Vitassay qPCR NG CiproRes reported 122 samples true positives for *N. gonorrhoeae*, 61 true positives for *N. gonorrhoeae* CIP resistant (point mutation detected in *gyrA* gene) and 55 true positives for *N. gonorrhoeae* CIP sensitive (wild type *gyrA* gene). Finally, 10 incongruent results were obtained with Vitassay kit: 2 false negatives for *N. gonorrhoeae*, 7 false negatives for the detection of *gyrA* mutation and one false negative for wild type *gyrA* gene detection. No false positive values were found. Statistical analysis was performed to obtain the sensitivity and specificity values indicated below (CI=95%):

Target	TP	TN	FP	FN	SE	SP
<i>N. gonorrhoeae</i>	122	0	0	2	0.98(0.94-0.99)	NA
<i>Point mutation gyrA gene</i>	61	56	0	7	0.89 (0.79-0.95)	1 (0.93-1)
WT <i>gyrA</i> gene	55	68	0	1	0.98 (0.90-1)	1 (0.94-1)

TP: true positive; TN: true negative; FP: false positive; FN: false negative; SE: sensitivity; SP: specificity.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of standards from the pathogens ranging from 10^7 to 10^1 copies/reaction. This assay has a detection limit of 4 DNA copies per reaction for *N. gonorrhoeae* ciprofloxacin resistance coding gene, 0.02 CFU per reaction for *N. gonorrhoeae*, and 1.28 CFU per reaction for *N. gonorrhoeae* ciprofloxacin sensitivity coding gene (positive rate $\geq 95\%$), on vaginal swab clinical samples.

Analytical specificity

The analytical specificity for *Neisseria gonorrhoeae* detection was tested within the panel of following microorganisms, where no cross-reactivity was observed with any of the species:

Cross-reactivity testing					
<i>Mycoplasma genitalium</i>	-	<i>Enterococcus faecalis</i>	-	<i>Acinetobacter baumannii</i>	-
<i>Mycoplasma hominis</i>	-	<i>Enterococcus faecium</i> serotype 11	-	<i>Aspergillus fumigatus</i>	-
<i>Bacteroides fragilis</i>	-	<i>Escherichia coli</i> strain 0.1285; O18:H7:K1	-	<i>Neisseria cinerea</i>	-
<i>Candida albicans</i>	-	<i>Gardnerella vaginalis</i>	-	<i>Neisseria meningitidis</i> serogroup A	-
<i>Candida dubliniensis</i>	-	<i>Haemophilus ducreyi</i> strain Class 1	-	<i>Proteus mirabilis</i>	-
<i>Candida glabrata</i>	-	<i>Haemophilus influenzae</i> strain Minn A	-	<i>Pseudomonas aeruginosa</i>	-
<i>Candida krusei</i> (<i>Issatchenkia orientalis</i>)	-	Hepatitis A Virus	-	<i>Serratia marcescens</i> subsp. <i>marcescens</i>	-
<i>Candida parapsilosis</i>	-	Herpes Simplex Virus 1 strain MacIntyre	-	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-
<i>Candida tropicalis</i>	-	Herpes Simplex Virus 2 strain MS	-	<i>Stenotrophomonas maltophilia</i>	-
<i>Chlamydia trachomatis</i> strain SW	-	Human Papillomavirus types 16 and 18	-	<i>Streptococcus agalactiae</i> strain Z019	-
<i>Chlamydia trachomatis</i> genovar F	-	<i>Klebsiella oxytoca</i>	-	<i>Streptococcus pneumoniae</i> strain Z022	-
<i>Chlamydia trachomatis</i> serovar E, I and K	-	<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i> serotype Capsular 2	-	<i>Treponema pallidum</i>	-
<i>Chlamydia trachomatis</i> strain LGV	-	<i>Listeria innocua</i> serotype 6a	-	<i>Trichomonas vaginalis</i>	-
Cytomegalovirus strain AD-169	-	<i>Listeria ivanovii</i> subsp. <i>ivanovii</i> serovar 5	-	<i>Ureaplasma parvum</i>	-
<i>Enterobacter aerogenes</i> serotype Cloaca B	-	<i>Listeria monocytogenes</i> serovar 4b	-	<i>Ureaplasma urealyticum</i>	-
<i>Enterobacter cloacae</i> serotype Cloaca A	-				

Analytical reactivity

The Vitassay qPCR NG CiproRes kit's reactivity for *N. gonorrhoeae* ciprofloxacin sensitive was evaluated against DNA extracted from *Neisseria gonorrhoeae* strain NCTC 8375 [B 5025], *Neisseria gonorrhoeae* strain 49226 and *Neisseria gonorrhoeae* strain Lvl Ng PorA as templates, showing positive results.

The Vitassay qPCR NG CiproRes kit's reactivity for *N. gonorrhoeae* ciprofloxacin resistant was evaluated against DNA extracted from clinical samples positive for *N. gonorrhoeae* ciprofloxacin resistant as templates, showing positive results.

Compatibles real-time PCR equipment

Vitassay qPCR NG CiproRes has been validated on the following equipment:

- 7500 Fast Real-Time PCR System (Applied Biosystems) ^I
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- CFX Opus 96 Touch™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTPPrime Real Time Detection Thermal Cycler (DNA-Technology)
- DTlite Real-Time PCR System (DNA-Technology)
- Rotor-Gene® Q (Qiagen) ^{II}
- Cobas z480 Analyzer (Roche)
- LightCycler® 480 Instrument II (Roche)
- MIC qPCR Cycler (Bio Molecular Systems) ^{II}
- NEOS-96 Real Time PCR System (Linear)

^I: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

^{II}: For Rotor-Gene® Q and MIC qPCR Cycler the product should be reconstituted following the procedure and transferred into specific tubes of each equipment.

Limitations

- All obtained results must be interpreted by a specialist in conjunction with available clinical information and laboratory findings.
- This assay has been validated with DNA extracted from urethral urine, urethral swabs, endocervical swabs and rectal swab samples already characterized as *N. gonorrhoeae* positive (using molecular assays). The use of other samples has not been established.
- The correct test performance depends on the sample's quality; proper DNA from clinical specimens must be extracted.

- This test does not provide quantitative values nor indicate the number of organisms present. This test is a qualitative assay.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination either by NG CiproRes Positive Control during its reconstitution, by samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.
- Detection may be affected by several factors and their combinations which may lead to false negative results, including a) inadequate specimen sampling, shipping, storage, handling; b) procedural errors (including DNA isolation); c) DNA degradation during specimen shipping, storage, and/or preparation; d) potential mutations of the target regions of the *Neisseria gonorrhoeae* CIP resistant or sensitive covered by this test which may result in DNA being undetectable e) pathogen load below the limit of detection for the assay; f) the presence of Real-Time amplification inhibitors or other types of interference (an interference study evaluating the effect of vaccines, antiviral therapeutics, antibiotics, chemotherapeutics, or immunosuppressant drugs used to prevent the infection or used during the treatment of the infection was not performed); g) failure to follow the manufacturer's instructions and procedures.
- Detection of DNA may not indicate the presence of viable and/or infectious bacteria or that *Neisseria gonorrhoeae* CIP resistant or sensitive is the causative agent for clinical symptoms.
- Negative results do not preclude *Neisseria gonorrhoeae* infection and should not be the sole basis of a patient treatment/management decision. Optimal specimen types and/or stage of infection most suitable for its collection have not been established yet. Consider the collection of multiple specimens from the same patient at different time points, which may increase the probability of detecting the bacteria.
- Fluorescence values may vary due to multiple factors such as PCR equipment, extraction system, sample type and its previous treatment, among others.
- If the patient clinical data, laboratory tests and epidemiological studies suggest possible *N. gonorrhoeae* CIP resistant or sensitive infection, and other sexually transmitted diseases have been discarded, a false negative result might not be discarded, and additional tests should be performed.

Attached I. Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 ⁽¹⁾
AriaDx Real-Time PCR System	Applied Biosystems
Applied Biosystems	7300 Real-Time PCR System ^{(3) (1)}
7500 Fast Real-Time PCR System ^{(1) (2)}	7500 Real-Time PCR System ⁽¹⁾
7500 Fast Dx Real-Time PCR System ^{(1) (2)}	7900 Real-Time PCR System ⁽³⁾
QuantStudio™ 12K Flex 96-well Fast	ABI PRISM 7000 ⁽³⁾
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7700 ⁽³⁾
QuantStudio™ 7 Flex 96-well Fast	QuantStudio™ 12K Flex 96-well
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
QuantStudio™ 5 Fast Real-Time PCR System	QuantStudio™ 7 Flex 96-well
QuantStudio™ 3 Fast Real-Time PCR System ⁽³⁾	QuantStudio™ 5 Fast Real-Time PCR System
StepOne Plus™ Real-Time PCR System ⁽³⁾	QuantStudio™ 5 Real-Time PCR System
StepOne™ Real-Time PCR System ⁽⁴⁾	QuantStudio™ 3 Real-Time PCR System ⁽³⁾
ViiA™ 7 Fast Real-Time PCR System	ViiA™ 7 Real-Time PCR System
Azure Biosystems	Agilent Technologies
Azure Cielo 3 ⁽⁵⁾	Mx3000P™ Real Time PCR System
Azure Cielo 6	Mx3005P™ Real Time PCR System
BIONEER	Analytik Jena
Exicycler™ 96	qTOWER
Bio-Rad	BIONEER
CFX96™ Real-Time PCR Detection System	Exicycler™ 96
Mini Opticon™ Real-Time PCR Detection System ⁽⁵⁾	BIOER
Roche	QuantGene 9600
LightCycler @480 Real-Time PCR System ^{(1) (6)}	Bio-Rad
LightCycler @96 Real-Time PCR System ⁽¹⁾	CFX96™ Deep Well Real-Time PCR
Cobas z480 Analyzer ^{(1) (6)}	iCycler iQ™ Real-Time PCR
	iCycler iQ™ 5 Real-Time PCR
	My iQ™ Real-Time PCR Detection System ⁽⁵⁾
	My iQ™ 2 Real-Time PCR Detection System ⁽⁵⁾
Special Formats ⁽⁷⁾	DNA-Technology
Bio Molecular Systems	DTlite Real-Time PCR System ⁽⁸⁾
Mic Real Time PCR Cyclers	DTprime Real-time Detection Thermal Cyclers ⁽⁸⁾
Cepheid	Eppendorf
SmartCycler®	Mastercycler™ ep <i>realplex</i>
Precision System Science Co., Ltd.	Qiagen
geneLEAD VIII System	QIAquant 96
Qiagen	
Rotor-Gene® Q	

- (1) A special grid is needed to fit these real-time PCR kits.
- (2) Select Ramp Speed "Standard".
- (3) No ROX caption.
- (4) No Cy5 caption.
- (5) Only FAM and HEX caption.
- (6) Specific compensation color is required.
- (7) The product must be reconstituted following the appropriate procedure (see Test procedure) and transferred to the specific tubes for Mic, SmartCycler®, Rotor-Gene® Q or geneLEAD VIII System.
- (8) See Attached III to configure exposure settings.

Attached II. Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR thermocyclers are specified in the following table:

Thermocycler	Vitassay Channel	Detection Channel	Observations
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Applied Biosystems ABI 7500	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Roche Cobas z 480	FAM	465/510	Colour Compensation required
	HEX	540/580	
	ROX	540/610	
	Cy5	610/670	
Cepheid Smartcycler®	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Agilent Technologies Mx3000P™ Mx 3005P™	FAM	FAM	Passive reference option ROX is none
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Bio Molecular Systems Mic Real Time PCR Cycler	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10).
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Agilent Technologies AriaMx	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Qiagen Rotor-Gene® Q	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
BIONEER Exicycler™ 96	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Attached III. Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
DTIite Real-Time PCR System (DNA-Technology)	FAM	500
	HEX	500
	ROX	500
	Cy5	500
DTprime Real-time Detection Thermal Cycler (DNA-Technology)	FAM	500*
	HEX	1000
	ROX	1000
	Cy5	1000

*If result in FAM channel is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.


Bibliography/Bibliografía

- Buder S., Schöfer H., Meyer T., Bremer V., Kohl P. K., Skaletz-Rorowski A., & Brockmeyer N. (2019). Bacterial sexually transmitted infections. *Journal der Deutschen Dermatologischen Gesellschaft = Journal of the German Society of Dermatology : JDDG*, 17(3), 287–315.
- LiverTox: Clinical and Research Information on Drug-Induced Liver Injury. (2012). National Institute of Diabetes and Digestive and Kidney Diseases. Available from: <https://www.ncbi.nlm.nih.gov/books/NBK548066/>
- Młynarczyk-Bonikowska, B., Majewska, A., Malejczyk, M., Młynarczyk, G., & Majewski, S. (2020). Multiresistant *Neisseria gonorrhoeae*: a new threat in second decade of the XXI century. *Medical microbiology and immunology*, 209(2), 95–108.
- Quillin, S. J., & Seifert, H. S. (2018). *Neisseria gonorrhoeae* host adaptation and pathogenesis. *Nature reviews. Microbiology*, 16(4), 226–240.
- Suay-García, B., & Pérez-Gracia, M. T. (2018). Future Prospects for *Neisseria gonorrhoeae* Treatment. *Antibiotics (Basel, Switzerland)*, 7(2), 49.
- Workowski, K. A., Bolan, G. A., & Centers for Disease Control and Prevention (2015). Sexually transmitted diseases treatment guidelines, 2015. *MMWR. Recommendations and reports : Morbidity and mortality weekly report. Recommendations and reports*, 64(RR-03), 1–137.

Trademarks

All trademarks that may appear in this package insert are the property of their respective owners.

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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