Real-Time PCR Kits

Vitassay qPCR Gram-negative Nosocomial Agents II

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- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTprime Real-time Detection Thermal Cycler (DNA-Technology)
- 7500 Fast Real-Time PCR System (Applied Biosystems)
- DTLite Real-Time Detection Thermal Cycler (DNA-Technology)
- Cobas LightCycler 480 System (Roche)
- Rotor-Gene® Q (Qiagen)
- CFX Opus 96 Touch™ Real-Time PCR System (Bio-Rad)
- NEOS-96 Real-Time PCR System (Linear)

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Termocicladores con bloque de bajo perfil	Termocicladores con bloque de alto perfil
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 ⁽¹⁾
AriaDx Real-Time PCR System	Agilent Technologies
Applied Biosystems	Mx3000P™ Real Time PCR System
7500 Fast Real-Time PCR System ^{(1) (2)}	Mx3005P™ Real Time PCR System
7500 Fast Dx Real-Time PCR System ^{(1) (2)}	Applied Biosystems
QuantStudio™ 12K Flex 96-well Fast	7300 Real-Time PCR System ^{(1) (4)}
QuantStudio™ 6 Flex 96-well Fast	7500 Real-Time PCR System ⁽¹⁾
QuantStudio™ 7 Flex 96-well Fast	7900 HT Real-Time PCR System ⁽³⁾
QuantStudio™ 5 Real-Time PCR System	ABI PRISM 7000 ⁽³⁾
QuantStudio™ 5 Fast Real-Time PCR System	ABI PRISM 7700 ⁽³⁾
QuantStudio™ 3 Fast Real-Time PCR System	QuantStudio™ 12K Flex 96-well
StepOne Plus™ Real-Time PCR System ⁽³⁾	QuantStudio™ 6 Flex 96-well
StepOne™ Real-Time PCR System ^{(3) (4)}	QuantStudio™ 7 Flex 96-well
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 5 Fast Real-Time PCR System
Azure Biosystems	QuantStudio™ 5 Real-Time PCR System
Azure Cielo 3 ⁽⁵⁾	QuantStudio™ 3 Real-Time PCR System
Azure Cielo 6	ViiA™ 7 Real-Time PCR System
BIONEER	Analytik Jena
Exicycler™ 96 Fast	qTOWER ⁽⁶⁾
Bio-Rad	BIONEER
CFX96™ Real-Time PCR Detection System	Exicycler™ 96
CFX96™ IVD Real-Time PCR Detection System	BIOER
Mini Opticon™ Real-Time PCR Detection System ⁽⁵⁾	QuantGene 9600
CFX Opus 96	Bio-Rad
Roche	CFX96™ Deep Well Real-Time PCR System
LightCycler®480 Real-Time PCR System ^{(1) (6)}	CFX96™ Deep Well IDV Real-Time PCR System
LightCycler®96 Real-Time PCR System	iCycler iQ™ Real-Time PCR Detection System
Cobas z480 Analyzer ^{(1) (6)}	iCycler iQ™5 Real-Time PCR Detection System
	My iQ™ Real-Time PCR Detection System ⁽⁵⁾
	My iQ™ 2 Real-Time PCR Detection System ⁽⁵⁾
Formatos especiales ⁽⁷⁾	DNA-Technology
Bio Molecular Systems	DTlite Real-Time PCR System ⁽⁸⁾
Mic Real Time PCR Cyclor	DTprime Real-time Detection Thermal Cyclor ⁽⁸⁾
Cepheid	Eppendorf
SmartCycler®	Mastercycler™ ep <i>realplex</i>
Qiagen	Qiagen
Rotor-Gene® Q	QIAquant 96 ⁽⁶⁾

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Intended use

Vitassay qPCR Gram-negative Nosocomial Agents II allows the qualitative detection and identification of DNA from *Pseudomonas aeruginosa*, *Klebsiella pneumoniae* and/or *Proteus mirabilis* by Real-Time PCR in blood culture and swab samples from patients with suspect of bacterial infection and/or multi-resistant infection or colonization and swabs specimens for epidemiological control by their healthcare professional. This product is intended to aid in the *Pseudomonas aeruginosa*, *Klebsiella pneumoniae* and/or *Proteus mirabilis* infections diagnosis, alongside the patient's clinical data, epidemiological risk factors, and other laboratory tests outcomes.

References

Vitassay qPCR Gram-negative Nosocomial Agents II 4x8-well strip, low profile	7041074
Vitassay qPCR Gram-negative Nosocomial Agents II 4x8-well strip, high profile	7042074

Reagents provided

In references 7041074 and 7042074:

Reference	Reagent/Material	Colour	Amount
7041S074/ 7042S074	Gram-negative II strips low/high profile	-	4 x 8-well strip
7C074	Gram-negative II Positive Control	red	1 vial
7002A	PCR grade water	white	1 vial x 1 mL
7007B	Resuspension buffer	green	1 vial x 1,8 mL
7008N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4 x 8-cap strip

Material and equipment required, not provided

- Collection and transport system.
- Laboratory freezers (- 30°C to - 10°C and/or ≤ -70°C).
- DNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Vortex
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated on the label.
- The resuspended positive control should be stored at -20°C. To avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

Summary

Pseudomonas aeruginosa, *Klebsiella pneumoniae* and *Proteus mirabilis* are gram-negative nosocomial multidrug-resistant pathogens. Their infections mainly affect patients in intensive care units and immunocompromised patients. The adaptation to the hospital environment together with their ability to easily acquire multiple genetic mobile elements harboring resistance and virulence genes are generating multidrug-resistant (MDR) strains of these pathogens, which is causing a major health problem worldwide.

Pseudomonas aeruginosa

Pseudomonas aeruginosa is a common, gram-negative bacterium found in the environment. *P. aeruginosa* may be a significant pathogenic factor of severe and often opportunistic infections in humans. Due to a huge arsenal of virulence factors that contribute to pathogenesis and antibiotic resistance determinants included in *P. aeruginosa*'s genome, these bacteria are capable of causing both acute and chronic infections. Some of its virulence factors are the outer membrane which has an asymmetric bilayer that limits the entry of harmful compounds, the biofilms formed by the bacterium which are highly resistant to antibiotics, disinfectants, and host defenses, the type IV pili which are essential structures for the initiation of the infection by mediating motility and adhesion, among others. But the most toxic *P. aeruginosa* virulence factor is the exotoxin A, which one inside is responsible for necrosing activity at the site of colonization.

It is often isolated from plants, fruits, soil, and water environments, such as rivers, lakes, and swimming pools. It infects airways and urinary tracts, causes blood infections, and is the most common cause of burn injury infections, hot-tub dermatitis, and outer ear infections. Because of its natural and acquired antibiotic resistance, treating *P. aeruginosa* infections can be problematic.

It is a major cause of healthcare-associated infections, being particularly problematic in intensive care units, being the most frequent colonizer of medical devices. Its infections are associated with high morbidity and mortality in many groups, including individuals with healthcare-associated pneumonia, chronic obstructive pulmonary disease, or cystic fibrosis, in which stage the pathogen is practically impossible to eradicate. It is included in the "critical" category of the WHO priority list of bacterial pathogens for which research and development of new antibiotics are urgently needed.

Klebsiella pneumoniae

Klebsiella pneumoniae is a gram-negative bacterium which belongs to the family *Enterobacteriaceae*, which includes the well-known genre *Salmonella* and *Escherichia*.

This bacterium can be located on the surface of mucosa in animals, but in humans is mostly found in the gastrointestinal tract, and a few in the nasopharynx, where the bacteria can enter the bloodstream or other tissues, and cause infection. Environmental niches include soil, water, a range of plant species, insects, birds, reptiles, and many different mammals in which this bacterium can be either a commensal organism or a potential pathogen.

Most of *K. pneumoniae* infections are opportunistic health-care-associated infections, sometimes referred to as "classical" *K. pneumoniae* infections. The most common manifestations are pneumonia, urinary tract, and wound infections, any of which can progress to bacteraemia. Vulnerable patient groups such as neonates, the elderly, those with inserted medical devices and the immunocompromised, are the most at danger. *K. pneumoniae* can also cause endophthalmitis, pneumonia, necrotizing fasciitis, non-hepatic abscess, meningitis, and pyogenic liver abscess in the absence of biliary tract disease, when it is outside the hospital. Some of its virulence factors are acquired siderophores, molecules that can competitively scavenge iron from host proteins or other sources, and surface receptors for internalization, as well as the *ybt* locus, which is found at 6-80% and 78-100% frequency in the worldwide MDR and hypervirulent clones, respectively, a genotoxic polyketide colibactin, which induces DNA damage in eukaryotic cells, among others.

The WHO recognizes extended-spectrum β -lactam producing and carbapenem-resistant *K. pneumoniae* as a critical public health threat. This pathogen, as well as *Pseudomonas aeruginosa*, are members of ESKAPE pathogens, which are described as the leading cause of resistant nosocomial infections. Both have inherent resistance to one or more antibiotics, and individual strains have accumulated resistance to a wide range of drugs.

Proteus mirabilis

Proteus mirabilis, a gram-negative bacterium, which belongs to the class *Gammaproteobacteria*. It has long been recognized as a member of the older *Enterobacteriales*, family *Eterobacteriaceae*, however in recent years it has been suggested that the genre *Proteus* is going to be placed within a new *Morganellaceae* family.

The bacterium can be found in a wide variety of environments, including soil, water sources, and sewage, but it is predominantly a commensal of the gastrointestinal tracts of humans (along with *Klebsiella species*, and *Escherichia coli*) and animals.

Proteus mirabilis is capable of causing a range of human infections, including wounds, and infection in eyes, gastrointestinal and urinary tracts. Urinary tract infections and catheter-associated urinary tract infections are sometimes aggravated by the production of bladder and kidney stones (urolithiasis), as well as persistent renal damage and can lead to bacteremia and sepsis.

The urinary tract infections are more common in people between the ages of 20 and 50 and most common in women. It is a common cause of infection in male infants of all ages, as well as elderly men and women and can lead to serious problems such as recurrent infection, pyelonephritis with sepsis, infant kidney injury, premature delivery, and antibiotic-related complications.

The mobility of *P. mirabilis* is due to its flagella, which not only aids colonization but it has also been linked to its ability to build biofilms and is thought to contribute to resistance to host defenses and certain antibiotics. Also, *P. mirabilis* produces two toxins, hemolysin

(HpmA) and *Proteus* toxigenin (Pta), both of which have been linked to tissue destruction and renal transmission, leading to acute pyelonephritis.

Test principle

Vitassay qPCR Gram-negative Nosocomial Agents II is based on real-time amplification of a conserved region of the *oprL* gene for *Pseudomonas aeruginosa*, *Hemolysin* gene for *Klebsiella pneumoniae* and *hns* gene for *Proteus mirabilis*. After DNA extraction, the presence of the viruses is detected by an increase in fluorescence observed during the reaction, following hydrolysis of the fluorescent probe.

The assay is based on 5' exonuclease activity using two primers and a fluorogenic hydrolysis probe to detect the accumulation of the amplified target sequence during the PCR reaction. When the polymerase begins to extend the primers, the probe is hydrolysed by its 5'-3' exonuclease activity resulting in spatial separation of the fluorophore and the quencher. The resulting increase in fluorescent signal is proportional to the amount of amplified product in the sample and is detected by real-time PCR equipment.

Vitassay qPCR Gram-negative Nosocomial Agents II is a ready-to-use test that contains in each well all the necessary reagents in a stabilised format to perform real-time PCR. In addition, an internal control allows the detection of a possible inhibition reaction. The wells contain the multiplex reaction mixture for the detection of *Pseudomonas aeruginosa*, *Klebsiella pneumoniae* and *Proteus mirabilis* genes as well as the Internal Control (IC). The detection channels of the target sequences are described below:

Target	Detection Channel
<i>Hns</i> gene	Cy5
<i>OprL</i> gene	FAM
<i>Hemolysin</i> gene	ROX
Internal Control	HEX, VIC, or JOE (depending on the equipment used)

Precautions

- For professional *in vitro* diagnostic use (use by qualified and trained clinical laboratory personnel).
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or is broken or if the foil has been broken or damaged.
- Do not remove desiccant from reagent pouches once is open.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles, and mask.
- Do not eat, drink, or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment, and reagents to the area in which the previous step was performed. Use separate areas for the patient samples and controls preparation.
- Avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile disposable aerosol resistant pipette tips or RNase/DNase positive displacement pipette tips is recommended.
- An appearance of the reaction mixture in stabilized format, normally found at the tube bottom, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or color different from whitish) does not alter the test functionality.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious and they must be managed according to the national safety legislation and national health waste legislation. Take necessary precautions during the collection, transport, storage, treatment, and disposal of samples.
- Use personal protective equipment (PPE) and biological safety cabinet for handling of potentially infectious samples and reagents according to current recommendations.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

Procedures

Specimen collection, transport, and conservation

For specimen collection, conservation, and transport, user-validated conditions must be followed. Vitassay qPCR Gram-negative Nosocomial Agents II has been tested on blood culture and swab samples such as biopsy, gluteal swabs, ulcer swabs, wound swabs, tongue swab, joint fluid swab, oral swabs, hair swabs, perineal swab, urethral swab, abscess swabs and ostomy swab from patients with suspect of bacterial infection and/or multi-resistant infection or colonization and swabs specimens for epidemiological control such as nasal, pharyngeal, rectal swabs, axillary swabs, inguinal swabs, epithelial swabs, colostomy swabs, one stoma swab, BAS, BAL and sputum specimens. Other sample types should be validated by the user.

Overall, clinical samples should be collected appropriately in clean containers with or without transport media (depending on sample type), correctly labelled and promptly processed to ensure the test's quality. Use of fresh samples for the test is also recommended. Transport must always be carried out in accordance with local and national regulations for the transport of biological samples. Specimens' transportation should be done either at room temperature for up to 2 hours. For long term transport (more than 2 hours), shipment at -80°C (or lower) is advisable. The samples can be stored frozen at -20°C or bellow (ideally -80°C) for conservation. Repeated freeze-thaw cycles should be avoided to prevent sample and nucleic acids degradation.

The clinical specimens must be collected, transported, and stored according to appropriate laboratory guidelines. For more information, refer to the following guidelines:

- IDSA (Miller, J. M., Binnicker, M. J., Campbell, S., ... & Pritt, B. S. (2018). A guide to utilization of the microbiology laboratory for diagnosis of infectious diseases: 2018 update by the Infectious Diseases Society of America and the American Society for Microbiology. *Clinical Infectious Diseases*, 67(6), e1-e94)
- García-Lechuz Moya JM, et al. (2017). Recogida, transporte y procesamiento general de las muestras en el Laboratorio de Microbiología. 2017. 1b. Procedimientos en Microbiología Clínica. Sociedad Española de Enfermedades Infecciosas y Microbiología Clínica (SEIMC).

DNA extraction

For nucleic acid isolation from clinical specimens, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. For the sample's pretreatment follow the instructions for use of the extraction kit used. The following extraction kits have been validated:

- MagDEA Dx SV kit, using the magLEAD® 12gC instrument (Precision System Science Co.).
- Invisorb® Spin Universal Kit (Invitex).

Positive control preparation

Reconstitute the lyophilized Gram-negative NA II Positive Control (red tube) with 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run for each assay.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of DNA sample, negative (yellow tube) or positive (red tube) controls and add them into the corresponding wells.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Polymerase activation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the FAM (*oprL* gene), Cy5 (*hns* gene), ROX (*Hemolysin* gene) and HEX, JOE, or VIC (Internal Control). Depending on the equipment used select the proper detection channel (Attached II).

Analysis and interpretation of results

The results analysis is done by the software itself of the used real-time PCR system following manufacturer's instructions. To verify the correct operation of the amplification mix, check the internal control (IC) signal emission.

Use the positive control amplification curve as a starting point during reaction validation (prior to interpretation of sample results), to ensure that the threshold falls within the exponential phase of the amplification curves and above any background noise signals. The threshold value may vary between different instruments due to different signal intensities. It is recommended to set the threshold values for each channel (target) independently by the end user.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive control used in each run must show an amplification curve (Ct ≤40) in FAM, ROX, and Cy5 channels, which validates the reaction.

Negative control

The negative control included in each run must show signal' absence (Ct >40 or no signal) in FAM, ROX, and Cy5, which validates the reaction.

The experiment seems to be failed if there is amplification signal in negative control or signal absence in the positive control for any channel. The assay should be repeated.

The Internal Control (IC) should show an amplification signal (Ct ≤40) in positive control and negative control wells.

The clinical samples test results assessment should be performed once controls' results have been validated. The result interpretation is summarized in the following table:

Vitassay qPCR Gram-negative Nosocomial Agents II					
<i>P. aeruginosa</i> (FAM) ³	<i>K. pneumoniae</i> (ROX) ⁴	<i>P. mirabilis</i> (Cy5) ⁴	Internal Control (HEX)	Interpretation	
+	+	-	+/- ¹	Valid	DNA of <i>P. aeruginosa</i> y <i>K. pneumoniae</i> detected
+	-	-	+/- ¹	Valid	DNA of <i>P. aeruginosa</i> detected
-	+	-	+/- ¹	Valid	DNA of <i>K. pneumoniae</i> detected
-	-	+	+/- ¹	Valid	DNA of <i>P. mirabilis</i> detected
+	-	+	+/- ¹	Valid	DNA of <i>P. aeruginosa</i> y <i>P. mirabilis</i> detected
+	+	+	+/- ¹	Valid	DNA of three pathogens detected
-	+	+	+/- ¹	Valid	DNA of <i>K. pneumoniae</i> y <i>P. mirabilis</i> detected
-	-	-	+ ²	Valid	Targets not detected
-	-	-	- ²	Invalid	Test failure

(+) Positive: Amplification signal

(-) Negative: No amplification signal

¹Sometimes, the Internal Control (IC) detection is not necessary since a high copies number of the target can cause a preferential amplification of target-specific nucleic acids. The Internal Control (IC) shows or not an amplification signal (Ct ≤40 or no signal).

²In the case of negative *Pseudomonas aeruginosa*, *Klebsiella pneumoniae* and/or *Proteus mirabilis* target genes detection, IC must show an amplification signal with Ct ≤35. If there is a signal' absence or Ct value > 35 of Internal Control, the result is considered 'Invalid', and retesting is required. It is recommended to repeat the RT-qPCR by diluting the DNA sample 1:10 and/or 1:100, or to re-extract and retest to check for possible failure in the extraction procedure and/or inhibition problems.

³For *oprL* gene from *P. aeruginosa* a cut-off of Ct 37 has been set. Therefore:

(+) Positive: Amplification signal (Ct ≤37)

(-) Negative: No amplification signal (Ct >37 or no signal)

⁴For *Hemolysin* gene from *K. pneumoniae* y *hns* gene from *P. mirabilis* a cut-off of Ct 35 has been set. Therefore:

(+) Positive: Amplification signal (Ct ≤35)

(-) Negative: No amplification signal (Ct >35 or no signal)

If the obtained result is confusing or doubtful, it is necessary to check that all the steps have been carried out correctly, to verify the correct performance of each RT-qPCR steps and to review all the parameters, the sigmoid shape of the curve and the fluorescence intensity. It is also recommended to repeat the assay, preferably in duplicate, depending on the available material (repeat RT-qPCR with the same isolated DNA sample, or re-extract and retest another aliquot of the same specimen or, obtain a new specimen and retest).

The test results must be evaluated by a health professional, together with the medical history, clinical symptoms and/or the results obtained in other diagnostic tests.

Quality Control

To confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Moreover, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

The clinical performance of Vitassay qPCR Gram-negative Nosocomial Agents II was evaluated in two comparative-retrospective studies, in which a total of 417 blood culture samples and 620 swab samples were used in first and second study, respectively. The samples were obtained every 15 days between March 2021 and March 2022 and from March to May 2021, respectively. Sample extraction in both studies was performed with the MagDEA DX SV kit using the automatic MagLEAD® 12 gC instrument (Precision System Science Co.) and run on DTprime Real-Time Detection Thermal Cycler (DNA-Technology).

In the first study, the results obtained with Vitassay qPCR Gram-negative Nosocomial Agents II were compared with the results obtained in culture. Discordant results were analysed with Magicplex™ Sepsis ID6 Kit (Seegene), Magicplex™ Sepsis ID4 Kit (Seegene) and by sequencing, confirming as true positive and true negative the results obtained with Vitassay assay. In addition, several random negative samples of each pathogen were also tested with Magicplex™ Sepsis ID6 Kit (Seegene) and Magicplex™ Sepsis ID4 Kit (Seegene). In the second study, the results obtained with Vitassay qPCR Gram-negative Nosocomial Agents II were compared with those obtained in culture, and discordant results were confirmed by sequencing.

After Vitassay qPCR Gram-negative Nosocomial Agents II analysis, it can be concluded that it represents an efficient diagnostic tool for the detection of *Pseudomonas aeruginosa*, *Klebsiella pneumoniae* and/or *Proteus mirabilis*, as the results show high agreement with comparator assays. The results are shown in the table below, and the global sensitivity and specificity values, and the values for each pathogen, were calculated from the study data (CI=95%).

Comparator assay	Target	TP	TN	FP	FN	SE	SP	PPV	NPV	
1	Culture	<i>P. aeruginosa</i>	5	412	0	0	1(0.47-1)	1(0.99-1)	1(0.56-1)	1(0.99-1)
		<i>K. pneumoniae</i>	22	395	0	0	1(0.84-1)	1(0.99-1)	1(0.84-1)	1(0.99-1)
		<i>P. mirabilis</i>	20	397	0	0	1(0.83-1)	1(0.99-1)	1(0.83-1)	1(0.99-1)
	MagicPlex™ Sepsis ID4/ID6 Real-Time Detection (Seegene)	<i>P. aeruginosa</i>	5	18	0	0	1(0.47-1)	1(0.81-1)	1(0.47-1)	1(0.81-1)
		<i>K. pneumoniae</i>	22	21	0	0	1(0.84-1)	1(0.83-1)	1(0.84-1)	1(0.83-1)
		<i>P. mirabilis</i>	20	22	0	0	1(0.83-1)	1(0.84-1)	1(0.83-1)	1(0.84-1)
2	Culture	<i>P. aeruginosa</i>	132	486	0	2	0.98(0.94-0.99)	1(0.99-1)	1(0.96-1)	0.99(0.98-0.99)
		<i>K. pneumoniae</i>	122	498	0	0	1(0.97-1)	1(0.99-1)	1(0.97-1)	1(0.99-1)
		<i>P. mirabilis</i>	35	585	0	0	1(0.90-1)	1(0.99-1)	1(0.90-1)	1(0.99-1)

TP: True Positive; TN: True Negative; FP: False Positive; FN: False Negative; SE: Sensitivity; SP: Specificity; PPV: Positive Predictive Values; NPV: Negative Predictive Values.

Analytical sensitivity

The analytical sensitivity was determined by the analysis of 10-fold dilution series of standards from the pathogens ranging from 10^7 to 10^1 copies/reaction. This assay has a detection limit of 0.8 CFU per reaction for *P. aeruginosa*, 0.01 CFU for *K. pneumoniae* y de 0.04 CFU for *P. mirabilis* (positive rate of 95%).

Analytical specificity

The analytical specificity for *P. aeruginosa*, *K. pneumoniae* and *P. mirabilis* detection was tested within the panel of following microorganisms, where no cross-reactivity was observed between any of the species.

Cross-reactivity testing					
<i>Human Adenovirus type 1,2, 3, 4, 5, 8, 15, 31, 40 and 41</i>	-	Dengue virus type 1 strain Hawaii A	-	<i>Pseudomonas aeruginosa</i> ATCC®10145™	+/-
<i>Aeromonas caviae</i>	-	Dengue virus type 2 strain New Guinea C	-	<i>Rickettsia conorii</i> strain Moroccan	-
<i>Aeromonas hydrophila</i> subsp. <i>Hydrophila</i>	-	Dengue virus type 3 strain H87	-	Rift Valley Fever Virus AR21229	-
<i>Anaplasma marginale</i>	-	Dengue virus type 4 strain H241	-	Rift Valley Fever Virus MP12	-
<i>Arcobacter butzleri</i>	-	<i>Dientamoeba fragilis</i> / <i>Blastocystis hominis</i>	-	<i>Salmonella bongori</i> serovar 66:z41	-
Astrovirus genotype I-VIII	-	<i>Entamoeba dispar</i> / <i>Blastocystis hominis</i>	-	<i>Salmonella enterica</i> subsp. <i>Enterica</i> serovar <i>enteritidis</i>	-
<i>Bacteroides fragilis</i>	-	<i>Entamoeba histolytica</i>	-	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>gallinarum</i>	-
<i>Bartonella henselae</i> strain Houston-1	-	<i>Enterococcus faecalis</i>	-	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>paratyphi A</i>	-
<i>Borrelia azfeli</i> strain P-Ko/1984	-	Enterohemorrhagic <i>Escherichia coli</i> serotype O157:H7	-	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>paratyphi B</i>	-

<i>Borrelia bavariensis</i>	-	Enteroinvasive <i>Escherichia coli</i>	-	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>pullorum</i>	-
<i>Borrelia bisetti</i>	-	Enteropathogenic <i>Escherichia coli</i>	-	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>typhi</i>	-
<i>Borrelia burgdorferi sensu stricto</i>	-	Enterotoxigenic <i>Escherichia coli</i> serotype O25:H42	-	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>typhimurium</i> . Serotype 4,5,12:i:1,2	-
<i>Borrelia burgdorferi</i> strain IRS	-	Enterovirus 68 and 71	-	Sapovirus	-
<i>Borrelia garinii</i>	-	Enterovirus <i>Coxsackievirus</i> A24, A9 and B3	-	<i>Serratia liquefaciens</i>	-
<i>Cryptosporidium parvum/hominis</i>	-	Enterovirus <i>Echovirus</i> 30	-	<i>Shigella dysenteriae</i> serotype 1	-
<i>Borrelia hermsii</i>	-	<i>Giardia intestinalis</i> strain WB clone C6	-	<i>Shigella flexneri</i>	-
<i>Borrelia japonica</i>	-	<i>Helicobacter heilmannii</i>	-	St Louis Encephalitis Virus	-
<i>Borrelia lusitanae</i>	-	<i>Helicobacter hepaticus</i>	-	<i>Staphylococcus aureus</i> subsp. <i>Aureus</i>	-
<i>Borrelia miyamotoi</i>	-	<i>Helicobacter pylori</i> J99	-	<i>Streptococcus pneumoniae</i>	-
<i>Borrelia spielmanii</i>	-	Human Rotavirus A	-	<i>Theileria annulata</i>	-
<i>Borrelia valaisiana</i>	-	Japanese encephalitis	-	Tick-Borne encephalitis virus strain Neudorfl	-
<i>Campylobacter coli</i>	-	Japanese encephalitis virus strain Nakayama	-	<i>Treponema phagedenis</i>	-
<i>Campylobacter fetus</i> subsp. <i>Fetus</i>	-	<i>Klebsiella oxytoca</i>	-	<i>Trypanosoma cruzi</i>	-
<i>Campylobacter hyointestinalis</i>	-	<i>Klebsiella pneumoniae</i> subsp. <i>Pneumoniae</i> serotype capsular 2 (CECT141)	-/+	Usutu Virus	-
<i>Campylobacter jejuni</i> subsp. <i>Jejuni</i>	-	<i>Klebsiella pneumoniae</i> subsp. <i>Pneumoniae</i>	-/+	<i>Vibrio parahaemolyticus</i> serotype O1:K1	-
<i>Campylobacter lari</i>	-	<i>Klebsiella pneumoniae</i> (NDM-1 resistant) ATCC®BAA-2146™	+/-	West Nile virus strain Heja	-
<i>Campylobacter upsaliensis</i>	-	<i>Klebsiella pneumoniae</i> (OXA-48 resistant) (ATCC® BAA-2524™)	+/-	West Nile virus strain NY99	-
<i>Candida albicans</i>	-	<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i> ATCC®10031™	+/-	West Nile virus strain Ug 37	-
Chikungunya Virus F24	-	<i>Leptospira</i>	-	Yellow Fever Virus strain 17D	-
Chikungunya Virus Martinique	-	<i>Listeria monocytogenes</i> serovar ½c	-	Yellow Fever Virus strain French Neurotropic	-
Chikungunya virus S27 Petersfield	-	Norovirus GI	-	<i>Yersinia enterocolitica</i> O:3	-
Chikungunya Virus WHO IS (R91064)	-	Norovirus GII. 4RNA	-	<i>Yersinia enterocolitica</i> O:9	-
<i>Citrobacter freundii</i>	-	<i>Plasmodium falciparum</i> strain 3D7	-	Zika Virus African strain	-
<i>Clostridium difficile</i>	-	<i>Proteus mirabilis</i> ATCC®25933™	+/-	Zika Virus Asian strain PF13/251013-18	-
<i>Clostridium difficile</i> O:27	-	<i>Proteus mirabilis</i> (CECT 170)	+/-	Zika Virus French Polynesian strain 11468/16	-
<i>Clostridium perfringens</i>	-	<i>Proteus vulgaris</i>	-	Zika Virus French Polynesian strain 11474/16	-
<i>Coxiella burnetii</i> strain Nine Mile Q	-	<i>Pseudomonas aeruginosa</i> (CECT 108)	+/-	Zika virus strain FB-GWUH-2016	-

Analytical reactivity

Vitassay qPCR Gram-negative Nosocomial Agents II reactivity for *Pseudomonas aeruginosa* was tested against DNA extracted from *Pseudomonas aeruginosa* (ATCC® 10145™) and *Pseudomonas aeruginosa* (CECT 108) (as templates), showing positive results.

Vitassay qPCR Gram-negative Nosocomial Agents II reactivity for *Klebsiella pneumoniae* was assessed against DNA extracted from *Klebsiella pneumoniae* subsp. *pneumoniae* (ATCC® 10031™), *Klebsiella pneumoniae* subsp. *pneumoniae* Serotype capsular 2 (CECT 141), *Klebsiella pneumoniae* NDM-1 reference strain (ATCC® BAA-2146™) and *Klebsiella pneumoniae* subsp. *pneumoniae* carbapenemase producer (OXA-48) (ATCC® BAA-2524™) (as templates), showing positive results.

Vitassay qPCR Gram-negative Nosocomial Agents II reactivity for *Proteus mirabilis* was assessed against DNA extracted from *Proteus mirabilis* (ATCC® 25933™) and *Proteus mirabilis* (CECT 170) (as templates), showing positive results.

Compatibles real-time PCR equipment

Vitassay qPCR Gram-negative Nosocomial Agents II has been validated on the following equipment:

- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTprime Real-time Detection Thermal Cycler (DNA-Technology)
- 7500 Fast Real-Time PCR System (Applied Biosystems)
- DTlite Real-Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)
- Cobas LightCycler 480 System (Roche)
- CFX Opus 96 Touch™ Real-Time PCR System (Bio-Rad)
- NEOS-96 Real-Time PCR System (Linear)

Limitations

- All obtained results should be interpreted by a specialist in conjunction with available clinical information and laboratory findings.
- This assay has been validated with DNA extracted from blood culture and samples from patients with suspect of bacterial infection and/or multi-resistant infection or colonization and swabs for epidemiological control such as biopsy, gluteal, ulcer, wound, tongue, joint fluid, oral, hair, perineal, urethral, abscess, ostomy, nasal, pharyngeal, rectal, axillary, inguinal epithelial, colostomy and stoma swabs; BAS, BAL and sputum. The use of other samples has not been established and must be validated by the user.
- The correct test performance depends on the sample's quality; DNA should be properly extracted from clinical samples.
- This test does not provide quantitative values nor indicate the number of organisms present. This test is a qualitative assay.
- Low copy numbers of the target DNA template may be detected below the detection limit, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination, either by Gram-negative II positive control during reconstitution, by samples containing high concentrations of target DNA template or by PCR products from previous reactions.
- Detection can be affected by several factors and their combinations that can lead to false negative results, including: (a) inadequate specimen sampling, shipping, storage and handling; (b) procedural errors (including DNA extraction); (c) DNA degradation during specimen shipping, storage and/or preparation; (d) pathogen load is below the detection limit for the assay; (e) presence of retrotranscription and/or real-time amplification inhibitors or other types of interference (an interference study evaluating the effect of antibiotics used to prevent infection or during treatment of infection was not performed); (f) mutations or polymorphisms in primer or probe binding regions may affect the detection of novel or unknown variants of *P. aeruginosa*, *K. pneumoniae* and/or *P. mirabilis*; (g) failure to follow the manufacturer's instructions and procedures.
- Detection of viral DNA may not indicate the presence of viable and/or infectious virus or that these viruses are the causative agents for clinical symptoms.
- Negative results do not preclude *P. aeruginosa*, *K. pneumoniae* and/or *P. mirabilis* infection and should not be the sole basis of a patient treatment/management decision. Optimal specimen types and/or stage of infection most suitable for its collection have not been established yet. Consider the collection of multiple specimens from the same patient at different time points, which may increase the probability of detecting the virus.
- If the patient clinical data, laboratory tests and epidemiological studies suggest possible *P. aeruginosa*, *K. pneumoniae* and/or *P. mirabilis* infection, and other gastrointestinal illnesses have been discarded, a false negative result might not be discarded, and additional tests should be performed.
- Fluorescence values may vary due to multiple factors such as PCR equipment, extraction system, sample type and its previous treatment, among others.

Attached I. Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following guidance table separated by tube type. Please refer to the table and check the equipment and its specifications before running the assay. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 ⁽¹⁾
AriaDx Real-Time PCR System	Agilent Technologies
Applied Biosystems	Mx3000P™ Real Time PCR System
7500 Fast Real-Time PCR System ^{(1) (2)}	Mx3005P™ Real Time PCR System
7500 Fast Dx Real-Time PCR System ^{(1) (2)}	Applied Biosystems
QuantStudio™ 12K Flex 96-well Fast	7300 Real-Time PCR System ^{(1) (4)}
QuantStudio™ 6 Flex 96-well Fast	7500 Real-Time PCR System ⁽¹⁾
QuantStudio™ 7 Flex 96-well Fast	7900 HT Real-Time PCR System ⁽³⁾
QuantStudio™ 5 Real-Time PCR System	ABI PRISM 7000 ⁽³⁾
QuantStudio™ 5 Fast Real-Time PCR System	ABI PRISM 7700 ⁽³⁾
QuantStudio™ 3 Fast Real-Time PCR System	QuantStudio™ 12K Flex 96-well
StepOne Plus™ Real-Time PCR System ⁽³⁾	QuantStudio™ 6 Flex 96-well
StepOne™ Real-Time PCR System ^{(3) (4)}	QuantStudio™ 7 Flex 96-well
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 5 Fast Real-Time PCR System
Azure Biosystems	QuantStudio™ 5 Real-Time PCR System
Azure Cielo 3 ⁽⁵⁾	QuantStudio™ 3 Real-Time PCR System
Azure Cielo 6	ViiA™ 7 Real-Time PCR System
BIONEER	Analytik Jena
Exicycler™ 96 Fast	qTOWER ⁽⁶⁾
Bio-Rad	BIONEER
CFX96™ Real-Time PCR Detection System	Exicycler™ 96
CFX96™ IVD Real-Time PCR Detection System	BIOER
Mini Opticon™ Real-Time PCR Detection System ⁽⁵⁾	QuantGene 9600
CFX Opus 96	Bio-Rad
Roche	CFX96™ Deep Well Real-Time PCR System
LightCycler®480 Real-Time PCR System ^{(1) (6)}	CFX96™ Deep Well IDV Real-Time PCR System
LightCycler®96 Real-Time PCR System	iCycler iQ™ Real-Time PCR Detection System
Cobas z480 Analyzer ^{(1) (6)}	iCycler iQ™5 Real-Time PCR Detection System
	My iQ™ Real-Time PCR Detection System ⁽⁵⁾
Special Formats ⁽⁷⁾	My iQ™ 2 Real-Time PCR Detection System ⁽⁵⁾
Bio Molecular Systems	DNA-Technology
Mic Real Time PCR Cyclor	DTlite Real-Time PCR System ⁽⁸⁾
Cepheid	DTprime Real-time Detection Thermal Cyclor ⁽⁸⁾
SmartCycler®	Eppendorf
Qiagen	Mastercycler™ ep <i>realplex</i>
Rotor-Gene® Q	Qiagen
	QIAquant 96 ⁽⁶⁾

(1) A special grid is needed to fit these real-time PCR kits.

(2) Select Ramp Speed "Standard" in New Experiment/Advanced Set-up/Experiment Properties.

(3) No Cy5 caption.

(4) No ROX caption.

(5) Only FAM and HEX caption.

(6) Specific compensation color is required.

(7) The product must be reconstituted following the appropriate procedure (see Test procedure) and transferred to the specific tubes for Mic, SmartCycler®, or Rotor-Gene® Q.

(8) See Attached III to configure exposure settings.

Attached II. Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real-Time PCR thermocyclers are specified in the following table:

Thermocycler	Vitassay Channel	Detection Channel	Observations
Bio-Rad CFX96™	FAM	FAM	During the first few cycles of an analysis, some wells may exhibit abnormal RFU values and show a non-sigmoid rising line. To correct this effect, select the Apply Fluorescence Drift Correction option in the Settings menu for the baseline configuration.
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Applied Biosystems ABI 7500	FAM	FAM	Passive reference option ROX is none. During the first few cycles of an assay, some wells may exhibit abnormal RFU values and show a non-sigmoidal rising line. To correct for this effect, select the Start cycle and End cycle values so that the baseline ends before the detection of significant fluorescence.
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Color Compensation required.
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Roche Cobas z 480	FAM	465/510	Color Compensation required.
	HEX	540/580	
	ROX	540/610	
	Cy5	610/670	
Cepheid Smartcycler®	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Agilent Technologies Mx3000P™ Mx3005P™	FAM	FAM	Passive reference option ROX is none.
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Bio Molecular Systems Mic Real Time PCR Cycler	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10).
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Agilent Technologies AriaMx	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Qiagen Rotor-Gene® Q	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range must be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
BIONEER Exicycler™ 96	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Attached III. Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
DTlite Real-Time PCR System (DNA-Technology)	FAM	500
	HEX	500
	ROX	500
	Cy5	500
DTprime Real-time Detection Thermal Cycler (DNA-Technology)	FAM	500*
	HEX	1000
	ROX	1000
	Cy5	1000

*If result in FAM channel is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.



Bibliography/Bibliografía

- Armbruster, C. E., Mobley, H., & Pearson, M. M. (2018). Pathogenesis of *Proteus mirabilis* Infection. *EcoSal Plus*, 8(1), 10.1128/ecosalplus.ESP-0009-2017.
- Davin-Regli, A., Lavigne, J. P., & Pagès, J. M. (2019). *Enterobacter* spp.: Update on Taxonomy, Clinical Aspects, and Emerging Antimicrobial Resistance. *Clinical microbiology reviews*, 32(4), e00002-19.
- Jurado-Martín, I., Sainz-Mejías, M., & McClean, S. (2021). *Pseudomonas aeruginosa*: An Audacious Pathogen with an Adaptable Arsenal of Virulence Factors. *International journal of molecular sciences*, 22(6), 3128.
- Mielko, K. A., Jabłoński, S. J., Milczewska, J., Sands, D., Łukaszewicz, M., & Młynarz, P. (2019). Metabolomic studies of *Pseudomonas aeruginosa*. *World journal of microbiology & biotechnology*, 35(11), 178.
- Jamil, R. T., Foris, L. A., & Snowden, J. (2022). *Proteus Mirabilis* Infections. In *StatPearls*. StatPearls Publishing.
- Wyres, K. L., Lam, M., & Holt, K. E. (2020). Population genomics of *Klebsiella pneumoniae*. *Nature reviews. Microbiology*, 18(6), 344–359
- Wang, G., Zhao, G., Chao, X., Xie, L., & Wang, H. (2020). The Characteristic of Virulence, Biofilm and Antibiotic Resistance of *Klebsiella pneumoniae*. *International journal of environmental research and public health*, 17(17), 6278.
- Wyres, K. L., & Holt, K. E. (2018). *Klebsiella pneumoniae* as a key trafficker of drug resistance genes from environmental to clinically important bacteria. *Current opinion in microbiology*, 45, 131–139.
- Yuan, F., Huang, Z., Yang, T., Wang, G., Li, P., Yang, B., & Li, J. (2021). Pathogenesis of *Proteus mirabilis* in Catheter-Associated Urinary Tract Infections. *Urologia internationalis*, 105(5-6), 354–361.

Trademarks

All trademarks that may appear in this package insert are the property of their respective owners.

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