
Vitassay qPCR

Entero + Parecho + Parotitis

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	Diana	TP	TN	FP	FN	SE	SP	PPV	NPV
1	Enterovirus	3	45	0	0	1 (0.29 – 1)	1 (0.92 – 1)	1 (0.29 – 1)	1 (0.92 – 1)
2	Virus de la Parotiditi	4	107	0	0	1 (0.39 – 1)	1 (0.96 – 1)	1 (0.39 – 1)	1 (0.96 – 1)
	Enterovirus	13	98	0	0	1 (0.75 – 1)	1 (0.96 – 1)	1 (0.75 – 1)	1 (0.96 – 1)

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- Cobas z480 Analyzer (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ¹
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTPRime Real Time Detection Thermal Cycler (DNA-Technology)
- DTIite Real-Time PCR System (DNA-Technology)
- Rotor-Gene® Q (Qiagen).
- QuantStudio™ 5 Real-Time PCR System (Applied Biosystems)

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Intended use

Vitassay qPCR Entero + Parecho + Parotitis allows the qualitative detection and identification of RNA from Mumps virus (Parotitis), Enterovirus and/or Parechovirus by Real-Time RT-PCR in clinical samples from patients with suspect of Mumps, Enterovirus and/or Parechovirus viral infection by their healthcare professional. This product is intended to aid in the Mumps virus, Enterovirus and/or Parechovirus infections diagnosis, alongside the patient's clinical data, epidemiological risk factors, and other laboratory tests outcomes.

References

Vitassay qPCR Entero + Parecho + Parotitis

4x8-well strip, low profile

7041075

Vitassay qPCR Entero + Parecho + Parotitis

4x8-well strip, high profile

7042075

Materials/reagents provided

Reference	Reagent/Material	Colour	Amount
7041S075/ 7042S075	Entero+Parecho+Parotitis strips low/high profile	-	4 x 8-well strip
7C075	Entero+Parecho+Parotitis Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4 x 8-cap strip

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated on the label.
- The resuspended positive control should be stored at -20°C. To avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

Additional equipment and material required

- Collection and transport system.
- Laboratory freezers (- 30°C to - 10°C and/or ≤ -70°C).
- RNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Vortex
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

Summary

Paramyxovirus (Mumps virus)

Mumps is a viral infection caused by a Paramyxovirus, which is a single-stranded RNA, linear non-segmented virus.

Humans are the only known hosts for the virus which causes mumps. The virus can be transmitted in saliva from 3 days before parotitis starts to 9 days after. Direct contact with respiratory secretions, saliva or fomites can spread mumps. That is why outbreaks usually occur among people in close contact, such as college dorms or sports teams. It usually affects young adults, likely related to imperfect vaccination effectiveness.

Fever, myalgias, anorexia, malaise, headache, and nonspecific upper respiratory symptoms are the symptoms that people with parotitis often manifest. But the most characteristic presentation is unilateral or bilateral swelling of the parotid glands.

The measles, mumps, and rubella (MMR) vaccine was introduced in 1967 and reported cases of mumps have decreased since then. Nowadays, it is recommended that children receive their first dose at age 12-15 months and the second at age 4-6 years. However, the number of outbreaks has increased over the last decade in New York, United Kingdom, and California, among other places.

The detection of mumps antibody alone is not enough to diagnose mumps, the use of RT-PCR can reveal more details about viral strains because it is an extremely sensitive tool.

Parechovirus and Enterovirus

Parechovirus and Enterovirus belong to the Picornaviridae family, which is one of the largest RNA virus families and contains an array of pathogens that infect both humans and animals. Their members are small, non-enveloped viruses containing a single-stranded ribonucleic acid (RNA). The genus Enterovirus (EV) is divided into 12 species (EV A-J and Rhinovirus A-C). EV A, B, C and D are found in human, and the rest of them

in animals. EV-A consist of 25 types, EV-B 63 of types, EV-C 23 of types and EV-D of 5 types. The genus Parechovirus is divided into 2 species (Parechovirus A and B). The species Parechovirus A consists of 16 types which infect humans, HPeV-1 to HPeV-16, whereas Parechovirus B consists of 4 types which infect rodents.

Transmission of Enterovirus occur through the fecal-oral, transplacental pathways, and respiratory droplets. Children, especially neonates and babies, are at high risk for developing aseptic meningitis due to Enterovirus and Parechovirus. While EV infections commonly affect older children and adults, HPeV infections in these populations are uncommon.

EVs can cause wide range of symptoms, from less serious respiratory and gastrointestinal conditions such as herpangina, and hand-foot-and-mouth disease, to more serious ones like pleurodynia, hepatitis, paralysis, among others.

For HPeV, HPeV-1 is the most prevalent type, and it mainly causes mild gastrointestinal and respiratory disease although more serious disease can occasionally be seen in young children. HPeV are the second most important cause of viral sepsis-like illness and meningitis in infants, both of which are mostly caused by HPeV-3 which is the most pathogenic HPeV type.

Nowadays, the gold standard for diagnosis of their infections is RT-qPCR, which has been demonstrated to be more accurate than viral culture in detection both EV and HPeV RNA.

Principle of the test

Vitassay qPCR Entero + Parecho + Parotitis is based on the real-time amplification of a conserved region of *F* (Fusion Protein) gene of Mumps virus, *5'UTR* gene of Enterovirus, and *5'UTR* gene of Parechovirus. After RNA isolation, the pathogens presence is detected by an increase in observed fluorescence during the reaction, after hydrolysis of the fluorescent probe.

The assay is based on 5' exonuclease activity using two primers and a fluorogenic hydrolysis probe to detect accumulation of the amplified target sequence during the PCR reaction. When the polymerase begins to spread the primers, the probe is hydrolyzed by its exonuclease 5'-3' activity causing the spatial separation of the fluorophore and the quencher. The increase in the resulting fluorescent signal is proportional to the amount of amplified product in the sample and is detected by means of real-time PCR equipment.

Vitassay qPCR Entero + Parecho + Parotitis is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible inhibition reaction. After the amplification reaction, Mumps virus is detected in FAM channel, Parechovirus is detected in ROX channel, Internal Control is detected in Cy5 channel, and the Enterovirus is detected in HEX, VIC, or JOE channels (depending on the equipment used).

Precautions

- For professional *in vitro* diagnostic use (use by qualified and trained clinical laboratory personnel).
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or is broken or if the foil has been broken or damaged.
- Do not remove desiccant from reagent pouches once is open.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles, and mask.
- Do not eat, drink, or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment, and reagents to the area in which the previous step was performed. Use separate areas for the patient samples and controls preparation.
- Avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile disposable aerosol resistant pipette tips or RNase/DNase positive displacement pipette tips is recommended.
- An appearance of the reaction mixture in stabilized format, normally found at the tube bottom, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or color different from whitish) does not alter the test functionality.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious and they must be managed according to the national safety legislation and national health waste legislation. Take necessary precautions during the collection, transport, storage, treatment, and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Use personal protective equipment (PPE) and biological safety cabinet for handling of potentially infectious samples and reagents according to current recommendations.

Procedures

Specimen collection, transport, and conservation

For specimen collection, conservation, and transport, user-validated conditions must be followed. The Vitassay qPCR Entero + Parecho + Parotitis has been tested swabs specimens (vesicles, cutaneous skin, mouth, saliva, mucosal ulcer lesions (for mumps virus and enterovirus detection, in nasopharyngeal swabs and plasma samples (for mumps virus detection), and in cerebrospinal fluid (CSF) (for enterovirus detection). Swabs specimens were collected in M4RT (Viral Transport Media) or in Cobas® PCR Media Dual Swab Sample Kit. Other sample types should be validated by the user.

Overall, clinical samples should be collected appropriately in clean containers with or without transport media (depending on sample type), correctly labelled and promptly processed to ensure the test's quality. It is recommended to use fresh specimens for the test. Transport must always be carried out in accordance with local and national regulations for the transport of biological samples. The specimens should be processed immediately or transported at -20°C or lower. The samples can be stored frozen at -20°C or lower (at -80°C ideally) for a long time. Repeated freeze-thaw cycles should be avoided to prevent sample and nucleic acids degradation.

The clinical specimens must be collected, transported, and stored according to appropriate laboratory guidelines. For more information, refer to the following guidelines:

- Miller JM, et al. (2018). A Guide to Utilization of the Microbiology Laboratory for Diagnosis of Infectious Diseases: 2018 Update by the Infectious Diseases Society of America and the American Society for Microbiology. *Clin Infect Dis*. Aug 31;67(6):e1-e94.
- García-Lechuz Moya, J.M., González López, J.J., Orta Mira, N., Sánchez Romero, M.I. (2017). Recogida, transporte y procesamiento general de las muestras de Laboratorio de Microbiología. 2017. 1b. Sánchez Romero, M.I., (coordinadora). Procedimientos en Microbiología Clínica. Cercenado Mansilla, E., Canton Moreno, R., (editores). *Sociedad Española de Enfermedades Infecciosas y Microbiología Clínica (SEIMC)*.

RNA extraction

For nucleic acid isolation from clinical samples, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. For the sample's pretreatment follow the instructions for use of the extraction kit used. The assay has been validated with the following extraction kits:

Invisorb® Spin Universal Kit (Invitex)

NucliSENS® eMAG® (Biomérieux).

MagDEA Dx SV kit, utilizando el magLEAD® 12gC instrument (Precision System Science Co.)

Positive control preparation

Reconstitute the lyophilized Entero+Parecho+Parotitis Positive Control (red tube) with 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend to open and manipulate it in a separate laboratory area away from the other components and samples.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run for each assay.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of RNA sample, negative (yellow tube) or positive (red tube) controls and add them into the corresponding wells.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Reverse transcription	45°C	15 min	1
Polymerase activation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the FAM (Mumps virus), ROX (Parechovirus), Cy5 (Internal Control), and HEX, JOE, or VIC (Enterovirus). Depending on the equipment used select the proper detection channel (Attached II).

Analysis and interpretation of results

The results analysis is done by the software itself of the used real-time PCR system following manufacturer's instructions. To verify the correct operation of the amplification mix, check the internal control (IC) signal emission.

Use the positive control amplification curve as a starting point during reaction validation (prior to interpretation of sample results), to ensure that the threshold falls within the exponential phase of the amplification curves and above any background noise signals. The threshold value may vary between different instruments due to different signal intensities. It is recommended to set the threshold values for each channel (target) independently by the end user.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive control used in each run must show an amplification curve ($Ct \leq 40$) in FAM, ROX, and HEX channels, which validates the reaction.

Negative control

The negative control included in each run must show signal' absence ($Ct > 40$ or no signal) in FAM, ROX, and HEX channels.

The experiment seems to be failed if there is amplification signal in negative control or signal absence in the positive control for any channel. The assay should be repeated.

The Internal Control (IC) should show an amplification signal ($Ct \leq 40$) in positive control and negative control wells.

The clinical samples test results assessment should be performed once controls' results have been validated. The result interpretation is summarized in the following table:

Entero + Parecho + Parotitis				Interpretation
Mumps virus (FAM)	Enterovirus (HEX)	Parechovirus(ROX)	Internal Control (Cy5)	
+	-	-	+/- ¹	Mumps virus RNA Detected
-	+	-	+/- ¹	Enterovirus RNA Detected
-	-	+	+/- ¹	Parechovirus RNA Detected
+	+	-	+/- ¹	Mumps virus and Enterovirus RNA Detected
+	-	+	+/- ¹	Mumps virus and Parechovirus RNA Detected
-	+	+	+/- ¹	Enterovirus and Parechovirus RNA Detected
+	+	+	+/- ¹	All targets' RNA Detected
-	-	-	+ ²	Targets not detected ²
-	-	-	- ²	Failed test

¹ Sometimes, the Internal Control (IC) detection is not necessary since a high copies number of the target can cause a preferential amplification of target-specific nucleic acids. The Internal Control (IC) shows or not an amplification signal (Ct ≤40 or no signal).

² In the case of negative pathogens target genes detection, IC must show an amplification signal with Ct ≤ 35. If there is a signal' absence or Ct value >35 of Internal Control, the result is considered 'Invalid', and retesting is required. It is recommended to repeat the qPCR by diluting the RNA sample 1:10 and/or 1:100, or to re-extract and retest to check for possible failure in the extraction procedure and/or inhibition problems.

If the obtained result is confusing or doubtful, it is necessary to check that all the steps have been carried out correctly, to verify the correct performance of each qPCR steps and to review all the parameters, the sigmoid shape of the curve and the fluorescence intensity. It is also recommended to repeat the assay, preferably in duplicate, depending on the available material (repeat qPCR with the same isolated RNA sample, or re-extract and retest another aliquot of the same specimen or, obtain a new specimen and retest).

The test results must be evaluated by a health professional, together with the medical history, clinical symptoms and/or the results obtained in other diagnostic tests.

Quality Control

To confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Moreover, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

To assess clinical sensitivity and specificity, two studies were conducted. In a first study (1), 48 cerebrospinal fluid (CSF) samples were analysed and in a second study (2) 111 swab samples (vesicles, cutaneous-skin, mouth, mucosal ulcer/lesions, mouth-saliva) were analysed.

In both studies the extraction was performed with the MagDEA Dx SV kit and magLEAD® 12 gC instrument (Precision System Science Co.) and analysed with the Vitassay qPCR Entero + Parecho + Parotitis kit with DTlite Real-Time PCR System (DNA-Technology) in the first study and with CFX96™ Real-Time PCR Detection System (Bio-Rad) in the second study.

The Vitassay qPCR Entero + Parecho + Parotitis kit's results were compared with the results obtained with the routine molecular methods used in the third level hospital where samples collected in the first study and Eurofins Biomnis and Cerba Xpert in the second study; the obtained results are shown in the following tables.

	Target	TP	TN	FP	FN	SE	SP	PPV	NPV
1	Enterovirus	3	45	0	0	1 (0.29 – 1)	1 (0.92 – 1)	1 (0.29 – 1)	1 (0.92 – 1)
2	Mumps Virus	4	107	0	0	1 (0.39 – 1)	1 (0.96 – 1)	1 (0.39 – 1)	1 (0.96 – 1)
	Enterovirus	13	98	0	0	1 (0.75 – 1)	1 (0.96 – 1)	1 (0.75 – 1)	1 (0.96 – 1)

TP = True Positive, TN = True Negative, FP = False Positive, FN = False Negative, SE = Sensibility, SP = Specificity, PPV = Positive Predictive Value, NPV = Negative Predictive Value.

In conclusion, the results show a high concordance to detect those pathogens using Vitassay qPCR Entero + Parecho + Parotitis.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of standards from the pathogens ranging from 10⁷ to 10¹ copies/reaction. This assay has a detection limit of 10 cDNA copies per reaction for Mumps virus, Enterovirus and Parechovirus (positive rate ≥ 95%).

Analytical specificity

The analytical specificity for the different pathogen detection associated to diseases in immunosuppressed patients, was tested within the panel of following microorganisms:

Cross-reactivity testing					
<i>Acinetobacter baumannii</i>	-	Epstein-Barr virus	-	<i>Plasmodium falciparum</i> 3D7	-
Adenovirus types 1, 2, 3, 4, 5, 8, 15,31,40 and 41	-	<i>Escherichia coli</i> 0.1285; O18:H7:K1	-	<i>Pneumocystis jirovecii</i> Type A1	-
<i>Aspergillus fumigatus</i>	-	<i>Haemophilus influenzae</i> MinnA	-	<i>Pneumocystis jirovecii</i> g885652	-
<i>Bacteroides fragilis</i>	-	Hepatitis A virus	-	<i>Pseudomonas aeruginosa</i>	-
BK Virus types Ib-2 and IV	-	Herpes simplex virus 1 (HSV-1) strain MacIntyre	-	St Louis Encephalitis Virus	-
<i>Candida albicans</i>	-	Herpes simplex virus 2 (HSV-2) MS	-	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-
Chikungunya S27 Petersfield	-	Human herpesvirus 6 (HHV6) Strain Z29	-	<i>Streptococcus agalactiae</i> Z019	-
Cytomegalovirus AD-169	-	Human herpesvirus 6 (HHV6) Type A	-	<i>Streptococcus pneumoniae</i> Z022	-
<i>Cryptococcus gatti</i> Z156	-	Human herpesvirus 6 (HHV6) Type B	-	Tick-Borne encephalitis virus strain Neudorf1	-
Dengue serotype 1 Hawaii A	-	JC Virus Type 1A	-	<i>Toxoplasma gondii</i> type II	-
Dengue serotype 2 New guinea C	-	JC Virus Type 2B	-	<i>Treponema pallidum</i>	-
Dengue serotype 3 H87	-	<i>Klebsiella oxytoca</i>	-	Varicella-Zoster Virus OKA	-
Dengue serotype 4 H241	-	<i>Listeria innocua</i> serotype 6a	-	Varicella-Zoster Virus 9/84	-
Japanese Encephalitis virus	-	<i>Listeria ivanovii</i> subsp. <i>Ivanovii</i> serovar 5	-	Varicella-Zoster Virus Ellen	-
<i>Enterobacter aerogenes</i> serotype Cloaca B	-	<i>Listeria monocytogenes</i> serovar 4b	-	West Nile virus Heja	-
<i>Enterobacter cloacae</i> serotype Cloaca A	-	<i>Listeria monocytogenes</i> serotype 1/2b	-	West Nile virus NY99	-
<i>Enterococcus durans</i>	-	Measles virus	-	West Nile virus Ug37	-
<i>Enterococcus faecalis</i>	-	<i>Neisseria meningitidis</i> serogroup A	-	Yellow Fever Virus strain 17D	-
<i>Enterococcus faecium</i> serotype 11	-	Parvovirus B19	-	<i>Yersinia enterocolitica</i> O:3	-

Analytical reactivity

The Vitassay qPCR Entero + Parecho + Parotitis kit's reactivity for Mumps virus was evaluated against extracted RNA from Mumps virus strain Enders (as template), showing positive results.

The Vitassay qPCR Entero + Parecho + Parotitis kit's reactivity for Enterovirus was evaluated against extracted RNA from the following strains (as template), showing positive results:

Enterovirus Echovirus Types 6, 7, 11, 18, 25, 30, Human Echovirus 4 strain Pesascek

Enterovirus types 68, D68 B3, 71 and A71

Enterovirus Coxsackievirus types A4, A9, A16, A21, A24, B3 and B4

Human Coxsackievirus B1 strain Conn-5

Enterovirus D68 strain US/MO/14-18949

The Vitassay qPCR Entero + Parecho + Parotitis kit's reactivity for Parechovirus was evaluated against extracted RNA from the following strains (as template), showing positive results:

Parechovirus types 1, 2, 3, 4, 6

Human Parechovirus 1 strain Harris

Compatible real-time PCR equipment

Vitassay qPCR Entero + Parecho + Parotitis has been validated on the following equipment:

- Cobas z480 Analyzer (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ¹
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTPPrime Real Time Detection Thermal Cycler (DNA-Technology)
- DTIite Real-Time PCR System (DNA-Technology)
- Rotor-Gene® Q (Qiagen)
- QuantStudio™ 5 Real-Time PCR System (Applied Biosystems)

¹: When using the Applied Biosystems 7500 Fast with strips it is recommended placing a plate holder for the tubes (Ref. PN 4388506).

Limitations

- All obtained results must be interpreted by a specialist in conjunction with available clinical information and laboratory findings.
- This assay has been validated with RNA extracted from swabs specimens (vesicles, cutaneous skin, mouth, saliva, mucosal ulcer/lesions (for mumps and enterovirus)), from nasopharyngeal swabs and plasma samples (for mumps), and from CSF (for enterovirus). The use of other samples has not been established and must be validated by the user.
- The correct test performance depends on the sample's quality; proper RNA from clinical specimens must be extracted.
- This test does not provide quantitative values nor indicate the number of organisms present. This is a qualitative test.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination either by Entero+Parecho+Parotitis Positive Control during its reconstitution, by samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.
- Detection may be affected by several factors and their combinations which may lead to false negative results, including a) inadequate specimen sampling, shipping, storage, handling; b) procedural errors (including RNA isolation); c) RNA degradation during specimen shipping, storage, and/or preparation; d) potential mutations of the target regions of the Mumps virus, Enterovirus and/or Parechovirus, covered by this test which may result in RNA being undetectable e) pathogen load below the limit of detection for the assay; f) the presence of Real-Time amplification inhibitors or other types of interference (an interference study evaluating the effect of vaccines, antiviral therapeutics, antibiotics, chemotherapeutics or immunosuppressant drugs used to prevent the infection or used during the treatment of the infection was not performed); g) failure to follow the manufacturer's instructions and procedures.
- Detection of RNA may not indicate the presence of viable and/or infectious bacteria or that these are the causative agents for clinical symptoms.
- Negative results do not preclude Mumps virus, Enterovirus and/or Parechovirus infection and should not be the sole basis of a patient treatment/management decision. Optimal specimen types and/or stage of infection most suitable for its collection have not been established yet. Consider the collection of multiple specimens from the same patient at different time points, which may increase the probability of detecting the bacteria.
- If the patient clinical data, laboratory tests and epidemiological studies suggest possible Mumps virus, Enterovirus and/or Parechovirus infection, and other

bacterial infections and/or multi-resistant infection or colonization have been discarded, a false negative result might not be discarded, and additional tests should be performed.

- Fluorescence values may vary due to multiple factors such as PCR equipment, extraction system, sample type and its previous treatment, among others.

Attached I. Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
Agilent Technologies	Abbott
AriaMx Real-Time PCR System	Abbott m2000 ⁽¹⁾
AriaDx Real-Time PCR System	Applied Biosystems
Applied Biosystems	7300 Real-Time PCR System ^{(3) (1)}
7500 Fast Real-Time PCR System ^{(1) (2)}	7500 Real-Time PCR System ⁽¹⁾
7500 Fast Dx Real-Time PCR System ^{(1) (2)}	7900 Real-Time PCR System ⁽⁴⁾
QuantStudio™ 12K Flex 96-well Fast	ABI PRISM 7000 ⁽⁴⁾
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7700 ⁽⁴⁾
QuantStudio™ 7 Flex 96-well Fast	QuantStudio™ 12K Flex 96-well
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
QuantStudio™ 5 Fast Real-Time PCR System	QuantStudio™ 7 Flex 96-well
QuantStudio™ 3 Fast Real-Time PCR System ⁽³⁾	QuantStudio™ 5 Fast Real-Time PCR System
StepOne Plus™ Real-Time PCR System ⁽³⁾	QuantStudio™ 5 Real-Time PCR System
StepOne™ Real-Time PCR System ^{(3) (4)}	QuantStudio™ 3 Real-Time PCR System ⁽³⁾
ViiA™ 7 Fast Real-Time PCR System	ViiA™ 7 Real-Time PCR System
Azure Biosystems	Agilent Technologies
Azure Cielo 3 ⁽⁵⁾	Mx3000P™ Real Time PCR System
Azure Cielo 6	Mx3005P™ Real Time PCR System
BIONEER	Analytik Jena
Exicycler™ 96	qTOWER ⁽⁶⁾
Bio-Rad	BIONEER
CFX96™ Real-Time PCR Detection System	Exicycler™ 96
Mini Opticon™ Real-Time PCR Detection System ⁽⁵⁾	BIOER
Roche	QuantGene 9600
LightCycler @480 Real-Time PCR System ^{(1) (6)}	Bio-Rad
LightCycler @96 Real-Time PCR System ⁽¹⁾	CFX96™ Deep Well Real-Time PCR
Cobas z480 Analyzer ^{(1) (6)}	iCycler iQ™ Real-Time PCR
	iCycler iQ™ 5 Real-Time PCR
	My iQ™ Real-Time PCR Detection System ⁽⁵⁾
	My iQ™ 2 Real-Time PCR Detection System ⁽⁵⁾
Special Formats ⁽⁷⁾	DNA-Technology
Bio Molecular Systems	DTlite Real-Time PCR System ⁽⁸⁾
Mic Real Time PCR Cyclers	DTprime Real-time Detection Thermal Cyclers ⁽⁸⁾
Cepheid	Eppendorf
SmartCycler®	Mastercycler™ ep <i>realplex</i>
Qiagen	Qiagen
Rotor-Gene® Q	QIAquant 96

- (1) A special grid is needed to fit these real-time PCR kits.
- (2) Select Ramp Speed "Standard".
- (3) No ROX caption.
- (4) No Cy5 caption.
- (5) Only FAM and HEX caption.
- (6) Specific compensation color is required.
- (7) The product must be reconstituted following the appropriate procedure (see Test procedure) and transferred to the specific tubes for Mic, SmartCycler®, Rotor-Gene® Q or geneLEAD VIII System.
- (8) See Attached III to configure exposure settings.

Attached II. Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR thermocyclers are specified in the following table:

Thermocycler	Vitassay Channel	Detection Channel	Observations
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Applied Biosystems ABI 7500	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Roche Cobas z 480	FAM	465/510	Colour Compensation required
	HEX	540/580	
	ROX	540/610	
	Cy5	610/670	
Cepheid Smartcycler®	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Agilent Technologies Mx3000P™ Mx 3005P™	FAM	FAM	Passive reference option ROX is none
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Bio Molecular Systems Mic Real Time PCR Cycler	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10).
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Agilent Technologies AriaMx	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Qiagen Rotor-Gene® Q	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range must be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
BIONEER Exicycler™ 96	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Attached III. Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
DTIite Real-Time PCR System (DNA-Technology)	FAM	500
	HEX	500
	ROX	500
	Cy5	500
DTprime Real-time Detection Thermal Cycler (DNA-Technology)	FAM	500*
	HEX	1000
	ROX	1000
	Cy5	1000

*If result in FAM channel is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.



Bibliography/Bibliografía

- de Crom, S. C., Rossen, J. W., van Furth, A. M., & Obihara, C. C. (2016). Enterovirus and parechovirus infection in children: a brief overview. *European journal of pediatrics*, 175(8), 1023–1029.
- Bockelman, C., Frawley, T. C., Long, B., & Koyfman, A. (2018). Mumps: An Emergency Medicine-Focused Update. *The Journal of emergency medicine*, 54(2), 207–214.

Trademarks

All trademarks that may appear in this package insert are the property of their respective owners.

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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